

Syddansk Universitet

High throughput formulation screening: Testing a novel two compartment microtiter plate approach

Jacobsen, Ann-Christin; Brandl, Martin; Bauer-Brandl, Annette

Publication date:
2018

Document version
Publisher's PDF, also known as Version of record

Document license
Unspecified

Citation for pulished version (APA):
Jacobsen, A-C., Brandl, M., & Bauer-Brandl, A. (2018). High throughput formulation screening: Testing a novel two compartment microtiter plate approach. Poster session presented at American Association of Pharmaceutical Scientists (AAPS) annual meeting - PharmSci 360 , Washington DC, United States.

General rights

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal ?

Take down policy

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

High throughput formulation screening: Testing a novel two compartment microtiter plate approach

Ann-Christin Jacobsen¹, Martin Brandl¹, Annette Bauer-Brandl¹

¹Drug Transport & Delivery Group, Department of Physics, Chemistry & Pharmacy,
University of Southern Denmark, Odense Denmark

CONTACT INFORMATION: annette.bauer@sdu.dk

PURPOSE

A novel two compartment microtiter plate-based **high throughput formulation screening (HTFS)** approach was tested for its suitability for formulation performance ranking of enabling formulations. Many new drug compounds have poor solubility and require complex formulation to enable oral administration. Identification of the optimum formulation design/composition can be a challenging and costly process, which may be bypassed by the present approach.

OBJECTIVES

The HTFS approach was tested using model formulations of **tadalafil (TDF)** with/without the amphiphilic polymer **Soluplus®** (crystalline and amorphous formulations). Formulations were evaluated with respect to both the amount of drug **dissolved** and the amount of drug **permeated**.

- Can the HTFS approach predict the formulation performance *in vivo*?
→ The *in vitro* HTFS results were compared to *in vivo* data from literature.
- How do screening parameters affect HTFS practicability and predictability?
→ The following test parameters were evaluated:
1) Incubation time, 2) acceptor media and 3) dispersion (i.e. donor) media

CONCLUSIONS

The two compartment microtiter plate approach was generally able to predict the formulation performance *in vivo*.

Varying the screening parameters of the HTFS approach indicated that the protocol can be simplified to achieve high practicability and simultaneously high predictability.

FUTURE OUTLOOK

Testing the novel approach with other formulations of poorly soluble drug compounds will reveal it's usefulness in early drug development.

RESULTS

In vitro in vivo correlation: Comparing the HTFS results to *in vivo* data from Krupa et al. 2016

Table 1 Coefficients of determination (R^2) when plotting the TDF concentration in either the donor compartment (dissolved TDF) or the acceptor compartment (permeated TDF) versus the *in vivo* AUC (Krupa et al. 2016). Incubation time: 6h, Acceptor medium: 1% VitE TPGS

Dispersion medium	R^2 (dissolved concentration vs AUC)	R^2 (permeated concentration vs AUC)
Phosphate buffer	0.8459	0.8499
SGF	0.826	0.8637
FaSSIF	0.8282	0.9391

A good *in vitro in vivo* correlation was achieved.

TDF concentration in the acceptor compartment (**permeated TDF**) generally yielded a **better prediction** at same time providing experimental ease (**no separation step**).




Figure 3 TDF concentration in A) the donor compartment (dissolved) or B) the acceptor compartment (permeated) vs *in vivo* AUC.

Influence of screening parameters

Acceptor media

Using surfactant solutions as acceptor medium increased TDF permeation as compared to neat buffer: *Sink conditions*

No difference between surfactants.




Figure 4 TDF concentration in the acceptor after 6h of incubation using different acceptor media. The dotted line indicates the TDF concentration when using neat buffer.

Dispersion media

FaSSIF as dispersion medium yielded the highest TDF concentration in the donor compartment followed by SGF and phosphate buffer.

The same clear difference in TDF concentration could not be seen in the acceptor compartment.




Figure 5 TDF concentration in the acceptor compartment after 1, 3 and 6h of incubation. Acceptor medium: 1% VitE TPGS.

Incubation time

Steady TDF transport over the dialysis membrane for at least 6h. Longer incubation increased analytical ease.

Figure 6 TDF concentration in A) the donor compartment or B) the acceptor compartment after 6 h of incubation using different media for dispersion of the formulations. Acceptor medium: 1% VitE TPGS




Figure 2 Flow diagram of the HTFS protocol.

METHODS

Figure 2 shows a flow diagram of the HTFS protocol used in all experiments.

- Four **TDF formulations** (i.e. suspensions) were prepared:
1) Crystalline TDF, 2) Crystalline TDF and Soluplus® (physical mixture), 3) freeze dried TDF and 4) 'co'-freeze dried TDF and Soluplus® (solid dispersion).

- Three **dispersion media** were used for the suspensions:
1) 50 mM phosphate buffer, 2) Simulated gastric fluid (SGF) and 3) Fasted state simulated intestinal fluid (FaSSIF).

- Five **acceptor media** containing surfactants were evaluated and compared to neat medium:
1) 1% Sodium dodecyl sulfate (SDS), 2) 0.2% d-α-Tocopheryl polyethylene glycol 1000 (VitE TPGS), 3) 1% VitE TPGS, 4) 0.2% Polysorbate 80 (PS 80) and 5) 1% PS 80.

300 μ L suspension and 200 μ L acceptor medium were transferred to the bottom-wells and top-wells, respectively, and incubated for 1, 3 or 6h. Figure 1 shows the assembled system. TDF concentration was determined via UHPLC.

REFERENCES / ACKNOWLEDGEMENTS

Krupa, A., et al., 2016. High-Energy Ball Milling as Green Process To Vitrify Tadalafil and Improve Bioavailability. Molecular pharmaceutics 13, 3891-3902
InnoME GmbH (D-Espelkamp) is acknowledged for co-financing this project.