

5-1995

The selective wavelength filtering characteristics of Tokai optical medical lenses and a discussion of current applications in the clinical setting.

David Dowe
Pacific University

Weon Jun
Pacific University

James Phelan
Pacific University

Recommended Citation

Dowe, David; Jun, Weon; and Phelan, James, "The selective wavelength filtering characteristics of Tokai optical medical lenses and a discussion of current applications in the clinical setting." (1995). *College of Optometry*. 1122.

<https://commons.pacificu.edu/opt/1122>

This Thesis is brought to you for free and open access by the Theses, Dissertations and Capstone Projects at CommonKnowledge. It has been accepted for inclusion in College of Optometry by an authorized administrator of CommonKnowledge. For more information, please contact CommonKnowledge@pacificu.edu.

The selective wavelength filtering characteristics of Tokai optical medical lenses and a discussion of current applications in the clinical setting.

Abstract

Tokai Optical Company in Okazaki, Japan, has developed a new line of CR-39 medical lenses designed to manage certain ocular and visual disorders. Before these lenses can be incorporated into the medical market, it is necessary to verify the selective filtration characteristics claimed by Tokai Optical Company. Using a modified Varian OMS 200 UV/Visible Spectrophotometer, Tokai Optical 'CC Filter' lenses were analyzed and compared to the transmission properties published by the manufacturer. It was concluded that the spectral transmission results differ by no more than 3.5 percent transmittance with the majority being less than 1.0 percent difference. The historical use of spectrally selective lenses in the treatment, management and prevention of ocular and visual disorders is also discussed.

Degree Type

Thesis

Degree Name

Master of Science in Vision Science

Committee Chair

Willard B. Bleything

Subject Categories

Optometry

Copyright and terms of use

If you have downloaded this document directly from the web or from CommonKnowledge, see the "Rights" section on the previous page for the terms of use.

If you have received this document through an interlibrary loan/document delivery service, the following terms of use apply:

Copyright in this work is held by the author(s). You may download or print any portion of this document for personal use only, or for any use that is allowed by fair use (Title 17, §107 U.S.C.). Except for personal or fair use, you or your borrowing library may not reproduce, remix, republish, post, transmit, or distribute this document, or any portion thereof, without the permission of the copyright owner. [Note: If this document is licensed under a Creative Commons license (see "Rights" on the previous page) which allows broader usage rights, your use is governed by the terms of that license.]

Inquiries regarding further use of these materials should be addressed to: CommonKnowledge Rights, Pacific University Library, 2043 College Way, Forest Grove, OR 97116, (503) 352-7209. Email inquiries may be directed to: copyright@pacificu.edu

**THE SELECTIVE WAVELENGTH FILTERING
CHARACTERISTICS OF TOKAI OPTICAL
MEDICAL LENSES AND A DISCUSSION
OF CURRENT APPLICATIONS IN
THE CLINICAL SETTING.**

By:

David Dowe

Weon Jun

James Phelan

**A thesis submitted to the faculty of the
College of Optometry
Pacific University
Forest Grove, OR
for the degree of
Doctor of Optometry
May 1995**


Adviser:

Willard B. Bleything, O.D., M.S.


Authors



David Dowe



Weon Jun



James Phelan

Advisor:



Willard B. Bleything O.D., M.S.

About the Authors

David Dowe was raised in Nanaimo, British Columbia, a small city on the east coast of beautiful Vancouver Island. After receiving his Bachelors of Science degree from the University of Victoria, he traveled to Guyana, South America where he participated on several medical and scientific missions into the remote South American interior. It was here that he decided, upon his return, to pursue a career in optometry. David recently completed his optometric preceptorship at Tyndall Air Force Base, Florida and plans to practice optometry in British Columbia after graduation.

Weon Jun was born in Jinju, Korea and moved to Fairbanks, Alaska in 1981. He completed his Visual Science degree at Pacific University and entered the Optometry School in 1991. Weon completed his preceptorship at the Eye Foundation of Utah. One of the most rewarding experiences in optometry school was being the president of AMIGOS and participating in eye care missions locally and internationally. He has been selected for the 1994 Dean's Award and was nominated for Who's Who Among Students in American Universities and Colleges. His future plans are to complete a residency program and practice in Oregon or Washington. He also plans to continue to support AMIGOS.

James Phelan is currently doing a preceptorship at Malmstrom Air Force Base in Great Falls, Montana. Raised in Helena, Montana, he studied Pre-Med., with a minor in chemistry at Carroll College. James was awarded the Thomas Sullivan Memorial Scholarship for Academic Achievement in Biology. He then shifted his focus from medical school to optometry school, choosing Pacific University because of its reputation and its ideal location. During his second year as an optometry student, James was awarded the Montana Optometric Association Scholarship for academic leadership and community involvement. James is planning to settle in Great Falls, Montana when he graduates in May.

Abstract

Tokai Optical Company in Okazaki, Japan, has developed a new line of CR-39 medical lenses designed to manage certain ocular and visual disorders. Before these lenses can be incorporated into the medical market, it is necessary to verify the selective filtration characteristics claimed by Tokai Optical Company. Using a modified Varian DMS 200 UV/Visible Spectrophotometer, Tokai Optical 'CC Filter' lenses were analyzed and compared to the transmission properties published by the manufacturer. It was concluded that the spectral transmission results differ by no more than 3.5 percent transmittance with the majority being less than 1.0 percent difference. The historical use of spectrally selective lenses in the treatment, management and prevention of ocular and visual disorders is also discussed.

Introduction

It is well known that illumination plays a crucial role in optimizing useful vision in patients with reduced visual function resulting from disease or the aging process. For this reason, practitioners prescribe neutral density or selective filtration (S.F.) lenses in an attempt to increase patient comfort and visual performance. It has also been proposed that by selectively filtering specific wavelengths from entering the eye, destructive pathological processes may be slowed.

Tokai Optical Company has developed a series of lenses specifically with this in mind. The lenses, presently marketed in Japan, are sold under the brand name of "CC Filter" lenses and include two types: CCP, and CCP-400, the latter allowing slightly more transmission of the visible spectrum. Both are manufactured from a CR-39 resin. Tokai medical lenses are available with six different transmission spectra and appear red, orange, yellow, light green, dark green, and brown in color. All are a fixed tint, and made of non photochromic material. This paper attempts to verify the claimed transmission characteristics of Tokai medical lenses through careful spectrophotometric analysis. Lenses with different transmission characteristics can be used to alleviate

ocular symptoms. Research on the clinical application of such lenses is reviewed and discussed.

Procedure

A total of 12 different lenses were verified using a modified Varian DMS 200 UV/Visible/IR Spectrophotometer. The substrate carrier was modified to allow measurement through plastic lenses as the spectrophotometer's primary use is to measure transmission characteristics of liquids. Modification entailed the simple removal of the carrier such that lenses could be positioned centrally in the path of the analyzer. Using the transmission through air as a 100% transmittance base line, each lens was scanned from 280 nm to 800 nm and a printout of the transmission curves was obtained. As well, transmissions measured *at specific wavelengths*, were compared to values reported by Tokai Optical Company. These wavelengths included 280 nm (harmful UV-C radiation), 315 nm (harmful UV-A radiation), 400nm (end of the UV spectrum), 450 nm (scotopic sensitivity maximum), 550 nm (photopic sensitivity maximum) and 850 nm (infrared). These two measuring approaches allowed for direct comparisons between the present study's values and those reported by Tokai Optical Company.

Results



Spectrophotometric analysis of the study lenses revealed only minimal deviations from Tokai Optical Company's reported transmission characteristics. The lens that differed the greatest was the Amber CCP400. *At 450 nm*, this lens showed 11.4% *greater* absorption than was reported by Tokai Optical Company (Table 1). Spectrophotometer plots from Tokai Optical Company as well as those from the present study can be found in Appendix 1.

Overall, the average deviation for each CCP lens varied by no more than 1.29% between the present study and Tokai Optical Company's values. The present study's results revealed both over and under-estimation of transmittances reported by Tokai (Figure 1).

Table 1.																		
Percent light transmission of CCP lenses at specific wavelengths																		
Colors:	Yellow			Orange			Red			Light Green			Dark Green			Brown		
	Total	PUCO	% Difference	Total	PUCO	% Difference	Total	PUCO	% Difference	Total	PUCO	% Difference	Total	PUCO	% Difference	Total	PUCO	% Difference
280 nm	0.16%	0.00%	0.16%	0.18%	0.00%	0.18%	0.18%	0.00%	0.18%	0.24%	0.00%	0.24%	0.19%	0.00%	0.19%	0.15%	0.00%	0.15%
315 nm	0.14%	0.00%	0.14%	0.14%	0.00%	0.14%	0.13%	0.00%	0.13%	0.16%	0.00%	0.16%	0.16%	0.00%	0.16%	0.14%	0.00%	0.14%
400 nm	0.05%	0.00%	0.05%	0.06%	0.00%	0.06%	0.06%	0.00%	0.06%	0.06%	0.00%	0.06%	0.06%	0.00%	0.06%	0.06%	0.00%	0.06%
450 nm	0.05%	0.00%	0.05%	0.05%	0.00%	0.05%	0.06%	0.00%	0.06%	0.06%	0.00%	0.06%	0.06%	0.00%	0.06%	0.09%	0.00%	0.09%
511 nm	2.97%	2.30%	0.67%	1.07%	4.60%	-3.53%	0.45%	0.02%	0.43%	8.98%	9.50%	-0.52%	4.84%	6.30%	-1.46%	5.39%	5.30%	0.09%
527 nm	14.20%	12.10%	2.10%	5.81%	11.10%	-5.29%	1.23%	0.08%	1.15%	21.10%	22.10%	-1.00%	13.40%	15.80%	-2.40%	11.40%	11.10%	0.30%
550 nm	49.20%	45.10%	4.10%	27.50%	29.10%	-1.60%	6.44%	4.90%	1.54%	37.40%	39.70%	-2.30%	25.40%	29.10%	-3.70%	24.70%	25.60%	-0.90%
850 nm	92.50%	90.70%	1.80%	92.10%	92.40%	-0.30%	92.80%	91.10%	1.70%	92.40%	90.00%	2.40%	92.20%	91.70%	0.50%	92.60%	92.30%	0.30%
		Avg. %	1.13%		Avg. %	-1.29%		Avg. %	0.66%		Avg. %	-0.11%		Avg. %	-0.82%		Avg. %	0.03%
Percent light transmission of CCP400 lenses at specific wavelengths																		
Colors:	Amber			Light Green			Yellow			Green			Olive			Evergreen		
	Total	PUCO	% Difference	Total	PUCO	% Difference	Total	PUCO	% Difference	Total	PUCO	% Difference	Total	PUCO	% Difference	Total	PUCO	% Difference
280 nm	0.12%	0.00%	0.12%	0.35%	0.00%	0.35%	0.21%	0.00%	0.21%	0.11%	0.00%	0.11%	0.34%	0.00%	0.34%	0.15%	0.00%	0.15%
315 nm	0.07%	0.00%	0.07%	0.12%	0.00%	0.12%	0.14%	0.00%	0.14%	0.25%	0.00%	0.25%	0.21%	0.00%	0.21%	0.10%	0.00%	0.10%
400 nm	9.31%	6.80%	2.51%	13.40%	10.80%	2.60%	3.02%	1.20%	1.82%	7.31%	3.40%	3.91%	0.08%	0.00%	0.08%	0.73%	0.03%	0.70%
450 nm	53.60%	42.20%	11.40%	59.20%	58.10%	1.10%	26.40%	19.70%	6.70%	40.90%	35.10%	5.80%	2.43%	2.00%	0.43%	9.11%	7.50%	1.61%
511 nm	79.70%	75.10%	4.60%	79.50%	79.50%	0.00%	43.10%	38.10%	5.00%	55.60%	55.90%	-0.30%	10.40%	10.10%	0.30%	17.30%	15.70%	1.60%
527 nm	83.70%	80.04%	3.66%	80.30%	81.40%	-1.10%	47.70%	43.20%	4.50%	55.50%	56.60%	-1.10%	13.70%	13.40%	0.30%	17.90%	16.60%	1.30%
550 nm	88.40%	86.30%	2.10%	78.60%	81.10%	-2.50%	55.80%	51.40%	4.40%	53.60%	53.30%	0.30%	19.90%	19.60%	0.30%	18.00%	17.30%	0.70%
850 nm	92.70%	91.10%	1.60%	92.40%	91.40%	1.00%	92.50%	90.90%	1.60%	92.40%	90.30%	2.10%	92.00%	91.60%	0.40%	92.60%	90.80%	1.80%
		Avg. %	3.26%		Avg. %	0.20%		Avg. %	3.05%		Avg. %	1.38%		Avg. %	0.30%		Avg. %	1.00%


PUCO = Pacific University College of Optometry

Figure 1. Transmission spectra of Tokai Optical CCP lenses.




PUCO = Pacific University College of Optometry

Red CCP




Light Green CCP




PUCO = Pacific University College of Optometry

Dark Green CCP



Brown CCP



PUCO = Pacific University College of Optometry

The CCP400 lenses average-percent-difference in transmittance did not exceed 3.26% (Table 1). Curiously, Tokai Optical Company consistently *overestimated* transmittance relative to the values found in the present study. Therefore, Tokai Optical's CCP400 Lenses consistently absorbed *more* radiation than was claimed by the Company (Figure 2).

Discussion


The present study supports Tokai Optical Companies claimed absorption characteristics of their medical lenses. Therefore, these lenses can be used for the purposes of their design, without hesitation, as they have been shown to perform to the levels claimed by the Company. The prescribing practitioner must, however, be aware of *how* patients will benefit from S.F. lenses as well as *who* will be best suited to them.

There are three proposed benefits from the use of tints and filters: (a) enhancement of visual performance in varying illuminations, (b) improved visual comfort; and (c) controlling the progression of ocular disease.¹ The latter is highly controversial and has proven difficult to assess. Patients with cataracts, retinitis pigmentosa, diabetic retinopathy, glaucoma, optic nerve atrophy, and maculopathies are among those who report improvements in visual functioning with S.F. lenses.²

Cataracts


Trokel³, in his investigations into the physical basis for transparency of the crystalline lens, concluded that the primary basis for visual distortion produced by a cataract was an increase in light scattering. This finding is supported in the clinic by observing patients with early cataracts or corneal scars. These patients can usually read one or two more lines on the Snellen chart if the exam room lights are turned down.⁴ This is because the room lights act as a source of glare thus decreasing the patients contrast sensitivity and therefore visual acuity. Neutral density filters will decrease the overall illumination entering the eye, reduce the scattering and hence increase contrast sensitivity and acuity; similar acuity

Figure 2. Transmission spectra of Tokai Optical CCP400 lenses.




PUCO = Pacific University College of Optometry

Yellow CCP400




Green CCP400




PUCO = Pacific University College of Optometry

Olive CCP400



Evergreen CCP400



PUCO = Pacific University College of Optometry

increases can be seen by reducing levels of ambient light. Although most of the scattering produced by a cataract or scarred cornea is not wavelength dependent, some Raleigh scattering is present.⁴ As shown in Figure 3., blue light Raleigh scatters approximately nine times more than the light we see at our maximum scotopic sensitivity level (555 nm, light green).

Figure 3. Raleigh scattering as a function of wavelength.

555nm, photopic sensitivity maximum	320 nm, blue light
$i=1/555^4$	$i=1/320^4$
$i=1.05 \times 10^{-11}$	$i=9.54 \times 10^{-11}$

Therefore, it would be logical that lenses which selectively absorb short wavelength light might further benefit the cataract patient. However, Bither⁵ demonstrated that patients with media opacities preferred neutral density filters over such S.F. lenses. This finding is supported by Chou and Cullen⁶ who suggested that shorter wavelengths only produce *slightly* more discomfort in situations of very intense illumination. Therefore, while Raleigh scattering does contribute to glare, more relief might result from reducing total illumination rather than selectively filtering high energy light.

Reducing scattered light, or overall illumination, is not the only means by which visual performance can be enhanced in patients with media opacities. Tupper⁴ states that S.F. lenses blocking the blue end of the spectrum will sharpen most objects in the real world. He states, for example, that yellow-green leaves, red-yellow plants as well as the gray-brown branches stand out more vividly against the blue sky. Therefore, regardless of light scattering, by using an optical filter that absorbs short light wavelengths, the world of colored objects is slightly distorted but its contrast can be much enhanced. Practitioners must also keep in mind the beneficial effects of reducing incident blue light to help slow the progression of cataracts.

Retinitis Pigmentosa

Retinitis Pigmentosa (R.P.) is a degenerative genetic condition resulting in debilitating vision loss. Initially, the R.P. patient loses night vision as a result of the initial loss of rod vision, also causing the patient to experience tunnel vision.⁷ This is followed eventually by the loss of cone vision. Patients with R.P. very often experience difficulties in bright conditions thought to be due to an abnormal recovery time after bleaching.² It has been proposed that extended exposure to light, especially high energy short wavelength light, might exacerbate pigmentary retinal degeneration's such as R.P. This phenomenon was demonstrated by Dowling and Sidman in 1962⁵, in certain pink eyed rats with inherited retinal dystrophy's. Berson⁸, however, showed through his conclusive studies in long term monocular occlusion, that the progression of R.P. *did not* slow in the absence of light. However, because of the known benefits to ocular tissue and the many subjective reports of increased comfort, practitioners should not hesitate to prescribe S.F. lenses to R.P. patients.^{9,7} When introduced in 1981, the Corning CPF 550 lens was touted as the first filter designed specifically for the R.P. patient.¹⁰ The CPF lens is photochromic, and appears red as it absorbs light in the blue end of the spectrum. Morrisette et al.¹¹ found that 26 out of 36 R.P. patients who tried this lens accepted it and reported such improvements as decreased dark adaptation time, fewer headaches and improved visual functioning and comfort when compared to a non tinted lens. Lynch and Brilliant¹², comparing the Corning CPF lens (S.F. lens) with Noir neutral density filters, found slight acuity increases in 23 out of 30 eyes tested, while no differences in contrast sensitivity or adaptation times could be elicited. Silver and Lyness¹³, in their study titled "Do retinitis pigmentosa patients prefer red *photochromic* lenses" concluded that most all R.P. patients subjectively prefer the use of a filtering lens, however there was no evidence that patients preferred a *photochromic* red lens over a fixed red tint. He continued that a significant number of patients *did* prefer the photochromic lenses but just as many subjects disliked them and therefore the consciences prescriber must have reservations about prescribing any tint without a trial.

Interestingly, it has been shown that neutral density filters actually result in a considerable *reduction* in visual performance (acuity measured under different conditions); no reduction in performance was demonstrated with S.F. lenses.¹

Macular degeneration

Age related macular degeneration (A.R.M.D.) accounts for 9500 new cases of legal blindness each year in the United States² and is the leading cause of blindness in the Western World.¹⁴ There are many reported functional changes associated with A.R.M.D. including: glare, reduced contrast sensitivity, reduced acuity and impaired color discrimination². It is well established that filters will decrease incident light thus decreasing amount of scatter and associated discomfort. Similarly, it is well established that in many conditions, S.F. lenses lead to a *subjective* increase in comfort and performance. Leat et al¹, studied low vision patients with 17 presenting conditions and objectively compared differences in performance through neutral density filters and S.F. lenses. Leat states that while subjective improvements were noted, patients suffering from macular or retinal degeneration's tend to benefit the least by using filters of any kind. Corning Optical¹⁵ however, states that A.R.M.D. patients are the primary users of their S.F. lenses, thus suggesting that visual performance or comfort, whether subjective or objective, is significantly enhanced through their use.

Leat¹ studied objective contrast sensitivities of patients with a wide variety of ocular conditions. Performance was evaluated under conditions of glare and under normal illumination with neutral density filters and a variety of S.F. lenses. Of the 44 subjects taking part, less than one-half showed improvement with any filter. She concluded that in the glare free conditions, people with posterior pole disease, such as maculopathies, retinopathies (diabetic, hypertensive) are more representative of the group gaining no improvement. Alternatively, anterior eye conditions (albinism, coloboma, aniridia, media opacities) are more numerous amongst the

group gaining improvement with S.F. lenses. Similarly under glare conditions she reports more improvement was noted with those with anterior eye conditions and more improvement was noted with S.F. lenses than with neutral density. While there may be a relationship between pathological etiology and the preference of a neutral density or a S.F. lens, there does not appear to be a correlation between etiology and the *absorptive properties* of the S.F. chosen. This was demonstrated in Nguyen and Hoeft's 1994 study on Corning blue blocker filters and related pathologies. Their investigation demonstrated that there was no correlation between location of the ocular pathology and the transmission cutoff point of the S.F. lens chosen.¹⁶

Conclusion:

Although the potential benefits of S.F. lenses, on comfort and visual performance, have long been realized, relatively few companies have invested in their production. The present paper supports Tokai Optical Companies reported transmittance curves of their medical lenses. Results show reported transmissions differed from the present study's, by no more than 3.5 percent transmittance with the majority differing by less than 1 percent. While some trends are evident, in general, research on the benefits derived from using such lenses remains uncertain. The reviewed literature suggests cataract patients prefer a more neutral density filter under normal conditions while in conditions of high glare, some selective filtering characteristics proved beneficial. Evidence suggests, to increase patient comfort, more emphasis should be placed on reducing the *total* amount of incident light. However, because blue light is thought to exacerbate cataract development, S.F. lenses would likely be of benefit.

Patients with retinitis pigmentosa experienced subjective improvements in visual comfort and performance with S.F. lenses relative to neutral density filters. An objective increase in acuity was demonstrated with a significant number of patients with S.F. lenses. Literature suggests there is no strong preference for *photochromic* lenses over a fixed tint. Visual performance may

actually be compromised with the use of a neutral density filter. Considering the literature reviewed, R.P. patients may greatly benefit from S.F. lenses, both in comfort and acuity. Patients suffering from macular degenerative conditions were shown to be among the group benefiting the least from *any* filter in normal lighting conditions however, under glare situations, S.F. lenses tended to be more beneficial than neutral density lenses. In general, people with age related macular degeneration see improvement in acuity and contrast sensitivity with increased illumination and thus a neutral density filter may be contraindicated.

It must be stressed that the literature reviewed reports *trends* , and in a clinical situation, individuals are likely to vary in their preferences when attempting to chose a medical lens. Unless patients are having difficulty meeting a desired acuity demand, the practitioner should rely more on the patients *subjective* responses rather than objective findings when choosing a medical lens. With new lenses entering the market this task may be made more difficult but the potential rewards to patients may be tremendous.


References:

1. Leat SJ, North RV, and Bryson H. Do long wavelength pass filters improve low vision performances? *Ophthal Physiol Optics* 1990; 10:219-224.
2. DeSylvia D. Low vision and aging. *Optom and Vision Sci* 1990; 67:(5) 319-322.
3. Trokel S. The physical basis for transparency of the crystalline lens. *Invest Ophthal* 1962; 1:493-500.
4. Tupper B, Miller D, and Miller R. The effect of a 550 nm cutoff filter on the vision of cataract patients. *Annals of Ophthal* 1985; 17(1): 67-72.
5. Bither PP, Hurt JD. CPF-550 vs. C-Lite: a comparison study. *Journal of the Amer Optom Assn* 1988; 59:(8) 623-628.
6. Chou BR, Cullen AP. Spectral characteristics of sports and occupational tinted lenses. *Canadian J of Optom* 1985; 47:77-88.
7. Yolton RL. Student's guide to visual information processing and perception. Pacific University College of Optometry. 14,168.
8. Berson EL. Light deprivation in retinitis pigmentosa. *Vis Res* 1980; 20:1179-84.
9. Woo GC, and Hackner AJ. An update on therapies for retinitis pigmentosa. *Canadian J of Optom* 1981; 43(1): 46-48.
10. Pensyl CD. Matching patient needs with the features of photochromic filters. *J of Vision Rehab* 1993; 7:(1) 10-12.
11. Morrisette DL, Mehr EB, Keswick CW and Lee PN. Users' and nonusers' evaluations of the CPF 550 lenses. *Am J of Optom and Physiol Optics* 1984; 61:(11) 704-710.
12. Lynch DM, Brilliant R. An evaluation of the Corning CPF550 Lens. *Optom Monthly* 1984; 75:36-42.
13. Silver JH and Lyness AL. Do retinitis pigmentosa patients prefer red photochromic lenses? *Ophthalmological Physiol Optics* 1985; 5(1): 87-88.
14. Kanski JJ. Clinical ophthalmology. London: Butterworth-Heinemann, 1990. 349-355.
15. Corning Medical Optics. Corning glare control lens manual 1991. Corning NY.
16. Nguyen TV and Hoeft WW. A study of Corning blue blocker filters and related pathologies. *J of Vision Rehab* 1994; (8):3 15-21.


Appendix 1. Spectrophotometer plots of percent transmission versus wavelength of Tokai Medical Lenses.

Pacific: Yellow CCP

Trace	A	Operator		Sample
Mode	%T	SBW (nm)	2.0	Smoothing (sec)
Ord	Max/Min	100.0	0.0	WL Max/Min (nm)
Speed (nm/min)	200	500.0	280.0	
5 Peaks, threshold 1.0				
Max	771.6 nm, 90.9	Max	706.9 nm, 90.0	
Max	571.3 nm, 89.5	Max	558.1 nm, 88.5	
Max	539.3 nm, 87.0			




Pacific: Orange CCP




Pacific: Light Green CCP

Trace			Operator		Sample
Mode	%T		SEW (nm)	2.0	Smoothing (sec)
Ord Max/Min	100.0	0.0	WL Max/Min (nm)	800.0	280.0
Speed (nm/min)	200				
Cycles	1		Time (min/cyc)	2.9	




Pacific: Dark Green CCP

Trace			Operator	Sample	
Mode	%T		SBW (nm)	2.0	Smoothing (sec)
Ord Max/Min	100.0	0.0	WL Max/Min (nm)	800.0	280.0
Speed (nm/min)	200				
Cycles	1		Time (min/cyc)	2.9	




Pacific: Brown CCP

Trace			Operator		Sample
Mode	%T		SBW (nm)	2.0	Smoothing (sec)
Ord Max/Min	100.0	0.0	WL Max/Min (nm)	800.0	280.0
Speed (nm/min)	200				
Cycles	1		Time (min/cyc)	2.9	




Pacific: Amber CCP400

Trace %T
 Mode %T
 Wavelength 100.0
 Speed (nm/min) 200
 Cycles 1
 Operator
 SBW (nm) 2.0
 WL Max/Min (nm) 800.0
 Smoothing (sec) 280.0
 Sample
 Time (min/cyc) 2.9




Pacific: Light Green CCP400




Pacific: Yellow CCP400

Trace				Operator		Sample
Mode	%T			SEW (nm)	2.0	Smoothing (sec)
Grd Max/Min	100.0	0.0		WL Max/Min (nm)	800.0	280.0
Speed (nm/min)	200					
Cycles	1			Time (min/cyc)	2.9	




Pacific: Green CCP400




Pacific: Olive CCP400

Trace %T
 Mode %T
 Snd Max/Min 100.0
 Speed (nm/min) 200
 Cycles 1
 Operator
 Sbw (nm) 2.0
 WL Max/Min (nm) 800.0
 Smoothing (sec) 1
 Sample
 Time (min/cyc) 2.9




Pacific: Evergreen CCP400

Trace				Operator		Sample	
Mode	%T			SBW (nm)	2.0	Smoothing (sec)	1
Ord Max/Min	100.0	0.0		WL Max/Min (nm)	800.0	280.0	
Speed (nm/min)	200						
Cycles	1			Time (min/cyc)	2.9		



94/10/31 11:06 AM

%T




Tokai: Yellow CCP

サンプル名 : CCP YL HARDCOATED

コメント :
スキャンスピード : 300(750) nm/min
スリット (近接外) : 自動制御
スリット (可視) : 6.00 nm
ホトマル電圧 : 自動制御
PbS 感度 : 2

94/10/31 11:13 AM



サンプル名 : CCP RO HARDCOATED

コメン

スキャンスピード : 300(750) nm/min


スリット(可視) : 6.00 nm

ホトマル電圧 : 自動制御

PbS 経度 : 2

Tokai: Orange CCP

94/10/31 11:25 AM



Tokai: Light Green CCP

サンプル名 : CCP YG HARDCOATED


コメント :

スキャンスピード : 300(750) nm/min

スリット(可換) : 6.00 nm

スリット(近郊外) : 自動制御 ホトマル電圧 : 自動制御 PbS感度 : 2


94/10/31 01:09 PM



Tokai: Dark Green CCP

サンプル名 : CCP UG HARDCOATED
コメント :
スキャンスピード : 600(1500) nm/min スリット(可視) : 6.00 nm
スリット(近紫外) : 自動制御 ホットマル電圧 : 自動制御 P b S 感度 : 2


94/10/31 01:19 PM



Tokai: Brown CCP

サンプル名 : CCP BR HARDCOATED
コメント :
スキャンスピード : 300(750) nm/min スリット(可視) : 6.00 nm
スリット(近紫外) : 自動制御 ホトマル電圧 : 自動制御 PbS感度 : 2


94/10/31 08:48 PM




サンプル名 : CCP400 AC HARDCOATED
コメント :
スキャンスピード : 300(750) nm/min スリット(可視) : 6.00 nm
スリット(近紫外) : 自動制御 ホトマル電圧 : 自動制御 P b S 感度 : 2

Tokai: Amber CCP400


Tokai: Light Green CCP400



Tokai: Yellow CCP400




Tokai: Green CCP400



サンプル名 : CCP400 NL HARDCOATED
コメント :
スキャンスピード : 300(750) nm/min スリット(可視) : 6.00 nm
スリット(近紫外) : 自動制御 ホトマル電圧 : 自動制御 P b S 感度 : 2

94/10/31 05:58 PM

2T



サンプル名 : CCP400 TR HARDCOATED


コメント :

スキャンスピード : 300(750) nm/min スリット(可変) : 6.00 nm

スリット(近赤外) : 自動制御 ホトマル電圧 : 自動制御 P b S 感度 : 2

Tokai: Olive CCP400

94/10/31 05:53 PM



CCP400

Evergreen

Tokai:

サンプル名 : CCP400 FR HARDCOATED


コメント :

スキャンスピード : 300(750) nm/min スリット(可視) : 6.00 nm

スリット(近紫外) : 自動制御 ホトマル電圧 : 自動制御 P b S 感度 : 2

94/10/31 11:11 AM

%T



サンプル名 : CCP OY HARDCOATED

コメント :

スキャンスピード : 300(750) nm/min

スリット(可変) : 6.00 nm

スリット(近接外) : 自動制御 ホトマル電圧 : 自動制御 P b S 腔壁 : 2

Tokai: Red CCP