

Pacific University

CommonKnowledge

---

College of Optometry

Theses, Dissertations and Capstone Projects

---

4-1986

## The effect of extended wear contact lenses on corneal thickness and the endothelial layer

David O. Marrs  
*Pacific University*

### Recommended Citation

Marrs, David O., "The effect of extended wear contact lenses on corneal thickness and the endothelial layer" (1986). *College of Optometry*. 1300.  
<https://commons.pacificu.edu/opt/1300>

This Dissertation is brought to you for free and open access by the Theses, Dissertations and Capstone Projects at CommonKnowledge. It has been accepted for inclusion in College of Optometry by an authorized administrator of CommonKnowledge. For more information, please contact [CommonKnowledge@pacificu.edu](mailto:CommonKnowledge@pacificu.edu).

---

# The effect of extended wear contact lenses on corneal thickness and the endothelial layer

## Abstract

Corneal thickness and the corneal endothelium were measured and observed during the extended wear of high (74 %) water content, medium (55%) water content, and low (38%) water content contact lenses. Thirty-nine patients were monitored for a five month period. Two different water content contact lenses were worn by each patient. The results of the study indicate no significant change in corneal thickness occurred with any of the lenses. Endothelial cell density was not effected by any of the different lenses, but morphological changes, increased polymegathism, did occur with all. The degree of successful wear of extended wear contact lenses will vary among individual subjects. Only nine (23%) of the subjects were recommended for continued extended wear at the end of this study.

## Degree Type

Dissertation

## Degree Name

Master of Science in Vision Science

## Committee Chair

Don C. West

## Subject Categories

Optometry

### Copyright and terms of use

If you have downloaded this document directly from the web or from CommonKnowledge, see the "Rights" section on the previous page for the terms of use.

**If you have received this document through an interlibrary loan/document delivery service, the following terms of use apply:**

Copyright in this work is held by the author(s). You may download or print any portion of this document for personal use only, or for any use that is allowed by fair use (Title 17, §107 U.S.C.). Except for personal or fair use, you or your borrowing library may not reproduce, remix, republish, post, transmit, or distribute this document, or any portion thereof, without the permission of the copyright owner. [Note: If this document is licensed under a Creative Commons license (see "Rights" on the previous page) which allows broader usage rights, your use is governed by the terms of that license.]

Inquiries regarding further use of these materials should be addressed to: CommonKnowledge Rights, Pacific University Library, 2043 College Way, Forest Grove, OR 97116, (503) 352-7209. Email inquiries may be directed to: [copyright@pacificu.edu](mailto:copyright@pacificu.edu)

PACIFIC UNIVERSITY LIBRARY  
FOREST GROVE, OREGON

THE EFFECT OF EXTENDED WEAR CONTACT LENSES  
ON CORNEAL THICKNESS AND THE  
ENDOTHELIAL LAYER

A Thesis  
Presented to  
the Faculty of  
Pacific University

In Partial Fulfillment of  
the Requirements for the  
Degree Master of Science in  
Clinical Optometry (Management Track)


Submitted by  
David O. Marrs, O.D.

April 1986


THE EFFECT OF EXTENDED WEAR CONTACT LENSES  
ON CORNEAL THICKNESS AND THE  
ENDOTHELIAL LAYER

Place: Pacific University

Approved:

  
Don C. West, B.S., O.D., Chairman


April 21, 1986  
Date

  
James E. Peterson, B.S., O.D.

April 21, 1986  
Date

  
A. Richard Reinke, B.S., O.D.

April 24, 1986  
Date

  
John E. Roggenkamp, B.S., O.D.

4/29/86  
Date

## TABLE OF CONTENTS

ABSTRACT . . . . .	i
LEGEND OF FIGURES . . . . .	ii
LEGEND OF TABLES . . . . .	iii
LEGEND OF GRAPHS . . . . .	iv
INTRODUCTION . . . . .	1
LITERATURE REVIEW . . . . .	3
STATEMENT OF THE PROBLEM . . . . .	20
METHODOLOGY . . . . .	21
TREATMENT OF DATA . . . . .	38
RESULTS . . . . .	41
DISCUSSION AND CONCLUSIONS . . . . .	51
REFERENCES . . . . .	54
APPENDICES . . . . .	59
A. CONTACT LENS INSTRUCTIONS . . . . .	59
B. INFORMED CONSENT FORM . . . . .	60
C. CONTACT LENS FORMS (SAM FORMS 1 THROUGH 8) . . . . .	61
D. PACHOMETER COMPUTER PROGRAM . . . . .	62
E. PACHOMETER INSTRUCTION MANUAL . . . . .	63
F. PACHOMETRY SUBJECT DATA FILES . . . . .	64
G. PACHOMETRY SPREADSHEET, THICKNESS AND PERCENTAGE CHANGES . . . . .	65
H. STATISTICAL DATA OF CORNEAL THICKNESS MEASUREMENT AT FIVE LOCATIONS OVER TIME . . . . .	66

I.	STATISTICAL DATA AND FREQUENCY DISTRIBUTION ON FIVE CORNEAL LOCATIONS PRIOR TO WEAR OF CONTACT LENSES . . . . .	67
J.	ENDOTHELIUM CELL COUNT AND STATISTICS - SUBJECTS AND CONTROL PATIENTS . . . . .	68
K.	KEYSTONE DIAGNOSTIC TESTING - AVIATORS UNIT . .	69

The opinions or assertions contained herein are the private views of the author and are not to be construed as official or as reflecting the views of the Department of the Army or the Department of Defense.



## ACKNOWLEDGEMENTS

My gratitude and appreciation goes to my advisor and teacher, Professor Don West, for his advice, guidance, interest and his ability to keep me motivated during the course of this study. I also wish to thank Doctors A. Richard "Dick" Reinke, James E. Peterson, John R. Roggenkamp, Niles Roth, Robert L. Yolton, Bradley M. Coffey and Eileen Brennan for their expert guidance and continual support. Thanks also to Laurel Gregory in her expert guidance of my tour of the library.

A special thank you to Barnes Hind, Inc., CooperVision, Inc. and Syntex Ophthalmics, Inc. for furnishing all of the contact lenses for this research project. Thanks also to Allergan Pharmaceuticals, Inc., American Optical, Inc. and Blairex Laboratories, Inc. for their donation of all the supplies and solutions used during the project and to Humphrey Instruments, Inc. for the loan of the Automatic Refractor and the Auto Keratometer.

To the United States Army, I extend my deepest gratitude for the opportunity to attend and the funding for this graduate program.

To Doctor Tim Allen and Doctor Gary Slater, thank you for your outstanding ability to cooperate, cajole and otherwise keep me and all those patients coming back for more.

To all of the patients, who put up with the frequent and lengthy examinations, thanks for the superb cooperation and steadfast support.

A very special thank you to Sharon, John and Ashley for their sacrifices, love and support while aiding me in the completion of this degree.

## ABSTRACT

Corneal thickness and the corneal endothelium were measured and observed during the extended wear of high (74%) water content, medium (55%) water content, and low (38%) water content contact lenses. Thirty-nine patients were monitored for a five month period. Two different water content contact lenses were worn by each patient. The results of the study indicate no significant change in corneal thickness occurred with any of the lenses. Endothelial cell density was not effected by any of the different lenses, but morphological changes, increased polymegathism, did occur with all. The degree of successful wear of extended wear contact lenses will vary among individual subjects. Only nine (23%) of the subjects were recommended for continued extended wear at the end of this study.

## LEGEND OF FIGURES

FIGURE	PAGE
1. PACHOMETER OPTICS . . . . .	10
2. TOUCH Vs.OVERLAP METHOD . . . . .	13
3. COUNTING GRIDS . . . . .	17
4. PHOTOGRAPH OF PACHOMETER . . . . .	31
5. PHOTOGRAPH OF ENDOTHELIAL CAMERA . . .	34
6. ARRANGEMENT OF CAMERA SYSTEM . . . . .	35
7. ENDOTHELIAL PHOTOGRAPH . . . . .	36

## LEGEND OF TABLES

TABLE	PAGE
1. DISTRIBUTION OF CONTACT LENS . . . . .	26
2. DISPOSITION OF PATIENTS NOT COMPLETING STUDY . . . . .	41
3. REASONS FOR REPLACEMENT OF CONTACT LENSES DURING RESEARCH PERIOD . . . . .	42
4. PATIENT RESPONSE TO QUESTIONARIE COMPARING LENSES . . . . .	43
5. RESULTS OF NON-CONTACT TONOMETER . . . . .	44
6. RESULTS OF KEYSTONE (AVIATORS UNIT). . . . .	45
7. GRADING OF POLYMEGATHISM FOUND IN RESEARCH SUBJECTS . . . . .	50

## LEGEND OF GRAPHS

GRAPH	PAGE
1. SPECTACLE Rx DISTRIBUTION . . . . .	22
2. DISTRIBUTION OF PATIENTS BY FLAT KERATOMETER READING . . . . .	23
3. PATIENT AGE DISTRIBUTION . . . . .	24
4. PERCENTAGE OF CENTRAL CORNEAL THICKNESS CHANGE WITH CSI-T LENSES . .	46
5. PERCENTAGE OF CENTRAL CORNEAL THICKNESS CHANGE WITH HYDROCURVE II LENSES . . . . .	47
6. PERCENTAGE OF CENTRAL CORNEAL THICKNESS CHANGE WITH PERMAFLEX LENSES . . . . .	48
7. PERCENTAGE OF CENTRAL CORNEAL THICKNESS CHANGES WITH HYDROCURVE II, PERMAFLEX, AND CSI-T LENSES . . . . .	49

## INTRODUCTION

Extended wear contact lenses are worn by approximately five million Americans and account for approximately 35 percent of the soft lens market with sales expanding rapidly.<sup>1-3</sup> The safety of extended wear lenses is of principal concern among the nation's eye care providers.<sup>2,4</sup>

Extended wear contact lenses are manufactured today by 18 United States manufacturers with several more presently undergoing Federal Drug Administration approval.<sup>4</sup> These hydrogel lenses range in water content from 32 to 79 percent.

The purpose of this research program was to study the effects of three types of extended wear contact lenses on the eye. The principle objective of the research is the effect that extended wear contact lenses have on corneal thickness and the endothelial layer of the cornea.

This research program was conducted in conjunction with two additional research programs in order to cover a broad scope of extended wear application. One of the additional programs investigated the changes in oxygen transmissibility of extended wear flexible contact lenses over time as a result of patient wear. The other program investigated the effects of extended wear contact lenses on contrast sensitivity function with time.

The Armed Forces of the United States are presently studying the use of extended wear contact lens on selected military personnel. The intent of this study is to evaluate

the overall safety of extended wear contact lenses and to provide procedural guidelines for the use of extended wear contact lenses.



## Literature Review

Leonardo da Vinci was credited with the concept of the first contact lens in 1508 and in 1888 A. Eugen Fick described the first contact lenses worn by a human subject.<sup>5</sup> The knowledge and research of contact lenses has grown by leaps and bounds since these early attempts. With the discovery in 1960 of hydroxyethyl methacrylate, or HEMA, in Czechoslovakia came the introduction of the soft contact lens. The rights to this material in the United States was sublicensed to Bausch & Lomb in 1966 and was approved by the Federal Food and Drug Administration in March of 1971.<sup>6</sup>

Approximately 120 million people (51% of total population) were using eyeglasses by 1982 and between 16.3 to 18 million (7%) were using contact lenses. It is estimated that 20 million Americans are wearing some form of contact lenses today.<sup>1</sup> The percentage using contact lenses has increased dramatically from 1.3% in 1971 to approximately 7% in 1982.<sup>7</sup>

Soft contact lenses have been used for the correction of refractive errors since their approval by the Food and Drug Administration in 1971. With the approval of extended wear contact lenses for correction of myopia in 1981, there has been an increasing demand for their use.<sup>8</sup> With the increasing popularity of contact lenses and the tremendous advancement in the technology of manufacturing, contact lenses may well overtake glasses in usage.<sup>7</sup> The use of

contact lenses has increased by about 1.5 million pairs per year since 1979.<sup>7</sup> It is estimated that approximately 17.5 million people, about 7 percent, were using contact lenses by 1982 and twenty million Americans are now wearing contact lenses.<sup>1,7</sup> With the large amount of advertisement by the optical companies and the optical chains, the popularity of extended wear contact lenses has exploded in today's market. The number of optometrists fitting extended wear contact lenses also has increased dramatically in the past four years. In a recent survey it is shown that approximately 90% of optometrists were fitting extended wear cosmetic and aphakic contact lenses.<sup>9</sup> In October, 1984, there were fourteen brands of extended wear lenses approved by the Food and Drug Administration for cosmetic wear while today there are eighteen brands available.<sup>10</sup> Additional lenses are currently under study by the FDA for approval.

There has been an abundance of studies to determine the safety and effects of extended wear contact lenses on the eye.<sup>11-16</sup> These studies have been accomplished by a variety of investigators with a wide assortment of controls. Some individuals are quite adamant in their belief that extended wear lenses are quite dangerous, and should not be fit for the routine cosmetic correction.<sup>17</sup> Others believe the lenses are highly successful.<sup>18,19</sup> It was the purpose of this study to evaluate the safety of extended wear lenses and to observe and measure the changes of the ocular system.

## CORNEAL THICKNESS

The first portion of this study dealt with changes in corneal thickness due to extended contact lens wear. Review of the literature shows various studies have been performed in the past few years with mixed results. These range in time from a few hours wear up to several months. Studies also vary in how the lenses were worn, with some wearing continuously and others a matter of a few days before removal.<sup>12-15,20</sup>

In one study of central corneal thickness with permanent (not removed for six months) wear, a thickness increase of 4-6% was found after the first 24 hour adaptation period.<sup>20</sup> The thickness increase of 4-6% gradually decreased to 2.5-3.1% after a six month period. Patients had as much as 9% increase during sleep, which took three to four hours to reverse after awakening.<sup>20</sup> Another study shows a mean increase of .010mm. after seven months of wearing extended wear lenses.<sup>14</sup> A study of extended wear lenses for thirteen weeks showed no statistically significant difference in pachometry, visual acuity, refraction, and keratometry between the beginning and the end of the study.<sup>15</sup> In another study of seven days continuous wear with three types of lenses, the mean overnight swelling ranged from 9.7% to 15.1%.<sup>13</sup> The mean swelling after the eyes were open for a period of twelve hours was from 1.6% to 5.8%, with the mean deswelling of the eyes over the twelve hour period to be

8.2%. Also included is a study of twenty (20) weeks of continual wear which resulted in no significant evidence of corneal swelling.<sup>12</sup> Larke summarizes his findings as follows<sup>21</sup>:

- a. Thick Hydrogel lenses provoke more stromal swelling than thin lenses.
- b. The extent of swelling declines with time.
- c. The extent of swelling rarely exceeds the maximum amount that can be provoked by anoxia alone (8%). In cases where this is not the case the posterior stroma shows evidence of gross change in the form of 'straie' and folds in Descemet's layer.
- d. The use of thin (0.1mm) lenses provokes on average less than one percent stromal swelling.

In reviewing these articles, there appears to be a wide variety of responses recorded.

It was the purview of this portion of the study to record the corneal thickness changes over a five month period on a regular basis and to evaluate the effects of extended wear contact lenses on the corneal thickness. The instrument utilized throughout this research was the Electronic Pachometer. The pachometer system consisted of a slitlamp, a fixation device, an optical doubling device, a potentiometer, and an Apple IIe computer. Further description of the pachometer may be found in the methodology section.

The pachometer is an extremely useful tool in evaluation of the cornea during contact lens wear. A French surgeon named Petit appears to be the first recorded (1723) individual to attempt measuring the thickness of the cornea of the human eye. His measurements were accomplished on

enucleated eyes of cadavers and were found to be about 0.4mm.<sup>22</sup> Physiologist, Blix in 1880 invented an instrument for a measurement in the living eye.<sup>23</sup> Blix' instrument had two identical microscope tubes, horizontally converging at an angle of 40° to a point in front of the tubes.

Gullstrand in 1924, also used the reflection from the two corneal surfaces for his measurements. His system was much more complicated but seemingly more accurate than that of Blix. Tscherning is said to have used the same principle.<sup>22</sup>

The invention of the corneal microscope and Gullstrand's slitlamp were utilized by Hartinger in 1921 to measure corneal thickness. He utilized the corneal microscope and a measurement drum devised by Ulbrich.<sup>22</sup> The corneal microscope was focused on the anterior surface and then on the posterior surface. The difference in the positions was read on the drum and gave the apparent thickness.

Koby and Juillierat also used the slit lamp in their attempt at measurement of corneal thickness. They measured the apparent length of an optical section of the slit lamp beam with an eyepiece micrometer. The real thickness of the cornea was then calculated by a formula that was based on the assumption that the corneal surfaces are concentric.<sup>23</sup> Goldmann in 1932 was apparently the first to use a split ocular to measure the optical section of the slit lamp beam as it passed through the cornea. He used a split ocular similar to the Lobeck eyepiece (a double prism in the plane

of the image of the ocular) and a diaphragm at the site of the exit pupil. In this manner two half images of the cornea are observed contiguous to each other. These are adjusted until they meet, but the hairlines must be kept clear and the refractive error of the observer must be fully corrected.<sup>24</sup>

Von Bahr, in 1948, developed a device based on the principle used by Blix. It allowed simultaneous adjustment of both observation and illumination system, which eliminated the problem of accurate fixation. Von Bahr's apparatus made use of two plane glass lamina which were symmetrically movable about vertical axes in front of the optical systems for both illumination and observation. Only the rays of the lower halves of the optical systems passed through the glass laminae. Thus, by proper rotation of these laminae, one could cause the specular reflection of the corneal surface to coincide with that of the endothelial surface. The rotation of the laminae was proportional to the thickness of the cornea, and from this angle could be calculated the real thickness of the cornea.<sup>23,24</sup>

Maurice and Giardini in 1951 used a modified von Bahr apparatus. They eliminated a number of the objections of the von Bahr instrument. Donaldson in 1965 proposed a split ocular with the Haag-Streit slitlamp to measure corneal thickness. He also concluded that the apparent corneal thickness measurement must be modified by a correction factor to obtain true corneal thickness.<sup>25</sup>

Mishima and Hedbys in 1968 developed a slightly different apparatus utilizing the Haag-Streit slitlamp and the standard Pachometer. The investigators developed the foundation for today's clinical pachometric measurements. They made a modification which allows the exact positioning of the slitlamp at the time of the measurement and an accurate determination of the corneal thickness to be made.<sup>26</sup>


The pachometer consists of a main attachment and split-image eyepiece. The main attachment contains two glass plates, the lower fixed and the upper rotatable, and these plates are placed in front of the right microscope. A thin metal diaphragm with a narrow vertical aperture extends from the main attachment. When the arm of the slitlamp and that of the microscope stand are set at an angle of  $35^{\circ}$  (the slitlamp at the observer's left of the microscope), the slitlamp beam passes through the vertical aperture of the diaphragm and makes an angle of  $40^{\circ}$  with the right half of the microscope.

Mishima and Hedbys modification consisted of extending the metal diaphragm approximately twice as long as the original and of attaching two small lamps, pinlites, on this extended diaphragm. The pinlites were placed about 10mm apart from each other on a vertical line.

The principle of the thickness measurement is similar to that described by Jaeger for the determination of the anterior chamber depth.<sup>26</sup> In Figure 1, the light beams (S)

meet the corneal surface at A and the endothelial surface at B. The direction of the beams is perpendicular to the corneal surface. The optical section of the cornea is observed at an angle of  $40^\circ$  from the direction of the slitlamp beams. The endothelium is seen on the extension of line M-C and therefore the apparent thickness of the corneal optical section corresponds to A-D.

Figure 1. The optical arrangement for the measurement of corneal thickness. Light beams meet corneal surface at A and the endothelial surface at B. Optical section of cornea observed at angle of  $40^\circ$  from direction of slit-lamp beams. The endothelium is seen on line extension M-C; thickness of corneal optical section corresponds to A-D (Mishima & Hedbys).






The light beams from the cornea pass through two glass plates and enter the microscope. The split-image eyepiece divides the visual field into upper and lower halves. The light beams which pass through the upper rotatable and the lower fixed plates form the images in the upper and lower visual field, respectively. When the upper glass plate is rotated, the light beam passing through it is displaced by refraction from the original light path, and therefore the upper half of the corneal image is displaced. When the endothelium of the upper half of the optical section is so displaced to be aligned with the epithelium of the lower half, the displacement of the light beam by the rotated glass plate is equal to the apparent thickness of the corneal optical section, A-D of Figure 1. The corneal thickness can, therefore, be calculated from the angle of rotation of the upper glass plate, using the thickness and refractive index of the glass plate, the curvature and refractive index of the cornea.

Mishima and Hedbys modification is used to secure the position of the slitlamp beam perpendicularly to the corneal surface.<sup>26</sup> This positioning of the slit lamp can be made by using images of two pinlites. In Figure 1, the image of pinlite (L) is made at I on the extension of line M-A, when the slitlamp beam is perpendicular to the corneal surface. Therefore, the images of the pinlites are seen on the epithelium of the optical section of the cornea.

Several variables in the measurement of the thickness play an active role in the reliability of the readings. Items such as subject fixation, slit beam width, slit focus, line voltage to the electronics in an electronic pachometer, ambient room illumination, and the subtleties of operator endpoint criterion can change from one moment to another.<sup>27</sup>

Of the above mentioned items, the endpoint of the measurement seems to have the largest influence. It has been shown that there is a statistically significant difference between the touch method and the overlap method (see Figure 2). The touch method is accomplished by having the endothelium of the upper section just touching the epithelium of the lower section.

Figure 2. The touch and overlap methods utilized in corneal thickness measurement during pachometry.



The overlapping measurements were accomplished by having the endothelium of the upper section overlap the epithelium of the lower section. In all cases the thickness indicated by the touch measurement was more than that indicated by the overlap method. It has been demonstrated that the touch technique provides the least variable method.<sup>25</sup> It is further shown that when using a pachometer for corneal studies it is important to realize that only changes greater than one to two percent will exceed the limitations (variability) in measurement and therefore demonstrate a true change in corneal thickness.

## ENDOTHELIUM

The next portion of the research deals with the change in the endothelial layer of the cornea during extended wear of contact lenses. The cornea is a transparent avascular structure in which there is a dynamic balance between salt and water, controlled by the endothelial pump, for the maintenance of corneal deturgescence (dehydration).<sup>21,28-30</sup> Several factors can impair the function of the endothelial layer of the cornea. Among these are normal aging, dystrophy and degeneration of the cornea, surgery and trauma. Considering the significance of the endothelial layer of the cornea, it is extremely important to be able to observe and to photograph this critical layer of the cornea.

The endothelial layer consists of a single layer of hexagonal cells having a thickness of approximately 5 micrometers ( $\mu\text{m}$ ).<sup>21,29</sup> The individual irregular hexagonal cell is normally about 20 micrometers in diameter. The endothelial membrane contains between 500,000 to 600,000 cells. The normal endothelium has between 2,500 - 4,000 cells/sq. mm in the young healthy eye and declines to approximately 1000 to 2000 per square millimeter in the elderly.<sup>21,31</sup> The average cell density for an individual of age 23 is approximately 3,060.<sup>32-34</sup> Capella suggests that between the ages of 15 and 50 the endothelial cell density remains fairly constant, and that pleomorphism occurs at about the age of 60 to 65.<sup>35</sup> Laing et al maintains that the

pleomorphism, due to aging, occurs gradually and progressively.<sup>36</sup>

The endothelial pump has been demonstrated to be a bicarbonate pump which controls the normal hydration, thickness, and transparency of the cornea.<sup>29,30</sup> The exact function of the endothelial pump is not fully understood.

The majority of articles report no mitosis of the endothelial cell except possibly in the very young, although a recent article by Olsen and Davanger has shown otherwise in studies on human eyes in vitro.<sup>37</sup>


Considering the significance of the endothelial layer it is important to observe and photograph this layer of the cornea. The two systems enabling the observer to accomplish this observation are the slitlamp and the specular microscope. The first of these, the slitlamp, the method used in this study is the more common and less expensive method. The slitlamp utilizes specular reflection from the endothelial cell layer to accomplish the observation. In 1970, Brown developed a macrophotographic system for anterior ocular segment photography which gave magnification of 10X.<sup>38</sup> Zantos and Holden and McMonnies and Zantos in 1977 and 1979 used a similar system and a process of rephotography to achieve higher magnification for the study of endothelial changes resulting from contact lens wear.<sup>39,40</sup> The other method for endothelial observation is the specular microscope. In this study, where the endothelium was

monitored, the use of topical anesthetics and the direct corneal contact that are necessary with the specular microscope were contraindicated. The slitlamp endothelial photographic system has the advantage of being noncontact and noninvasive, and is capable of repeatable observation of specific areas.<sup>41</sup> Slit lamp photographs will provide approximately  $0.08 \text{ mm}^2$  for study while the specular microscope provides approximately  $0.035 \text{ mm}^2$ .<sup>31</sup>

After obtaining photographs of the corneal endothelial layer, there are several methods of measurement and quantification. There have been four methods recorded for measurement of the endothelial mosaic.<sup>42</sup> The four are: A rectangle, planimeter, digitizer, and cell sizer. In the rectangle method, also called fixed frame analysis, the observer counts the number of cells in a rectangle and reports a cell density. The Planimeter is used to trace the outline of a single cell or group of cells and measure their area directly. When using the Digitizer, a computer compiles the measurements. A digitizer records the X-Y coordinates of cell apices directly when the photograph is placed on an electronic grid and apices of the cell outline are touched with a pen. The fourth method is the cell sizer, which works by simply comparing the endothelial photograph with a picture of cells of known size. The results from the rectangle, planimeter, and digitizer agree within 10% of one another.<sup>42</sup>

When using slitlamp photography, the most common measurement technique is to project slides of the endothelium onto a counting grid as shown in figure 3.

Figure 3. Counting grid utilized in determination of endothelial cell density.



The cells in six rectangles of equal area are counted and averaged. Each rectangle is approximately the size of the area that would be photographed by a specular microscope. Cells that touch the superior and left edge of the rectangle are counted and those touching the inferior and right side are not counted. A measurement scale photographed at the same magnification as the endothelial cells is projected onto the counting grid to determine the area being photographed. The average number of cells per rectangle is then multiplied by the magnification factor as determined with the measurement scale slide. The result is the cell density

(cells/ sq. mm.).<sup>34</sup> The instrument used in this study was the "Nikon Non-Contact Endothelial Camera".<sup>43</sup>

It has been demonstrated that the endothelium undergoes pleomorphism with age.<sup>36</sup> The mean endothelial cell area is known to increase with increasing age. The mean area of the cells is approximately twice as large in people in their 80's as found in individuals in their 20's.<sup>36</sup> An enlargement of cells is found also in corneas suffering injury or those that have undergone surgery, i.e. cataract extraction.<sup>29</sup> Several investigators have documented transient changes in the endothelium due to contact lens wear.<sup>39,41,44-46</sup> In the majority of new contact lens patients the appearance of endothelial blebs is noted within ten minutes and seems to peak at approximately 30 minutes and then to gradually reduce to normal. The blebs also rapidly disappear upon removal of the lenses. Investigators have further stated that a decrease in number of endothelial cells occurs with the wearing of hard contact lenses. Schoessler and Woloschak state there is a decrease in typical hexagonal cells and increased polymegathism.<sup>47</sup> A study, by Caldwell et al reported that prolonged hard contact lens wear lowers the endothelial cell count.<sup>48</sup> Stocker and Schoessler found that prolonged wear of PMMA contact lenses caused polymegathism (a large variation in cell size) but could not confirm the cell loss reported by Caldwell et al.<sup>48,49</sup>

The purview of this portion of the research was to



monitor the endothelial layer in extended wear contact lens patients and to measure changes in the endothelial cells over a five month period. Both male and female patients were examined to provide a fair representation of the United States Military Forces as well as the population in general. A number of researchers have attempted to correlate ocular changes and the female menstrual cycle.<sup>50-52</sup> El Hage and Beaulne found in their study on this subject, that there was no apparent change in corneal thickness during the cycle.<sup>51</sup>

This research makes a comparison of three types of contact lenses that are manufactured from three different materials and have different water contents. Patients were fit with a different lens on the right and left eye. In recent research in this area, when corneal thickness, corneal oxygen consumption, and the endothelial bleb response were monitored, no significant contralateral effect was observed.<sup>53</sup>

## STATEMENT OF THE PROBLEM

The problem addressed in this study is that of the "acceptable changes" found in the endothelial layer of the cornea and that of corneal thickness during wear of extended wear contact lenses. The first hypothesis which this study will address, assumes changes in corneal thickness will be proportional to the differences in contact lens water content (of various manufacturers) with the use of extended wear contact lenses. The second hypothesis assumes that with differences in water content and therefore oxygen transmissibility there will be morphological changes of the endothelial layer.

The sample used was a group of thirty-nine subjects (seventy-eight eyes) composed of twenty-five females and fourteen males. Age of the subjects ranged between 18 and 35 years of age to parallel the primary military age. The study dealt with spherical lens fittings only of refractive errors between -1.00D to -8.00D. The following conditions were considered absolute contraindications for extended wear contact lenses: Anterior segment infections, blepharitis, dry eye (tear break up time of less than 10 seconds), poor personal hygiene, history of past ocular surgery, trauma or pathology, pregnancy.

## METHODOLOGY

The investigative procedures of this study were conducted at the clinic facilities of the Pacific University, College Of Optometry.

Following a comprehensive visual examination, patients meeting the predetermined criteria were selected for participation in the study. The patients were selected from Pacific University, Forest Grove and surrounding communities.

The criteria established for patient selection included the following:


1. Age - between 18 and 35 years of age (to parallel the primary military age).
2. Refractive Error - Between -1.00 Diopter to -8.00 Diopters (again to stay within military requirements).
3. Tear Break Up Time (TBUT) - Greater than 10 seconds.
4. Subjects may not exhibit any of the following:
  - a. Anterior segment infections.
  - b. Blepharitis.
  - c. Poor personal hygiene.
  - d. History of past ocular surgery, trauma or pathology.
  - c. Pregnancy.

Approximately 100 individuals were screened for participation in this study. Twenty-five (25) females and fourteen (14) males were selected for the study. Patients who had previously worn soft contact lenses were accepted for the study. Abstension from wearing their prior contact lenses for a minimum of fourteen days was required before fitting procedures began. Refraction and Keratometer readings were monitored to insure corneal stability prior to fitting

on those who had been wearing lenses immediately prior to the study. Eleven of the patients selected had worn soft contact lenses at some time prior to the study with only six patients wearing contact lenses at the time of selection for the program.


The following graph illustrates the distribution of spectacle correction of the patients involved in this study. Seventy-eight eyes are represented.

Graph 1. Distribution of spectacle correction of the patients in the extended wear study.



The following graph illustrates the distribution of patients by their flattest keratometer readings. Seventy-eight eyes are represented.

Graph 2. Distribution of Keratometer (Flat) Readings of the patients who participated in the extended wear study.



The following graph illustrates the ages of the patients selected for this study, thirty-nine (39) in number. The mean age was 27.1 within a range of 18 to 35 years.

Graph 3. Ages of the patients who participated in the extended wear study.



The research made a comparison of three types of lenses of different water content. They were the Permaflex (74%) by CooperVision, the Hydrocurve II (55%) by Barnes Hind and the CSI-T (38.5%) by Syntex.

Permaflex™ hydrophilic contact lenses are manufactured by CooperVision Optics, Ltd. using the material Surfilcon A. Permaflex™ has a diameter of 13.5 to 15.0mm. The center thickness for minus lenses ranges from 0.08 to 0.22mm. It is available in 8.70 and 8.90 base curves and from -0.25D to -10.00D in minus power. The lens material, surfilcon A, is a copolymer of methyl methacrylate, vinyl pyrrolidone, and other methacrylates. CooperVision specifies the physical

properties of the Permaflex™ contact lens as follows:

- a. Refractive Index - 1.385 at 25°C
- b. Light Transmittance - >96 percent
- c. Water Content - 74 percent
- d. Oxygen Permeability -  $34.0 \times 10^{-11}$  (cm<sup>2</sup>/sec) (ml O<sub>2</sub>/ml x mm Hg) at 25°C

Hydrocurve II® hydrophilic contact lenses are manufactured by Barnes-Hind, Inc. using the material Bufileon A. Hydrocurve II® has a diameter of 14.0 to 16.0mm. The center thickness for low minus lenses is 0.05mm. It is available in 8.50, 8.80 and 9.10 base curves and from -0.25D to -12.00D in minus power. The lens material, bufileon A, is a hydrophilic random copolymer of 2-hydroxyethyl methacrylate, N-(1,1-dimethyl-3-oxobutyl)-acrylamide, and methacrylic acid. Barnes-Hind indicates that the physical properties of the Hydrocurve II® contact lens are as follows:

- a. Refractive Index - 1.41
- b. Light Transmittance - 90 percent
- c. Water Content - 55 percent
- d. Oxygen transmissibility - Dk/L  $22.7 \times 10^{-9}$

CSI®T hydrophilic contact lenses are manufactured by Syntex Ophthalmics, Inc. using the material Crofileon A. CSI®T has a diameter of 13.0 to 15.0mm. The center thickness for minus lenses ranges from 0.02 to 0.40mm. It is available in 8.30, 8.60, 8.90 and 9.35 base curves and from plano to -7.00D in minus power. The lens material, crofileon A, is a 61.5% 2,3-dihydroxypropyl methacrylate and methyl methacrylate cross-linked with ethylene glycol dimethacrylate. Syntex Ophthalmics, Inc. specifies the

physical properties of the CSI<sup>®</sup>T contact lens as follows:

- a. Refractive Index - 1.44
- b. Light Transmittance - >95 percent
- c. Water Content - 38.5 percent
- d. Oxygen Permeability -  $12.0 \times 10^{-11}$

Patients were supplied with a different water content contact lens on each eye. The study was single blind since the patients were not told that they were wearing two different brands of lenses. All patients, with one exception, were fit with the Hydrocurve II (55%) on one eye and either the Permaflex (74%) or the CSI-T (38.5%) on the other eye. The one exception, due to fitting difficulty, was fit with a Permaflex on one eye and the CSI-T on the other. This subject was unable to adapt to extended wear and was dropped from the study during the first week and that subject's data is not included in the results.

The following table shows the distribution of contact lenses by brand, base curve-diameter and number of eyes wearing that lens.

Table 1. Distribution of contact lenses by brand, base curve-diameter and number of eyes wearing the lens.

Brand	Base Curve-Diameter	Number Of Eyes
Permaflex <sup>™</sup>	8.70 - 14.4	20
Permaflex <sup>™</sup>	8.90 - 14.4	1
Hydrocurve II <sup>®</sup>	8.50 - 14.0	3
Hydrocurve II <sup>®</sup>	8.80 - 14.5	29
Hydrocurve II <sup>®</sup>	9.10 - 14.5	6
CSI <sup>®</sup> T	8.60 - 13.8	11
CSI <sup>®</sup> T	8.60 - 14.8	3
CSI <sup>®</sup> T	8.90 - 14.8	4
CSI <sup>®</sup> T	9.35 - 14.8	1



An extensive trial lens set was made available by the three companies so that optimal fittings could be accomplished. At least two researchers, usually three, were present at each fitting to insure the subjects received the best fitting lenses. At the time of dispensing, each subject was instructed in insertion, removal and care by the same researcher to insure that each received the same instructions. Each patient also received written instructions on insertion, removal and care which are attached as appendix A. All subjects signed the Informed Consent Form required by the Institutional Review Board. A copy of the Informed Consent Form is attached as appendix B.

All patients were given the same wearing schedule as follows:

- Day One - Six (6) hours of wear and removed overnight.
- Day Two - Ten (10) hours of wear and removed overnight.
- Day Three - Worn all waking hours and removed overnight.
- Day Four - Extended wear.

After initial adaptation, subjects were instructed to wear lenses for a period of one week and then to remove overnight for cleaning. They were further instructed that contact lens removal at any time for cleaning with Pliagel® was permitted if the contact lenses were returned to the eye. The American Optical Septicon system was provided to all patients for their weekly cleaning regimen. The Blairex System II, consisting of the Blairex Deionizer and salt

tablets was provided to each subject to prepare their normal saline. Patients were told that they could use a drop of Clerz<sup>®</sup> 2 for lubrication as needed.

The patients reported to the clinic for follow-up examinations according to the following schedule: Day 1 of extended wear, day 4, day 7, day 14, day 21, day 28, day 60, day 90, day 120 and day 150.

During each visit the following were completed as part of the examination:

- a. History/Discussion.
- b. Visual acuity with a Clason Projector.
- c. Contrast Sensitivity Function with contact lenses.
- d. Automatic refraction with the Humphrey's Automatic Refractor with contact lenses.
- e. Keratometry with the Humphrey's Auto Keratometer with contact lenses.
- f. Slit Lamp examination.
- g. Lens permeability and thickness measurement.
- h. Contrast Sensitivity Function without contact lenses.
- i. Endothelial photography.
- j. Pachometry.
- k. Automatic refraction with the Humphrey's Automatic Refractor without contact lenses.
- l. Keratometry with the Humphrey's Auto Keratometer without contact lenses.
- m. Keratometry with the Bausch & Lomb Keratometer.
- n. Corneal photography of both corneas.
- o. Intra Ocular Pressure with the American Optical Non-Contact Tonometer.

The initial fitting of contact lenses was recorded on the "Contact Lens Fitting Form" (SAM FORM 1) and each follow-up examination was recorded on the "Contact Lens Followup Form" (SAM FORM 2). A performance depth perception test was taken prior to starting the study and at the conclusion of the study. The Keystone diagnostic series (aviators unit) was

used as a visual performance determination. Baseline testing (pre-contact lens wear) was done with best spectacle correction prior to start of contact lens wear. Final testing (post-contact lens wear) was preformed (a) while wearing contact lenses and (b) with best spectacle correction at the conclusion of the study. The results of this testing was recorded on "Performance Depth Perception Test" (SAM FORM 3).

Meticulous records of the endothelial photographs were maintained in the Photograph Log Book on the "Photo Record" (SAM FORM 4).

Four control patients were utilized during the course of this research. Two female and two male patients were selected that met the general criteria utilized for subjects. These control patients had not previously worn contact lenses and agreed not to wear contact lenses during the duration of the study. The results of the testing accomplished on these patients was recorded on the "Control Patient Form" (SAM FORM 5).

All patients were fit with a single brand of contact lens at the conclusion of the research period. A questionnaire was completed by each contact lens wearer at the completion of the 150 day examination. "Extended Wear Contact Lens Patient Questionnaire" (SAM FORM 7) was utilized for this purpose. One of the questions asked the patients was: "Do you have a preference between your right and left lenses during the study with:

- (a) Vision - Right\_\_\_ Left\_\_\_ No Difference\_\_\_
- (b) Comfort - Right\_\_\_ Left\_\_\_ No Difference\_\_\_
- (c) Ease of Handling - Right \_\_\_ Left\_\_\_ No Difference\_\_\_"?

This and other information was utilized during the post-study refitting procedure. The researchers' judgement was always of paramount importance during this procedure, but in cases of equally fitting research lenses, the patients' preference would be considered.

The promise of a new pair of contact lenses upon completion of the program as well as the exceptional followup care was used as an incentive during the study. Upon dispensing of the new pair of contacts the patients were asked to read and sign a statement enforcing the need for routine care and follow-up examinations. A copy of this form (SAM FORM 8) and their contact lens specifications (SAM FORM 6) were given to the patients. A copy of these eight "SAM" forms are included as appendix C.

The two instruments of primary importance to this research were the Pachometer and the Endothelial Camera. The purpose of the Pachometer is to provide the mechanical to electronic digital conversion, data analysis, operator interaction and the output functions necessary in order to accurately measure apparent corneal thickness. The total measurement system consists of a slitlamp, a fixation device, an optical doubling device (with eyepiece), a potentiometer, and an Apple IIe computer.

Figure 4. Photograph of Pachometer



The fixation device allows fixation of the patient to be changed at the discretion of the operator in eight (8) meridians and seventeen (17) locations. The fixation device contains nineteen (19) L.E.D.s (Light Emitting Diodes). There are seventeen (17) red L.E.D.s arranged in eight (8) meridians and two (2) yellow L.E.D.s that allow for alignment in the vertical plane to insure that readings are taken

perpendicular to the cornea. Any one of the 17 red L.E.D.s may be illuminated at random from the Apple IIe keyboard. The two yellow L.E.D.s remain illuminated at all times to insure proper alignment.

The potentiometer assembly attaches to the slitlamp by slipping it on the mounting post used for a Goldman Applanation Tonometer. The potentiometer is slid over the large pin and is locked into place with the small short pin. The potentiometer is a device for measuring electromotive force or potential difference by comparison with a known voltage.

The bi-prism eyepiece is inserted into the right ocular of the slitlamp with the eyepiece slit in a horizontal orientation. This allows for vernier alignment of the upper and lower slits.

The Apple IIe computer along with the "Contact Lens Research Software" performs all of the analysis and hardware control functions. A copy of the "Contact Lens Research Software" program is attached as appendix D. The computer receives all of the input from the potentiometer, performs the analysis, controls the fixation lights, and provides an interface to the Apple Imagewriter Printer. The computer program computes the mean, standard deviation, and range of the current set of thickness samples. These items are recomputed after each measurement or after deleting a measurement. The computer and program also calibrates each measurement to compensate for operator biases in the

measurement process. All pachometer readings were accomplished by this researcher in order to control the endpoint criterion, touch versus overlap, as much as possible. The ambient room illumination was the same at each reading and the slit beam adjustment knob was glued in place in order to maintain the same width throughout the study.

Pachometry was accomplished at five locations on both corneas at each of the eleven visits. Central, superior, inferior, nasal and temporal reading were performed. The superior, inferior, nasal and temporal readings were taken at points 3.55mm from the center of the cornea. The central reading is represented by the point "Q" in the patients data files (see Appendix F). Point "I" is the superior position, and point "M" represents the inferior position. Point "K" is the nasal reading on the right eye and the temporal reading for the left eye. Point "O" is the temporal reading on the right eye and the nasal reading for the left eye. The "Pachometer Instruction Manual" compiled by this researcher is included in appendix E.

The other instrument of primary importance to this research was the endothelial camera. The Nikon Non-Contact Endothelial Camera was utilized through out the research period. The camera body used was the Nikon F3 with the MD4 motor drive. The Nikon DW-4 6X viewfinder was utilized for the critical high magnification photography. The objective lens used for all photography was the Nikon 5X in conjunction with a Nikon TC-201 2x Teleconverter. The F-number of the

objective lens was 2.6 for all exposures. Total observation magnification of the system was approximately 60X. Ektachrome 100 ASA slide film was used for all photographs. The unit is powered by the Nikon model T1 power supply. The illuminating light (observation light) was maintained at the medium setting. The flash (strobe output) was maintained at the high or number 5 position.


Figure 5. Photograph of Endothelial Camera





An angle of  $60^\circ$  was maintained between the objective arm and the illumination source arm, being approximately equally split on either side of the midline. Figure 6 illustrates the arrangement of the camera system.

Figure 6. Illustration of the Objective and Illumination Source arrangement of the Endothelial Camera




In order to obtain photographs of the same area of the cornea, a system was devised based upon one reported by Barr and Schoessler.<sup>41</sup> A sighting telescope and the camera assembly base were mounted in fixed positions on the instrument stand. The patient was positioned in the same manner at each examination so that the crosshairs of the telescope coincided with the pupil of the left eye while observing a fixed fixation target. Utilizing a protractor that was attached to the head rest of the instrument, the fixation target was set in the same position at each setting

as recorded for that particular subject.

All photographs were taken by this researcher to eliminate as much operator variability as possible. A minimum of three exposures were taken of each eye on every visit. The slides were reviewed after development to ascertain the highest quality photograph from that series. The slides were then printed at the largest possible enlargement on photographic paper of 2 1/2 inch by 10 inch strips. The resulting print gave an overall magnification factor of 400X. This was determined by photographing a high resolution photographically reduced scale from a measuring reticle marked in tenths of a millimeter. Figure 7 is an example of the photographs employed during this study.

Figure 7. Photograph of Endothelial Cells.



In order to determine the cell density from the photograph, a template was used to inscribe two rectangles and two squares of equal area on the print. The cells were counted in each of the four areas and then averaged. Partial cells were counted that touched the left and top edges and not counted if they touched the bottom and right edges. The cell average was then multiplied by a factor to give the number of cells per square millimeter.

## TREATMENT OF DATA

This section will be divided into two separate divisions, one dealing with the data collected from pachometry and the other dealing with the endothelium photography data.

### PACHOMETRY

Microsoft File, a database manager used with the Macintosh™ Computer, was used to store all of the thickness measurements that were collected during the course of the research. The completed pachometry subject data files are contained in Appendix F. This information was then transferred into Microsoft Excel, which is a fully integrated program consisting of a spreadsheet, database and charting package. A spreadsheet showing all patient data regarding pachometry is located in Appendix G. This spreadsheet also shows the percentage change for each thickness measurement from day one to day 150. The Excel program was used to manipulate the data so that it could be loaded directly into a StatView™ statistical program.

There were eleven visits for each of the patients and therefore eleven variables within each group of patients. The analysis of variance and the corresponding test of significance based upon the "F" distribution allows testing differences among all of the means at the same time. Edwards describes the test as follows: "The analysis of variance, as the name indicates, deals with variances rather than with standard deviations and standard errors."<sup>54</sup> The rationale of

the analysis of variance is that the total sum of squares of a set of measurements composed of several groups can be analyzed or broken down into specific parts, each part identifiable with a given source of variation. In the simplest case, the total sum of squares is analyzed into two parts: a sum of squares based upon variation within the several groups, and a sum of squares based upon the variation between the group means. Then, from these two sums of squares, independent estimates of the population variance are computed".

Analysis of Variance by lens type and by position of measurement (central, superior, etc.) are located in Appendix H. There were twelve (12) eyes represented in the Permaflex™ group, fourteen (14) in the Syntex CSI®T group, and twenty-six (26) in the Hydrocurve II® group. Also included is the study of the control group of eight (8) eyes.

All of the Analysis of Variance tests proved statistically insignificant. Results of the one way analysis of variance and the corresponding test of significance based on the "F" distribution demonstrated all five points tested to be below 1.174 with  $p > .25$ . The mean for the central measurement was .525 with a standard deviation of .039. There was a range of .152 with the maximum being .613 and the minimum being .461. This, therefore, supports the statistical null hypothesis and disproves the working hypothesis.

Located in Appendix I are statistical data that were formulated from the five different points measured during this study. These pages show the various statistical data on 52 eyes that finished the study. The data represents the findings that were taken prior to dispensing of contact lenses.

#### ENDOTHELIUM

Counting of the endothelial cells was accomplished as described in the Methodology Section. Results of the counting at day zero and day 150 were recorded on an Excel Worksheet along with the results of the polymegathism grading. The Polymegathism Grading System was utilized as described in the February 1986, Contactology.<sup>55</sup> A copy of this worksheet for the 52 subjects as well as the control patients is located in Appendix J. These data were transferred from Excel into the Stat View Statistical Program for analysis. No statistical significance was found between beginning and final cell count. Results of the one way analysis of variance and the corresponding test of significance based on the "F" distribution was .079 with  $p > .25$ . The mean was 2900 cells/mm<sup>2</sup> with a standard deviation of 244.147. There was a range of 900 with a maximum of 3300 cells/mm<sup>2</sup> and a minimum of 2400 cells/mm<sup>2</sup>. The statistical data are located in Appendix J.

## RESULTS

This section on results will be separated into three different sections subtitled: General Comments, Pachometry and Endothelium.

## GENERAL COMMENTS

Twenty-six (26) of the thirty-nine (39) patients who started the program completed the full 150 days of the study. The following table describes what happened to the thirteen patients who did not complete the study.

Table 2. Disposition of 13 patients not completing the full 150 days of the extended wear contact lens study.

Patient Number	Visit Completed	Reason
101	28	Joined Army
104	120	EKC like symptoms interrupted wear
106	60	EKC like symptoms interrupted wear
110	90	Moved from area
114	28	Pregnancy
116	00	Inability to adapt to extended wear
117	90	EKC like symptoms interrupted wear
124	60	EKC like symptoms interrupted wear
125	60	EKC like symptoms interrupted wear
127	00	Pregnancy
128	00	Unable to meet study requirements
131	07	Inability to adapt to extended wear
138	60	EKC like symptoms interrupted wear

The six patients with Epidemic Keratoconjunctivitis (EKC) like symptoms were able to continue wearing contact lenses after a wearing break. These six cases were not typical EKC, but did have several similar symptoms. Mild epithelial punctate keratitis with small corneal infiltrations were

present in most cases. Preauricular adenopathy was present only in a few of the cases. Symptoms such as photophobia, lacrimation, discharge, or swelling were rarely described by any of the patients effected. Symptoms affected five female and one male patient. The six cases involved five eyes that were wearing Permaflex™ and one that was wearing CSI®T. All cases cleared in approximately three to four weeks.

Thirty-four (43.5%) of the seventy-eight contact lenses that were originally ordered for the patients completed the entire study (150 days). The following table gives the breakdown on the twenty-eight lenses that had to be replaced during the study.

Table 3. Reasons that contact lenses were replaced during the extended wear study.

Number replaced	Explanation
5	Lost lens
2	Fitting Change
2	Power Change
8	Crack or Split in Lens
8	Torn Lens
2	Chip in Edge
1	Growth on Lens

All three lens types gave similar results as to replacement and to percentage finishing the program. At the termination of the research program, thirty-five patients were refit with new contact lenses: eighteen received



Hydrocure II<sup>®</sup>, eleven received Permaflex<sup>™</sup> and six received Syntex CSI<sup>®</sup>T.

The following table shows the response to question six from the patient questionnaire that was taken at the termination of the program.

Table 4. Patient response to questionnaire comparing Hydrocure II (HC), Permaflex (PF), and CSI-T (CS) contact lenses in three categories: Vision, Comfort and Handling. Numbers indicated show preferred lens. ND indicates patient could not determine a difference in the lenses.

	VISION	COMFORT	HANDLING
<u>GROUP 1</u>	HC - 2	HC - 3	HC - 4
Hydrocure &	PF - 4	PF - 6	PF - 2
Permaflex	ND - 13	ND - 10	ND - 13
<u>GROUP 2</u>	HC - 2	HC - 1	HC - 6
Hydrocure &	CS - 5	CS - 4	CS - 2
Syntex CSI-T	ND - 7	ND - 9	ND - 6

Non-Contact Tonometry was completed on each visit. Table 5 is a compilation of the results of the tonometry readings. Previous research has shown that the Non-Contact Tonometer has a  $\pm 3$ mm Hg range at the 95% accuracy level.<sup>56,57</sup> No significant clinical difference was noted on any of the patients during the course of the study. Each entry on table 5 represents the measurements that were taken on the various visits.

Table 5. Results of Non-Contact Tonometer. The right eye measurement is recorded first followed by a slash (/) and the left eye measurement. The triple X (XXX) notation indicates that no measurement was accomplished on that occasion.

PT. #	DAY 0	DAY 1	DAY 4	DAY 7	14	21	28	60	90	120	150
101	17/17	14/13	16/12	16/13	18/18	17/15	15/14	XXX	XXX	XXX	XXX
102	11/14	11/11	09/15	10/10	12/11	09/10	11/15	10/13	11/10	11/08	09/09
103	09/10	11/11	XXX	12/11	10/09	10/11	10/08	12/10	14/12	15/12	11/11
104	14/15	12/10	15/14	15/16	14/12	12/14	11/09	10/10	13/11	16/12	XXX
105	09/13	07/09	07/10	09/09	07/09	08/09	10/11	08/11	07/09	10/12	08/10
106	16/14	12/13	18/14	14/12	11/12	14/14	15/11	10/13	XXX	XXX	XXX
107	20/20	16/16	19/18	19/20	18/20	20/20	16/19	16/17	19/20	06/18	19/19
108	12/13	11/11	11/10	10/10	11/12	11/13	14/13	11/14	13/10	10/08	11/09
109	11/15	12/14	12/15	13/14	11/14	12/16	13/15	10/12	14/16	12/12	12/13
110	12/13	13/12	15/15	16/15	13/11	13/13	15/12	12/12	16/16	XXX	XXX
111	14/13	09/09	09/10	11/09	12/13	09/11	10/11	08/09	11/12	10/09	10/08
112	14/14	17/15	15/15	14/17	12/14	17/19	17/19	14/15	12/15	12/14	12/15
113	09/09	10/10	11/11	14/11	11/09	12/11	11/11	10/10	10/10	13/11	11/11
114	15/17	15/14	11/12	14/15	13/13	11/13	13/11	XXX	XXX	XXX	XXX
115	10/13	10/12	09/11	08/10	10/12	10/12	09/12	11/12	08/10	08/11	11/12
116	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX
117	11/12	13/13	13/10	09/09	14/10	12/11	13/09	12/09	14/10	XXX	XXX
118	13/11	09/11	11/08	11/13	13/15	11/15	12/12	14/12	13/13	13/14	11/10
119	14/15	15/11	13/12	14/13	12/12	17/13	13/15	13/13	16/16	13/13	13/12
120	11/12	12/13	12/12	14/14	13/14	15/14	15/13	15/12	12/12	12/11	17/15
121	12/09	13/10	11/14	12/12	09/08	13/10	11/10	10/11	09/10	12/09	09/10
122	17/16	14/15	13/13	14/15	16/19	11/14	12/15	15/16	14/16	15/13	18/16
123	09/11	11/11	11/10	12/11	14/13	09/11	13/12	12/08	12/13	12/09	09/08
124	20/20	22/20	18/17	18/20	17/18	19/13	17/16	20/18	XXX	XXX	XXX
125	18/19	15/12	14/13	10/14	15/19	14/12	12/12	09/11	XXX	XXX	XXX
126	12/12	12/11	11/11	12/13	11/11	11/11	13/13	11/11	12/12	12/09	13/13
127	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX
128	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX
129	14/13	12/10	12/15	10/13	13/12	10/09	11/11	09/10	11/09	12/14	11/12
130	09/09	08/11	11/09	09/08	08/10	09/10	08/10	08/11	10/11	08/08	06/09
131	17/16	14/15	16/15	16/17	XXX	XXX	XXX	XXX	XXX	XXX	XXX
132	14/15	15/15	16/13	14/15	15/13	15/12	14/12	14/13	12/16	12/12	15/13
133	10/12	11/10	09/11	08/10	08/08	09/09	08/10	08/09	11/08	09/09	10/10
134	19/15	14/13	16/12	16/13	12/12	13/11	14/12	14/13	17/13	15/14	15/12
135	11/09	11/10	11/10	12/11	11/09	11/10	12/12	09/11	10/11	11/10	09/09
136	13/12	09/10	12/12	13/13	11/13	12/15	14/14	14/12	10/12	15/12	12/13
137	11/09	09/10	10/09	11/12	09/08	11/11	12/08	11/08	10/07	09/10	09/09
138	15/15	12/14	11/14	11/12	12/12	11/10	14/17	13/14	XXX	XXX	XXX
139	10/13	10/12	13/12	12/12	10/11	10/09	10/11	10/10	11/12	10/12	11/09

Results of the visual performance test, Keystone Aviators Unit, is shown in Table 6.

Table 6. Results of the Keystone Diagnostic Series Aviators Unit. Test was administered prior to and after contact lens wear. Results shown are for best spectacle Rx and best contact lens Rx at the end of the research period.

	IMPROVED	REDUCED	NO CHANGE
BEST SPECTACLE Rx	21	3	9
Best Contact Lens	15	5	13


The full results are located in chart form in Appendix K.

#### PACHOMETRY


Relative to the first hypothesis, data analysis performed with the analysis of variance and the corresponding test of significance are located in Appendix H. The statistical analysis of corneal thickness change, at all five locations, measured over a period of 150 days was found to be insignificant. The following graphs (4,5,6,7) illustrate the percentage of change that the central corneal thickness experienced with the three types of contact lenses.

Analysis performed on the five locations measured on the control patients also proved to be insignificant with the maximum "F" test being .683 with  $p > .25$ .


Graph 4. Percentage of central corneal thickness change with Syntex CSI-T Contact Lenses over a period of 150 days.



Graph 5. Percentage of central corneal thickness change with Hydrocurve II Contact Lenses over a period of 150 days.




Graph 6. Percentage of central corneal thickness change with Permafex Contact Lenses over a period of 150 days.



Graph 7 illustrates simultaneously the percentage of change that the central corneal thickness experienced with the three types of contact lenses. The thickness increase experienced at day sixty by all three lens types cannot be explained by the contextual contents of this research design.

Graph 7. Percentage of central corneal thickness change with Hydrocurve II (HC), Permaflex (PF), and Snytex CSI-T (CS) contact lenses over a period of 150 days.



## ENDOTHELIUM

Relative to the second hypothesis, data analysis performed with the analysis of variance and the corresponding test of significance are located in Appendix J. The statistical analysis of the change in endothelial cell count between day zero and day 150 was not significant (.079 with  $p > .25$ ). The mean number of cells found was 2900/mm<sup>2</sup>. While there was not a significant change noted in cell density, there was an increase in polymegathism. 65.4% of the patients showed a small change in the amount of polymegathism demonstrated. No significant change in the amount of polymegathism could be found in the remaining 34.6% of the subjects. The following table shows the breakdown of the grading of the polymegathism, or lack of, in the patient population.

Table 7. Representation of the amount of polymegathism demonstrated by subjects in the extended wear research study. 52 eyes are represented.

NORMAL (NO POLYMEGATHISM)	6
GREATER THAN NORMAL - LESS THAN MODERATE	18
MODERATE	24
GREATER THAN MODERATE - LESS THAN SEVERE	4

No change was noted in either cell density or cell morphology among any of the control patients.



## DISCUSSION AND CONCLUSION

The objective of this research project was to investigate the correlation between extended wear contact lenses and their effect on corneal thickness and the endothelial layer of the cornea.

The research questions were as follows: (1) Will changes in corneal thickness be proportional to the differences in contact lens water content during extended wear? (2) What effect will different (water content and oxygen transmissibility) contact lenses make on the morphological changes of the endothelial layer.

To examine these questions, measurements of corneal thickness at five locations on the cornea were taken before the wearing of contact lenses and on ten occasions, over a five month period, during extended wear. Endothelial photographs were also taken at all of the aforementioned occasions.

Based on the results of the analysis of variance and the corresponding test of significance, using the "F" distribution, the changes in corneal thickness throughout the study were not significant. None of the corneal thickness changes, at any of the five locations measured, proved to be statistically significant with any of the three different contact lenses.

The photographs of the corneal endothelium taken over the five month research period revealed no significant change in the endothelial cell density. 65.4 percent of the patients

demonstrated morphological changes of their endothelium. This was exhibited in the form of small increases in the amount of polymegathism.

The majority of the subjects did experience the transient appearance of endothelial blebs. One subject experienced significant temporary swelling of the cornea. These changes were easily observed with the slitlamp as folds in Descemet's Membrane prior to measurement by pachometry. This patient was dropped from the extended wear program, but continues to wear hydrophilic daily wear contact lenses in a successful manner.

None of the control patients exhibited either a change in cell density or cell morphology. Neither did any of the control patients exhibit a significant change in corneal thickness in any of the measured locations.

Only nine (23%) of the original 39 subjects were released from the research program wearing their lenses in an extended wear mode. These patients were advised to remove their lenses for overnight cleaning twice a week (approximately every three to four days). They were also told of the importance of continued followup examinations.

The remaining subjects were advised to wear their lenses on a daily wear basis. Continual engorgement of the corneoscleral limbal plexus was the primary professional concern in advising these subjects against further extended wear. Clearly, one of the most significant findings of this study was the fact that only nine of the original 39 patients

were recommended for continued extended wear. This is exceedingly important from the standpoint that these patients were carefully selected from approximately 100 candidates for the program and that the compliance rendered by these patients was of an exceptional nature.

The conclusions are as follows:

(1) Corneal thickness was not significantly affected by the wearing of high, medium, or low water content contact lenses on an extended wear basis.

(2) Corneal endothelial cell density is not affected by the wearing of high, medium, or low water content contact lenses on an extended wear basis.

(3) Corneal endothelial morphological changes (increased polymegathism) does occur with wear of high, medium, and low water content contact lenses on an extended wear basis.

(4) Determination of successful wear of any extended wear contact lenses can only be done on an individual basis. Extensive followup examinations done on a timely basis, for the duration of the period that the patient wears extended wear contact lenses, is of utmost importance.

## REFERENCES

1. Heltzel EE Users Blamed For Problems With Contact Lenses. The Oregonian 1985 Dec 15:13(col 1).
2. Nazario SL Some Eye Doctors Begin to Doubt Safety of Extended-Wear Contacts. The Wall Street Journal 1985 Nov 8: 35(col 5).
3. Stuller J What's New In Contact Lenses? Reader's Digest 1985 July:109-112.
4. Anonymous. Extended-Wear Lens Safety Issue Explored At Symposium. American Optometric Association News 1985 Dec 15 24:1 and 11.
5. Mandell RB Contact Lens Practice (3rd Ed.). Illinois Charles C. Thomas 1981.
6. Miller D, White PF Complications of Contact Lenses International Ophthalmology Clinics. Boston: Little, Brown and Co.,1981.
7. Anonymous. Oklahoma Optometric Association Newsletter 1985 Jan (col 1).
8. Houde WL. Rubin ML. Extended-Wear Lenses. An Update. Survey Of Ophthalmology 1981;26:103-105.
9. Pavlichko JL Times Change For Extended Wear. Review of Optometry 1984;(Apr):99-100.
10. White P, Scott C Contact Lenses and Solutions Summary. Manual used in Continuing Education Seminar Pacific University College of Optometry 1984;(Oct):10-11.
11. Cordrey P The Eye in Extended Contact Lens Wear. American Journal of Optometry & Physiological Optics 1980; 57:113-117.
12. Hirji NK, Larke JR Corneal Thickness In Extended Wear Of Soft Contact Lenses. British Journal of Ophthalmology 1979;63:274-276.
13. Holden BA Mertz GW, McNally JJ Corneal Swelling Response to Contact Lenses Worn Under Extended Wear Conditions. Investigative Ophthalmology and Visual Science, 1983;24:218-226.
14. Miller JP, Coon LJ, Meier RF Extended Wear Of Hydrocurve II<sub>55</sub> Soft Contact Lenses. Journal of the American Optometric Association 1980;51:225-233.

15. Millodot M Clinical Evaluation of an Extended Wear Lens. International Contact Lens Clinic 1984;11:16-23.
16. Snyder AC, Schoessler JP Corneal Thickness Changes Associated With Daily and Extended Contact Lens Wear. American Journal of Optometry & Physiological Optics 1983; 60:830-838.
17. Weissman BA Danger: EXW. Optometric Monthly 1983; (Jan):21-22.
18. Binder PS Myopic Extended Wear with the Hydrocurve II Soft Contact Lens. Ophthalmology 1983;90:623-626.
19. Rengstorff RH, Nilsson KT Long-Term Effects of Extended Wear Lenses: Changes in Refraction, Corneal Curvature, and Visual Acuity. American Journal of Optometry & Physiological Optics 1985;62:66-68.
20. Sorensen T, Corydon, L Changes in Central Corneal Thickness During Permanent Lens Wear. Contacto 1979;23: 28-30.
21. Larke JR The Eye In Contact Lens Wear. Boston: Butterworth, 1985.
22. Von Bahr G Corneal Thickness. American Journal of Ophthalmology 1956;42:251-263.
23. Donaldson DD A New Instrument For The Measurement Of Corneal Thickness. Archives of Ophthalmology 1966;76:25-31.
24. Borish IM Clinical Refraction(3rd Ed.). Illinois: Professional Press, 1975.
25. Molinari JF, Bonds T Pachometry: A Comparison Between Touch and Overlap Measurement Method. American Journal of Optometry & Physiological Optics 1983;60:61-66.
26. Mishima S, Hedbys BO Measurement of Corneal Thickness With the Haag-Streit Pachometer. Archives of Ophthalmology, 1968;80:710-713.
27. Snyder AC Optical Pachometry Measurements: Reliability and Variability. American Journal of Optometry & Physiological Optics 1984;61:408-413.
28. Newell FW Ophthalmology Principles and Concepts (5th ed). St. Louis The C.V. Mosby Company, 1982.
29. Mishima S Clinical Investigations on the Corneal Endothelium. American Academy of Ophthalmology 1982;89: 525-530.

30. Hodson S, Miller F The Bicarbonate Ion Pump in the Endothelium Which Regulates the Hydration of Rabbit Cornea. *Journal of Physiology* 1976;263:563-577.
31. Barr JT Corneal Endothelial Observation Techniques. *American Journal of Optometry & Physiological Optics* 1981; 58:618-620.
32. Bourne WM, Kaufman HE Specular Microscopy of Human Corneal Endothelium In Vivo. *American Journal of Ophthalmology* 1976;81:319-323.
33. Laule A, Cable MK, Haffman CE, et al Endothelial Cell Population Changes of Human Cornea During Life. *Archives of Ophthalmology* 1973;96:2031-2035.
34. Barr JT, Schoessler JP Measurement of Endothelial Cell Density and Quantification of Endothelial Mosaic Change. *International Contact Lens Clinic* 1979;(Mar-Apr):36-39.
35. Capella JA The Pathology of Corneal Endothelium. *Ann. Ophthalmology* 1971;3:397-400.
36. Laing RA, Sandstrom MM, Berrospi AR, et al Changes in the Corneal Endothelium as a Function of Age. *Experimental Eye Research* 1976;22:587-594.
37. Olsen EG, Davanger M The Healing of Human Corneal Endothelium. *ACTA Ophthalmologica* 1984;62:885-892.
38. Brown N Macrophotography of the Anterior Segment of the Eye. *British Journal of Ophthalmology* 1970;54:697-701.
39. Zantos SG, Holden BA Transient Endothelial Changes Soon After Wearing Soft Contact Lenses. *American Journal of Optometry & Physiological Optics* 1977;54:856-858.
40. McMonnies CW, Zantos SG Endothelial Bedewing of the Cornea in Association With Contact Lens Wear. *British Journal of Ophthalmology* 1979;63:478-481.
41. Barr JT, Schoessler JP Corneal Endothelial Response to Rigid Contact Lenses. *American Journal of Optometry & Physiological Optics* 1980;57:267-274.
42. Waring GO, Krohn MA, Ford GE, et al Four Methods of Measuring Human Corneal Endothelial Cells From Specular Photomicrographs. *Archives of Ophthalmology* 1980;98: 848-855.
43. Anonymous. Nikon Non-Contact Endothelial Camera Instructions. Nippon Kogaku K.K., Tokyo, Japan.

44. Antti V, Jukka M, Jukka S, et al Contact Lens Induced Transient Changes in Corneal Endothelium. ACTA Ophthalmologica 1981;59:552-559.
45. Kamiya C A Study of Corneal Endothelial Response To Contact Lenses. Contact Lens 1982;8:92-95.
46. Schoessler JP, Woloschak MJ, Mauger TF Transient Endothelial Changes Produced By Hydrophilic Contact Lenses. American Journal of Optometry & Physiological Optics 1982; 59:764-5.
47. Schoessler JP, Woloschak MJ Corneal Endothelium in Veteran PMMA Contact Lens Wearers. International Contact Lens Clinic 1981;(Nov-Dec):19-25.
48. Caldwell DR, Kastl PR, Dabezies OH, et al The Effect of Long-Term Hard Lens Wear On Corneal Endothelium. Contact Lens. 1982;8:87-91.
49. Stocker EG, Schoessler JP Corneal Endothelial Polymegathism Induced by PMMA Contact Lens Wear. Investigative Ophthalmology & Visual Science 1985;26: 857-863.
50. Kiely PM, Carney LG, Smith G Menstrual Cycle Variations of Corneal Topography and Thickness. American Journal of Optometry & Physiological Optics 1983;60:822-829.
51. El Hage SG, Beaulne C Changes In Central And Peripheral Corneal Thickness With Menstrual Cycle. American Journal Of Optometry and Archives Of American Academy Of Optometry 1973;(Nov):863-871.
52. Feldman F, Bain J, Matuk A.R Daily Assessment of Ocular and Hormonal Variables Throughout the Menstrual Cycle. Archives of Ophthalmology 1978;96:1835-1836.
53. Efron N, Kotow M, Martin DK, et al Physiological Response of the Contralateral Cornea to Monocular Hydrogel Contact Lens Wear. American Journal of Optometry & Physiological Optics 1984;61:517-522.
54. Edwards AL Statistical Methods For The Behavioral Sciences. New York: Holt, Rhinehart, and Winston, 1964.
55. Anonymous. Polymegathism May Indicate Contact Lens Ill Effect. Contactology International 1986;2:12.
56. Grolman B A New Tonometer System. American Journal of Optometry and Archives of American Academy of Optometry 1972;49:646-660.

57. Myers KJ, Scott CA The Non-Contact ("Air Puff")  
Tonometer: Variability and Corneal Staining. American  
Journal of Optometry and Physiological Optics 1975;52:  
36-46.



APPENDIX A  
CONTACT LENS INSTRUCTIONS

CARING FOR YOUR EXTENDED WEAR CONTACT LENSES

We want you to follow the steps listed below, always starting with the right lens first:

1. Set up all of the items you will need to clean and disinfect your extended wear lenses (Pliagel, extenzyme cleaner, Septicon disinfection system, and Blairex deionizer).
2. Wash your hands with a pure soap that does not contain additives. Rinse thoroughly.
3. Remove your lenses and clean them with Pliagel (see directions supplies with the package). This removes inorganic material from the lens surface.
4. Rinse your lenses off with non-thimerosal saline (Blairex system for making normal saline) to completely remove cleaner.
5. Place the lenses in the Septicon lens baskets and soak in the extenzyme cleaner for 15-20 minutes using the Lensept cup #1. This removes organic material (protein complex from your tears etc.) from the lens.
6. Replace the extenzyme cleaner with fresh saline and shake for 15 seconds.
7. Replace saline with Lensept solution and soak lenses for 10 minutes.
8. Fill Rinse cup #2 with saline solution. Remove the lens basket from cup #1 and place into cup #2. The lenses must be left in cup #2 for 6 hours to neutralize the Lensept disinfection solution.
9. After the 6 hour neutralizing cycle, remove the Septicon lens baskets from the Rinse cup #2. Replace the solution with fresh saline solution and allow the lenses to soak in it for another 3 minutes before placing the lenses back on the eyes.

FOLLOW-UP CARE FOR EXTENDED WEAR CONTACT LENS PATIENTS

Patients wearing contact lenses for extended periods of time require additional professional follow-up care to more closely monitor how their lenses are fitting. For this reason and for the purposes of this research project, the following schedule of follow-up evaluations will be followed:

- \*After 24 hours of extended wear,
- \*After 4 days of extended wear,
- \*After 1 week of extended wear and then weekly for three weeks,
- \*After each month of extended wear for the remainder of the 5 months of the study.

Your eyes should always "FEEL GOOD, LOOK GOOD AND SEE GOOD." If your eyes bother you at anytime, or if you have any questions or concerns, we want to be the first to know. You can always reach one of our research staff 24 hours a day, 7 days a week. Office 357-6151 ext. 2453, Dr. Tim Allen 357-3953, Dr. David Marrs 357-4026, Dr. Gary Slater 357-8378.

## Extended Wear Lens Care and Handling

Proper care is necessary for successful wear, normal lens life, and good eye health. You will be provided with products to clean, disinfect and store your extended wear lenses. Use them as instructed.

Your lens cleaner is \_\_\_\_\_

Your lens disinfectant is \_\_\_\_\_

Your soaking solution is \_\_\_\_\_

Your rinsing solution is \_\_\_\_\_

Eyedrops to use before sleep and upon waking \_\_\_\_\_

**NOTE:** These products have been prescribed specifically for your lenses and eyes. Do not change or substitute brands unless you check with us first. Use of improper solutions may result in lens damage or eye irritation.

**SPECIAL INSTRUCTIONS:** \_\_\_\_\_

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

Please note that although certain brands of lenses may be FDA approved for 7, 14, or even 30 days of wear, the adaptability of your eyes is the key factor in determining wearing time. Trust us to recommend a schedule suited to your individual needs. And remember, like any medical device contact lenses must be monitored on a regular basis. Professional follow-up care is the most important element in successful long term lens wear.

### IN THE BEGINNING IT IS NORMAL IF:

- 1) Your lenses itch or feel funny.
- 2) One lens is more noticeable than the other.
- 3) Your vision seems fuzzier than with glasses.
- 4) One eye sees better than the other.
- 5) You have trouble handling your lenses.

### REMOVE AND DO NOT SLEEP IN YOUR LENSES IF:

- 1) You develop unusual pain or redness.
- 2) You develop unusually cloudy or foggy vision.
- 3) You experience a decrease in vision that does not clear up.
- 4) You suspect something is wrong.

### WEARING SCHEDULE

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

Next appointment: \_\_\_\_\_

Date \_\_\_\_\_

Dispenser \_\_\_\_\_

Patient \_\_\_\_\_

APPENDIX B  
INFORMED CONSENT FORM

## INFORMED CONSENT FORM

1. INSTITUTION

- A. TITLE OF PROJECT:
- a. The Effect of Extended Wear Contact Lenses on Corneal Thickness and the Endothelial Layer.
  - b. Changes in Oxygen Transmissibility (DK/L) in Extended Wear Contact Lenses with Time.
  - c. The Effect of Extended Wear Contact Lenses on Contrast Sensitivity Function with Time.
- B. PRINCIPAL INVESTIGATORS:
- a. David Marrs 357-4026
  - b. Tim Allen 357-3953
  - c. Gary Slater 357-8378
- C. ADVISOR: Dr. Don West 357-9036
- D. LOCATION: Pacific University College of Optometry  
Forest Grove, Oregon
- E. DATE: 1985 - 1986

2. DESCRIPTION OF PROJECT

This project is designed to observe, monitor and record (1) the effects of extended wear contact lenses on corneal thickness and the endothelial layer. (2) The changes in oxygen transmissibility (DK/L) of extended wear contact lenses with time. (3) The effect of extended wear contact lenses on contrast sensitivity function with time. This will be accomplished by fitting extended wear contact lenses on human subjects and monitoring and measuring the above items.

3. DESCRIPTION OF RISKS

Participants will be at no greater risk than if fit with contact lenses in a normal clinical situation. As with all contact lens wearers, participants will experience a normal adaptation period during which symptoms may occur. Several clinical instruments will be used in close proximity to the eye, presenting a minor risk of contact injury. These instruments are routinely used by optometrists, and the investigators are extensively trained in their use.

4. DESCRIPTION OF BENEFITS

This study will serve to increase the basic understanding of the response of the eye to extended wear contact lenses. It will further allow study of the changes in oxygen transmission of the lenses as a function of the wearing time. The study will also determine the effect of extended wear contact lenses on visual performance with time. This study will help to determine the efficacy of the extended wear lens and be useful in establishment of the optimum fitting and wearing programs for the extended wear contact lens. All examination fees will be waived for all participants, and those completing the study will receive an additional pair of extended wear contact lenses. Since the study involves careful monitoring of changes, participants will receive care that is equal to or more intensive than that received by routine contact lens patients.

5. COMPENSATION AND MEDICAL CARE

If you are injured in this experiment it is possible that compensation or medical care will not be available from Pacific University, the experimenters, or any organization associated with the experiment. All reasonable care will be used to prevent injury.

6. ALTERNATIVES ADVANTAGEOUS TO SUBJECTS

Participants may find it advantageous to pursue fitting with an alternative brand or type of contact lens or spectacles. If it is determined that such an alternative would be advantageous, a referral will be made to the general clinic of Pacific University College of Optometry.

7. OFFER TO ANSWER ANY INQUIRES

The experimenter will be pleased to answer any questions that might arise at any time during the course of this study. If you are not satisfied with the answers you receive, please call Dr. James Peterson at 357-0442. During participation in this research project you are not considered a clinic patient. All questions should be directed to the researchers and/or the faculty advisor who will be solely responsible for any treatment, except in the case of an emergency.

8. FREEDOM TO WITHDRAW

You are free to withdraw your consent and to discontinue participation in this project or activity at any time without prejudice.

I have read and understand the above. I am 18 years of age or over (or this form is signed for me by my parent or guardian).

PRINTED NAME\_\_\_\_\_

SIGNED\_\_\_\_\_DATE\_\_\_\_\_

ADDRESS\_\_\_\_\_PHONE\_\_\_\_\_

CITY\_\_\_\_\_STATE/ZIP\_\_\_\_\_

NAME AND ADDRESS OF A PERSON NOT LIVING WITH YOU WHO WILL  
ALWAYS KNOW YOUR ADDRESS.\_\_\_\_\_

\_\_\_\_\_

## APPENDIX C

CONTACT LENS FORMS  
SAM (SLATER, ALLEN, MARRS) FORMS  
1 THROUGH 8



## CONTACT LENS FITTING FORM

NAME \_\_\_\_\_ I.D. NUMBER \_\_\_\_\_

DATE \_\_\_\_\_

SUBJECTIVE:

OD

05

K's OD

05

PUPIL SIZE \_\_\_\_\_

CORNEAL DIAMETER \_\_\_\_\_

PALPEBRAL FISSURE \_\_\_\_\_

TEAR BUT \_\_\_\_\_

## DIAGNOSTIC FITTING

TRIAL 1 LENS - OD \_\_\_\_\_ TRIAL 2 LENS - OD \_\_\_\_\_  
OS \_\_\_\_\_ OS \_\_\_\_\_

	COVERAGE	CENTER	MOVEMENT	OVER K	OVER RET	OVER REF
TRIAL 1 OD	_____	_____	_____	_____	_____	_____
OS	_____	_____	_____	_____	_____	_____
TRIAL 2 OD	_____	_____	_____	_____	_____	_____
OS	_____	_____	_____	_____	_____	_____

SLIT LAMP \_\_\_\_\_

### FINAL CONTACT LENS PRESCRIPTION

	BASE CURVE	DIAMETER	POWER	MANUFACTURE
OD	_____	_____	_____	_____
OS	_____	_____	_____	_____

PATIENT NAME \_\_\_\_\_ I.D. NUMBER \_\_\_\_\_

DATE \_\_\_\_\_ VISIT \_\_\_\_\_ DAYS

DISCUSSION \_\_\_\_\_

\_\_\_\_ CLASON VA - C.L. - OD 20/ OS 20/ OU 20/Spec. - OD 20/ OS 20/ OU 20/

OVER REFRACTION - OD \_\_\_\_\_ OS \_\_\_\_\_

\_\_\_\_ CSF - SEE ATTACHED

\_\_\_\_ AUTO REFRACTION (WITH LENS) - SEE ATTACHED NUMBER \_\_\_\_\_

\_\_\_\_ AUTO "K" (WITH LENS) - SEE ATTACHED NUMBER \_\_\_\_\_

\_\_\_\_ SLIT LAMP - Injection \_\_\_\_\_ Neo-Vascularization \_\_\_\_\_ Iritis \_\_\_\_\_ Edema \_\_\_\_\_ Staining \_\_\_\_\_

\_\_\_\_ OD - DK/L \_\_\_\_\_ L \_\_\_\_\_ DK \_\_\_\_\_ HUM \_\_\_\_\_ TEMP \_\_\_\_\_

OS - DK/L \_\_\_\_\_ L \_\_\_\_\_ DK \_\_\_\_\_ HUM \_\_\_\_\_ TEMP \_\_\_\_\_

\_\_\_\_ ENDO PHOTO - \_\_\_\_\_ PER SQ MM - FIXATION OD \_\_\_\_\_° OS \_\_\_\_\_°

\_\_\_\_ PACHOMETRY

\_\_\_\_ AUTO REFRACTION (WITHOUT LENS) - SEE ATTACHED NUMBER \_\_\_\_\_

\_\_\_\_ AUTO "K" (WITHOUT LENS) - SEE ATTACHED NUMBER \_\_\_\_\_

\_\_\_\_ CORNEASCOPE - PHOTO ATTACHED

\_\_\_\_ B&amp;L "K" - OD \_\_\_\_\_@\_\_\_\_\_, \_\_\_\_\_@\_\_\_\_\_

OS \_\_\_\_\_@\_\_\_\_\_, \_\_\_\_\_@\_\_\_\_\_

\_\_\_\_ IOP - NCT - OD \_\_\_\_\_ OS \_\_\_\_\_

NOTES - \_\_\_\_\_

NEXT APPOINTMENT \_\_\_\_\_

PERFORMANCE DEPTH PERCEPTION TEST  
KEYSTONE DIAGNOSTIC SERIES - AVIATORS UNIT

NAME\_\_\_\_\_ I.D. NUMBER\_\_\_\_\_

BASELINE (PRE-CONTACT LENS WEAR)

SLIDE #\_\_\_\_\_

PERCENT STEREOPSIS\_\_\_\_\_

SECONDS OF ARC\_\_\_\_\_

EXAMINER\_\_\_\_\_

DATE\_\_\_\_\_

FINAL (POST-CONTACT LENS WEAR)

SLIDE #\_\_\_\_\_

PERCENT STEREOPSIS\_\_\_\_\_

SECONDS OF ARC\_\_\_\_\_

EXAMINER\_\_\_\_\_

DATE\_\_\_\_\_

COMMENTS\_\_\_\_\_

PHOTO RECORD

61D

<u>PHOTO</u>	<u>PATIENT'S NAME</u>	<u>ID #</u>	<u>EYE</u>	<u>NOTE</u>	<u>FRAME</u>
*					
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					
16					
17					
18					
19					
20					
21					
22					
23					
24					
25					

ROLL NUMBER \_\_\_\_\_

DATE TO LAB \_\_\_\_\_

RECEIVED FROM LAB \_\_\_\_\_

SAM FORM 4, 5 AUG 85

-

PATIENT NAME \_\_\_\_\_ I.D. NUMBER \_\_\_\_\_

DATE \_\_\_\_\_ VISIT \_\_\_\_\_ DAYS

DISCUSSION \_\_\_\_\_  
\_\_\_\_\_\_\_\_ CLASON VA - Spec. - OD 20/ OS 20/ OU 20/

REFRACTION - OD \_\_\_\_\_ OS \_\_\_\_\_

\_\_\_ CSF - SEE ATTACHED

\_\_\_ ENDO PHOTO - \_\_\_\_\_ PER SQ MM - FIXATION OD \_\_\_\_\_° OS \_\_\_\_\_°

\_\_\_ PACHOMETRY

\_\_\_ AUTO REFRACTION (WITHOUT LENS)

\_\_\_ AUTO "K" (WITHOUT LENS)

\_\_\_ CORNEASCOPE - PHOTO ATTACHED

\_\_\_ B&amp;L "K" - OD \_\_\_\_\_ @ \_\_\_\_\_ , \_\_\_\_\_ @ \_\_\_\_\_

OS \_\_\_\_\_ @ \_\_\_\_\_ , \_\_\_\_\_ @ \_\_\_\_\_

\_\_\_ IOP - NCT - OD \_\_\_\_\_ OS \_\_\_\_\_

NOTES - \_\_\_\_\_

DEPTH PERCEPTION TEST \_\_\_\_\_

NEXT APPOINTMENT \_\_\_\_\_

**PACIFIC UNIVERSITY  
COLLEGE OF OPTOMETRY**

2043 COLLEGE WAY  
FOREST GROVE, OR 97116

**EXTENDED WEAR CONTACT LENS RESEARCH  
CONTACT LENS INFORMATION**

PATIENT \_\_\_\_\_

Original "K"      O.D. \_\_\_\_\_ Latest "K"      O.D. \_\_\_\_\_

Date \_\_\_\_\_ O.S. \_\_\_\_\_ Date \_\_\_\_\_ O.S. \_\_\_\_\_

Spectacle Rx      O.D. \_\_\_\_\_ Latest Refraction      O.D. \_\_\_\_\_

Date \_\_\_\_\_ O.S. \_\_\_\_\_ Date \_\_\_\_\_ O.S. \_\_\_\_\_

**CURRENT SOFT LENS SPECIFICATIONS**

Date Prescribed \_\_\_\_\_ Manufacturer \_\_\_\_\_

	Base	Power	Diameter	VA
O.D.				20/
O.S.				20/

**PRESCRIPTION IS VALID FOR SIX MONTHS FROM DATE PRESCRIBED**

## EXTENDED WEAR CONTACT LENS PATIENT QUESTIONNAIRE

1. Do you feel that you received adequate and professional care during the research period? YES\_\_\_\_ NO\_\_\_\_ COMMENT \_\_\_\_  
\_\_\_\_\_
2. Are there any changes you would suggest in future research projects of this nature? YES\_\_\_\_ NO\_\_\_\_ COMMENT\_\_\_\_\_  
\_\_\_\_\_
3. Do you plan to remain an extended wear contact lens wearer? YES\_\_\_\_ NO\_\_\_\_ COMMENT\_\_\_\_\_  
\_\_\_\_\_
4. Would you be interested in participating in another research project? YES\_\_\_\_ NO\_\_\_\_ COMMENT\_\_\_\_\_  
\_\_\_\_\_
5. Did you have any problems in following the instructions given during the research period? YES\_\_\_\_ NO\_\_\_\_ COMMENT\_\_\_\_\_  
\_\_\_\_\_
6. Did you have a preference between your right and left lenses during the study with ...
  - a. Vision - Right\_\_\_\_ Left\_\_\_\_ No Difference\_\_\_\_
  - b. Comfort - Right\_\_\_\_ Left\_\_\_\_ No Difference\_\_\_\_
  - c. Ease of Handling - Right\_\_\_\_ Left\_\_\_\_ No Difference\_\_\_\_Comments\_\_\_\_\_  
\_\_\_\_\_

**PACIFIC UNIVERSITY**  
**COLLEGE OF OPTOMETRY**  
2043 COLLEGE WAY  
FOREST GROVE, OR 97116

**EXTENDED WEAR CONTACT LENS RESEARCH**  
**CONTACT LENS INFORMATION**

PATIENT \_\_\_\_\_

1. Recommend you continue to follow the contact lens care procedures that were prescribed at the beginning of this study.
2. Recommend you have continuing follow-up examinations at 90 day (3 months) to 180 day (6 months) intervals.
3. Recommend you wear your contact lenses in the following manner:

\_\_\_\_\_

I have read and understand the above. I have received a new pair of contact lenses and the prescription for these lenses.

\_\_\_\_\_  
Signature

\_\_\_\_\_  
Date



APPENDIX D  
PACHOMETER COMPUTER PROGRAM

# JLIST

```

10  REM *****
11  REM CONTACT LENS RESEARCH HELLO PROGRAM
12  REM CREATED FOR PACIFIC UNIVERSITY
13  REM COPYRIGHT (C) 1985 BY KEVIN SPREIER
14  REM VERSION 1.1  09 SEP 85
15  REM *****
100 D$ = CHR$(4);G$ = CHR$(7); TEXT : HOME : PRINT D$"NOMONC,1,0": POKE
    - 17814,8: POKE - 17807,8
110  FOR I = 1 TO 39: PRINT "*";: NEXT : PRINT
120  PRINT "*      Contact Lens Research Software      *"
125  PRINT "* Copyright (c) 1985 by Kevin Spreier  *"
130  FOR I = 1 TO 39: PRINT "*";: NEXT : PRINT
140  PRINT : PRINT
150  HTAB 7: PRINT "Options:": PRINT
160  READ N
165  FOR I = 1 TO N
170  READ N$(I),ADDRESS(I)
175  HTAB 10: PRINT I") "N$(I)
180  NEXT I
200  HTAB 10: PRINT I") Exit to BASIC"
330  PRINT : HTAB 7: PRINT "Choose one -->";
340  GET A$:A = ASC (A$) - 48: IF A < 1 OR A > I THEN PRINT G$;: GOTO
    340
350  PRINT A$
360  IF A = 1 THEN HOME : END
380  VTAB 9 + A: INVERSE : HTAB 13: PRINT N$(A);: NORMAL : PRINT
385  POKE 47722,ADDRESS(A): POKE 47729,ADDRESS(A)
390  PRINT D$"RUN"N$(A): END
1000 DATA 3
1005 DATA "PERMEOMETER",8
1010 DATA "PACHOMETER",8
1015 DATA "PACH.REPORT",8

```

1222  
JLIST

```
10 REM *****
11 REM PACHOMETER CREATED FOR PACIFIC UNIVERSITY
12 REM COPYRIGHT (C) 1985 BY KEVIN SPREIER
13 REM VERSION 1.4 03 SEP 85
15 REM *****
18 X$ = "": REM Dummy variable for special INPUT$ routine
19 INPT = 940: REM Special input routine to allow commas
20 D$ = CHR$ (4):G$ = CHR$ (7): ONERR GOTO 9005
25 PRINT D$"BLOADINPUT$,A940"
28 CAL = 3000: REM CALIBRATION FACTOR (1 MM=3000 A/D POINTS)
29 BIAS = 810: REM .270 MM SYSTEM BIAS TO ADD TO ZERO POINT
30 EYE = 0: REM INITIAL EYE=0D
35 LED$ = "A":LED = 1:LD = 1: REM DEFAULT STARTING LIGHT
40 S1 = 4: REM SLOT OF LED INTERFACE CARD
41 POKE 49280 + S1 * 16,7: REM TURN ON "A" LED INITIALLY
45 S2 = 3: REM SLOT OF A/D CARD
46 AI13 = - 16256 + S2 * 16: REM A/D CARD I/O PORT
47 POKE AI13,18: REM DO A DUMMY A/D CONVERSION
50 DIM L(16): REM LED POSITION TRANSLATION TABLE
55 DATA 7,0,1,2,3,4,5,6,9,10,11,12,13,14,15,8
60 FOR I = 1 TO 16: READ L(I): NEXT : REM READ IN LED TABLE
61 DIM X%(17),Y%(17): REM LED SCREEN COORDINATES
62 DATA 20,24,25,24,20,16,15,16,20,22,22,22,20,18,18,18,20
63 DATA 9,10,12,14,15,14,12,10,11,11,12,13,13,13,12,11,12
64 FOR I = 1 TO 17: READ X%(I): NEXT : FOR I = 1 TO 17: READ Y%(I): NEXT

75 DIM DT%(1,17,12): REM POTENTIOMETER A/D RAW DATA ARRAY
80 DIM DT(1,17,3): REM CALC. DATA
85 DR = 2: REM DATA DISK DRIVE #
99 REM ***MAIN MENU
100 TEXT : HOME
105 FOR I = 1 TO 39: PRINT "*";: NEXT : PRINT
110 HTAB 11: PRINT "Pachometer Program"
115 HTAB 3: PRINT "Copyright (c) 1985 by Kevin Spreier"
116 FOR I = 1 TO 39: PRINT "*";: NEXT : PRINT
117 PRINT : PRINT
118 POKE 34,4: REM SET TOP OF SCROLL WINDOW
119 IF FIRST < > 0 THEN 130: REM SKIP CALIBRATION IF NOT FIRST TIME TH
ROUGH
120 FIRST = 1
121 PRINT "New zero mark (Y/N) N";
122 HTAB 20: INPUT A$: IF A$ = "" THEN A$ = "N"
124 IF A$ < > "Y" AND A$ < > "y" AND A$ < > "N" AND A$ < > "n" THEN
PRINT G$;: GOTO 122
125 VTAB PEEK (37): HTAB 21: PRINT A$
126 IF A$ = "N" OR A$ = "n" THEN PRINT D$"OPENPACHOMETER.ZERO": PRINT
D$"READPACHOMETER.ZERO": INPUT ZERO: PRINT D$"CLOSEPACHOMETER.ZERO":
GOTO 130
127 PRINT : INPUT "Zero Pachometer head and press <RETURN>:";A$
128 POKE AI13,18:ZERO = PEEK (AI13 + 1) * 256 + PEEK (AI13) + BIAS: REM
ADD BIAS TO NEW ZERO VALUE
129 PRINT D$"OPENPACHOMETER ZERO": PRINT D$"WRITEPACHOMETER.ZERO": PRINT
ZERO: PRINT D$"CLOSEPACHOMETER.ZERO"
130 PRINT "Subject name:?" :SUB$;
133 HTAB 15: CALL INPT:A$ = MID$ (X$,1): REM Input string with commas
134 IF A$ = "" THEN VTAB PEEK (37): HTAB 15: PRINT SUB$: GOTO 139
136 SUB$ = A$
139 PRINT "Id #: " :ID$;
```

```

139 PRINT "Id #: ";ID$;
142 HTAB 6: INPUT A$: IF A$ = "" THEN VTAB PEEK (37): HTAB 7: PRINT I
D$: GOTO 148
145 ID$ = A$
148 PRINT "Date: ";DT$;
151 HTAB 6: INPUT A$: IF A$ = "" THEN VTAB PEEK (37): HTAB 7: PRINT D
T$: GOTO 157
154 DT$ = A$
157 PRINT : PRINT "Eye (1=OD, 2=OS): ";
160 IF EYE = 0 THEN PRINT "OD";
163 IF EYE = 1 THEN PRINT "OS";
166 HTAB 18: INPUT A$: IF A$ = "" THEN 172
169 EYE = VAL (A$) - 1
172 VTAB PEEK (37): HTAB 19: IF EYE = 0 THEN PRINT "OD"
175 IF EYE = 1 THEN PRINT "OS"
180 HOME : VTAB 9
181 HTAB 15: PRINT "      A      "
182 HTAB 15: PRINT " H          B "
183 HTAB 15: PRINT "      P I J      "
184 HTAB 15: PRINT "G  O Q K  C"
185 HTAB 15: PRINT "      N M L      "
186 HTAB 15: PRINT " F          D "
187 HTAB 15: PRINT "      E          "
188 POKE 34,16: REM SET TOP OF SCROLL WINDOW
189 KEY = 1:OKEY = 1: REM INITIAL KEYPRESS VALUE
199 REM ***MAIN PROGRAM LOOP
200 IF KEY < 0 THEN 210: REM SKIP LIGHT DISPLAY IF NON-LIGHT KEY PRESSE
D
201 IF OKEY > 0 AND KEY < 17 THEN VTAB Y%(OKEY): HTAB X%(OKEY): PRINT
CHR$ (OKEY + 64): REM TURN OFF PREVIOUS SCREEN CHAR.
202 INVERSE
204 VTAB Y%(KEY): HTAB X%(KEY): IF KEY = 17 THEN IF C < 16 THEN NORMAL

205 IF KEY = 17 THEN PRINT "Q": GOTO 210
206 PRINT LED$: REM LIGHT SCREEN CHAR.
210 NORMAL : VTAB 6: PRINT " "; IF EYE = 0 THEN PRINT "OD";
211 IF EYE = 1 THEN PRINT "OS";
213 PRINT " Light #: ";
214 IF C > 15 THEN PRINT "Q";: GOTO 217
215 PRINT LED$;
217 HTAB 26: PRINT "<ESC>=New eye"
218 HTAB 25: PRINT "^D=Delete point"
219 REM **KEYBOARD HANDLER
220 IF KEY > 0 AND KEY < 17 THEN OKEY = KEY
221 VTAB 6: HTAB 14
222 IF PEEK (49249) > 127 THEN PRINT G$: GOTO 275: REM PUSHBUTTON PRE
SSED
223 IF PEEK (49152) < 128 THEN 222: REM LOOP UNTIL KEY PRESSED
224 KEY = PEEK (49152) - 192: POKE 49168,0:A$ = CHR$ (KEY + 64): IF KE
Y = - 51 THEN VTAB PEEK (37) + 1: HTAB 10: PRINT : GOTO 270: REM
RETURN KEY PRESSED
225 IF KEY = - 32 THEN PRINT : GOTO 275: REM SPACE KEY PRESSED
230 IF KEY = - 37 THEN PRINT : GOTO 440: REM ESC KEY PRESSED
231 IF KEY = - 60 THEN PRINT : GOTO 1000
235 IF KEY < 1 OR KEY > 17 THEN PRINT G$;: GOTO 221
240 IF KEY = 17 THEN C = C + 16: IF C > 16 THEN C = 0
245 IF KEY < > 17 THEN LED = KEY:LED$ = A$
270 POKE 49280 + S1 * 16,L(LED) + C: REM TURN ON LED
271 GOTO 300: REM SKIP A/D CODE
274 REM READ PACHOMETER HEAD POSITION
275 POKE A113,18:N = DT%(EYE,LD,0) + 1
276 IF N > 12 THEN N = 12: VTAB 24: HTAB 7: PRINT G$G$"Maximum # of sam
ples = 12": GOTO 200
280 DT%(EYE,LD,N) = PEEK (A113 + 1) * 256 + PEEK (A113) - ZERO:DT%(EYE
,LD,0) = N

```

```

301 IF C < 10 THEN LD = 17: PRINT "Q"; GOTO 304
302 PRINT LED$;
304 PRINT "Data (mm):"
305 N = DT%(EYE,LD,0): REM # OF SAMPLES
310 K = 0: FOR I = 1 TO N STEP 4
312 FOR J = 0 TO 3
315 IF K >= N THEN 320
316 HTAB 2 + J * 10: IF K > 8 THEN HTAB PEEK (36)
317 PRINT K + 1;";"; INT (DT%(EYE,LD,K + 1) / CAL * 1000) / 1000;;K = K
+ 1
320 NEXT J
325 PRINT : NEXT I
349 REM CALCULATE STATISTICS
350 MN = 0:SD = 0:RG = 0:SQ = 0
353 IF N = 0 THEN 390
354 HI% = DT%(EYE,LD,1):LO% = HI%
355 FOR I = 1 TO N
360 D% = DT%(EYE,LD,1):MN = MN + D%:SQ = SQ + D% * D%
365 IF D% < LO% THEN LO% = D%
370 IF D% > HI% THEN HI% = D%
375 NEXT I
380 SD = SQR (SQ - (MN * MN / N)) / CAL
381 MN = MN / N / CAL
382 RG = (HI% - LO%) / CAL
390 PRINT : PRINT " Mean:"; INT (MN * 1000 + .5) / 1000;; HTAB 13: PRINT
" S.D.:"; INT (SD * 1000 + .5) / 1000;; HTAB 26: PRINT " Range:";
INT (RG * 1000 + .5) / 1000
395 DT(EYE,LD,1) = MN:DT(EYE,LD,2) = SD:DT(EYE,LD,3) = RG
400 GOTO 200
415 HTAB 7: INPUT "Press <RETURN> to continue:";A$
440 POKE 34,7: HOME
442 S$ = "N": IF EYE = 1 THEN S$ = "Y"
444 EYE = EYE + 1: IF EYE > 1 THEN EYE = 0
450 PRINT "Save data (Y/N): "S$;
455 HTAB 17: INPUT A$: IF A$ = "" THEN 470
460 IF A$ < > "Y" AND A$ < > "y" AND A$ < > "N" AND A$ < > "n" THEN
PRINT G$;; GOTO 450
465 S$ = A$
470 IF S$ = "Y" OR S$ = "y" THEN 500
490 GOTO 100
500 FI$ = LEFT$ (SUB$,1) + ID$ + ".PACH." + DT$
505 PRINT "Filename: "FI$;
510 HTAB 10: INPUT A$: IF A$ = "" THEN VTAB PEEK (37): HTAB 11: PRINT
FI$: GOTO 520
515 FI$ = A$
520 PRINT "Drive #: "DR;
521 HTAB 9: INPUT A$: IF A$ = "" THEN VTAB PEEK (37): HTAB 10: PRINT
DR: GOTO 529
525 A = VAL (A$): IF A < 1 OR A > 2 THEN VTAB PEEK (37): PRINT G$;; GOTO
521
526 DR = A
529 PRINT "Saving data..."
530 PRINT D$"OPEN"FI$,D"DR: PRINT D$"WRITE"FI$
535 PRINT SUB$: PRINT ID$: PRINT DT$
538 PRINT CAL: PRINT ZERO
540 PRINT "OD": FOR I = 1 TO 17: FOR J = 1 TO 3: PRINT DT(0,I,J): NEXT
J: NEXT I
545 PRINT "OS": FOR I = 1 TO 17: FOR J = 1 TO 3: PRINT DT(1,I,J): NEXT
J: NEXT I
550 PRINT "OD": FOR I = 1 TO 17: FOR J = 0 TO DT%(0,I,0): PRINT DT%(0,I
,J): NEXT J: NEXT I
555 PRINT "OS": FOR I = 1 TO 17: FOR J = 0 TO DT%(1,I,0): PRINT DT%(1,I
,J): NEXT J: NEXT I
560 PRINT D$"CLOSE"FI$
565 VTAB PEEK (37): HTAB 15: PRINT "finished."
570 PRINT : PRINT "Exit program (Y/N): "S$;
575 HTAB 22: INPUT " ";A$

```

```

580 IF A$ < > "Y" AND A$ < > "y" AND A$ < > "N" AND A$ < > "n" THEN
    PRINT G$;: GOTO 575
585 S$ = A$
590 IF S$ = "N" OR S$ = "n" THEN 100
595 PRINT D$"RUNHELLO,D1"
999 REM DELETE A DATA POINT
1000 PRINT "Delete data #: ";
1001 HTAB 16: INPUT A$: IF A$ = "" THEN 1020
1003 PT = VAL (A$): IF PT < 1 OR PT > 12 THEN VTAB PEEK (37): PRINT G
    $;: GOTO 1001
1004 IF PT = N THEN 1015
1005 IF PT > N THEN 1020
1009 FOR I = PT TO 11
1010 DT%(EYE,LD,I) = DT%(EYE,LD,I + 1): REM SHIFT CA DOWN BY 1
1014 NEXT I
1015 N = N - 1:DT%(EYE,LD,0) = N
1020 VTAB PEEK (37): HTAB 1: PRINT " ": REM ERASE LI
    NE
1025 GOTO 300: REM REDRAW DATA
9000 REM General purpose error handler
9005 TEXT :I = PEEK (222): PRINT G$"Error: ";: IF I = 0 OR I > 15 THEN
    J = 53856 + I: GOTO 9015
9010 J = 43377 + PEEK (43583 + I)
9015 IF I = 254 THEN PRINT "Bad input";
9020 IF I = 255 THEN PRINT "Cntrl C attempted.": GOTO 9040
9025 K = PEEK (J): PRINT CHR$ (K);: IF K < 192 THEN J = J + 1: GOTO 90
    25
9030 PRINT " at line "; PEEK (218) + 256 * PEEK (219): IF I = 10 THEN
    PRINT D$"UNLOCK"FI$: PRINT "File has been unlocked. Try again."
9035 IF I = 4 OR I = 8 OR I = 9 THEN PRINT "Insert another disk and tr
    y again.": GOSUB 9050: RESUME
9040 POKE 216,0: GOSUB 9050: GOTO 100
9050 HTAB 7: INPUT "Press <RETURN> to continue: ";A$: RETURN

```

# LIST

```

10 REM *****
15 REM PACH.REPORT Created for Pacific University
17 REM Copyright (c) 1985 by Kevin Spreier
20 REM Version 1.02 09 Sep 85
25 REM *****
30 X$ = "": REM Dummy variable for special INPUT$ routine
35 INPT = 940: REM Special input routine to allow commas
50 D$ = CHR$ (4):G$ = CHR$ (7):E$ = CHR$ (27)
55 PRINT D$"BLOADINPUT$,A940"
58 PR = 1: REM Default printer interface slot
60 DR = 2: REM Default data drive
65 DIM DT(1,17,3): REM Calculated data array
90 PRINT D$"PR#3": REM Enable 80-column card
100 REM Menu selection
105 POKE 49168,0: REM Clear keyboard strobe
110 PRINT CHR$ (17): REM Display 40 columns
115 TEXT : HOME : GOSUB 985: HTAB 7: PRINT "Pachometer Report Generator
": HTAB 3: PRINT "Copyright (c) 1985 by Kevin Spreier": GOSUB 985: VTAB
8: ONERR GOTO 1005
119 POKE 34,5: REM Don't scroll over top 4 screen lines
120 PRINT
125 IF FIRST < > 0 THEN 200: REM Skip questions if not first time thro
ugh
130 P = 0: PRINT "Do you want a printout? Y";
131 HTAB 25: INPUT "":A$: IF A$ = "" THEN A$ = "Y"
132 IF A$ = "Y" OR A$ = "y" THEN P = 1: GOTO 135
133 IF A$ < > "N" AND A$ < > "n" THEN PRINT G$;: GOTO 131
135 VTAB PEEK (37): HTAB 25: PRINT A$
160 PRINT : PRINT "Insert data disk in drive "DR" and": PRINT " press
<RETURN>: ";
161 INPUT "":A$
195 FIRST = 1
200 HOME :T = PEEK (34): PRINT
205 PRINT D$"CATALOG,D"DR
210 FOR V = PEEK (34) TO PEEK (37) + 1: VTAB V: HTAB 5
215 CH = PEEK (( PEEK (40) + PEEK (41) * 256) + PEEK (36)): IF CH = 1
60 THEN T = T + 1: NEXT V: GOTO 225
220 HTAB 4: PRINT "-": CHR$ (V - T + 65): "-": NEXT V
225 POKE 49168,0: PRINT "Enter letter (1,2=Dr; <RET>=Done):";
230 GET A$:A = ASC (A$): IF A = 49 OR A = 50 THEN DR = A - 48: GOTO 20
0
231 IF A = 13 THEN 900: REM Return key pressed
232 IF A > 96 THEN A = A - 32: REM Convert lower case letters to upper
case
233 IF A - 64 > (V - T) OR A - 65 < 0 THEN PRINT G$;: GOTO 230
235 VTAB A + T - 66: HTAB 2:CH = PEEK (( PEEK (40) + PEEK (41) * 256)
+ PEEK (36)): IF CH < > 212 THEN PRINT G$;: GOTO 230
240 HTAB 8:FI$ = "":LCHAR = 0:SPCNT = 0
241 FOR X = PEEK (36) TO 39
242 CH = PEEK (( PEEK (40) + PEEK (41) * 256) + X): IF CH = 160 THEN IF
LCHAR = 160 THEN SPCNT = SPCNT + 1
244 IF SPCNT = 6 THEN FI$ = LEFT$ (FI$, LEN (FI$) - 2): GOTO 250
245 FI$ = FI$ + CHR$ (CH):LCHAR = CH: NEXT X
250 HTAB 8: INVERSE : PRINT FI$: NORMAL
299 REM Read data from file
310 PRINT D$"OPEN"FI$

```

```

310 PRINT D$ "READ FI$
320 CALL INPT:SUB$ = MID$(X$,1): REM Input string with commas
322 INPUT ID$: INPUT DT$
324 INPUT CAL: INPUT ZERO
325 INPUT OD$: FOR I = 1 TO 17: FOR J = 1 TO 3: INPUT DT(0,I,J): NEXT J
: NEXT I
330 INPUT OS$: FOR I = 1 TO 17: FOR J = 1 TO 3: INPUT DT(1,I,J): NEXT J
: NEXT I
335 PRINT D$"CLOSE"FI$
395 PRINT : REM Fix enhanced 80-col. bug
399 REM Print report
400 PRINT CHR$(18): POKE 34,0: HOME : REM Display 80 columns and clear
screen
401 IF P = 0 THEN 410
405 PRINT D$"PR#1"PR: REM Turn printer on
410 PRINT "Name: "SUB$
415 PRINT "Id #: "ID$
420 PRINT "Date: "DT$
421 PRINT
422 PRINT SPC(13);"-OS (mm)-"; SPC(31);"-OD (mm)-"
423 PRINT
425 PRINT SPC(6);"Mean"; SPC(5);"Range"; SPC(6);"S.D."; SPC(16);"M
ean"; SPC(5);"Range"; SPC(6);"S.D."
426 PRINT SPC(6);"----"; SPC(5);"----"; SPC(6);"----"; SPC(16);"
----"; SPC(5);"----"; SPC(6);"----"
428 FOR I = 1 TO 17: REM Print LED data
429 REM OS data
430 D = DT(1,I,1): IF D = 0 THEN IF DT(0,I,1) = 0 THEN 475: REM Don't p
rint a line if both OS, OD values are 0
435 PRINT CHR$(I + 64)": "; REM Label LED #
436 IF D = 0 THEN PRINT SPC(35); GOTO 455: REM No data for that spe
cific OS LED
440 GOSUB 800: REM Print data value
445 FOR J = 2 TO 3:D = DT(1,I,J): GOSUB 800: NEXT J
449 REM OD data
450 PRINT SPC(5); REM Space between eyes
455 D = DT(0,I,1): PRINT CHR$(I + 64)": "; REM Label LED #
457 IF D = 0 THEN 470: REM No data for that specific OD LED
460 GOSUB 800: REM Print data value
465 FOR J = 2 TO 3:D = DT(0,I,J): GOSUB 800: NEXT J
470 PRINT
475 NEXT I
480 IF P = 0 THEN GOSUB 995
481 IF P = 1 THEN PRINT CHR$(12): REM Send form-feed char. to printe
r
485 PRINT D$"PR#3": REM Turn printer off and clear screen
490 GOTO 100
799 REM Print the variable D in the format "0.000"
800 D = INT(D * 1000 + .5): IF D > 1000 THEN D = D - 1000: PRINT "1.";
: GOTO 810
805 PRINT "0.";
810 IF D < 100 THEN PRINT "0"; IF D < 10 THEN PRINT "0";
815 PRINT D; SPC(5);
820 RETURN
899 REM Go back to main menu
900 PRINT : PRINT : PRINT "Exit program (Y/N): Y";
905 HTAB 21: INPUT "":A$: IF A$ = "" THEN A$ = "Y"
910 IF A$ < > "Y" AND A$ < > "y" AND A$ < > "N" AND A$ < > "n" THEN
PRINT G$; GOTO 905
915 IF A$ = "N" OR A$ = "n" THEN FIRST = 0: GOTO 100
916 VTAB PEEK(37): HTAB 21: PRINT A$
920 PRINT D$"RUNHELLO,D1"
985 FOR I = 1 TO 39: PRINT "*"; NEXT : PRINT : RETURN
995 PRINT : HTAB 7: PRINT "Press any key to continue."; GET A$: RETURN

1000 REM General purpose error handler

```



```

1005 POKE 216,0: TEXT :I = PEEK (222): PRINT G$"Error: ";: IF I = 0 OR
    I > 15 THEN J = 53856 + I: GOTO 1015
1010 J = 43377 + PEEK (43583 + I)
1015 IF I = 254 THEN PRINT "Bad input";
1020 IF I = 255 THEN PRINT "Cntrl C attempted.": GOTO 1040
1025 K = PEEK (J): PRINT CHR$ (K);: IF K < 192 THEN J = J + 1: GOTO 10
    25
1030 PRINT " at line "; PEEK (218) + 256 * PEEK (219): IF I = 10 THEN
    PRINT D$"UNLOCK"FI$: PRINT "File has been unlocked. Try again."
1035 IF I = 4 OR I = 8 OR I = 9 THEN PRINT "Insert another disk and tr
    y again.": GOSUB 995: RESUME
1040 GOSUB 995: GOTO 100

```

]

APPENDIX E  
PACHOMETER INSTRUCTION MANUAL

## PACHOMETER INSTRUCTION MANUAL

### GENERAL DESCRIPTION

The purpose of the Pachometer is to provide the mechanical to electronic digital conversion, data analysis, operator interaction and the output functions necessary in order to accurately measure corneal thickness. The total measurement system consists of a slitlamp, a fixation device, an optical doubling device (with eyepiece), a potentiometer, and an Apple IIe (enhanced) computer.

The fixation device allows fixation of the patient to be changed at the discretion of the operator in eight (8) meridians and seventeen (17) locations. The fixation device contains nineteen (19) L.E.D.s (Light Emitting Diodes). There are seventeen (17) red L.E.D.s arranged in eight (8) meridians and two (2) yellow L.E.D.s that allow for fixation in the vertical plane to insure that readings are taken perpendicular to the cornea. Any one of the 17 red L.E.D.s may be illuminated at random from the Apple IIe keyboard. The two yellow L.E.D.s remain illuminated at all times to insure proper alignment.

The potentiometer assembly attaches to the slitlamp by slipping it on the mounting post used with a Goldman Applanation Tonometer. The potentiometer is slid over the large pin and is locked into place with the small short pin. The potentiometer is a device for measuring electromotive force or potential difference by comparison with a known voltage.

The bi-prism eyepiece is inserted into the right ocular of the slitlamp with the eyepiece slit in the horizontal orientation. This allows for vernier alignment of the upper and lower slits.

The Apple IIe computer along with the "Contact Lens Research Software" (Program written by Kevin Spreier) performs all of the analysis and hardware control functions. The computer receives all of the input from the potentiometer, performs the analysis, controls the fixation device, and provides the interface to the Apple Imagewriter Printer. The computer program computes the mean, standard deviation, and range of the current set of thickness samples. These items are recomputed after each measurement or after deleting a measurement. The computer and program also calibrates each measurement to compensate for operator biases in the measurement process.

#### OPERATING INSTRUCTION

Insert the program (floppy disk) titled Contact Lens Research Software into Drive 1 of the Apple IIe computer. Next insert a formatted disk into Drive 2 for data storage. Turn on the computer using the switch on the system saver on the left side of the computer. The computer screen will display the following:

```
*****
*               Contact Lens Research Software               *
*               Copyright (c) 1985 by Kevin Spreier           *
*****
```

#### **Options:**

1. PERMEOMETER
2. PACHOMETER
3. PACH. REPORT
4. Exit to BASIC

**Choose one -->**

Press the "2" key on the keyboard and you will get the following on the computer screen:

\*\*\*\*\*

**Pachometer Program**  
**Copyright (c) 1985 by Kevin Spreier**

\*\*\*\*\*

**New zero mark (Y/N) ? N**

It is advisable to zero the instrument each time a new operator uses the instrument. If you wish to zero the instrument press the "Y" key and then the "Return" key. The computer screen will then display the following:

\*\*\*\*\*

**Pachometer Program**  
**Copyright (c) 1985 by Kevin Spreier**

\*\*\*\*\*

**New zero mark (Y/N) ?**

**Zero Pachometer head and press <RETURN>:**

At this time turn on the slitlamp and insert the test (focus) bar into the slitlamp. Lightly grasp the adjustment rod from the potentiometer and while observing the narrow slit of light that is focused on the test bar, move the adjustment rod until the slit of light becomes straight and single. Figure 1 shows the proper alignment picture.




Figure 1

When the observer has established proper alignment of the slit then he/she should press the "RETURN" key. This is also the proper time to focus the eyepiece. The computer screen will then display the following:

**Subject name: ?**

Type in the name (Last name first) and press "RETURN".  
The screen will then display the following:

**Id #: ?**

You may then enter up to a three (3) digit number to identify your patient. Again press the "RETURN" key and the screen will then display the following:

**Date: ?**

Enter the date and other information that may be desired i.e., Day 7 (of wearing extended wear contact lenses).

\* The date line can only be one line long on the screen.\*

Press the "RETURN" key and the screen will then appear as follows:

```
*****
                          Pachometer Program
                          Copyright (c) 1985 by Kevin Spreier
*****
```

**Subject name: ? Doe, John D.**

**Id #: ? 123**

**Date: ? 15 Dec 99 Day 7**

**Eye (1=OD, 2=OS): ? OD**

The computer program has preset the right eye to measure first, and you do not have to change this unless you do not need to measure the right eye or if you choose to measure the left first. If you are ready to measure the right eye simply press the "RETURN" key. The screen will then appear as follows:

\*\*\*\*\*

**Pachometer Program**  
**Copyright (c) 1985 by Kevin Spreier**

\*\*\*\*\*

OD Light #: A

<ESC> = New eye  
 /\ D = Delete Point

```

      A
    H   B
  P I J
G  O Q K  C
    N M L
  F       D
      E
  
```

You may then pick the fixation point that you want to measure by pressing the appropriate key. If you choose to measure central thickness first the press the "Q" key and the screen will appear as follows:

\*\*\*\*\*

**Pachometer Program**  
**Copyright (c) 1985 by Kevin Spreier**

\*\*\*\*\*

OD Light #: Q

<ESC> = New eye  
 /\ D = Delete Point

```

      A
    H   B
  P I J
G  O Q K  C
    N M L
  F       D
      E
  
```

'Q' Data (mm):

Mean: Ø

S.D.: Ø

Range: Ø

Up to twelve (12) values may be recorded in any of the fixation locations. If more than twelve (12) are attempted, the following will appear on the screen:

**Maximum # of samples = 12**

Due to the wiring of the fixation device, the central fixation point (Q) can only be illuminated in conjunction with another point. It is suggested that a peripheral point be selected and then the "Q". When finished with the central fixation point (Q), it is necessary to press the "Q" key to release that point and then make the next selection.

When all of the measurements have been completed for the right eye, press <ESC> and the screen will appear as follows:

```
*****
                        Pachometer Program
                Copyright (c) 1985 by Kevin Spreier
*****
```

**OD Light #: Q**

**<ESC> = New eye**

**/\ D = Delete Point**

**Save data (Y/N): ? N**

If both eyes are to be measured then simply press "RETURN" and continue. If you only want to record the right eye then change the "N" to "Y" and press "RETURN". When continuing to record the left eye measurements the "Subject Name" line will appear with the name already recorded. Again press the "RETURN" key and the "Id #" will appear and so on until you again have the fixation select pattern. Simply follow the procedures outlined for the right eye in accomplishing the measurements for the left.

If a value is recorded that is significantly different than the others, it may be deleted. Press the "CONTROL" key simultaneously with the "D" key and then select the number of the value to be deleted and press the corresponding key followed by the "RETURN" key. The point will be deleted and the mean, S.D., and range values will be recalculated.

After measurements for both eyes have been completed, the values should be recorded on the data disk. Press the "ESC" key and the following will appear:

**Save data (Y/N): ? Y**



If you do want to save the data, simply press the "RETURN" key and the following appears:

**Filename: ? D123 .PACH. 15DEC99 DAY 7**

This represents the file for John Doe, Id # 123, taken on 15 Dec 99, and the information sequence of Day 7. Again press the "RETURN" key and the following will appear:

**Drive #: ? 2**

Since the Data Storage Disk is in drive 2, this is correct and you need only press the "RETURN" key. the following appears:

**Saving data. . .**

followed by the following when finished:

**Saving data. . . finished**

**Exit program (Y/N): ? Y**

If finished, press the "RETURN" key and you will return to the main menu. If another patient is to be measured, then change the "Y" to "N" and press "RETURN" and the sequence starts again.

### PACHOMETER REPORT

In order to obtain a hard copy of the measurements taken, you must again view the main menu that appears as below:

```

*****
*               Contact Lens Research Software               *
*               Copyright (c) 1985 by Kevin Spreier           *
*****

```

**Options:**

1. PERMEOMETER
2. PACHOMETER
3. PACH. REPORT
4. Exit to BASIC

**Choose one - ->**

Press The "3" key and the following screen will appear:

```

*****
*               Pachometer Report Generator                   *
*               Copyright (c) 1985 by Kevin Spreier           *
*****

```

**Do you want a printout? Y**

If you want a printout then press the "RETURN" key and the following appears:

**Insert data disk in drive 2 and press <RETURN>:**

Once you have pressed the "RETURN" key the following will appear:

```

*****
*               Pachometer Report Generator                   *
*               Copyright (c) 1985 by Kevin Spreier           *
*****
T 005 D123 .PACH. 15DEC99 DAY 7
T 005 M321 .PACH. 15DEC99 DAY 14
T 005 J456 .PACH. 15DEC99 DAY 1
T 005 C999 .PACH. 16DEC99 DAY 30

```

The above is an example of the Filenames that could appear on the data storage disk. The data disk will store approximately eighty (80) files. It will display eighteen (18) filenames at one time. If the one you are wishing to obtain is not on the screen, then press any key and the

following eighteen files will appear. Once you have found the file you need, Press the "CONTROL" key simultaneously with the "C" key and the screen will appear as below:

\*\*\*\*\*

**Pachometer Report Generator  
Copyright (c) 1985 by Kevin Spreier**

\*\*\*\*\*

T -A- D123 .PACH. 15DEC99 DAY 7  
T -B- M321 .PACH. 15DEC99 DAY 14  
T -C- J456 .PACH. 15DEC99 DAY 1  
T -D- C999 .PACH. 16DEC99 DAY 30  
ENTER LETTER (1,2 = DR; <RET> = DONE):

Again the screen will display all files that were previously seen, up to 18. Press the key corresponding to the letter between the dashes in front of the file you wish to copy. Be sure the printer is turned on. If the file is located on a disk where there are not yet eighteen (18) files, you will not have to press the "CONTROL" and the "C" keys as the computer will do it automatically. Also if the file is located in the last group of files displayed, the computer will accomplish this step. The following is an example of the report that will be printed.

**Name: Doe, John**

**Id # : 123**

**Date : 15DEC99 day 7**

	-OS (mm)-				-OD (mm)-		
	Mean	Range	S.D.		Mean	Range	S.D.
	----	-----	----		----	-----	----
I:	0.526	0.010	0.013	I:	0.539	0.008	0.009
K:	0.521	0.008	0.010	K:	0.531	0.009	0.010
M:	0.515	0.009	0.012	M:	0.545	0.010	0.013
O:	0.533	0.012	0.012	O:	0.514	0.013	0.014
Q:	0.499	0.011	0.013	Q:	0.495	0.012	0.013

If when first entering PACH. REPORT and the question "Do you want a printout ? Y" appears and you do not want a printout, then change the "Y" to "N" and press "RETURN". Follow the same instructions as above and the report will appear on the screen in the same format as above.

When you are ready to exit the program, you must have the statement "Enter letter (1,2 = Dr; <RET> = Done):" at the bottom of the screen. Press the "RETURN" key and "Exit program (Y/N): Y" will appear. If you are ready

to exit the program, simply press the "RETURN" key and you will return to the main menu. If not ready to exit, change the "Y" to "N" and press the "RETURN" and the sequence will start again.

Always return to the main menu before quitting. Once the main menu is on the screen, remove both disks and then turn off the computer.

### MEASUREMENT PROCEDURE

After preparing the slitlamp and choosing the fixation point desired, the observer is ready to take the thickness measurement. The fixation board has two functions. First, the two yellow L.E.D.s assure that the readings are taken perpendicular to the cornea. Second, the fixation board insures that the readings are taken at the same location on the cornea of each patient even on multiple sittings. The computer allows the observer to have immediate statistical feedback as to the accuracy of the readings being taken.

### PATIENT INSTRUCTIONS

1. Inform the patient. Tell them that you are going to measure the thickness of the cornea but will not touch the eye.
2. Point out the Red L.E.D. that is illuminated (usually start with the center one). Ask the patient to maintain steady fixation on the L.E.D. that is illuminated.
3. Remind them to stay as still as possible with their forehead against the headrest.
4. Inform the patient that they may blink whenever they need to.
5. Inform the patient that they will hear a beep from the computer each time you take a reading.

### TAKING THE MEASUREMENT

1. The observer's left hand should be on the slitlamp joy stick. The right hand on the potentiometer adjustment rod. Keeping the right elbow on the table is helpful in making a more steady adjustment during the measurement.

2. The observer may use either the right or the left eye when looking through the bi-prism eyepiece. The angle between the light beam and the microscope should be set at forty (40°) degrees, with the light beam to the left of the microscope.

3. Up to twelve (12) values may be recorded for each fixation point. Readjustments in alignment will be necessary throughout the collection of the values. The following is necessary before each value is collected:

a. The slitlamp parallelepiped should be in the best possible focus, obtained by moving the slitlamp joy stick forward and backward. The observer will have to make minor adjustments to keep the sharpest focus at all times. The observer should first focus the inside surface of the upper image and then focus the inside surface of the lower image. One will always be slightly better than the other, so that the "best" focus is a compromise between the two.

b. The upper and lower yellow L.E.D.s should be an equal distance above and below the horizontal line of the bi-prism eyepiece as shown in figure 2. The yellow L.E.D.s should be superimposed by the split parallelepiped.

c. The red fixation L.E.D. should be positioned just above the horizontal line created by the bi-prism eyepiece when taking the central measurement. In most other positions, this is not possible and the red fixation L.E.D. should be as near as possible. It is assured that the readings being taken are perpendicular to the cornea as long as the yellow L.E.D.s are properly positioned.

d. To make an accurate reading, the inside surface of the upper and the lower parallelepiped images should be just touching, but not overlapping. Start the reading with the two surfaces well separated and bring them together slowly until the observer sees first touch as shown in figure 2. If the surfaces do not touch or if they overlap, the observer will record an incorrect value (see figure 3).

e. Occlusion of the eye not being tested is helpful in maintaining accurate fixation. \*The scale that is physically located on the potentiometer is not accurate when utilized with this computer setup.\*




FIGURE 2




FIGURE 3

4. When all of the above (3., a.-d.) have been met, press the record button (located on top of the joy stick or the space bar on the computer) to record that value. After recording all of the values desired for that fixation point, choose the next fixation point and repeat the process.

APPENDIX F  
PACHOMETRY SUBJECT DATA FILES  
AND CONTROL PATIENT  
DATA FILES

## INFORMATION MISSING FROM PATIENT FILES

<u>SUBJECT #</u>	<u>EYE</u>	<u>DAY</u>	<u>POINTS MISSING</u>
101	OS	7	I
102	OU	7	ALL
103	OS	0	I, K
103	OU	4	ALL
104	OU	1	ALL
107	OS	14	I
109	OU	21	ALL
113	OU	4	ALL
119	OD	21	Q
119	OS	21	I
132	OD	118	I
133	OS	1	Q, O, K, M
133	OU	7	ALL
136	OU	21	ALL
137	OU	14	ALL
139	OU	7	ALL
139	OS	14	I, K

ABOVE DATA POINTS MISSING DUE TO MALFUNCTION OF COMPUTER  
CAUSING DATA STORAGE DISK TO BE OVER WRITTEN.



SUBJECT # 101

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.542	0.577	0.569	0.579	0.585
1	0.542	0.595	0.575	0.607	0.613
4	0.6	0.62	0.615	0.631	0.627
7	0.587	0.631	0.589	0.607	0.587
14	0.593	0.618	0.578	0.617	0.586
20	0.57	0.607	0.58	0.61	0.603
28	0.567	0.575	0.549	0.574	0.573

JOINED ARMY - DROPPED

SUBJECT # 101

EYE: OS

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.56	0.574	0.579	0.567	0.584
1	0.534	0.561	0.59	0.531	0.577
4	0.585	0.569	0.606	0.581	0.598
7	0.531		0.561	0.554	0.554
14	0.555	0.599	0.601	0.571	0.58
20	0.58	0.589	0.6	0.568	0.57
28	0.546	0.576	0.565	0.566	0.552

JOINED ARMY - DROPPED

SUBJECT # 102

EYE: OD

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.495	0.54	0.541	0.545	0.545
1	0.482	0.494	0.503	0.516	0.525
4	0.511	0.512	0.498	0.521	0.539
7					
14	0.477	0.527	0.506	0.518	0.533
21	0.468	0.519	0.507	0.523	0.533
28	0.478	0.52	0.494	0.524	0.517
62	0.516	0.54	0.548	0.542	0.555
91	0.51	0.567	0.557	0.557	0.562
120	0.49	0.53	0.516	0.489	0.528
149	0.483	0.514	0.505	0.519	0.521

SUBJECT # 102

EYE: OS

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.489	0.499	0.516	0.518	0.531
1	0.479	0.499	0.534	0.503	0.513
4	0.498	0.518	0.541	0.508	0.541
7					
14	0.492	0.525	0.512	0.49	0.508
21	0.483	0.5	0.525	0.504	0.512
28	0.468	0.505	0.53	0.51	0.528
62	0.504	0.532	0.561	0.534	0.545
91	0.523	0.568	0.568	0.548	0.56
120	0.461	0.499	0.507	0.478	0.505
149	0.495	0.527	0.521	0.495	0.532

SUBJECT # 103

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.488	0.512	0.513	0.531	0.528
1	0.525	0.573	0.579	0.583	0.585
4					
11	0.546	0.561	0.539	0.565	0.571
15	0.52	0.516	0.529	0.538	0.539
21	0.495	0.537	0.51	0.535	0.538
28	0.488	0.515	0.517	0.54	0.541
60	0.52	0.55	0.518	0.552	0.554
110	0.517	0.534	0.534	0.555	0.554
123	0.515	0.543	0.534	0.553	0.564
152	0.514	0.547	0.537	0.563	0.567

SUBJECT # 103

EYE: OS

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.519		0.56		0.568
1	0.539	0.559	0.557	0.528	0.561
4					
11	0.547	0.576	0.567	0.536	0.561
15	0.508	0.522	0.552	0.509	0.529
21	0.502	0.52	0.525	0.481	0.533
28	0.503	0.532	0.53	0.512	0.547
60	0.539	0.548	0.573	0.512	0.571
110	0.531	0.552	0.567	0.512	0.558
123	0.523	0.518	0.56	0.528	0.554
152	0.535	0.562	0.58	0.532	0.555

SUBJECT # 104

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.534	0.556	0.551	0.563	0.558
1					
4	0.512	0.549	0.523	0.529	0.545
7	0.497	0.563	0.517	0.55	0.54
14	0.557	0.584	0.546	0.595	0.583
24	0.53	0.56	0.538	0.541	0.548
30	0.511	0.529	0.543	0.565	0.54
60	0.524	0.56	0.558	0.568	0.561
90	0.502	0.543	0.528	0.541	0.521
117	0.513	0.552	0.53	0.538	0.535

COMPLETED DAY 120 ONLY

SUBJECT # 104

EYE: OS

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.534	0.568	0.574	0.547	0.557
1					
4	0.516	0.526	0.545	0.518	0.53
7	0.52	0.53	0.527	0.503	0.54
14	0.574	0.586	0.599	0.57	0.581
24	0.538	0.568	0.567	0.518	0.548
30	0.534	0.56	0.567	0.537	0.56
60	0.534	0.566	0.567	0.545	0.56
90	0.517	0.535	0.546	0.528	0.542
117	0.541	0.545	0.545	0.518	0.541

COMPLETED DAY 120 ONLY

SUBJECT # 105

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.505	0.517	0.514	0.566	0.55
1	0.489	0.541	0.533	0.534	0.558
4	0.53	0.559	0.56	0.587	0.582
7	0.529	0.54	0.557	0.575	0.553
14	0.521	0.504	0.51	0.547	0.545
21	0.532	0.546	0.543	0.585	0.589
31	0.56	0.568	0.568	0.588	0.602
60	0.54	0.551	0.55	0.565	0.597
90	0.522	0.56	0.561	0.586	0.573
125	0.538	0.554	0.556	0.58	0.596
151	0.54	0.548	0.56	0.585	0.59

SUBJECT # 105

EYE: OS

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.523	0.534	0.55	0.518	0.547
1	0.536	0.522	0.566	0.542	0.541
4	0.576	0.576	0.585	0.58	0.579
7	0.55	0.57	0.589	0.566	0.591
14	0.523	0.539	0.539	0.526	0.547
21	0.554	0.554	0.589	0.559	0.576
31	0.567	0.58	0.613	0.578	0.613
60	0.557	0.562	0.59	0.577	0.591
90	0.546	0.564	0.589	0.563	0.574
125	0.564	0.563	0.597	0.556	0.586
151	0.536	0.551	0.592	0.581	0.589

SUBJECT # 106

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.538	0.568	0.517	0.534	0.554
1	0.525	0.553	0.527	0.563	0.55
4	0.54	0.576	0.554	0.57	0.565
7	0.544	0.564	0.545	0.578	0.578
14	0.535	0.564	0.561	0.574	0.582
21	0.551	0.556	0.544	0.581	0.575
28	0.554	0.556	0.557	0.574	0.583
61	0.552	0.579	0.588	0.602	0.596

COMPLETED DAY 61 ONLY

SUBJECT # 106

EYE: OS

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.519	0.551	0.54	0.52	0.55
1	0.521	0.536	0.53	0.51	0.533
4	0.546	0.538	0.561	0.551	0.564
7	0.542	0.564	0.56	0.533	0.55
14	0.524	0.55	0.575	0.549	0.566
21	0.547	0.558	0.57	0.552	0.572
28	0.547	0.57	0.569	0.559	0.582
61	0.578	0.579	0.594	0.597	0.61

COMPLETED DAY 61 ONLY

SUBJECT # 107

EYE: OD

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.562	0.577	0.581	0.589	0.594
1	0.56	0.596	0.599	0.606	0.614
4	0.604	0.612	0.615	0.629	0.646
7	0.529	0.561	0.555	0.586	0.577
14	0.548	0.577	0.547	0.571	0.557
22	0.557	0.557	0.537	0.554	0.552
29	0.537	0.556	0.551	0.56	0.542
60	0.564	0.585	0.565	0.582	0.591
88	0.565	0.586	0.568	0.575	0.588
120	0.533	0.551	0.544	0.558	0.548
150	0.579	0.617	0.599	0.604	0.593

SUBJECT # 107

EYE: OS

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.542	0.58	0.561	0.556	0.593
1	0.585	0.607	0.62	0.582	0.601
4	0.619	0.62	0.618	0.608	0.621
7	0.54	0.573	0.565	0.554	0.553
14	0.543		0.555	0.547	0.558
22	0.539	0.528	0.553	0.549	0.563
29	0.549	0.563	0.559	0.526	0.553
60	0.559	0.572	0.586	0.566	0.593
88	0.56	0.563	0.564	0.548	0.555
120	0.543	0.554	0.565	0.548	0.555
150	0.593	0.598	0.6	0.589	0.596

SUBJECT # 108

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.545	0.559	0.553	0.562	0.573
1	0.589	0.574	0.595	0.603	0.605
4	0.588	0.605	0.629	0.611	0.618
7	0.579	0.602	0.605	0.61	0.61
13	0.611	0.619	0.647	0.642	0.647
20	0.6	0.606	0.625	0.623	0.62
33	0.617	0.648	0.645	0.654	0.637
61	0.589	0.606	0.614	0.624	0.611
90	0.571	0.596	0.595	0.585	0.614
121	0.534	0.549	0.557	0.555	0.583
149	0.59	0.607	0.601	0.605	0.615

SUBJECT # 108

EYE: OS

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.538	0.574	0.563	0.564	0.571
1	0.571	0.592	0.58	0.568	0.595
4	0.578	0.593	0.608	0.6	0.619
7	0.581	0.613	0.593	0.595	0.605
13	0.602	0.609	0.615	0.609	0.62
20	0.6	0.613	0.612	0.628	0.634
33	0.603	0.634	0.632	0.643	0.641
61	0.606	0.619	0.621	0.624	0.617
90	0.565	0.587	0.597	0.587	0.59
121	0.548	0.572	0.566	0.552	0.568
149	0.6	0.613	0.605	0.583	0.624



SUBJECT # 109

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.568	0.58	0.591	0.591	0.581
1	0.522	0.536	0.531	0.552	0.569
4	0.53	0.57	0.56	0.596	0.581
7	0.522	0.57	0.556	0.586	0.575
14	0.528	0.552	0.542	0.55	0.557
21					
28	0.493	0.534	0.526	0.549	0.538
60	0.558	0.58	0.574	0.566	0.597
90	0.543	0.573	0.573	0.592	0.603
125	0.514	0.556	0.54	0.56	0.57
150	0.53	0.549	0.549	0.559	0.566

SUBJECT # 109

EYE: OS

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.563	0.581	0.587	0.571	0.598
1	0.562	0.569	0.552	0.548	0.552
4	0.563	0.582	0.588	0.566	0.552
7	0.563	0.591	0.599	0.552	0.573
14	0.527	0.556	0.565	0.543	0.565
21					
28	0.521	0.533	0.532	0.521	0.534
60	0.596	0.588	0.601	0.59	0.572
90	0.584	0.569	0.589	0.569	0.604
125	0.551	0.568	0.568	0.557	0.572
150	0.529	0.556	0.554	0.548	0.547

SUBJECT # 110

EYE: OD

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.472	0.518	0.493	0.528	0.532
1	0.499	0.529	0.534	0.53	0.551
4	0.522	0.557	0.551	0.595	0.604
7	0.51	0.554	0.56	0.569	0.584
14	0.517	0.534	0.543	0.559	0.539
21	0.545	0.553	0.556	0.579	0.567
28	0.551	0.568	0.563	0.575	0.592
62	0.54	0.567	0.563	0.573	0.577
89	0.495	0.532	0.52	0.525	0.542

COMPLETED DAY 89 ONLY

SUBJECT # 110

EYE: OS

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.498	0.509	0.493	0.507	0.52
1	0.525	0.547	0.555	0.545	0.539
4	0.559	0.551	0.552	0.556	0.573
7	0.543	0.579	0.566	0.55	0.573
14	0.544	0.55	0.535	0.536	0.538
21	0.553	0.573	0.574	0.548	0.573
28	0.557	0.575	0.59	0.556	0.583
62	0.561	0.576	0.59	0.561	0.565
89	0.504	0.529	0.541	0.518	0.541

COMPLETED DAY 89 ONLY

SUBJECT # 111

EYE: OD

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.512	0.542	0.521	0.546	0.566
1	0.536	0.554	0.557	0.6	0.587
4	0.555	0.583	0.584	0.611	0.604
7	0.546	0.579	0.57	0.589	0.605
13	0.597	0.634	0.603	0.648	0.65
21	0.588	0.609	0.609	0.623	0.639
28	0.57	0.594	0.571	0.593	0.609
61	0.576	0.611	0.588	0.613	0.619
90	0.555	0.562	0.563	0.554	0.565
124	0.577	0.589	0.61	0.619	0.606
158	0.533	0.545	0.54	0.554	0.556

SUBJECT # 111

EYE: OS

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.522	0.532	0.535	0.541	0.551
1	0.553	0.563	0.579	0.561	0.572
4	0.565	0.594	0.607	0.575	0.604
7	0.574	0.593	0.609	0.585	0.603
13	0.615	0.648	0.657	0.638	0.651
21	0.593	0.622	0.628	0.614	0.64
28	0.59	0.606	0.616	0.597	0.609
61	0.604	0.604	0.638	0.617	0.639
90	0.56	0.586	0.584	0.571	0.59
124	0.599	0.604	0.605	0.596	0.61
158	0.533	0.542	0.559	0.559	0.547

SUBJECT # 112

EYE: OD

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.613	0.621	0.615	0.637	0.62
1	0.578	0.608	0.595	0.615	0.635
4	0.637	0.675	0.668	0.684	0.706
8	0.585	0.638	0.609	0.645	0.654
14	0.572	0.592	0.59	0.593	0.6
21	0.588	0.604	0.597	0.61	0.617
28	0.603	0.631	0.638	0.639	0.651
61	0.583	0.602	0.601	0.61	0.598
90	0.561	0.575	0.583	0.582	0.591
123	0.559	0.589	0.603	0.586	0.605
151	0.577	0.595	0.61	0.617	0.619

SUBJECT # 112

EYE: OS

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.596	0.633	0.621	0.606	0.627
1	0.607	0.62	0.615	0.614	0.606
4	0.657	0.683	0.706	0.678	0.688
8	0.613	0.636	0.63	0.629	0.629
14	0.58	0.599	0.61	0.593	0.594
21	0.599	0.602	0.62	0.604	0.612
28	0.629	0.634	0.653	0.636	0.657
61	0.602	0.61	0.626	0.627	0.621
90	0.558	0.589	0.61	0.583	0.607
123	0.586	0.619	0.608	0.603	0.62
151	0.588	0.609	0.612	0.597	0.624

SUBJECT # 113

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.468	0.487	0.502	0.505	0.514
1	0.521	0.537	0.551	0.541	0.548
4					
7	0.511	0.514	0.528	0.551	0.547
14	0.517	0.509	0.515	0.52	0.513
21	0.482	0.516	0.513	0.512	0.509
29	0.511	0.527	0.518	0.534	0.534
60	0.56	0.557	0.554	0.574	0.572
90	0.542	0.566	0.548	0.558	0.555
118	0.504	0.525	0.54	0.543	0.53
151	0.505	0.52	0.53	0.526	0.525

SUBJECT # 113

EYE: OS

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.507	0.495	0.515	0.501	0.534
1	0.525	0.528	0.546	0.536	0.553
4					
7	0.524	0.54	0.551	0.536	0.549
14	0.506	0.533	0.532	0.51	0.53
21	0.471	0.478	0.482	0.483	0.492
29	0.511	0.529	0.532	0.541	0.537
60	0.56	0.573	0.563	0.567	0.569
90	0.545	0.549	0.563	0.556	0.56
118	0.535	0.536	0.546	0.527	0.552
151	0.506	0.524	0.535	0.502	0.524

SUBJECT # 114

EYE: OD

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.565	0.575	0.554	0.591	0.587
1	0.546	0.584	0.575	0.596	0.605
4	0.6	0.601	0.602	0.622	0.631
7	0.582	0.61	0.61	0.617	0.634
14	0.612	0.619	0.61	0.654	0.631
21	0.59	0.606	0.598	0.617	0.632
28	0.564	0.569	0.581	0.589	0.608

DROPPED - DUE TO PREGNANCY - DAY 28

SUBJECT # 114

EYE: OS

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.561	0.572	0.591	0.584	0.576
1	0.584	0.59	0.601	0.589	0.602
4	0.597	0.602	0.633	0.618	0.621
7	0.606	0.605	0.626	0.605	0.634
14	0.603	0.625	0.627	0.611	0.629
21	0.601	0.608	0.634	0.612	0.631
28	0.582	0.593	0.615	0.579	0.607

DROPPED - DUE TO PREGNANCY - DAY 28

SUBJECT # 115

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.461	0.487	0.482	0.498	0.511
1	0.446	0.487	0.482	0.512	0.498
4	0.47	0.495	0.503	0.519	0.529
7	0.504	0.496	0.487	0.501	0.52
14	0.452	0.459	0.467	0.497	0.484
21	0.437	0.491	0.491	0.506	0.507
28	0.479	0.51	0.489	0.519	0.508
58	0.476	0.501	0.508	0.515	0.541
86	0.502	0.527	0.518	0.519	0.54
120	0.488	0.512	0.506	0.523	0.525
150	0.483	0.503	0.481	0.5	0.519

SUBJECT # 115

EYE: OS

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.494	0.522	0.523	0.504	0.525
1	0.501	0.484	0.516	0.504	0.533
4	0.498	0.533	0.549	0.519	0.53
7	0.506	0.521	0.529	0.518	0.53
14	0.494	0.497	0.519	0.515	0.524
21	0.484	0.504	0.533	0.501	0.522
28	0.499	0.526	0.52	0.514	0.516
58	0.569	0.527	0.549	0.512	0.549
86	0.524	0.536	0.541	0.524	0.534
120	0.496	0.52	0.539	0.501	0.541
150	0.505	0.512	0.504	0.492	0.508

SUBJECT # 116

EYE: OD

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:

DROPPED - INABILITY TO WEAR EW

SUBJECT # 116

EYE: OS

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:

DROPPED - INABILITY TO WEAR EW



SUBJECT # 117

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.536	0.571	0.551	0.589	0.588
1	0.544	0.606	0.577	0.593	0.586
4	0.574	0.608	0.591	0.623	0.635
7	0.559	0.581	0.588	0.576	0.607
14	0.555	0.59	0.562	0.593	0.577
21	0.544	0.571	0.543	0.586	0.578
28	0.539	0.586	0.553	0.573	0.566
60	0.52	0.558	0.549	0.552	0.568
88	0.54	0.561	0.557	0.566	0.552

COMPLETED DAY 88 ONLY

SUBJECT # 117

EYE: OS

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.561	0.573	0.584	0.541	0.573
1	0.568	0.595	0.593	0.582	0.615
4	0.581	0.618	0.63	0.6	0.621
7	0.569	0.601	0.605	0.564	0.606
14	0.559	0.578	0.587	0.57	0.575
21	0.553	0.568	0.596	0.55	0.569
28	0.546	0.572	0.607	0.558	0.575
60	0.541	0.563	0.576	0.537	0.55
88	0.521	0.565	0.576	0.535	0.557

COMPLETED DAY 88 ONLY

SUBJECT # 118

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.549	0.556	0.557	0.573	0.582
1	0.554	0.586	0.564	0.588	0.615
3	0.541	0.601	0.583	0.612	0.621
7	0.558	0.579	0.561	0.596	0.59
14	0.553	0.617	0.573	0.623	0.6
21	0.574	0.614	0.596	0.638	0.62
28	0.579	0.622	0.586	0.615	0.622
62	0.62	0.639	0.606	0.622	0.648
91	0.561	0.585	0.583	0.602	0.603
126	0.554	0.591	0.561	0.602	0.586
150	0.573	0.595	0.586	0.604	0.614

SUBJECT # 118

EYE: OS

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.538	0.543	0.553	0.537	0.557
1	0.58	0.58	0.62	0.576	0.594
3	0.574	0.586	0.621	0.571	0.606
7	0.563	0.583	0.602	0.555	0.581
14	0.576	0.606	0.636	0.567	0.602
21	0.599	0.619	0.637	0.587	0.631
28	0.609	0.622	0.63	0.593	0.615
62	0.602	0.627	0.626	0.591	0.625
91	0.555	0.585	0.587	0.572	0.577
126	0.551	0.58	0.587	0.549	0.606
150	0.58	0.615	0.626	0.598	0.607

SUBJECT # 119

EYE: OD

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.566	0.601	0.56	0.571	0.607
1	0.55	0.597	0.573	0.574	0.627
4	0.525	0.586	0.553	0.573	0.595
7	0.552	0.565	0.563	0.591	0.583
14	0.515	0.564	0.565	0.575	0.584
21		0.544	0.501	0.527	0.535
28	0.544	0.548	0.551	0.559	0.572
60	0.539	0.556	0.554	0.567	0.566
90	0.554	0.592	0.575	0.574	0.575
118	0.602	0.624	0.622	0.643	0.641
151	0.518	0.539	0.528	0.535	0.565

SUBJECT # 119

EYE: OS

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.559	0.578	0.614	0.558	0.581
1	0.567	0.576	0.623	0.577	0.627
4	0.551	0.568	0.583	0.551	0.595
7	0.56	0.586	0.599	0.559	0.605
14	0.554	0.571	0.577	0.569	0.572
21	0.496		0.521	0.52	0.537
28	0.541	0.567	0.574	0.561	0.577
60	0.548	0.564	0.582	0.566	0.591
90	0.575	0.584	0.6	0.567	0.603
118	0.616	0.633	0.645	0.63	0.634
151	0.539	0.541	0.554	0.53	0.547

SUBJECT # 120

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.57	0.585	0.564	0.604	0.587
1	0.516	0.541	0.533	0.597	0.579
4	0.524	0.554	0.526	0.575	0.565
7	0.532	0.538	0.523	0.56	0.555
14	0.544	0.546	0.54	0.578	0.539
21	0.542	0.571	0.538	0.581	0.585
28	0.527	0.545	0.535	0.586	0.565
59	0.542	0.552	0.547	0.571	0.56
88	0.522	0.521	0.539	0.557	0.558
118	0.525	0.557	0.55	0.556	0.567
146	0.523	0.547	0.545	0.548	0.554

SUBJECT # 120

EYE: OS

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.57	0.6	0.614	0.57	0.596
1	0.524	0.523	0.564	0.516	0.533
4	0.55	0.539	0.558	0.517	0.555
7	0.512	0.516	0.571	0.529	0.538
14	0.534	0.537	0.552	0.524	0.546
21	0.547	0.565	0.6	0.536	0.569
28	0.552	0.547	0.579	0.511	0.565
59	0.552	0.541	0.554	0.537	0.567
88	0.531	0.565	0.592	0.551	0.571
118	0.523	0.553	0.571	0.549	0.568
146	0.542	0.555	0.567	0.519	0.551

SUBJECT # 121

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.556	0.546	0.549	0.575	0.575
1	0.545	0.583	0.57	0.589	0.596
4	0.534	0.538	0.545	0.587	0.598
7	0.516	0.521	0.524	0.578	0.545
14	0.525	0.547	0.513	0.553	0.587
21	0.515	0.533	0.539	0.547	0.559
28	0.545	0.552	0.558	0.562	0.591
60	0.54	0.552	0.557	0.563	0.6
87	0.551	0.562	0.561	0.591	0.597
122	0.491	0.505	0.515	0.537	0.518
150	0.54	0.556	0.551	0.572	0.589

SUBJECT # 121

EYE: OS

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.551	0.557	0.564	0.55	0.588
1	0.552	0.586	0.563	0.58	0.63
4	0.551	0.56	0.588	0.563	0.598
7	0.526	0.523	0.54	0.527	0.55
14	0.543	0.553	0.553	0.543	0.571
21	0.53	0.54	0.557	0.538	0.565
28	0.547	0.552	0.566	0.563	0.584
60	0.55	0.547	0.579	0.562	0.601
87	0.561	0.557	0.567	0.562	0.594
122	0.504	0.501	0.524	0.519	0.52
150	0.561	0.544	0.561	0.562	0.558

SUBJECT # 122

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.581	0.603	0.614	0.628	0.628
1	0.579	0.604	0.592	0.627	0.652
4	0.569	0.619	0.598	0.62	0.643
7	0.559	0.596	0.598	0.606	0.614
14	0.616	0.65	0.641	0.661	0.672
22	0.61	0.652	0.62	0.636	0.651
28	0.596	0.629	0.612	0.647	0.665
62	0.623	0.656	0.65	0.668	0.674
88	0.616	0.648	0.647	0.65	0.654
120	0.614	0.631	0.646	0.656	0.637
142	0.579	0.628	0.604	0.621	0.642

SUBJECT # 122

EYE: OS

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.584	0.588	0.617	0.595	0.607
1	0.595	0.612	0.631	0.6	0.611
4	0.592	0.617	0.612	0.59	0.605
7	0.595	0.6	0.608	0.572	0.608
14	0.651	0.652	0.666	0.651	0.673
22	0.623	0.645	0.666	0.627	0.652
28	0.621	0.632	0.644	0.613	0.651
62	0.639	0.658	0.67	0.656	0.673
88	0.629	0.648	0.659	0.645	0.659
120	0.628	0.655	0.659	0.634	0.647
142	0.604	0.621	0.65	0.608	0.637

SUBJECT # 123

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.484	0.531	0.502	0.53	0.52
1	0.495	0.539	0.514	0.531	0.545
4	0.486	0.526	0.508	0.538	0.564
7	0.507	0.549	0.523	0.539	0.563
15	0.499	0.532	0.513	0.54	0.555
22	0.514	0.536	0.539	0.578	0.573
29	0.504	0.55	0.521	0.521	0.547
60	0.509	0.532	0.522	0.529	0.546
92	0.525	0.55	0.523	0.522	0.546
120	0.482	0.513	0.507	0.511	0.508
150	0.507	0.529	0.515	0.537	0.521

SUBJECT # 123

EYE: OS

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.48	0.511	0.531	0.506	0.515
1	0.499	0.526	0.533	0.48	0.515
4	0.501	0.52	0.523	0.495	0.519
7	0.496	0.531	0.55	0.505	0.535
15	0.522	0.545	0.563	0.511	0.547
22	0.556	0.552	0.573	0.539	0.573
29	0.52	0.519	0.546	0.506	0.526
60	0.524	0.539	0.552	0.522	0.551
92	0.516	0.521	0.546	0.518	0.542
120	0.478	0.506	0.513	0.474	0.485
150	0.489	0.515	0.528	0.482	0.535

SUBJECT # 124

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.567	0.612	0.622	0.614	0.634
1	0.561	0.622	0.596	0.585	0.602
4	0.553	0.569	0.559	0.579	0.569
7	0.537	0.565	0.548	0.56	0.542
14	0.578	0.598	0.58	0.595	0.599
21	0.558	0.557	0.558	0.565	0.569
28	0.554	0.569	0.588	0.586	0.591
59	0.573	0.586	0.582	0.592	0.597

COMPLETED DAY 59 ONLY

SUBJECT # 124

EYE: OS

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.595	0.606	0.604	0.62	0.622
1	0.571	0.576	0.593	0.581	0.585
4	0.57	0.56	0.558	0.563	0.559
7	0.559	0.545	0.552	0.566	0.567
14	0.566	0.591	0.586	0.574	0.612
21	0.539	0.549	0.561	0.552	0.57
28	0.569	0.581	0.595	0.59	0.59
59	0.588	0.583	0.588	0.598	0.593

COMPLETED DAY 59 ONLY



SUBJECT # 125

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.531	0.533	0.526	0.536	0.552
1	0.484	0.524	0.49	0.527	0.52
3	0.48	0.509	0.523	0.523	0.532
7	0.504	0.541	0.522	0.546	0.553
14	0.484	0.514	0.507	0.522	0.546
21	0.516	0.535	0.524	0.544	0.543
28	0.57	0.595	0.575	0.587	0.604
61	0.517	0.552	0.543	0.561	0.558

COMPLETED DAY 61 ONLY

SUBJECT # 125

EYE: OS

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.537	0.555	0.551	0.529	0.583
1	0.521	0.496	0.516	0.509	0.513
3	0.522	0.55	0.566	0.509	0.544
7	0.539	0.553	0.576	0.535	0.568
14	0.51	0.506	0.526	0.508	0.53
21	0.537	0.555	0.551	0.525	0.55
28	0.563	0.581	0.609	0.569	0.584
61	0.577	0.555	0.567	0.55	0.57

COMPLETED DAY 61 ONLY

SUBJECT # 126

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.495	0.499	0.479	0.506	0.52
1	0.518	0.529	0.513	0.54	0.512
4	0.48	0.517	0.502	0.511	0.519
8	0.511	0.547	0.522	0.559	0.539
14	0.502	0.514	0.495	0.518	0.517
21	0.497	0.519	0.507	0.527	0.523
28	0.492	0.518	0.502	0.528	0.509
59	0.514	0.553	0.528	0.553	0.551
91	0.524	0.566	0.542	0.563	0.541
115	0.483	0.541	0.503	0.526	0.552
154	0.503	0.546	0.507	0.532	0.522

SUBJECT # 126

EYE: OS

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.469	0.498	0.519	0.481	0.5092
1	0.528	0.529	0.547	0.535	0.551
4	0.497	0.501	0.523	0.502	0.515
8	0.511	0.521	0.536	0.533	0.541
14	0.514	0.512	0.524	0.524	0.527
21	0.517	0.51	0.523	0.489	0.509
28	0.494	0.521	0.533	0.501	0.506
59	0.521	0.538	0.554	0.536	0.552
91	0.533	0.569	0.582	0.548	0.581
115	0.518	0.53	0.541	0.506	0.541
154	0.521	0.539	0.551	0.524	0.547

SUBJECT # 127

EYE: OD

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:

DROPPED FROM STUDY DUE TO PREGNANCY

SUBJECT # 127

EYE: OS

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:

DROPPED FROM STUDY DUE TO PREGNANCY

SUBJECT # 128

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:

DROPPED - UNABLE TO CONFORM TO GUIDELINES OF STUDY

SUBJECT # 128

EYE: OS

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:

DROPPED - UNABLE TO CONFORM TO GUIDELINES OF STUDY

SUBJECT # 129

EYE: OD

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.52	0.541	0.525	0.537	0.54
1	0.502	0.538	0.541	0.562	0.558
4	0.538	0.541	0.544	0.553	0.568
7	0.55	0.571	0.551	0.566	0.568
14	0.536	0.563	0.547	0.571	0.578
21	0.566	0.592	0.573	0.605	0.606
28	0.546	0.574	0.572	0.604	0.59
62	0.511	0.532	0.517	0.545	0.551
90	0.52	0.542	0.543	0.549	0.549
120	0.527	0.553	0.549	0.551	0.561
147	0.519	0.53	0.522	0.535	0.535

SUBJECT # 129

EYE: OS

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.505	0.537	0.536	0.517	0.525
1	0.534	0.543	0.557	0.538	0.551
4	0.514	0.552	0.561	0.555	0.556
7	0.546	0.574	0.572	0.569	0.569
14	0.53	0.555	0.573	0.537	0.573
21	0.559	0.57	0.593	0.583	0.6
28	0.562	0.561	0.58	0.589	0.591
62	0.515	0.524	0.55	0.529	0.548
90	0.528	0.54	0.563	0.539	0.548
120	0.533	0.549	0.555	0.547	0.568
147	0.52	0.515	0.552	0.517	0.554

SUBJECT # 130

EYE: OD

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.544	0.494	0.548	0.575	0.563
1	0.513	0.514	0.548	0.554	0.59
4	0.534	0.534	0.566	0.588	0.574
7	0.54	0.564	0.579	0.585	0.598
14	0.49	0.496	0.517	0.545	0.55
21	0.462	0.49	0.503	0.508	0.532
29	0.521	0.524	0.535	0.548	0.535
60	0.547	0.554	0.568	0.586	0.586
88	0.556	0.574	0.584	0.608	0.58
118	0.527	0.545	0.544	0.558	0.545
151	0.498	0.5	0.534	0.531	0.524

SUBJECT # 130

EYE: OS

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.526	0.584	0.591	0.558	0.585
1	0.513	0.54	0.576	0.569	0.554
4	0.542	0.566	0.57	0.592	0.589
7	0.541	0.564	0.56	0.567	0.571
14	0.513	0.508	0.53	0.509	0.52
21	0.497	0.471	0.5	0.503	0.53
29	0.509	0.512	0.547	0.546	0.564
60	0.564	0.564	0.581	0.579	0.573
88	0.576	0.585	0.593	0.592	0.575
118	0.521	0.552	0.558	0.527	0.56
151	0.515	0.522	0.552	0.523	0.545

SUBJECT # 131

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.52	0.525	0.529	0.545	0.533
1	0.576	0.604	0.585	0.603	0.589
4	0.61	0.633	0.609	0.626	0.616
7	0.597	0.622	0.623	0.619	0.631

DROPPED - UNABLE TO WEAR EW

SUBJECT # 131

EYE: OS

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.529	0.524	0.544	0.515	0.545
1	0.571	0.581	0.578	0.562	0.585
4	0.576	0.585	0.582	0.586	0.589
7	0.565	0.573	0.581	0.576	0.589

DROPPED - UNABLE TO WEAR EW

SUBJECT # 132

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.544	0.574	0.576	0.58	0.583
1	0.556	0.593	0.607	0.608	0.61
4	0.577	0.597	0.59	0.603	0.593
7	0.561	0.567	0.572	0.585	0.571
14	0.587	0.621	0.58	0.604	0.595
21	0.57	0.607	0.599	0.616	0.609
28	0.579	0.603	0.58	0.591	0.593
60	0.532	0.56	0.561	0.548	0.549
90	0.569	0.605	0.593	0.611	0.593
118	0.569		0.577	0.608	0.568
147	0.553	0.589	0.575	0.579	0.584

SUBJECT # 132

EYE: OS

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.532	0.568	0.576	0.567	0.58
1	0.586	0.596	0.599	0.596	0.605
4	0.567	0.591	0.6	0.576	0.57
7	0.565	0.575	0.591	0.564	0.574
14	0.574	0.592	0.619	0.572	0.589
21	0.597	0.61	0.615	0.605	0.594
28	0.594	0.59	0.61	0.594	0.595
60	0.542	0.554	0.565	0.554	0.529
90	0.563	0.594	0.611	0.59	0.58
118	0.547	0.556	0.585	0.58	0.578
147	0.557	0.569	0.594	0.566	0.568



SUBJECT # 133

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.537	0.555	0.575	0.598	0.571
1	0.519	0.562	0.548	0.583	0.601
4	0.538	0.547	0.521	0.571	0.571
7					
14	0.535	0.573	0.561	0.599	0.578
21	0.532	0.568	0.524	0.569	0.547
27	0.53	0.56	0.562	0.587	0.576
60	0.535	0.565	0.564	0.57	0.593
94	0.5	0.546	0.547	0.548	0.553
118	0.483	0.514	0.51	0.526	0.531
147	0.531	0.554	0.545	0.578	0.567

SUBJECT # 133

EYE: OS

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.548	0.561	0.571	0.541	0.59
1		0.543			
4	0.523	0.548	0.539	0.555	0.56
7					
14	0.544	0.571	0.564	0.559	0.576
21	0.54	0.544	0.543	0.532	0.528
27	0.541	0.586	0.564	0.561	0.556
60	0.553	0.582	0.577	0.555	0.591
94	0.53	0.548	0.545	0.531	0.576
118	0.515	0.54	0.533	0.527	0.54
147	0.535	0.558	0.569	0.534	0.567

SUBJECT # 134

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.465	0.495	0.491	0.519	0.502
1	0.472	0.502	0.455	0.506	0.497
4	0.469	0.524	0.516	0.547	0.55
7	0.468	0.493	0.476	0.515	0.503
14	0.581	0.623	0.619	0.594	0.565
21	0.492	0.527	0.491	0.518	0.523
28	0.482	0.514	0.494	0.53	0.539
60	0.494	0.535	0.522	0.53	0.534
102	0.462	0.501	0.496	0.527	0.506
123	0.474	0.508	0.493	0.508	0.494
145	0.477	0.49	0.489	0.486	0.489

SUBJECT # 134

EYE: OS

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.464	0.499	0.503	0.465	0.5
1	0.458	0.489	0.511	0.471	0.496
4	0.49	0.527	0.559	0.5	0.533
7	0.467	0.519	0.518	0.489	0.517
14	0.558	0.592	0.628	0.582	0.598
21	0.5	0.504	0.519	0.504	0.516
28	0.494	0.513	0.518	0.511	0.515
60	0.512	0.517	0.536	0.504	0.527
102	0.471	0.507	0.517	0.478	0.495
123	0.469	0.503	0.509	0.477	0.505
145	0.465	0.496	0.507	0.47	0.511

SUBJECT # 135

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.493	0.543	0.551	0.537	0.547
1	0.526	0.577	0.56	0.575	0.587
4	0.522	0.543	0.53	0.548	0.544
7	0.515	0.552	0.561	0.53	0.573
14	0.476	0.513	0.502	0.521	0.534
21	0.527	0.543	0.554	0.555	0.579
28	0.572	0.568	0.552	0.592	0.611
60	0.559	0.599	0.584	0.62	0.625
90	0.561	0.577	0.576	0.591	0.6
120	0.52	0.535	0.548	0.537	0.56
151	0.525	0.543	0.555	0.563	0.562

SUBJECT # 135

EYE: OS

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.518	0.536	0.564	0.535	0.565
1	0.539	0.583	0.599	0.551	0.611
4	0.526	0.54	0.537	0.534	0.565
7	0.545	0.543	0.567	0.546	0.572
14	0.5	0.5	0.511	0.496	0.53
21	0.526	0.55	0.567	0.55	0.576
28	0.572	0.596	0.604	0.575	0.613
60	0.58	0.56	0.595	0.588	0.598
90	0.566	0.574	0.587	0.576	0.594
120	0.527	0.529	0.566	0.546	0.568
151	0.513	0.525	0.555	0.521	0.551

SUBJECT # 136

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.561	0.59	0.586	0.584	0.601
1	0.528	0.535	0.523	0.532	0.559
4	0.527	0.573	0.55	0.602	0.588
8	0.527	0.538	0.549	0.589	0.562
14	0.551	0.561	0.546	0.583	0.576
21					
28	0.542	0.553	0.563	0.597	0.583
61	0.566	0.587	0.57	0.581	0.587
89	0.533	0.564	0.551	0.559	0.56
121	0.539	0.554	0.547	0.568	0.584
149	0.566	0.576	0.567	0.582	0.579

SUBJECT # 136

EYE: OS

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.569	0.59	0.61	0.575	0.611
1	0.517	0.554	0.568	0.533	0.552
4	0.528	0.555	0.562	0.554	0.572
8	0.537	0.558	0.578	0.532	0.556
14	0.544	0.562	0.58	0.541	0.569
21					
28	0.547	0.553	0.588	0.558	0.571
61	0.548	0.588	0.607	0.557	0.599
89	0.542	0.558	0.574	0.552	0.568
121	0.527	0.556	0.565	0.539	0.56
149	0.547	0.566	0.586	0.558	0.581

SUBJECT # 137

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.465	0.544	0.512	0.512	0.498
1	0.449	0.505	0.493	0.515	0.518
4	0.507	0.54	0.537	0.563	0.56
7	0.458	0.507	0.494	0.533	0.494
14					
21	0.444	0.473	0.488	0.49	0.476
28	0.45	0.501	0.483	0.504	0.494
61	0.513	0.544	0.529	0.562	0.547
91	0.484	0.505	0.503	0.53	0.523
119	0.486	0.516	0.485	0.52	0.492
149	0.507	0.537	0.517	0.523	0.53

SUBJECT # 137

EYE: OS

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.485	0.52	0.51	0.49	0.494
1	0.52	0.55	0.568	0.539	0.55
4	0.567	0.587	0.58	0.553	0.534
7	0.508	0.51	0.542	0.515	0.539
14					
21	0.493	0.514	0.525	0.502	0.513
28	0.474	0.491	0.507	0.458	0.499
61	0.53	0.544	0.563	0.529	0.529
91	0.496	0.519	0.537	0.491	0.506
119	0.48	0.497	0.537	0.496	0.527
149	0.513	0.53	0.543	0.5	0.54

SUBJECT # 138

EYE: OD

LENS TYPE PF

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.524	0.572	0.547	0.583	0.578
1	0.557	0.57	0.555	0.569	0.574
4	0.544	0.568	0.551	0.576	0.57
6	0.565	0.582	0.559	0.582	0.582
13	0.591	0.614	0.608	0.645	0.637
20	0.587	0.598	0.587	0.603	0.599
27	0.605	0.624	0.591	0.632	0.624
60	0.564	0.59	0.559	0.574	0.583

COMPLETED DAY 60 ONLY

SUBJECT # 138

EYE: OS

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.549	0.573	0.537	0.55	0.55
1	0.55	0.582	0.587	0.562	0.57
4	0.559	0.581	0.596	0.554	0.574
6	0.564	0.605	0.604	0.576	0.574
13	0.608	0.655	0.625	0.619	0.628
20	0.579	0.6	0.602	0.612	0.593
27	0.595	0.616	0.616	0.595	0.609
60	0.573	0.587	0.593	0.57	0.588

COMPLETED DAY 60 ONLY

SUBJECT # 139

EYE: OD

LENS TYPE HC II

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.483	0.513	0.507	0.525	0.511
1	0.516	0.528	0.517	0.549	0.56
4	0.504	0.539	0.498	0.541	0.542
7					
14	0.473	0.534	0.499	0.538	0.529
21	0.512	0.53	0.507	0.537	0.534
28	0.523	0.543	0.532	0.56	0.556
59	0.525	0.54	0.521	0.548	0.548
87	0.52	0.539	0.524	0.54	0.532
120	0.503	0.54	0.53	0.548	0.545
151	0.464	0.483	0.483	0.508	0.513

SUBJECT # 139

EYE: OS

LENS TYPE CSI-T

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.476	0.507	0.518	0.505	0.512
1	0.506	0.563	0.539	0.519	0.529
4	0.509	0.525	0.536	0.503	0.522
7					
14	0.491		0.54		0.523
21	0.509	0.527	0.544	0.522	0.532
28	0.54	0.556	0.556	0.516	0.54
59	0.511	0.524	0.553	0.525	0.546
87	0.512	0.531	0.539	0.512	0.54
120	0.505	0.542	0.554	0.51	0.53
151	0.462	0.479	0.514	0.453	0.5

PATIENT NO: 301

EYE: OD

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.488	0.5	0.486	0.493	0.51
14	0.481	0.507	0.499	0.519	0.538
30	0.458	0.488	0.486	0.508	0.518
60	0.465	0.489	0.488	0.506	0.525
90	0.47	0.49	0.478	0.508	0.522
120	0.482	0.484	0.49	0.495	0.508

PATIENT NO: 301

EYE: OS

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.488	0.508	0.528	0.493	0.512
14	0.494	0.492	0.534	0.519	0.546
30	0.489	0.488	0.523	0.488	0.525
60	0.504	0.509	0.524	0.497	0.521
90	0.488	0.485	0.524	0.493	0.524
120	0.463	0.489	0.528	0.468	0.511



PATIENT NO: 303

EYE: OD

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.497	0.524	0.489	0.543	0.512
14	0.496	0.524	0.48	0.537	0.531
28	0.507	0.536	0.506	0.544	0.515
60	0.487	0.521	0.515	0.543	0.539
90	0.451	0.487	0.467	0.515	0.475
120	0.488	0.509	0.493	0.521	0.494

PATIENT NO: 303

EYE: OS

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.508	0.534	0.522	0.512	0.515
14	0.51	0.542	0.56	0.509	0.528
28	0.509	0.528	0.534	0.523	0.533
60	0.51	0.544	0.531	0.53	0.516
90	0.46	0.49	0.497	0.469	0.502
120	0.498	0.523	0.524	0.501	0.511

PATIENT NO: 304

EYE: OD

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.542	0.568	0.546	0.571	0.57
14	0.539	0.565	0.566	0.571	0.575
30	0.552	0.575	0.561	0.56	0.57
60	0.491	0.532	0.522	0.53	0.532
90	0.52	0.564	0.547	0.559	0.558
120	0.536	0.571	0.554	0.554	0.56

PATIENT NO: 304

EYE: OS

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.547	0.57	0.549	0.558	0.568
14	0.563	0.581	0.592	0.554	0.576
30	0.562	0.59	0.562	0.575	0.57
60	0.527	0.536	0.533	0.533	0.545
90	0.559	0.548	0.568	0.555	0.554
120	0.55	0.568	0.576	0.535	0.558

PATIENT NO: 305

EYE: OD

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.558	0.589	0.551	0.572	0.564
14	0.53	0.566	0.543	0.572	0.558
30					
60	0.516	0.546	0.542	0.564	0.551
90	0.512	0.552	0.533	0.565	0.577
120	0.522	0.569	0.546	0.563	0.569

PATIENT NO: 305

EYE: OS

DAY:	POINT Q:	POINT I:	POINT O:	POINT K:	POINT M:
0	0.557	0.584	0.594	0.571	0.599
14	0.55	0.567	0.574	0.544	0.575
30					
60	0.541	0.549	0.577	0.54	0.558
90	0.54	0.567	0.562	0.543	0.566
120	0.533	0.574	0.591	0.544	0.571

APPENDIX G  
PACHOMETRY SPREADSHEET, THICKNESS  
AND PERCENTAGE CHANGES





AW										AZ										BA										BB										BC										BD										BE										BF										BG										BH										BI										BJ										BK										BL										BM										BN										BO										BP										BQ										BR										BS										BT										BU										BV										BW										BX										BY										BZ										CA										CB										CC										CD										CE										CF										CG										CH										CI										CJ										CK										CL										CM										CN										CO																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																									
M90										M91										M92										M93										M94										M95										M96										M97										M98										M99										M00										M01										M02										M03										M04										M05										M06										M07										M08										M09										M10										M11										M12										M13										M14										M15										M16										M17										M18										M19										M20										M21										M22										M23										M24										M25										M26										M27										M28										M29										M30										M31										M32										M33										M34										M35										M36										M37										M38										M39										M40										M41										M42										M43										M44										M45										M46										M47										M48										M49										M50										M51										M52										M53										M54										M55										M56										M57										M58										M59										M60										M61										M62										M63										M64										M65										M66										M67										M68										M69										M70										M71										M72										M73										M74										M75										M76										M77										M78										M79										M80										M81										M82										M83										M84										M85										M86										M87										M88										M89										M90										M91										M92										M93										M94										M95										M96										M97										M98										M99										M00										M01										M02										M03										M04										M05										M06										M07										M08										M09										M10										M11										M12										M13										M14										M15										M16										M17										M18										M19										M20										M21										M22										M23										M24										M25										M26										M27										M28										M29										M30										M31										M32										M33										M34										M35										M36										M37										M38										M39										M40										M41										M42										M43										M44										M45										M46										M47										M48										M49										M50										M51										M52										M53										M54										M55										M56										M57										M58										M59										M60										M61										M62										M63										M64										M65										M66										M67										M68										M69										M70										M71										M72										M73										M74										M75										M76										M77										M78										M79										M80										M81										M82										M83										M84										M85										M86										M87										M88										M89										M90										M91										M92										M93										M94										M95										M96										M97										M98										M99										M00										M01										M02										M03										M04										M05										M06										M07										M08										M09										M10										M11										M12										M13										M14										M15										M16										M17										M18										M19										M20										M21										M22										M23										M24										M25										M26										M27										M28										M29										M30										M31										M32										M33										M34										M35										M36										M37										M38										M39										M40										M41										M42										M43										M44										M45										M46										M47										M48										M49										M50										M51										M52										M53										M54										M55										M56										M57										M58										M59										M60										M61										M62										M63										M64										M65										M66										M67										M68										M69										M70										M71										M72										M73										M74										M75										M76										M77										M78										M79										M80										M81										M82										M83										M84										M85										M86										M87										M88										M89										M90										M91										M92										M93										M94										M95										M96										M97										M98										M99										M00										M01										M02										M03										M04										M05										M06										M07										M08										M09										M10										M11										M12										M13										M14										M15										M16										M17										M18										M19										M20										M21										M22										M23										M24										M25										M26										M27										M28										M29										M30										M31										M32										M33										M34										M35										M36										M37										M38										M39										M40										M41										M42										M43										M44										M45										M46										M47										M48										M49										M50										M51										M52										M53										M54										M55										M56										M57										M58										M59										M60										M61										M62										M63										M64										M65										M66										M67										M68										M69										M70										M71										M72										M73										M74										M75										M76										M77										M78										M79										M80										M81										M82										M83										M84										M85										M86										M87										M88										M89										M90										M91										M92										M93										M94										M95										M96										M97										M98										M99										M00										M01										M02										M03										M04										M05										M06										M07										M08										M09										M10										M11										M12										M13										M14										M15										M16										M17										M18										M19										M20										M21										M22										M23										M24										M25										M26										M27										M28										M29										M30										M31										M32										M33										M34										M35										M36										M37										M38										M39										M40										M41										M42										M43										M44										M45										M46										M47										M48										M49										M50										M51										M52										M53										M54										M55										M56										M57										M58										M59										M60										M61										M62										M63										M64										M65										M66										M67										M68										M69										M70										M71										M72										M73										M74										M75										M76										M77										M78										M79										M80										M81										M82										M83										M84										M85										M86										M87										M88										M89										M90										M91										M92										M93										M94										M95										M96										M97										M98										M99										M00										M01										M02										M03										M04										M05										M06										M07										M08										M09										M10										M11										M12										M13										M14										M15										M16										M17										M18										M19										M20										M21										M22										M23										M24										M25										M26										M27										M28										M29										M30										M31										M32										M33										M34										M35										M36										M37										M38										M39										M40										M41										M42										M43										M44										M45										M46										M47										M48										M49										M50										M51										M52										M53										M54										M55										M56										M57										M58										M59										M60										M61										M62										M63										M64										M65										M66										M67										M68										M69										M70										M71										M72										M73										M74										M75										M76										M77										M78										M79										M80										M81										M82										M83										M84										M85										M86										M87										M88										M89										M90										M91										M92										M93										M94										M95										M96										M97										M98										M99										M00										M01										M02										M03										M04										M05										M06										M07										M08										M09										M10										M11										M12										M13										M14										M15										M16										M17										M18										M19										M20										M21										M22										M23										M24										M25										M26										M27										M28										M29										M30										M31										M32										M33										M34										M35										M36										M37										M38										M39										M40										M41										M42										M43										M44										M45										M46										M47										M48										M49										M50										M51										M52										M53										M54										M55										M56										M57										M58										M59										M60										M61										M62										M63										M64										M65										M66										M67										M68										M69										M70										M71										M72										M73										M74										M75										M76										M77										M78										M79										M80										M81										M82										M83										M84										M85										M86										M87										M88										M89										M90										M91										M92										M93										M94										M95										M96										M97										M98										M99										M00										M01										M02										M03										M04										M05										M06										M07										M08										M09										M10										M11										M12										M13										M14										M15										M16										M17										M18										M19										M20										M21										M22										M23										M24										M25										M26										M27										M28										M29										M30										M31										M32										M33										M34										M35										M36										M37										M38										M39										M40										M41										M42										M43										M44										M45										M46										M47										M48										M49										M50										M51										M52										M53										M54										M55										M56										M57										M58										M59										M60										M61										M62										M63										M64										M65										M66										M67										M68										M69										M70										M71										M72										M73										M74										M75										M76										M77										M78										M79										M80										M81										M82										M83										M84										M85										M86										M87										M88										M89										M90										M91										M92										M93										M94										M95										M96										M97										M98										M99										M00										M01										M02										M03										M04										M05										M06										M07										M08										M09										M10										M11										M12										M13										M14										M15										M16										M17										M18										M19										M20										M21										M22										M23										M24										M25										M26										M27										M28										M29										M30										M31										M32										M33										M34										M35										M36										M37										M38										M39										M40										M41										M42										M43										M44										M45										M46										M47										M48										M49										M50										M51										M52										M53										M54										M55										M56										M57										M58										M59										M60										M61										M62										M63										M64										M65										M66										M67										M68										M69										M70										M71										M72										M73										M74										M75										M76										M77										M78										M79										M80										M81										M82										M83										M84										M85										M86										M87										M88										M89										M90										M91										M92										M93										M94										M95										M96										M97										M98										M99										M00										M01										M02										M03										M04										M05										M06										M07										M08										M09										M10										M11										M12										M13										M14										M15										M16										M17										M18										M19										M20										M21										M22										M23										M24										M25										M26										M27										M28										M29										M30										M31										M32										M33										M34										M35										M36										M37										M38										M39										M40										M41										M42										M43										M44										M45										M46										M47										M48										M49										M50										M51										M52										M53										M54										M55										M56										M57										M58										M59										M60										M61										M62										M63										M64										M65										M66										M67										M68										M69										M70										M71										M72										M73										M74										M75										M76										M77										M78										M79										M80										M81										M82										M83										M84									

CN	CO	CP	CQ	CR	CS	CT	CU	CV	CW	CX	CY	CZ	DA	DB	DC	DD	DE	DF	DG	DH	DI	DJ	DK	DL	DM	DN	DO	DP	DQ	
0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	0% D21	
5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188	5.188		
3.571	2.613	3.627	0.178	-2.397	-2.500	0.348	-2.415	-0.176	2.478	3.854	0.000	1.294	-0.550	1.835	3.030	5.000	2.957	2.202	3.119	-1.010	1.852	-4.621	-10.275	-3.119	-2.424	-4.815	-8.654	-4.771	-4.404	
-4.455	-3.889	-6.285	-4.037	-2.202	-3.434	-3.704	-8.688	-3.853	-5.138	3.854	0.000	1.294	-0.550	1.835	3.030	5.000	2.957	2.202	3.119	-1.010	1.852	-4.621	-10.275	-3.119	-2.424	-4.815	-8.654	-4.771	-4.404	
-1.227	0.220	1.744	-2.703	-3.578	-4.284	1.202	2.713	-1.544	-0.585	3.067	6.813	8.721	3.089	2.637	6.953	13.828	10.078	5.792	5.481	5.776	0.000	-1.744	-7.722	-4.896	1.227	8.611	0.969	-4.440	0.160	
1.434	4.889	-0.585	0.753	-1.884	0.000	0.586	0.780	1.695	2.462	8.557	7.422	0.975	3.855	4.924	5.943	4.297	4.095	4.520	4.924	8.300	6.055	4.095	4.143	6.818	5.328	8.208	4.678	7.386		
-3.276	-6.250	-6.250	-6.250	-6.182	-3.083	-5.357	-3.697	3.854	2.321	0.528	2.312	1.250	0.528	2.312	1.250	0.528	2.312	1.250	0.528	2.312	1.250	0.528	2.312	1.250	0.528	2.312	1.250	0.528	2.312	
-0.749	0.719	-2.359	-3.908	-1.792	-4.307	-4.856	-1.452	0.355	-3.226	-1.873	0.719	1.270	0.808	0.538	-5.993	-2.338	-4.174	-3.908	-6.631	-3.933	-0.719	-3.811	-4.440	-4.122	-2.424	-4.815	-8.654	-4.771	-4.404	
0.749	0.000	-1.220	-5.302	-1.616	-6.000	-1.408	-1.220	-1.828	0.539	0.000	-0.352	-1.220	-0.366	0.539	-3.184	-5.810	-4.878	-3.473	-2.693	1.311	-0.409	-5.052	-5.302	-2.873	-2.424	-4.815	-8.654	-4.771	-4.404	
5.927	3.745	7.091	7.915	5.302	4.435	8.614	11.455	11.505	12.065	6.501	5.243	7.273	11.390	8.044	4.398	5.618	7.091	8.687	4.938	7.839	5.431	8.545	7.338	7.130	2.486	3.184	5.998	8.948	3.357	7.273
2.418	-2.113	5.222	8.801	3.791	2.974	-2.113	7.737	7.441	8.254	2.602	1.937	13.733	12.794	7.581																
5.395	1.270	5.556	6.154	4.000	5.395	3.448	8.970	7.500	5.818	11.368	5.082	10.000	14.808	10.909																
-0.890	-3.466	-7.573	-6.945	-7.071	-4.448	-6.540	-8.164	-4.804	-6.754	0.356	1.386	-2.754	-1.168	-0.505	0.534	1.562	-2.238	-2.377	-1.018	-5.160	-4.508	-4.365	-5.263	7.744	3.025	8.932	3.098	2.547	-0.168	
-0.854	-8.968	-1.428	-1.253	-5.059	1.232	-2.931	-2.937	-3.824	-6.745	3.137	-1.379	4.456	1.799	0.000	3.321	-2.931	0.535	-1.438	-5.408	0.185	-4.483	0.213	-1.439	-6.408	8.419	3.103	6.954	5.935	0.506	
10.092	8.408	13.023	10.854	8.202	13.211	15.921	16.681	16.370	11.182	8.073	8.408	11.031	11.032	6.832	4.771	8.619	7.585	4.083	7.155	-2.018	-1.780	0.723	-1.748	1.743	8.257	8.587	8.680	7.651	7.390	
11.524	6.794	8.722	11.344	11.033	12.082	10.453	12.256	14.507	12.259	12.038	7.840	10.302	10.638	8.058	5.019	2.268	6.038	4.078	3.327	1.859	-0.348	0.653	-2.128	-0.525	11.524	6.794	7.460	3.369	9.287	
					-13.204	-7.931	-10.999	-7.107	-7.401	-1.761	0.000	-2.678	-4.230	2.754	-4.401	-1.207	-3.048	0.169	3.787	-9.507	-4.138	-8.823	-5.245	-1.897	-6.690	-5.345	-7.107	-5.415	-2.587	
					-7.460	-8.282	-9.370	-8.757	-10.702	5.861	1.205	2.381	3.327	-4.948	3.730	-2.065	0.041	-0.350	1.003	-2.131	-2.238	-0.213	-2.452	-4.348	-6.033	-4.303	-5.622	-4.028	-8.520	
15.488	6.757	12.774	9.659	6.579	10.737	9.653	14.199	8.902	11.278	14.407	9.454	14.199	8.523	8.459	4.873	2.703	5.477	-0.568	1.880											
11.044	12.574	16.430	8.087	10.192	11.847	12.967	19.675	9.655	12.115	12.651	13.165	19.675	10.651	8.654	1.205	3.929	9.738	2.170	4.030											
14.844	12.982	9.891	14.103	12.898	11.328	8.594	8.597	8.608	7.597	12.503	12.731	12.860	12.271	9.365	8.398	8.890	8.061	1.485	-0.177	12.695	11.672	17.835	13.970	7.067	4.105	0.554	3.647	1.465	-1.767	
13.606	19.197	17.383	13.494	16.152	13.027	13.910	15.140	10.351	10.528	15.709	13.534	19.253	14.048	11.971	7.280	10.150	9.159	5.545	0.708	14.751	13.594	13.894	10.166	10.706	2.107	1.880	4.486	3.327	-0.728	
-4.078	-2.738	-2.927	-4.239	-0.484	-1.631	1.610	3.740	0.914	5.000	-4.894	-3.095	-2.276	-4.239	-3.545	-8.483	-7.407	-5.203	-8.634	-4.677	-8.800	-5.153	-1.847	-8.008	-2.419	-5.272	-4.187	-0.813	-3.140	-0.151	
0.503	-4.897	-0.161	-3.309	-2.892	5.537	0.158	5.159	4.950	4.785	1.007	-3.630	0.825	3.465	-0.957	-8.376	-6.951	-1.771	-3.798	-3.190	-1.670	-2.212	-2.093	-0.495	-1.115	-1.344	-3.791	-1.445	-0.478	-0.161	
2.981	5.955	2.181	1.305	-0.973	9.188	8.214	3.187	5.743	9.891	19.558	14.374	10.353	13.663	11.284	15.812	16.222	9.163	10.498	9.777	7.692	7.670	7.350	7.525	3.117	7.908	6.776	5.578	4.158	2.140	
-7.101	-3.434	-6.408	-3.393	-7.865	0.789	6.669	3.304	7.984	0.565	10.454	16.754	9.320	13.174	6.554	7.495	10.909	9.320	10.978	4.869	5.523	8.207	6.019	5.190	3.371	-0.197	5.858	3.883	0.200	-1.873	
4.423	5.391	7.942	4.360	7.668	-0.177	-1.043	4.874	-0.338	3.578																					
7.130	6.294	7.276	4.794	9.549	3.743	3.671	4.061	-0.856	5.382																					
-5.206	0.821	1.867	1.606	-0.783	3.905	4.723	1.452	4.217	-0.587	3.254	2.874	5.394	3.414	5.871	8.894	8.214	7.469	4.217	5.675	5.857	5.133	4.870	5.020	2.740	4.772	2.885	-0.207	0.407	1.566	
-2.024	-3.448	1.912	-0.595	-0.571	1.012	0.766	-0.574	1.984	-1.714	15.182	0.958	4.871	1.587	4.571	6.073	2.688	3.442	3.966	1.714	0.405	-0.389	3.043	-0.595	3.048	2.227	-1.818	-1.850	2.381	-3.238	
1.493	0.000	-1.452	-0.509	-1.701	0.560	2.627	8.970	-2.718	-3.741	-2.985	-2.277	-0.363	-6.282	-3.401	0.746	-1.751	1.089	-3.905	-6.122											
-1.426	-0.873	2.055	1.864	-0.698	-2.674	-0.175	8.970	3.142	0.258	-3.565	-1.745	-1.370	-0.739	-4.014	-7.130	-1.398	-1.970	-1.709	-2.708											
4.554	10.432	7.002	11.344	6.529	5.464	11.871	5.200	7.330	6.873	12.823	16.823	12.787	8.551	11.340	2.186	5.218	4.868	5.091	3.608	0.911	8.987	8.713	5.081	8.687	4.377	7.014	8.278	8.410	5.498	
11.339	9.996	15.190	9.311	13.285	13.197	14.540	18.054	10.428	10.413	15.808	15.474	13.301	10.058	12.988	3.160	7.795	8.148	6.516	3.591	2.418	6.814	6.148	2.235	8.792	2.607	13.260	13.201	11.598	8.977	
11.339	10.496	-10.596	-7.708	-11.862	-3.887	-8.819	-1.807	-2.102	-5.786	-4.770	-7.488	-1.071	-0.701	-6.755	-2.120	-1.498	2.679	0.625	-2.722	6.360	8.414	6.148	12.609	5.601	-8.481	-10.318	-5.714	-6.305	-6.919	
11.270	-15.147	-6.810	-7.973	-5.232	-2.201	-1.903	-8.511	0.538	-0.688	-1.188	-1.435	-5.414	1.434	1.721	2.622	1.038	-2.260	1.613	3.787	10.197	11.693	5.044	12.903	9.122	-3.578	-6.401	-9.772	-5.018	-5.852	
-4.912	-2.393	-4.610	-3.808	-0.341	7.544	-6.838	-8.747	-2.980	-3.748	-4.915	-6.841	-3.014	-5.464	-4.600	-8.421	-10.94	-4.433	-7.781	-4.940	-7.895	-7.789	-2.483	-7.947	-3.407	-8.246	-6.496	-3.369	-9.272	-5.622	
-1.070	-5.335	-0.500	-0.800	-0.840	8.970	8.970	8.970	-10.188	-6.301	-3.188	-8.288	-8.778	-8.708	-4.888	-8.888	-2.588	-2.588	-2.588	-2.588	-4.188	-8.288	-7.850	-7.070	-3.684	-4.998	-4.912	-7.500	-7.655	-8.947	-7.550
-2.734	-2.381	-1.821	-1.870	-2.783	-1.978	1.099	1.688	-2.261	2.783	-2.678	1.000	1.487	-2.087	4.348	-2.899	2.830	2.186	2.783	3.826	-11.691	10.654	-6.193	-6.609	-9.913	-2.878	1.832	0.364	-0.522	2.432	
-8.811	-3.052	-1.241	-2.182	-0.912	-0.776	-0.898	0.385	2.964	-0.680	-0.141	-1.735	2.960	2.182	2.211	1.815	0.000	0.532	2.182	1.020	-8.530	-10.054	-7.090	-5.636	-11.565	1.815	-2.934	-0.532	2.182	-5.105	
4.991	8.126	0.977	1.274	3.669	2.582	4.312	-3.226	3.025	5.832	7.259	8.783	8.955	6.969	7.325	6.024	7.463	5.375	3.503	4.140	5.680	6.404	5.074	4.453	1.433	-0.344	4.146	-1.629	-1.115	2.229	
6.678	9.694	7.942	5.778	7.414	3.656	7.493	4.376	3.025	7.249	8.918	11.905	10.252	10.763	7.705	10.2															

APPENDIX H  
STATISTICAL DATA OF CORNEAL THICKNESS  
MEASUREMENT AT FIVE LOCATIONS  
OVER TIME



# CENTRAL THICKNESS CHANGES - PERMAFLEX

66A

## One Way ANOVA 11 Groups

### Analysis of Variance Table

Source	DF:	Sum Squares:	Mean Square:	F-test:
Between groups	10	.009	.001	.531
Within groups	121	.198	.002	p > .25
Total	131	.206		

1

Model II estimate of between component variance = -.00006392

Group	Count:	Mean:
Column 1	12	.531
Column 2	12	.535
Column 3	12	.545
Column 4	12	.546
Column 5	12	.537

2

Group:	Count:	Mean:
Column 6	12	.536
Column 7	12	.544
Column 8	12	.561
Column 9	12	.545
Column 10	12	.532

3

Group:	Count:	Mean:
Column 11	12	.537

4

# SUPERIOR THICKNESS CHANGES - PERMAFLEX

66B

## One Way ANOVA 11 Groups

### Analysis of Variance Table

Source	DF	Sum Squares	Mean Square	F-test
Between groups	10	.007	.001	.394
Within groups	121	.214	.002	p > .25
Total	131	.221		

Model II estimate of between component variance = -.00008929

1

Group:	Count:	Mean:
Column 1	12	.552
Column 2	12	.549
Column 3	12	.56
Column 4	12	.564
Column 5	12	.551

2

Group:	Count:	Mean:
Column 6	12	.554
Column 7	12	.569
Column 8	12	.568
Column 9	12	.568
Column 10	12	.552

3

Group:	Count:	Mean:
Column 11	12	.554

4

# NASAL THICKNESS CHANGES - PERMAFLEX

66C

## One Way ANOVA 11 Groups

### Analysis of Variance Table

Source	DF:	Sum Squares:	Mean Square:	F-test:
Between groups	10	.008	.001	.413
Within groups	121	.242	.002	p > .25
Total	131	.25		

Model II estimate of between component variance = -.00009777

1

Group:	Count:	Mean:
Column 1	12	.545
Column 2	12	.551
Column 3	12	.563
Column 4	12	.558
Column 5	12	.551

2

Group:	Count:	Mean:
Column 6	12	.549
Column 7	12	.561
Column 8	12	.566
Column 9	12	.561
Column 10	12	.547

3

Group:	Count:	Mean:
Column 11	12	.541

4

# TEMPORAL THICKNESS CHANGES - PERMAFLEX

66D

## One Way ANOVA 11 Groups

### Analysis of Variance Table

Source	DF	Sum Squares	Mean Square	F-test
Between groups	10	.006	.001	.367
Within groups	121	.212	.002	p > .25
Total	131	.218		

Model II estimate of between component variance = -.0000923

1

Group:	Count:	Mean:
Column 1	12	.565
Column 2	12	.567
Column 3	12	.571
Column 4	12	.572
Column 5	12	.561

2

Group:	Count:	Mean:
Column 6	12	.563
Column 7	12	.575
Column 8	12	.582
Column 9	12	.576
Column 10	12	.559

3

Group:	Count:	Mean:
Column 11	12	.561

4

# INFERIOR THICKNESS CHANGES - PERMAFLEX

66E

## One Way ANOVA 11 Groups

### Analysis of Variance Table

Source	DF	Sum Squares	Mean Square	F-test
Between groups	10	.006	.001	.327
Within groups	121	.212	.002	p > .25
Total	131	.218		

Model II estimate of between component variance = -.00009835

1

Group:	Count:	Mean:
Column 1	12	.564
Column 2	12	.57
Column 3	12	.576
Column 4	12	.573
Column 5	12	.564

2

Group:	Count:	Mean:
Column 6	12	.564
Column 7	12	.574
Column 8	12	.583
Column 9	12	.577
Column 10	12	.564

3

Group:	Count:	Mean:
Column 11	12	.562

4

# CENTRAL THICKNESS CHANGES - SYNTAX C01-T

66F

## One Way ANOVA 11 Groups

### Analysis of Variance Table

Source	DF:	Sum Squares:	Mean Square:	F-test:
Between groups	10	.009	.001	.724
Within groups	143	.168	.001	p > .25
Total	153	.176		

Model II estimate of between component variance = -.00002312

1

Group:	Count:	Mean:
Column 1	14	.526
Column 2	14	.534
Column 3	14	.546
Column 4	14	.534
Column 5	14	.537

2

Group:	Count:	Mean:
Column 6	14	.541
Column 7	14	.542
Column 8	14	.55
Column 9	14	.541
Column 10	14	.53

3

Group:	Count:	Mean:
Column 11	14	.526

4

# SUPERIOR THICKNESS CHANGES - SYNTEX CSI-T

66G

## One Way ANOVA 11 Groups

### Analysis of Variance Table

Source	DF	Sum Squares	Mean Square	F-test
Between groups	10	.009	.001	.756
Within groups	143	.163	.001	p > .25
Total	153	.171		

Model II estimate of between component variance = -.00001988

1

Group:	Count:	Mean:
Column 1	14	.545
Column 2	14	.557
Column 3	14	.566
Column 4	14	.556
Column 5	14	.566

2

Group:	Count:	Mean:
Column 6	14	.552
Column 7	14	.552
Column 8	14	.564
Column 9	14	.558
Column 10	14	.549

3

Group:	Count:	Mean:
Column 11	14	.544

4

# NASAL THICKNESS CHANGES - SYNTAX CSI-T

66H

## One Way ANOVA 11 Groups

### Analysis of Variance Table

Source	DF:	Sum Squares:	Mean Square:	F-test:
Between groups	10	.01	.001	.905
Within groups	143	.153	.001	p > .25
Total	153	.163		

Model II estimate of between component variance = -.000007278

1

Group:	Count:	Mean:
Column 1	14	.556
Column 2	14	.568
Column 3	14	.582
Column 4	14	.574
Column 5	14	.578

2

Group:	Count:	Mean:
Column 6	14	.564
Column 7	14	.566
Column 8	14	.58
Column 9	14	.569
Column 10	14	.565

3

Group:	Count:	Mean:
Column 11	14	.559

4



# TEMPORAL THICKNESS CHANGES - SYNTAX CSI-T

661

## One Way ANOVA 11 Groups

### Analysis of Variance Table

Source	DF:	Sum Squares:	Mean Square:	F-test:
Between groups	10	.006	.001	.505
Within groups	143	.182	.001	$p > .25$
Total	153	.189		

Model II estimate of between component variance = -.00004505

1

Group:	Count:	Mean:
Column 1	14	.535
Column 2	14	.547
Column 3	14	.556
Column 4	14	.544
Column 5	14	.549

2

Group:	Count:	Mean:
Column 6	14	.542
Column 7	14	.543
Column 8	14	.556
Column 9	14	.55
Column 10	14	.541

3

Group:	Count:	Mean:
Column 11	14	.537

4

# INFERIOR THICKNESS CHANGES - SYNTAX CSI-T

66J

## One Way ANOVA 11 Groups

### Analysis of Variance Table

Source	DF:	Sum Squares:	Mean Square:	F-test:
Between groups	10	.006	.001	.425
Within groups	143	.185	.001	p > .25
Total	153	.191		

Model II estimate of between component variance = -.00005325

1

Group:	Count:	Mean:
Column 1	14	.559
Column 2	14	.57
Column 3	14	.572
Column 4	14	.568
Column 5	14	.571

2

Group:	Count:	Mean:
Column 6	14	.563
Column 7	14	.561
Column 8	14	.574
Column 9	14	.563
Column 10	14	.558

3

Group:	Count:	Mean:
Column 11	14	.555

4

# CENTRAL THICKNESS CHANGES - HYDROCURVE 11

66K

## One Way ANOVA 11 Groups

### Analysis of Variance Table

Source	DF:	Sum Squares:	Mean Square:	F-test:
Between groups	10	.013	.001	.79
Within groups	275	.447	.002	p > .25
Total	285	.46		

Model II estimate of between component variance = -.00001312

1

Group:	Count:	Mean:
Column 1	26	.521
Column 2	26	.527
Column 3	26	.535
Column 4	26	.534
Column 5	26	.535

2

Group:	Count:	Mean:
Column 6	26	.527
Column 7	26	.535
Column 8	26	.547
Column 9	26	.539
Column 10	26	.526

3

Group:	Count:	Mean:
Column 11	26	.53

4

# SUPERIOR THICKNESS CHANGES - HYDROCURVE 11

66L

## One Way ANOVA 11 Groups

### Analysis of Variance Table

Source	DF	Sum Squares	Mean Square	F-test
Between groups	10	.011	.001	.732
Within groups	275	.418	.002	p > .25
Total	285	.429		

Model II estimate of between component variance = -.00001564

1

Group:	Count:	Mean:
Column 1	26	.547
Column 2	26	.553
Column 3	26	.564
Column 4	26	.556
Column 5	26	.555

2

Group:	Count:	Mean:
Column 6	26	.551
Column 7	26	.558
Column 8	26	.566
Column 9	26	.563
Column 10	26	.55

3

Group:	Count:	Mean:
Column 11	26	.55

4

# NASAL THICKNESS CHANGES - HYDROCURVE 11

66M

## One Way ANOVA 11 Groups

### Analysis of Variance Table

Source	DF:	Sum Squares:	Mean Square:	F-test:
Between groups	10	.014	.001	.942
Within groups	275	.408	.001	p > .25
Total	285	.422		

1

Model II estimate of between component variance = -.00000332

Group:	Count:	Mean:
Column 1	26	.558
Column 2	26	.568
Column 3	26	.58
Column 4	26	.57
Column 5	26	.568

2

Group:	Count:	Mean:
Column 6	26	.562
Column 7	26	.572
Column 8	26	.576
Column 9	26	.572
Column 10	26	.56

3

Group:	Count:	Mean:
Column 11	26	.558

4

# TEMPORAL THICKNESS CHANGES - HYDROCURVE II

66N

## One Way ANOVA 11 Groups

### Analysis of Variance Table

Source	DF:	Sum Squares:	Mean Square:	F-test:
Between groups	10	.011	.001	.701
Within groups	275	.438	.002	p > .25
Total	285	.449		

1

Model II estimate of between component variance = -.0000183

Group:	Count:	Mean:
Column 1	26	.541
Column 2	26	.546
Column 3	26	.555
Column 4	26	.549
Column 5	26	.544

2

Group:	Count:	Mean:
Column 6	26	.544
Column 7	26	.55
Column 8	26	.561
Column 9	26	.556
Column 10	26	.544

3

Group:	Count:	Mean:
Column 11	26	.543

4

# INFERIOR THICKNESS CHANGES - HYDROCURVE II

660

## One Way ANOVA 11 Groups

### Analysis of Variance Table

Source	DF	Sum Squares	Mean Square	F-test
Between groups	10	.018	.002	1.174
Within groups	275	.424	.002	p > .25
Total	285	.442		

1

Model II estimate of between component variance = .0000103

Group:	Count:	Mean:
Column 1	26	.557
Column 2	26	.57
Column 3	26	.582
Column 4	26	.568
Column 5	26	.564

2

Group:	Count:	Mean:
Column 6	26	.564
Column 7	26	.573
Column 8	26	.582
Column 9	26	.571
Column 10	26	.56

3

Group:	Count:	Mean:
Column 11	26	.561

4

# CENTRAL THICKNESS CHANGES - CONTROLS

66P

## One Way ANOVA 6 Groups

### Analysis of Variance Table

Source	DF:	Sum Squares:	Mean Square:	F test:
Between groups	5	.003	.001	.637
Within groups	42	.042	.001	p > .25
Total	47	.045		

1

Model II estimate of between component variance = -.00004491

Group	Count:	Mean:
Column 1	8	.523
Column 2	8	.52
Column 3	8	.513
Column 4	8	.505
Column 5	8	.5

2

Group	Count:	Mean:
Column 6	8	.509

3



# SUPERIOR THICKNESS CHANGES - CONTROLS

## One Way ANOVA 6 Groups

66Q

### Analysis of Variance Table

Source	DF	Sum Squares	Mean Square	F-test
Between groups	5	.003	.001	.561
Within groups	42	.049	.001	p > .25
Total	47	.052		

1

Model II estimate of between component variance = -.00006348

Group:	Count:	Mean:
Column 1	8	.547
Column 2	8	.543
Column 3	8	.534
Column 4	8	.528
Column 5	8	.523

2

Group:	Count:	Mean:
Column 6	8	.536

3

# NASAL THICKNESS CHANGES - CONTROLS

66R

## One Way ANOVA 6 Groups

### Analysis of Variance Table

Source	DF:	Sum Squares:	Mean Square:	F-test:
Between groups	5	.003	.001	.683
Within groups	42	.031	.001	p > .25
Total	47	.033		

1

Model II estimate of between component variance = -.000029

Group:	Count:	Mean:
Column 1	8	.546
Column 2	8	.557
Column 3	8	.537
Column 4	8	.538
Column 5	8	.537

2

Group:	Count:	Mean:
Column 6	8	.544

3

# TEMPORAL THICKNESS CHANGES - CONTROLS

66S

## One Way ANOVA 6 Groups

### Analysis of Variance Table

Source	DF:	Sum Squares:	Mean Square:	F-test:
Between groups	5	.002	.0003226	.327
Within groups	42	.041	.001	p > .25
Total	47	.043		

1

Model II estimate of between component variance = -.00008299

Group.	Count:	Mean:
Column 1	8	.526
Column 2	8	.527
Column 3	8	.525
Column 4	8	.521
Column 5	8	.511

2

Group	Count:	Mean:
Column 6	8	.516

3

# INFERIOR THICKNESS CHANGES - CONTROLS

## One Way ANOVA 6 Groups

66T

Analysis of Variance Table

Source	DF	Sum Squares	Mean Square	F-test
Between groups	5	.002	.0004185	.547
Within groups	42	.032	.001	$p > .25$
Total	47	.034		

1

Model II estimate of between component variance = -.00004328

Group:	Count:	Mean:
Column 1	8	.544
Column 2	8	.553
Column 3	8	.539
Column 4	8	.536
Column 5	8	.535

2

Group:	Count:	Mean:
Column 6	8	.535

3

## APPENDIX I

STATISTICAL DATA AND FREQUENCY DISTRIBUTION ON  
FIVE CORNEAL LOCATIONS PRIOR TO WEAR OF  
CONTACT LENSES

# **CENTRAL MEASUREMENT** **52 EYES**

## **STATISTICAL DATA**

Column 1					
Mean:	Std. Dev.:	Std. Error:	Variance:	Coef. Var.:	Count:
.525	.039	.005	.002	7.431	52
Minimum:	Maximum:	Range:	Sum:	Sum Squared:	* Missing:
.461	.613	.152	27.297	14.407	0
Median:	Mode:	Geo. Mean:	Har. Mean:	Kurtosis:	Skewness:
.525	•	.524	.522	-.945	.07

## **FREQUENCY DISTRIBUTION**

Column 1				
Bar:	From: (>)	To: (<)	Count:	Percent:
1	.461	.476	7	13.462
2	.476	.492	6	11.538
3	.492	.507	6	11.538
4	.507	.522	6	11.538
5	.522	.538	4	7.692
6	.538	.553	9	17.308
7	.553	.568	7	13.462
8	.568	.583	4	7.692
9	.583	.599	2	3.846
10	.599	.614	1	1.923

-Mode

## SUPERIOR MEASUREMENT 52 EYES

### STATISTICAL DATA

Column 1					
Mean:	Std. Dev.:	Std. Error:	Variance:	Coef. Var.:	Count:
.548	.037	.005	.001	6.81	52
Minimum:	Maximum:	Range:	Sum:	Sum Squared:	* Missing:
.487	.633	.146	28.471	15.659	0
Median:	Mode:	Geo. Mean:	Har. Mean:	Kurtosis:	Skewness:
.543	.499	.546	.545	-.833	.149

### FREQUENCY DISTRIBUTION

Column 1				
Bar:	From: (≥)	To: (<)	Count:	Percent:
1	.487	.502	9	17.308
2	.502	.516	4	7.692
3	.516	.531	4	7.692
4	.531	.546	10	19.231
5	.546	.561	6	11.538
6	.561	.575	4	7.692
7	.575	.59	8	15.385
8	.59	.605	5	9.615
9	.605	.619	0	0
10	.619	.634	2	3.846

-Mode

## NASAL MEASUREMENT 52 EYES

### STATISTICAL DATA

Column 1					
Mean:	Std. Dev.:	Std. Error:	Variance:	Coef. Var.:	Count:
.554	.037	.005	.001	6.669	52
Minimum:	Maximum:	Range:	Sum:	Sum Squared:	* Missing:
.481	.628	.147	28.822	16.045	0
Median:	Mode:	Geo. Mean:	Har. Mean:	Kurtosis:	Skewness:
.551	●	.553	.552	-.885	.112

### FREQUENCY DISTRIBUTION

Column 1				
Bar:	From: (≥)	To: (<)	Count:	Percent:
1	.481	.496	1	1.923
2	.496	.511	7	13.462
3	.511	.525	5	9.615
4	.525	.54	8	15.385
5	.54	.555	6	11.538
6	.555	.57	5	9.615
7	.57	.585	8	15.385
8	.585	.599	6	11.538
9	.599	.614	3	5.769
10	.614	.629	3	5.769



## TEMPORAL MEASUREMENT 52 EYES

### STATISTICAL DATA

Column 1					
Mean:	Std. Dev.:	Std. Error:	Variance:	Coef. Var.:	Count:
.545	.038	.005	.001	6.993	52
Minimum:	Maximum:	Range:	Sum:	Sum Squared:	# Missing:
.465	.637	.172	28.35	15.531	0
Median:	Mode:	Geo. Mean:	Har. Mean:	Kurtosis:	Skewness:
.55	•	.544	.543	-.529	.069

### FREQUENCY DISTRIBUTION

Column 1				
Bar:	From: (≥)	To: (<)	Count:	Percent:
1	.465	.482	3	5.769
2	.482	.499	2	3.846
3	.499	.517	9	17.308
4	.517	.534	5	9.615
5	.534	.551	9	17.308
6	.551	.568	10	19.231
7	.568	.585	6	11.538
8	.585	.603	4	7.692
9	.603	.62	3	5.769
10	.62	.637	1	1.923

-Mode

## INFERIOR MEASUREMENT 52 EYES

### STATISTICAL DATA

Column 1					
Mean:	Std. Dev.:	Std. Error:	Variance:	Coef. Var.:	Count:
.559	.037	.005	.001	6.664	52
Minimum:	Maximum:	Range:	Sum:	Sum Squared:	* Missing:
.494	.628	.134	29.077	16.33	0
Median:	Mode:	Geo. Mean:	Har. Mean:	Kurtosis:	Skewness:
.566	•	.558	.557	-1.096	-.067

### FREQUENCY DISTRIBUTION

Column 1				
Bar:	From: (>)	To: (<)	Count:	Percent:
	.4	.507	4	7.692
2	.507	.521	8	15.385
3	.521	.534	5	9.615
4	.534	.548	4	7.692
5	.548	.561	3	5.769
6	.561	.574	7	13.462
7	.574	.588	8	15.385
8	.588	.601	7	13.462
9	.601	.615	3	5.769
10	.615	.628	3	5.769

## APPENDIX J

ENDOTHELIUM CELL COUNT AND STATISTICS  
SUBJECTS AND CONTROL PATIENTS

PT. NO.	EYE	DAY 0 # CELLS	DAY 150 # CELLS	NOTES CELLS RATED - CHANGE NOTED
102	OD	3200	3200	MODERATE POLYMEGATHISM - NO SIGNIFICANT CHANGE
102	OS	3100	3100	MODERATE POLYMEGATHISM - NO SIGNIFICANT CHANGE
103	OD	3200	3200	MODERATE POLYMEGATHISM - NO SIGNIFICANT CHANGE
103	OS	3200	3200	MODERATE POLYMEGATHISM - NO SIGNIFICANT CHANGE
105	OD	3200	3100	NORMAL - NO SIGNIFICANT CHANGE
105	OS	3100	3200	NORMAL - NO SIGNIFICANT CHANGE
107	OD	3100	3100	NORMAL TO SLIGHT POLYMEGATHISM - SMALL INCREASE
107	OS	3100	3100	NORMAL TO SLIGHT POLYMEGATHISM - SMALL INCREASE
108	OD	3200	3200	NORMAL - NO SIGNIFICANT CHANGE
108	OS	3200	3300	NORMAL - NO SIGNIFICANT CHANGE
109	OD	2800	2800	MODERATE POLYMEGATHISM - NO SIGNIFICANT CHANGE
109	OS	2800	2800	MODERATE POLYMEGATHISM - NO SIGNIFICANT CHANGE
111	OD	3100	3200	NORMAL - NO SIGNIFICANT CHANGE
111	OS	3300	3300	NORMAL - NO SIGNIFICANT CHANGE
112	OD	2800	3000	SLIGHT POLYMEGATHISM - SMALL INCREASE
112	OS	2600	2900	SLIGHT POLYMEGATHISM - SMALL INCREASE
113	OD	2800	2800	MODERATE POLYMEGATHISM - NO SIGNIFICANT CHANGE
113	OS	2800	2800	MODERATE POLYMEGATHISM - NO SIGNIFICANT CHANGE
115	OD	2700	2700	MODERATE POLYMEGATHISM - SLIGHT INCREASE
115	OS	2700	2700	MODERATE POLYMEGATHISM - SLIGHT INCREASE
118	OD	3200	3200	SLIGHT POLYMEGATHISM - SMALL INCREASE
118	OS	3200	3200	SLIGHT POLYMEGATHISM - SMALL INCREASE
119	OD	2800	2800	SLIGHT TO MODERATE POLYMEGATHISM - SMALL INCREASE
119	OS	2800	2800	SLIGHT TO MODERATE POLYMEGATHISM - SMALL INCREASE
120	OD	3000	3000	MODERATE TO SEVERE - SMALL INCREASE
120	OS	3000	3000	MODERATE TO SEVERE - SMALL INCREASE
121	OD	2600	2600	MODERATE POLYMEGATHISM - SMALL INCREASE
121	OS	2600	2600	MODERATE POLYMEGATHISM - SMALL INCREASE
122	OD	2500	2500	NORMAL TO SLIGHT - NO SIGNIFICANT CHANGE
122	OS	2500	2500	NORMAL TO SLIGHT - NO SIGNIFICANT CHANGE
123	OD	2600	2600	MODERATE POLYMEGATHISM - SLIGHT INCREASE
123	OS	2600	2600	MODERATE POLYMEGATHISM - SLIGHT INCREASE
126	OD	2600	2600	NORMAL TO SLIGHT - NO SIGNIFICANT CHANGE
126	OS	2400	2500	NORMAL TO SLIGHT - NO SIGNIFICANT CHANGE
129	OD	2800	2800	MODERATE POLYMEGATHISM - SLIGHT INCREASE
129	OS	2800	2800	MODERATE POLYMEGATHISM - SLIGHT INCREASE
130	OD	3100	3100	SLIGHT TO MODERATE POLYMEGATHISM - SMALL INCREASE
130	OS	3100	3100	SLIGHT TO MODERATE POLYMEGATHISM - SMALL INCREASE
132	OD	2400	2400	MODERATE POLYMEGATHISM - SLIGHT INCREASE
132	OS	2400	2400	MODERATE POLYMEGATHISM - SLIGHT INCREASE
133	OD	2900	2800	MODERATE POLYMEGATHISM - SLIGHT INCREASE
133	OS	2800	2800	MODERATE POLYMEGATHISM - SLIGHT INCREASE
134	OD	2900	2900	SLIGHT POLYMEGATHISM - SMALL INCREASE
134	OS	2900	2900	SLIGHT POLYMEGATHISM - SMALL INCREASE
135	OD	2900	2900	SLIGHT POLYMEGATHISM - SMALL INCREASE
135	OS	2900	2900	SLIGHT POLYMEGATHISM - SMALL INCREASE
136	OD	3100	3100	MODERATE POLYMEGATHISM - SLIGHT INCREASE

136	OS	3100	3100	MODERATE POLYMEGATHISM - SLIGHT INCREASE
137	OD	3000	3000	MODERATE POLYMEGATHISM - SLIGHT INCREASE
137	OS	3000	3000	MODERATE POLYMEGATHISM - SLIGHT INCREASE
139	OD	2800	2800	MODERATE TO SEVERE - SMALL INCREASE
139	OS	2800	2800	MODERATE TO SEVERE - SMALL INCREASE

PT .	EYE	DAY 0	DAY 120	NOTES
NO.		#CELLS	#CELLS	CELLS RATED - CHANGE NOTED
301	OD	3000	3000	MODERATE POLYMEGATHISM - NO CHANGE
301	OS	3000	3000	MODERATE POLYMEGATHISM - NO CHANGE
303	OD	3100	3100	NORMAL - NO CHANGE
303	OS	3100	3100	NORMAL - NO CHANGE
304	OD	3100	3100	NORMAL - NO CHANGE
304	OS	3100	3100	NORMAL - NO CHANGE
305	OD	2800	2800	SLIGHT TO MODERATE POLYMEGATHISM - NO CHANGE
305	OS	2800	2800	SLIGHT TO MODERATE POLYMEGATHISM - NO CHANGE

# ENDOTHELIAL COUNT

## STATISTICAL DATA DAY 0 VS. DAY 150

### One Way ANOVA 2 Groups

#### Analysis of Variance Table

Source	DF:	Sum Squares:	Mean Square:	F-test:
Between groups	1	4711.538	4711.538	.079
Within groups	102	6120576.923	60005.656	p > .25
Total	103	6125288.462		

Model II estimate of between component variance = -1063.348

Group:	Count:	Mean:
Column 1	52	2886.538
Column 2	52	2900

# ENDOTHELIAL COUNT

## STATISTICAL DATA DAY 150

### Column 2

Mean:	Std. Dev.:	Std. Error:	Variance:	Coef. Var.:	Count:
2900	244.147	33.857	59607.843	8.419	52
Minimum:	Maximum:	Range:	Sum:	Sum Squared:	# Missing:
2400	3300	900	150800	440360000	0
Median:	Mode:	Geo. Mean:	Har. Mean:	Kurtosis:	Skewness:
2900	2800	2889.714	2879.231	-.884	-.261

## FREQUENCY DISTRIBUTION

### Column 2

Bar:	From: ( $\geq$ )	To: ( $<$ )	Count:	Percent:	
1	2400	2490.1	2	3.846	
2	2490.1	2580.2	3	5.769	
3	2580.2	2670.3	5	9.615	
4	2670.3	2760.4	2	3.846	
5	2760.4	2850.5	12	23.077	-Mode
6	2850.5	2940.6	5	9.615	
7	2940.6	3030.7	5	9.615	
8	3030.7	3120.8	8	15.385	
9	3120.8	3210.9	8	15.385	
10	3210.9	3301	2	3.846	



# **ENDOTHELIAL COUNT CONTROL PATIENTS**

## **STATISTICAL DATA**

Column 1					
Mean:	Std. Dev.:	Std. Error:	Variance:	Coef. Var.:	Count:
3000	130.931	46.291	17142.857	4.364	8
Minimum:	Maximum:	Range:	Sum:	Sum Squared:	# Missing:
2800	3100	300	24000	72120000	0
Median:	Mode:	Geo. Mean:	Har. Mean:	Kurtosis:	Skewness:
3050	3100	2997.441	2994.825	-1	-.816

## **FREQUENCY DISTRIBUTION**

Column 1				
Bar:	From: ( $\geq$ )	To: ( $<$ )	Count:	Percent:
1	2800	2830.1	2	25
2	2830.1	2860.2	0	0
3	2860.2	2890.3	0	0
4	2890.3	2920.4	0	0
5	2920.4	2950.5	0	0
6	2950.5	2980.6	0	0
7	2980.6	3010.7	2	25
8	3010.7	3040.8	0	0
9	3040.8	3070.9	0	0
10	3070.9	3101	4	50

-Mode

## APPENDIX K

KEYSTONE DIAGNOSTIC TESTING - AVIATORS UNIT

## EXPLANATION OF ABBREVIATIONS ON FOLLOWING PAGE

PT. # = PATIENT NUMBER

PRE S # = PRETEST HIGH SCORE CARD NUMBER

% S = PERCENT STEREOPSIS

SEC = SECONDS OF ARC

POST S # = POSTTEST (WITH SPECTACLE Rx) HIGH SCORE CARD  
NUMBER

P % S = POSTTEST (WITH SPECTACLE Rx) PERCENT STEREOPSIS

POST SEC = POSTTEST (WITH SPECTACLE Rx) SECONDS OF ARC

CL S # = POSTTEST (WITH CONTACT LENS) HIGH SCORE CARD  
NUMBER

CL % S = POSTTEST (WITH CONTACT LENS) PERCENT STEREOPSIS

CL SEC = POSTTEST (WITH CONTACT LENS) SECONDS OF ARC

XXX = NOT TAKEN

PT. #	PRE S #	% S	SEC	POST S #	P % S	POST SEC	CL S #	CL % S	SCL SEC
101	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX
102	48	85	32	49	90	26	48	85	32
103	38	35	175	43	60	75	42	55	80
104	44	65	65	46	75	45	45	70	55
105	42	55	80	44	65	65	41	50	85
106	48	85	32	48	85	32	47	80	39
107	45	70	55	47	80	39	46	75	45
108	47	80	39	47	80	39	48	85	32
109	48	85	32	49	90	26	49	90	26
110	46	75	45	XXX	XXX	XXX	XXX	XXX	XXX
111	46	75	45	45	70	55	45	70	55
112	48	85	32	48	85	32	48	85	32
113	43	60	75	42	55	80	42	55	80
114	44	65	65	44	65	65	47	80	39
115	47	80	39	48	85	32	48	85	32
116	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX
117	46	75	45	48	85	32	46	75	45
118	46	75	45	46	75	45	46	75	45
119	48	85	32	49	90	26	46	75	45
120	47	80	39	48	85	32	47	80	39
121	45	70	55	46	75	45	45	70	55
122	44	65	65	46	75	45	44	65	65
123	47	80	39	51	100	19	48	85	32
124	47	80	39	47	80	39	47	80	39
125	40	45	115	48	85	32	49	90	26
126	46	75	45	46	75	45	46	75	45
127	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX
128	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX
129	45	70	55	44	65	65	46	75	45
130	39	40	140	42	55	80	44	65	65
131	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX	XXX
132	44	65	65	44	65	65	44	65	65
133	43	60	75	47	80	39	45	70	55
134	39	40	140	39	40	140	40	45	115
135	41	50	85	42	55	80	44	65	65
136	48	85	32	49	90	26	48	85	32
137	44	65	65	46	75	45	46	75	45
138	38	35	175	39	40	140	38	35	175
139	37	30	210	38	35	175	37	30	210