

**Scientific Institute
of Public Health**

Relative potency estimates for PAHs by means of the CALUX assay

Cappuyns V., Van Overmeire I, Windal I,
Hellebosch L, Carboneille S, Van Wouwe N,
Hanot V, Degroodt JM, Goeyens L

Introduction

CALUX bioassay : binding and activation of the aryl hydrocarbon receptor (AhR) by dioxinlike compounds

















- → dioxins
- → PCBs
- → polycyclic aromatic hydrocarbons (PAHs)

Polycyclic aromatic hydrocarbons (PAHs)



- > 100 compounds with 2 or more aromatic rings
- mutagenic and carcinogenic effects
- **EPA:** 16 priority PAHs
 - several other PAHs that are not routinely analyzed are responsible for mutagenic and carcinogenic effects. (e.g. 'heavy PAHs')

US EPA PAHS

















16 US EPA priority PAHs	MW	Structure
<i>naphtalene</i>	128	
<i>acenaphthylene</i>	152	
<i>acenaphtene</i>	154	
<i>fluorene</i>	166	
<i>anthracene</i>	178	
<i>phenanthrene</i>	178	
<i>fluoranthene</i>	202	
<i>pyrene</i>	202	
<i>benzo(a)anthracene</i>	228	
<i>chrysene</i>	228	
<i>benzo(b)fluoranthene</i>	252	
<i>benzo(k)fluoranthene</i>	252	
<i>benzo(a)pyrene</i>	252	
<i>Indeno(1,2,3-c,d)pyrene</i>	276	
<i>benzo(g,h,i)perylene</i>	276	
<i>dibenzo(a,h)anthracene</i>	278	

Polycyclic aromatic hydrocarbons (PAHs)



- > 100 compounds with 2 or more aromatic rings
- mutagenic and carcinogenic effects
- **US-EPA:** 16 priority PAHs
 - several other PAHs that are not routinely analyzed are responsible for mutagenic and carcinogenic effects. (e.g. 'heavy PAHs')
- ◆ **European Union (2002)**
 - list of 15 PAHs (+ BcF) that are of major concern for human health due to their toxic properties
 - monitor those 15 priority PAHs (+BcF) in food and the environment to enable long-term exposure assessments

EU priority PAHs

Name	MW	Structure
<i>benzo(c)fluorene</i>	216	
<i>cyclopenta(c,d)pyrene</i>	226	
<i>benzo(a)anthracene</i>	28	
<i>chrysene</i>	228	
<i>5-methylchrysene</i>	242	
<i>benzo(b)fluoranthene</i>	252	
<i>benzo(k)fluoranthene</i>	252	
<i>benzo(a)pyrene</i>	252	
<i>benzo(j)fluoranthene</i>	252	
<i>Indeno(1,2,3-c,d)pyrene</i>	276	
<i>benzo(g,h,i)perylene</i>	276	
<i>dibenzo(a,h)anthracene</i>	278	
<i>dibenzo(a,e)pyrene</i>	302	
<i>dibenzo(a,h)pyrene</i>	302	
<i>dibenzo(a,i)pyrene</i>	302	
<i>dibenzo(a,l)pyrene</i>	302	


CALUX bioassay




- mouse H1L6.1 cell line from XDS
- 24 h exposure time
- maximum concentrations tested in the *in vitro* bioassays:
 - 50-400 mg/L (in DMSO)

Dose response curves

- Normal response



Benzo(k)fluoranthene





24 h exposure

Dose response curves

- Low response


Chrysene




24 h exposure

Dose response curves

- ♦ *High response*




Indeno(1,2,3-c,d)pyrene




24 h exposure

Dose response curves

- Decrease of response at high concentrations



Benzo(a)pyrene



24 h exposure

REP values



- ♦ *Non-parallel dose response curves → responses from EC_{20} and EC_{80} to derive REP_{20-80} range (Villeneuve et al., 2000)*

24 h exposure


REP-values



*relative to
TCDD,
24 h exposure*

	CALUX REP weight	REP range EC ₂₀ -EC ₈₀	CALUX REP molar	Efficacy % of TCDD
B[c]F	NI	NI	NI	9
CPP	NI	NI	NI	15
B[a]A	1 E-5 ± 2 E-6	8 E-5 - 3 E-6	7 E-6 ± 1 E-6	52
Chry	8 E-5 ± 6 E-6	4 E-4 - 8 E-5	6 E-5 ± 4 E-6	51
5MeChry	1 E-4 ± 1 E-5	3 E-4 - 1 E-4	1 E-4 ± 1 E-5	64
B[a]P	6 E-5 ± 9 E-6	6 E-5 - 6 E-5	5 E-5 ± 7 E-6	88
B[b]F	2 E-4 ± 2 E-5	4 E-4 - 1 E-4	1 E-4 ± 1 E-5	95
B[k]F	1 E-3 ± 2 E-4	4 E-3 - 1 E-3	1 E-3 ± 2 E-4	110
B[j]F	8 E-5 ± 1 E-5	2 E-4 - 6 E-5	6 E-5 ± 1 E-5	120
I[123-cd]P	3 E-5 ± 3 E-6	9 E-5 - 8 E-6	2 E-5 ± 3 E-6	143
B[ghi]Pe	NI	NI	NI	23
DB[ah]A	2 E-5 ± 3 E-6	4 E-4 - 2 E-6	2 E-5 ± 3 E-6	125
DB[ae]P	3 E-4 ± 3 E-5	5 E-4 - 4 E-4	3 E-4 ± 3 E-5	210
DB[ah]P	NC	NC	NC	117
DB[ai]P	6 E-4 ± 7 E-5	6 E-3 - 6 E-4	6 E-4 ± 6 E-5	116
DB[al]P	NI	NI	NI	24

REP values



Relative to TCDD, 24 h exposure; weight derived

Comparison with Toxicity equivalent Factors (TEF) (US-EPA, 1993)

Substance	TEF	REP
Benzo(a)pyrene	1	1
Dibenzo(a,h)anthracene	1	*
benzo(a)anthracene	0.1	0.2
benzo(b)fluoranthene	0.1	3
benzo(k)fluoranthene	0.1	22
Indeno(1,2,3-cd)pyrene	0.1	0.4
chrysene	0,01	1.4

Relative to TCDD, 24 h exposure

Comparison with literature data

PAH	This study	Ziccardi et al.(2000)	Ziccardi et al. (2002)	Clemons et al. (1998)
<i>cell line</i>	<i>H1L6.1</i>	<i>H1L1.1c2</i>	<i>H1L1.1c2</i>	<i>Hepa1c1c7</i>
B[c]F	NI	-	-	-
CPP	NI	-	-	-
B[a]A	7 E-6	2 E-3	NI	1 E-5
Chry	6 E-5	-	1 E-4	1 E-2
5Me-chry	1 E-4	-	-	-
B[a]P	5 E-5	4 E-3	8 E-5	1 E-5
B[b]F	1 E-4	-	2 E-4	-
B[k]F	1 E-3	4 E-1	2 E-2	5 E-2
B[j]F	6 E-5	-	-	-
I[123-cd]P	2 E-5	-	3 E-3	-
B[ghi]Pe	NI	-	NI	-
DB[ah]A	2 E-5	-	2 E-3	5 E-2
DB[ae]P	3 E-4	-		-
DB[ah]P	NC	-		-
DB[ai]P	6 E-4	-		-
DB[al]P	NI	-		-

Mouse cell lines

Conclusion (1)



- ♦ PAHs can significantly contribute to the TEQ-value determined by means of the 24 h CALUX bioassay.
- ♦ Nevertheless, the ranking of REP values for different PAHs obtained in the present study is comparable with the results from other researchers using an *in vitro* luciferase bioassay.
- ♦ Additionally, the CALUX response was investigated for the EU priority PAHs B[c]F and B[j]F, for which no REP-values were known until now.

European regulation on PAHs in food


- ♦ European regulation (EC 208/2005) on food contaminants:
 - Maximum levels for benzo(a)pyrene (BaP) in food
 - Use of BaP as a marker substance for all PAHs
- better to express PAH toxicity in BaP-TEQ instead of TCDD-TEQ

Use of the CALUX assay for PAHs




- ◆ Literature:
 - exposure times between 3h and 24h
 - Reference compound:
 - TCDD
 - Benzo(a)pyrene
 - ◆ Maximum induction of luciferase activity after 6-12 h (Postlind (1993), Machala et al. (2001); Villeneuve et al. (2002); Masunaga et al. (2004))
- Metabolic decomposition of PAHs


Influence of exposure time



24 h exposure time, TCDD





6 h exposure time, BaP



Dose response curves, 6 h exposure

Indeno(1,2,3-c,d)pyrene




not always a sigmoidal shape

Conclusion (2)



- ♦ Response of PAH = $f(\text{exposure time})$
- ♦ 6 h exposure: higher variability in response compared to 24 h incubation
- ♦ not always a sigmoidal shape of the curve (6 h exposure)
- ♦ RLU: RSD = 6-15%
TEQ: much higher RSD because of 'flat' calibration curve
=> problem for quality control


Influence of incubation temperature



37 °C
31 °C



24 h exposure

Influence of incubation temperature



6 h exposure

Antagonistic effects?



1: TCDD: 25 000 fg/wel
BaP: 93 000 pg/well

24 h exposure