

Membrane Enhanced Peptide Synthesis

Sheung So,^a Ludmila G. Peeva,^a Edward W. Tate,^b Robin J. Leatherbarrow,^b Andrew G. Livingston*^a

Received (in XXX, XXX) Xth XXXXXXXXX 200X, Accepted Xth XXXXXXXXX 200X


First published on the web Xth XXXXXXXXX 200X

5 DOI: 10.1039/b000000x

This paper reports a new technology platform that advantageously combines organic solvent nanofiltration (a newly emerging technology capable of molecular separations in organic solvents) with solution phase peptide synthesis - Membrane Enhanced Peptide Synthesis (MEPS).

In the last decade the market for peptide based pharmaceuticals has been growing rapidly: as of 2005 an estimated 40 peptide-based pharmaceutical products were for sale, with approximately 270 more in different phases of clinical trials¹. Large-scale peptide manufacture is essential to bring these drugs to market. Solid-phase peptide synthesis (SPPS) is the most widely used technology, since it neatly solves the critical purification problems encountered at each stage in solution phase synthesis. However it faces serious challenges including mass transfer, steric hindrance, and resin handling^{2,3,4}. The concept of membrane separation coupled to solution phase synthesis offers major advantages over SPPS by combining the advantages of “classical” solution phase synthesis with the ease of purification of the solid phase method². Compared to SPPS, reactions in solution phase provide faster rates, and are less affected by steric hindrance due to peptide folding, or reactions within confined space which result in transepeptidation, and are not limited by intraparticle diffusional mass transfer phenomena^{3,4}. This enables reduction in the large reactant excesses, in some cases of up-to 5-fold^{5,6}, that are employed in SPPS to compensate for mass transfer limitations. Solution phase reactions are also easier to scale-up due to the absence of swelling effects and cake formation within solid resins, which require more complex reactor design. Prior applications of membrane separation in peptide chemistry are restricted to re-concentration of peptides⁷, amino acid recovery^{8,9}, and a single report on membrane separation for purification between reaction cycles during peptide synthesis^{2,10}. In the latter case, peptides built on poly(ethylene glycol) (PEG) were separated from impurities by ultrafiltration. However, due to a lack of organic solvent-compatible membranes, this required evaporation of organic solvent after each coupling and each deprotection step, neutralisation after deprotection, and uptake in water prior to ultrafiltration. Water was then removed by evaporation and/or azeotropic distillation before re-dissolving the PEG-peptide into organic solvent for the next coupling step. This complex process limited the synthesis to volatile organic solvents, and the lengthy solvent switching system made the separation process unsuitable for large-scale production. Molecular separation in organic solvents via nanofiltration (Organic Solvent Nanofiltration – OSN) is an emerging technology¹¹, which should be an ideal separation

method for in-cycle purification during peptide synthesis.



55 **Figure 1:** Schematic representation of membrane enhanced peptide synthesis (MEPS). Peptide chain assembly was performed following this scheme using the apparatus presented in **Figure 4**.

. Here we report for the first time Membrane Enhanced Peptide Synthesis (MEPS), a new technology platform that advantageously combines OSN with solution phase peptide synthesis.

The MEPS concept is illustrated in **Figure 1**. Peptide chain assembly occurs via: 1) amide coupling; 2) a washing step for removal of excess reagents *via* constant volume diafiltration; 3) deprotection; 4) a washing step for removal of deprotection by-products and excess reagents again *via* diafiltration. The cycle is repeated as many times as necessary, adding a further amino acid each cycle, until the desired peptide sequence is obtained. In contrast to the previously reported studies^{2,10} washing is carried out immediately after the coupling and deprotection steps using the reaction solvent, and does not require any solvent exchanges.

The peptide is assembled on a soluble polymeric support, methoxy-amino-PEG with molecular weight (MW) 5,000 g.mol⁻¹ (MeO-PEG-NH₂), to increase retention by the membrane³. Since PEGylated peptides have received much recent attention due to their enhanced therapeutic and pharmacokinetic potential¹²⁻¹⁴, in some cases a bioactive peptide made by MEPS might also be used directly without cleavage from PEG.

Fmoc (9-fluorenylmethoxycarbonyl) peptide synthesis chemistry was chosen due to its widespread application and mild deprotection conditions (piperidine/DMF). Hydroxymethylphenylacetic acid (HMPA) was used as linker to attach the first amino acid in the peptide sequence to the MeO-PEG-NH₂ enabling facile cleavage of the final peptide via acidolysis. Peptide coupling proceeds optimally in polar aprotic solvents such as dimethylformamide (DMF) or N-methyl pyrrolidinone (NMP), and DMF was used in both reaction and diafiltration steps. The use of these solvents in combination with OSN membranes has only recently become possible with the development of membranes¹⁵ which possess good stability in these aprotic solvents and make this approach viable.

The most important consideration for successful realisation of this approach was the choice of membrane. This must possess excellent long term stability in the reaction solvent (DMF) and high selectivity between MeO-PEG-peptide, and side reaction products and excess reagents, including unreacted amino acids, activators and deprotection reagents. Membrane performance should not be affected by frequent switching of the reaction media between DMF solution in the coupling step and 20% piperidine/DMF solution in the deprotection step.

After screening a range of commercial and developmental membranes we identified a ceramic OSN membrane that met these requirements: the Inopor ZrO₂ coated membrane with 3nm pore size and hydrophobic surface modification (Innocermic, Germany). Based on the membrane characterisation data it was estimated that 10-12 solvent wash volumes should be sufficient for removal of all excess reagents from around 1.0 to less than 0.01 equivalents.

To prove the MEPS concept we started by producing a model peptide H-Tyr-Ala-Tyr-Ala-Tyr-OH. This was chosen as it includes one of the largest protected amino acids Fmoc-Tyr(Bu)-OH, and one of the smallest hydrophobic amino acids, Fmoc-Ala-OH, thus providing information on the performance of the MEPS process with respect to different molecular sizes and properties of amino acids. PyBOP was chosen for the coupling reaction as it is one of the largest of the commercial activators available, and its successful removal presents a challenge for the MEPS process. DIC was used for the esterification linking the first amino acid onto MeO-PEG-HMPA. Thus the synthesis also provided insight into the behaviour of the post-reaction species derived from both activators during diafiltration. Finally, this first experiment also sought to establish membrane stability at high concentrations of organic base (piperidine) during the deprotection step (see Table S1, Supporting Information).

The model peptide produced contained no peptide by-products as confirmed by RP-HPLC and MALDI-TOF analysis. This absence of any detectable peptide impurities indicates that membrane purification is efficient at removing un-reacted protected amino acids throughout the synthesis. Our calculations suggested that traces of impurities would still remain in the system after each purification step even after 10 volumes of washing solvent, but this did not affect the final peptide purity. Apparently the level of impurities that can be




Figure 2: HPLC chromatograms of peptide TP-5 produced by MEPS and SPPS processes, and TP-5 standard purchased from Sigma-Aldrich (UK). The target TP-5 peptide was eluted at 10.3 minutes. Both syntheses (MEPS and SPPS) were performed under the same reaction conditions of 2 equivalents of reagents per 1 equivalent peptide, and single reaction cycles. Peptide purity was determined as a ratio between the target peptide TP-5 peak area and the total area of the peaks corresponding to peptide sequences in the solution. The purity of TP-5 produced by MEPS was determined as ~94% (two impurities eluted at 10.0 minutes and 10.4 minutes) while TP-5 produced by SPPS was ~77% pure (one impurity eluted at 10.5 minutes). The large peaks eluted between 19 - 23 minutes were PEGylated wastes such as MeO-PEG-HMPA and the peak eluted at 13 minutes was not of peptide origin as confirmed by MALDI-TOF analysis and were not taken into account for the purity calculations.

tolerated in the system without provoking side reactions is higher than anticipated, and so there is potential to reduce the washing volume.

Encouraged by this first success, Thymopentin (H-Arg-Lys-Asp-Val-Tyr-OH) was synthesized as a second demonstration of the MEPS process. Thymopentin (TP-5) is a derivative of the naturally occurring hormone thymopietin, with potential for the treatment of rheumatoid arthritis, AIDS and other primary immunodeficiencies¹⁶. Besides being a potential active pharmaceutical, this peptide includes a range of amino acids from aromatic (Tyr), acidic (Asp) and basic (Lys and Arg) to hydrophobic (Val). It also contains the largest Fmoc/Boc protected amino acid (Fmoc-Arg(Boc)₂, M.W. 597 g.mol⁻¹, and so this synthesis represents a significant further challenge for MEPS. RP-HPLC analyses of TP-5 produced by both MEPS and SPPS are illustrated in **Figure 2**.

The purity of the MEPS product (as a percentage of total (TP-5 plus peptidic by-products)) was estimated at 94%. MALDI-TOF analysis (**Figure 3**) confirmed the target product molecular weight of MH+ 680 and identified the two impurities as formed by deletion of Asp, MH+ 564, and Lys, MH+ 550. TP-5 produced by SPPS under the same reagent excess (2 equivalents) was only 77% pure and the main impurity was identified as deletion of Arg, MH+ 524. This result demonstrates the key advantage of liquid phase synthesis over SPPS – a higher purity was obtained using the same excess of reagents. Typically more equivalents are needed for SPPS⁵. The overall yield of TP-5 produced by the MEPS process was estimated to be 92%, with respect to the starting MeO-PEG-NH₂ material.




Figure 3: MALDI-TOF mass spectrum of the different peptides produced by MEPS and SPPS TP-5 synthesis and isolated using semi-preparative HPLC. Spectrum B) – from the MEPS process and D) – from the SPPS process correspond to the peak eluted at 10.3 minutes. It was identified as TP-5 and showed the target molecular mass, MH^+ of 680 Da. Spectrum A) and C) – from the MEPS process correspond to the two impurities eluted at 10.0 and 10.4 minutes respectively and were identified as deletion of Lys, MH^+ 550 Da and Asp, MH^+ 564 Da. Spectrum E) from the SPPS process corresponds to the impurity eluted at 10.5 minutes and was identified as deletion of Arg, MH^+ 524 Da.





Figure 4: Experimental set-up used for peptide chain assembly: Both coupling and deprotection reactions were performed in the Reaction Vessel where mixing was provided via the Circulation Pump. Upon completion of each reaction, the system was pressurised to 7 barg using nitrogen gas. Fresh DMF solvent was pumped via an HPLC pump from the Solvent Reservoir into the system to replace the permeated solvent and maintain constant liquid volume within the Reaction Vessel; (a) Inopor ZrO_2 ceramic membrane – front view; (b) SEM image of the Inopor ZrO_2 ceramic membrane - edge view at magnification 370x. The line bar corresponds to 50 μm .

For this investigation batch size of 0.9 mmol TP-5 was produced which yielded ~0.6 g of product respectively. With the current laboratory set-up (Figure 4) > 20 mmol batches of peptide can be produced by simply increasing the feed volume and using identical operating conditions. Further scale-up to kilogram or ton scale may be possible by increasing the size of the equipment.

The MEPS process proposed in this work integrates the advantages of performing peptide synthesis in solution with a direct membrane purification of the post-reaction mixture. The process is less constrained by mass-transfer limitations, and requires a smaller excess of reagents, than SPPS, yet demonstrates excellent purity and yield of the final peptide. We anticipate that further optimisation of the separation step and wash solvent volume will result in solvent savings and improved process economics. Thus, we conjecture that the MEPS process offers an important alternative technology platform for peptide and PEGylated peptide production at industrial scale.

This work was funded by the U.K. Engineering and Physical Sciences Research Council (EPSRC) Grant NoEP/C50982X/1; The authors would like to acknowledge also the financial support from the FP7 IAPP EC Grant NoPIAP-GA-2008-218068.

Notes and references

^a Department of Chemical Engineering, Imperial College London, London SW7 2AZ, UK. Fax: +44 (0)20 75945629; Tel: +44 (0)20 75945582; E-mail: a.livingston@imperial.ac.uk

^b Department of Chemistry, Imperial College London, London SW7 2AZ, UK. E-mail e.tate@imperial.ac.uk

† Electronic Supplementary Information (ESI) available: Experimental Methods and Procedures; Peptide Synthesis; Analysis; Peptide purity and yield. See DOI:

- 1 B. Bray, *Nature Review Drug Discovery*, 2003, **2**, 587 – 593.
- 2 E. Bayer, M. Mutter, *Nature*, 1972, **237**, 512 – 513.
- 3 E. Bayer, M. Mutter, R. Uhmann, J. Polster, H. Mauser, *J. Am. Chem. Soc.*, 1974, **96**, 7333-7336.
- 4 V. Gut , J. Rudinger in *Peptides 1968*, (Eds.: E. Bricas) Amsterdam, 1968, 185-188.
- 5 M. Amblard, J. Fehrentz, J. Martinez, G. Subra, *Mol. Biotech.*, 2006, **33**, 239 – 254.
- 6 W. Chan, P. White, *Fmoc solid phase peptide synthesis: a practical approach*, Oxford University Press, **2000**.
- 7 T. Tsuru, S. Nakao, S. Kimura, T. Shutou, *Sep. Sci. Technol.*, 1994, **29**, 971 – 984.
- 8 S-L. Li, C. Li , Y-S. Lui, X-L. Wang, Z-A. Cao, *J. Membr. Sci.*, 2003, **222**, 191 – 201.
- 9 X. Wang, A. Ying, W. Wang, *J. Membr. Sci.*, 2002, **196**, 59 – 67.
- 10 E. Bayer, D. Gillessen, H. Kuenzi, M. Mutter, R. Studer (F. Hoffmann-La Roche & Co.), DE-F2047413, **1980**.
- 11 P. Vandezande, L.E.M. Gevers, I.F.J. Vankelecom, *Chem. Soc. Rev.*, 2008, **37**, 365-405.
- 12 M. Hamidi, A. Azadi, P. Rafiei, *Drug Delivery.*, 2006, **13**, 399 – 409.
- 13 S.-H. Lee, S. Lee, Y. Youn, D. Na, S. Chae, Y. Byun, K. Lee, *Bioconjugate Chem.*, 2005, **16**, 377 – 382.
- 14 M. Roberts, M. Bentley, J. Harris, *Adv. Drug Delivery Rev.*, 2002, **54**, 459 – 476.
- 15 Y. See Toh, F. Lim, A. Livingston, *J. Membr. Sci.*, 2007, **301**, 3–10.
- 16 E. Ochoa, S. Medrano, M. de Carlin, A. Dilonardo, *Cell. Mol. Neurobiol.*, 1988, **8**, 325 – 331.