Data in Brief 3 (2015) 234-239



Contents lists available at ScienceDirect

Data in Brief

journal homepage: www.elsevier.com/locate/dib

Data Article

# Quantitative proteomics of rat and human pancreatic beta cells



# B. Brackeva<sup>a,b</sup>, G. Kramer<sup>c</sup>, J.P.C. Vissers<sup>d</sup>, G.A. Martens<sup>a,b,\*</sup>

<sup>a</sup> B-Probe, Diabetes Research Center, Vrije Universiteit Brussel (VUB), Belgium

<sup>b</sup> Department of Clinical Chemistry and Radio-immunology, Universitair Ziekenhuis Brussel, Brussels, Belgium

<sup>c</sup> Department of Medical Biochemistry, Academisch Medisch Centrum, Amsterdam, The Netherlands

<sup>d</sup> Waters Corporation, Wilmslow, United Kingdom

# ARTICLE INFO

Article history: Received 13 February 2015 Received in revised form 23 February 2015 Accepted 23 February 2015 Available online 10 March 2015

#### ABSTRACT

Data set description: This data set is composed by label-free alternate-scanning LC-MS/MS proteomics analysis human and Wistar rat pancreatic islet endocrine cells.

The mass spectrometry data of the human and rat pancreatic beta cells and the resulting proteome search output from ProteinLynx GlobalSERVER (PLGS) have been deposited to the ProteomeXchange Consortium [1] via the PRIDE partner repository with the dataset identifiers PXD001539 (human) and PXD001816 (rat). From these mass spectrometry data, 'relative molar amount units' between cell types and across species were calculated.

Biological relevance: These data provide a quantitative view on the unfractionated proteomes of human and rat beta and alpha cells. It is likely biased towards the proteins with higher molar abundance, relating to core functional pathways, but also includes several proteins with an islet-enriched expression. The quality of the cell preps is state-of-the-art, and the label-free quantitation is both precise and accurate, allowing detailed quantitative analysis.

© 2015 The Authors. Published by Elsevier Inc. This is an open access article under the CC BY license (http://creativecommons.org/licenses/hu/4.0()

(http://creativecommons.org/licenses/by/4.0/).

DOI of original article: http://dx.doi.org/10.1016/j.jprot.2015.01.009

\* Corresponding author. Tel.: + 32 2 476 37 81.

http://dx.doi.org/10.1016/j.dib.2015.02.019

2352-3409/© 2015 The Authors. Published by Elsevier Inc. This is an open access article under the CC BY license (http://creativecommons.org/licenses/by/4.0/).

E-mail address: geert.martens@uzbrussel.be (G.A. Martens).

# Specifications table

Subject area	Biology
More specific subject area	Rat and human pancreatic beta cell proteomics
Type of data	Excel tables
How data was	Label-free alternate-scanning LC-MS proteomics; Nanoscale LC separation of tryptic peptides was
acquired	performed with
	a nanoACQUITY system (Waters Corporation), equipped with a
	Symmetry C18 5 mm, 2 cm 6180 mm precolumn and an Atlantis
	C18 3 mm, 25 cm 675 mm or a Atlantis C18 3 mm, 15 cm
	675 mm analytical reversed phase column (Waters Corporation).
Data format	Raw/Analyzed
Experimental factors	Proteomes of rat alpha and beta cells were studied immediately after cell isolation, without in vitro culture or treatment. Proteins were typically extracted within 4 h of animal killing. Human beta cells
	were prepared from islet-enriched pancreas fractions, cultured for 6–14 days after human organ
	procurement. Selected fractions cleared for research were then trypsin-dissociated into single cell
	suspensions and FACS-sorted according to cell size, granularity and zinc content to an average
	insulin-purity of 60%. Proteomes of these cells was studied immediately after cell isolation without
	additional treatment or further culture.
Experimental	Proteomes were recorded using label-free LC-MS/MS. Accurate mass precursor and fragment ion
features	LC-MS data were collected in data independent, alternate scanning (LC-MSE) mode of acquisition
	[2,3]. Biological interpretation of the "relative molar amount units" was done after geometric
	normalization to a network of stable reference proteins [4–8]
Data source location	Rat data: Waters Corp., Manchester, UK, Analyst: J.P.C. Vissers
	Human data: Department of Medical Biochemistry, Academisch Medisch Centrum, Amsterdam,
	the Netherlands, Analyst: G. Kramer
Data accessibility	Raw data uploaded to ProteomeXchange Consortium via the PRIDE partner repository with the
	dataset identifiers
	PXD001539: human data
	Project DOI: 10.6019/PXD001539
	Reviewer account details:
	Username: reviewer99564@ebi.ac.uk
	Password: pJzPiClj
	PXD001816: rat data
	Project DUI: 10.6019/PXD001816
	Username: reviewer43213@eDi.ac.uk
	Passworu: yzbodylw

# Value of the data

- The methods for FACS-purification of rat alpha and beta cells, and FACS-enrichment of human islet endocrine cells were developed by our lab, and results in state-of-the-art endocrine purities. This results in minimization of exocrine contamination, which is a problem affecting other proteomics studies on these cell types
- Human and rat beta cells were analyzed using exactly the same methods of protein extraction, trypsinization, LC-MS/MS analysis and downstream data analysis by the same people (Vissers and Kramer). This makes it possible for the first time to make a detailed quantitative comparison of the core proteomes of this cell type in two species.
- We applied rigorous statistical criteria for peptide/protein identification and quantification, thereby reducing false positive identifications perhaps at the expense of claimed coverage.
- We showed that alternate-scanning LC-MS/MS generated highly reproducible protein quantifications (with overall technical plus biological imprecision < 20% and high accuracies as judged by correct measurement of multi-enzyme stoichiometry [4]).
- A limitation is our depth of coverage: with an average of 400–500 protein identifications/ quantifications per cell type, our analysis is likely biased to the core proteome and proteins with higher molar abundance.

### 1. Data, full experimental design, materials and methods

Description of the analyzed islet endocrine cell preparations:

- (i) FACS-purified rat insulin-producing beta cells (90% insulin+, 3% glucagon+, 1% somatostatin+ and 2% pancreatic polypeptide+ cells); n=3 biological replicates, 3 technical replicate injections each:
- (ii) FACS-enriched rat glucagon-producing alpha cells (2% insulin+, 94% glucagon+, 1% somatostatin+ and 2% pancreatic polypeptide+ cells); n=3 biological replicates, 3 technical replicates each. Rat alpha and beta cells were FACS-purified as described [9,10]. Cells were processed for proteomics immediately after cell isolation without additional culture or treatment.
- (iii) FACS-enriched human beta cell preparations (n=4 different donors, 3 technical replicates each) and
- (iv) human pancreatic exocrine/duct cells (n=1 biological replicate, 4 technical replicates). Human beta cells were FACS-purified from pancreatic islet-enriched fractions that were prepared for clinical islet transplantation in non-uremic type 1 diabetic patients. Sporadic cell preparations not suitable for islet transplantation were cleared for research use. These cell preparations were typically cultured for 6–14 days in HamF10 supplemented with 6 mM glucose, 2 mM glutamine, 0.5% BSA, penicillin/streptomycin, with refreshing of the medium and removal of dead cells every 3 days. Islet-enriched fractions were dissociated into single cell suspensions, and viable endocrine cells FACS-isolated according to cell size, granularity and zinc content, to an average insulin-positive purity of  $60 \pm 6\%$  (range: 53–68%). These cells were then immediately processed for proteomics, exactly as done for the rat cells, without further treatment or culture.

Protein extraction After isolation, cell preparations were washed and soluble protein was extracted with 0.5% (w/v) RapiGest detergent in 50 mM ammonium bicarbonate in the presence of Complete Protease Inhibitor Cocktail and bovine DNAse II solution, followed by centrifugation to remove cellular debris. One human beta cell pellet was subjected to an additional re-extraction in 1% (w/v) Rapigest to verify if this resulted in higher recovery of hard-to-extract/membrane-bound proteins (HP2010, Table 1). To remove protease inhibitors and most of reduced insulin molecules, protein extract was reduced with 100 mM dithiothreitol (DTT). Proteins were denatured by heating at 80 °C for 15 min, followed by 30 min at 60 °C after addition of 2.5  $\mu$ l 200 mM iodoacetamide. Trypsinization was carried out overnight at 37 °C (1:25 w/w trypsin ratio) in final volume of 100  $\mu$ l. Finally, RapiGest detergent was removed by acidifying digest to pH=2 with trifluoroacetic acid and incubation for 15 min at 37 °C.

LC-MS configuration Nanoscale LC separation of the tryptic peptides was performed with a NanoAcquity system (Waters Corporation). Samples were loaded on to a Symmetry C18 5 µm,  $2 \text{ cm} \times 180 \text{ }\mu\text{m}$  trap column (Waters) at a flow rate of 5  $\mu$ l/min prior to separation on a Bridged Ethyl Hybrid C18 1.7  $\mu$ m, 25 cm  $\times$  75  $\mu$ m analytical reversed-phase column (Waters) by application of a 90 min gradient from 1% ACN and 0.1% formic acid to 40% ACN and 0.1% formic acid at a column flow rate of 0.25  $\mu$ l/min. The column temperature was maintained at 35 °C. Analysis of the eluted tryptic peptides was performed using a Synapt G2 Q-TOF (quadrupole time-of-flight) mass spectrometer (Waters Corporation) equipped with a nanolockspray source (Waters Corporation) fitted with a pico-tip emitter (New Objective) operated at a capillary voltage of approximately 3 kV. For all measurements, the mass spectrometer was operated in v-mode with a typical resolution of at least 20,000 full width at half maximum. All analyses were performed in positive mode ESI. The time-offlight analyzer of the mass spectrometer was externally calibrated with a NaI mixture from m/z 50 to 1990. The collision gas used was argon, maintained at a constant pressure of  $2.0 \times 10-3$  mbar in the collision cell. The lock mass, [Glu1]-fFibrinopeptide B, was delivered from the auxiliary pump of the NanoAcquity system with a concentration of 100 fmol/ $\mu$ l at 0.5  $\mu$ l/min to the reference sprayer of the nanolockspray source. The data were post-acquisition lock-mass corrected using the monoisotopic mass of the doubly charged precursor of [Glu1]-Fibrinopeptide B, delivered through the reference

#### Table 1

Overview of all technical replicate LC-MS/MS injections of the indicated independent biological isolates. First column indicated the sample name, second column listed the codes of individual technical replicates, as used in the supplementary raw MZml data files.

Human pancreas number	Code tech. Replicate		% insulin+	% dead	Rapigest
HP2010 3121	20111012_s4b 20111103_s4a 20111103_s4a	human beta cells, donor 1	68	4	0.50%
HP2064 1421	20111103_s4c 20111103_s5a 20111103_s5b	human beta cells, donor 2	62	3	0.50%
HP2785 2721	20111103_s5c 20111103_s6a 20111103_s6b	human beta cells, donor 3	58	2	0.50%
HP2789 1821	20111103_s6c 20111103_s7a 20111103_s7b	human beta cells, donor 4	53	4	0.50%
HP2703 2	20111103_s7c 20111103_s8a 20111103_s8b	human duct/exocrine	2	24	0.50%
HP2010 3121 1% Rapigest	20111103_s8c 20111103_s16a 20111103_s16b 20111103_s16c	re-extraction of HP2010 at 1% Rapigest (data not further used)	68	4	1.00%
Rat beta cells	Code tech. Replicate		% insulin+	% glucagon	Rapigest
beta 001	011108_002 011108_003	biological replicate 1 (pool of 25 animals)	90	3 (average)	0.50%
beta 002	011108_004 011108_010 011108_011	biological replicate 2 (pool of 25 animals)			
beta 003	011108_012 011108_019 011108_020 011108_021	biological replicate 3 (pool of 25 animals)			
Rat alpha cells	Code tech. Replicate		% insulin+	% glucagon	Rapigest
alpha 001 alpha 002	011108_006 011108_007 011108_008 011108_014	biological replicate 1 (pool of 25 animals) co-isolated with beta 001 biological replicate 2	2	94 (average)	0.50%
alpha 003	011108_015 011108_016 011108_022 011108_023 011108_024	(pool of 25 animals) co-isolated with beta 002 biological replicate 3 (pool of 25 animals) co-isolated with beta 003			

sprayer, which was sampled every 120 s. Accurate mass precursor and fragment ion LC-MS data were collected in data independent, alternate scanning (LC-MSE) mode of acquisition [2,3]

LC-MS data processing and protein identification Continuum LC-MS data were processed and searched using ProteinLynx GlobalSERVER v2.5 (Waters Corporation). Protein identifications were obtained by searching databases of Rattus norvegicus databases (v15.12, 7449 entries) and Homo Sapiens release (2011\_11, 20,335 entries). Sequence information of Alcohol dehydrogenase Saccharomyces cerevisiae was added to the databases to afford the ability to normalize the data sets and to estimate amounts and concentration and that of known contaminant proteins (e.g. serum albumin Bos taurus and trypsin Sus scrofa). A decoy was generated on the fly with every database search experiment conducted to estimate the protein false positive rate of identification. Data independent scanning protein identifications were accepted when more than three fragment ions per peptide, seven fragment ions per protein and more than 2 peptides per protein were

#### Table 2

Overview of total proteins identified and quantified by LC-MS/MS in the various cell preparations, along criteria specified in Materials & Methods.

		rat		human	
		beta	alpha	beta	exocrine
number of proteins identified	$n \ge 1$	506	417	462	372
number of proteins quantified	$n \ge 2$	465	353	413	300
GEO-normalized molar amount	highest	3.61	8.22	8.95	14.74
	lowest	0.02	0.02	0.01	0.06
	dynamic range	226	548	746	246
biological replicates (x)		3	3	4	1
technical replicates (n)		3	3	3	3

identified, in at least one technical replicate per sample. Protein quantitation was only reported when the protein was detected in at least 2 out of 3 technical replicate of at least one biological replicate. Typical search criteria used for protein identification included automatic peptide and fragment ion tolerance settings (approximately 10 and 25 ppm, respectively), 1 allowed missed cleavage, fixed carbamidomethyl-cysteine modification and variable methionine oxidation. Raw data were expressed as 'relative molar amount units' calculated by dividing the determined molar amount for a given protein by the summed determined amount for all identified proteins as this accounts for both technical and biological variations [3,4,11]. To enable reuse of the continuum LC-MS data by other proteome search engines, MS1 and reconstructed MSn mass spectrometric data generated by PLGS were exported in mzML format.

Geometric normalization between cell types and across species: to compare these relative molar amounts between cell types and across species, data were biologically interpreted after geometric normalization to a set of stable housekeeping proteins, as described in the associated 4 Pubmed-indexed studies thus far: analysis of the rat beta cell proteome [4], identification of doublecortin [5], protein phosphatase 1-inhibitor 1 [6] and ubiquitin thioesterase-L1 [8] as candidate real-time biomarkers for beta cell destruction (Table 2).

#### Acknowledgment

This study was supported by research grants from the Research Foundation Flanders (FWO G.0492.12 project grant and Senior Clinical Investigator career support grant to G.M.), from the European Foundation for the Study of Diabetes (EFSD/JDRF/Roche Young Investigator Award to G.M.), by the Agentschap voor Innovatie door Wetenschap en Technologie (IWT to G.M. and O.C.), by the Wetenschappelijk Fonds Willy Gepts from the Universitair Ziekenhuis Brussel (to G.M and K.V.) and by the European Union, 7th FP grant (HEALTH-F5-2009-241883) to the Center for Beta Cell Therapy in Diabetes (www.betacelltherapy.org). Funding organizations did not influence data collection and interpretation.

### References

- [1] J.A. Vizcaino, E.W. Deutsch, R. Wang, A. Csordas, F. Reisinger, et al., ProteomeXchange provides globally coordinated proteomics data submission and dissemination, Nat. Biotechnol. 32 (2014) 223–226.
- [2] J.C. Silva, R. Denny, C.A. Dorschel, M. Gorenstein, I.J. Kass, et al., Quantitative proteomic analysis by accurate mass retention time pairs, Anal. Chem. 77 (2005) 2187–2200.
- [3] S.J. Geromanos, J.P. Vissers, J.C. Silva, C.A. Dorschel, G.Z. Li, et al., The detection, correlation, and comparison of peptide precursor and product ions from data independent LC-MS with data dependant LC-MS/MS, Proteomics 9 (2009) 1683–1695.
- [4] G.A. Martens, L. Jiang, K. Verhaeghen, J.B. Connolly, S.G. Geromanos, et al., Protein markers for insulin-producing beta cells with higher glucose sensitivity, PLoS ONE 5 (2010) e14214.
- [5] L. Jiang, B. Brackeva, G. Stange, K. Verhaeghen, O. Costa, et al., LC-MS/MS identification of doublecortin as abundant beta cell-selective protein discharged by damaged beta cells in vitro, J. Proteomics 80 (2013) 268–280.

- [6] L. Jiang, B. Brackeva, Z. Ling, G. Kramer, J.M. Aerts, et al., Potential of protein phosphatase inhibitor 1 as biomarker of pancreatic beta-cell injury in vitro and in vivo, Diabetes 62 (2013) 2683–2688.
- [7] G.A. Martens, E. Motte, G. Kramer, L.W. Gaarn, G. Stange, et al., Functional characteristics of neonatal rat beta cells with distinct markers, J. Mol. Endocrinol. 52 (1) (2013) 11–28.
- [8] B. Brackeva, P. De, V. G. Kramer, O. Costa, K. Verhaeghen, et al., Potential of UCHL1 as biomarker for destruction of pancreatic beta cells, J. Proteomics (2015).
- [9] G. Stangé, M. Van de Casteele, H. Heimberg, Purification of rat pancreatic B-cells by fluorescence-activated cell sorting, Methods Mol. Med. 83 (2003) 15–22.
- [10] M. Van De Winkel, E. Maes, D. Pipeleers, Islet cell analysis and purification by light scatter and autofluorescence, Biochem. Biophys. Res. Commun. 107 (1982) 525–532.
- [11] J.C. Silva, M.V. Gorenstein, G.Z. Li, J.P. Vissers, S.J. Geromanos, Absolute quantification of proteins by LCMSE: a virtue of parallel MS acquisition, Mol. Cell. Proteomics 5 (2006) 144–156.