



E1660
JACC March 12, 2013
Volume 61, Issue 10

TCT@ACC-i2: Invasive and Interventional Cardiology

COMPARISON OF IN VIVO LONGITUDINAL STRENGTH AND CONFORMABILITY FOLLOWING STENT IMPLANTATION IN RABBIT ILIAC ARTERY

Poster Contributions

Poster Sessions, Expo North

Saturday, March 09, 2013, 10:00 a.m.-10:45 a.m.

Session Title: Coronary Stents

Abstract Category: 47. TCT@ACC-i2: Coronary Intervention, Devices

Presentation Number: 2101-238

Authors: Sho Torii, Gaku Nakazawa, Takeshi Ijichi, Ayako Yoshikawa, Masayoshi Tokunaga, Yuji Ikari, Tokai University School of Medicine, Isehara, Japan

Background: Recently, longitudinal coronary stent deformation has been highlighted and some bench tests have demonstrated difference in longitudinal strength among the stent platforms, however, this has not been investigated in an in vivo setting. This is of interest because there may be a trade-off between high longitudinal strength and high conformability. The aim of the current study was to determine the in vivo longitudinal strength and conformability of the various stent platforms following stent implantation in rabbit iliac arteries.

Methods: We evaluated 4 types of commercially available stents; MultiLink8 (Abbott Vascular), Omega (Boston Scientific), Integrity (Medtronic) and Nobori (Terumo Corporation). To investigate the longitudinal strength, constant axial force was applied to the stent edge by a guiding catheter after deployment in the rabbit iliac artery and the compression rate was calculated by measuring stent length. In order to evaluate conformability, stents were deployed crossing over the iliac bifurcation and the bifurcation angles were measured before and after the stent implantation. If the change in the angle was small, the stent was considered to be more conformable.

Results: The Omega demonstrated significantly greater longitudinal compression compared with other stents (Omega $17.4 \pm 9.3\%$, MultiLink8 $2.8 \pm 2.3\%$, Integrity $2.8 \pm 1.4\%$, and Nobori $3.8 \pm 3.2\%$, $p=0.0036$), but Omega showed the greatest conformability evidenced by the least change in the bifurcation angle (Omega $12.7 \pm 0.8^\circ$, Integrity $25.7 \pm 2.4^\circ$, Integrity $28.3 \pm 1.1^\circ$, and Nobori $28.1 \pm 6.8^\circ$, $p=0.0347$).

Conclusions: In this rabbit model, the Omega stent of Element platform, showed the least longitudinal strength but the greatest conformability compared with the other stent platforms.