

Data in Brief 6 (2016) 630-633



Contents lists available at ScienceDirect

Data in Brief





Data Article

Data in support of effect of blue LED irradiation in human lymphoma cells



Phil-Sun Oh, Hyosook Hwang, Hwan-Seok Jeong, Jeongil Kwon, Hyun-Soo Kim, Minjoo Kim, SeokTae Lim, Myung-Hee Sohn, Hwan-leong leong*

Department of Nuclear Medicine, Molecular Imaging & Therapeutic Medicine Research Center, Cyclotron Research Center, Institute for Medical Sciences, Biomedical Research Institute, Chonbuk National University Medical School and Hospital, Jeonju, Jeonbuk 561-712, Republic of Korea

ARTICLE INFO

Article history: Received 7 December 2015 Received in revised form 11 January 2016 Accepted 11 January 2016 Available online 15 January 2016

ABSTRACT

As a new and preferred light source for phototherapy, blue light emitting diodes (LEDs) with wavelengths of 400-500 nm have been used to treat hyperbilirubinaemia in infantile jaundice [1]. Recent studies report that blue LED irradiation induces apoptosis by stimulating a mitochondrial pathway and reduces the early growth rate of melanoma cells in mice [2]. Here, we detected the induction of apoptotic cell death and formation of autophagosome in human B lymphoma cells after irradiation with blue LED. This paper provides data in support of the research article entitled "Blue light emitting diode induces apoptosis in lymphoid cells by stimulating autophagy" [3].

© 2016 The Authors. Published by Elsevier Inc. This is an open access article under the CC BY license

(http://creativecommons.org/licenses/by/4.0/).

Specifications table

Subject area **Biology** More specific sub-**Apoptosis** ject area

DOI of original article: http://dx.doi.org/10.1016/j.biocel.2015.11.004

* Corresponding author.

E-mail address: jayjeong@jbnu.ac.kr (H.-J. Jeong).

How data was Microscope, flow cytometry

acquired
Data format Raw data, analyzed
Experimental Blue LED irradiation

Experimental Apoptosis of blue LED-irradiated cells was detected by TUNEL staining for features measuring DNA fragments and annexin V/PI staining to distinguish apoptotic

and necrotic cells

Data source Jeonju, Republic of Korea

location

factors

Data accessibility Data are provided in this article

Value of the data

• The data provide information about the effect of blue LED irradiation in B lymphoma cells.

• The data can be used to identify the interaction between autophagy and apoptosis.

• The data inform future study for application of the blue LED to kill cancer cells including lymphoma.

1. Data

Here, we evaluated the increase of fragmented DNA, levels of intracellular superoxide anion $(O_2^{\bullet-})$, LC3 conversion and caspase activation in RAMOS cells. Irradiation with blue LED induced apoptotic cell death through autophagosome activation in RAMOS cells.

2. Experimental design, materials and methods

2.1. Induction of apoptotic cell death under blue LED irradiation in RAMOS cells

The human B cell lymphoma RAMOS cell line was grown in IMDM supplemented with 10% FBS, $10 \,\mu\text{g/ml}$ gentamicin, and $0.25 \,\mu\text{g/ml}$ amphotericin B. We used a blue LED with a wavelength of 450 nm at a power of $6.3 \, \text{mW/cm}^2$ for the experiment. RAMOS cells ($1 \times 10^6 \, \text{cells/well}$) were exposed to blue LED for 4 h. TUNEL staining to detect DNA fragmentation carried out according to the manufacturer's instructions (in situ Apoptosis Detection Kit, Takara, Japan). As shown in Fig. 1, numerous cells were stained green, indicating that apoptotic cell death occurred in RAMOS cells exposed to blue LED. The experiments were performed three times with similar results.

2.2. Autophagosome activation by blue LED irradiation in RAMOS cells

Intracellular $(O_2^{\bullet-})$ levels were measured using an oxidation-sensitive fluorescent probe dye, dihydroethidium (DHE, Ex/Em=518 nm/605 nm), as previously described [3]. The $(O_2^{\bullet-})$ levels were increased up to 349.2% in blue LED-irradiated RAMOS cells at 2 h as compared to the control. Also, the band of cleaved caspase 3 showed weakly after 4 h exposure to blue LED (Fig. 2A). Protein preparation and western blotting were carried out as previously described [2,3]. All experiments were performed at least three times and the results of the immunoblot assays were calculated as relative intensity using Image J software. LC3 conversion levels (LC3-I to LC3-II) were detected. LC3-II level which is positively correlated with the number of autophagosomes [4] increased by 139.3% and 221.3% after 3 h and 4 h of exposure to blue LED (Fig. 2B). To assess interaction between autophagy and apoptosis, 3-methyladenine (3-MA) which blocks the initiation stage of autophagy was pretreated. The cells were suspended in 500 μ l annexin V working solution containing 5 μ l of annexin V-FITC and 10 μ l of




Fig. 1. TUNEL staining of RAMOS cells. Cells were treated with the TdT enzyme and stained with dUTP-fluorescein isothiocyanate using a TUNEL staining kit (Takara Inc., Japan).

propidium iodide (PI) for 15 min, according to the manufacturer's instructions. Fluorescence was measured by flow cytometer (Becton Dickinson, San Jose, CA) and analyzed with Cell Quest software (Becton Dickinson). Moreover, the population of annexin V positive cells was decreased in 3-MA (autophagy inhibitor)-pretreated RAMOS cells as compared to the blue LED group. After blue LED irradiation, the viable cells were 70.7%, the primary necrotic cells were 12.6%, the early apoptotic cells were 4.7%, and the late apoptotic/secondary necrotic cells were 11.9% in 3-MA-pretreated RAMOS cells, respectively.

3. Statistical analysis

All data represent the mean of at least three independent experiments. The results were analyzed using the unpaired two-tailed Student's t-test. P < 0.05 was designated as the level of significance.

Acknowledgments

This work was supported by National Research Foundation of Korea funded by the Ministry of Science, ICT and Future Planning (Nos. 2012M2A2A7014020 and 2015M2A2A6A04044884), and the National R&D Program for Cancer Control, Ministry of Health, Welfare and Family Affairs (0620220).




Fig. 2. Apoptotic effect of blue LED via autophagy in RAMOS cells. (A) Intracellular $O_2^{\bullet-}$ production using DHE staining was analyzed by flow cytometry. (B) Expression of caspase-3 and LC3-I/II under LED irradiation was determined via Western blot analysis. (C) The cells were exposed to blue LED in the presence of 3-MA for 3 h. FACS analysis indicated that cells with apparent early apoptosis (lower right quadrant) and late apoptosis (upper right quadrant) show positive annexin V-FITC staining. *p < 0.05 compared with the control group.

Appendix A. Supplementary material

Supplementary data associated with this article can be found in the online version at http://dx.doi.org/10.1016/j.dib.2016.01.018.

References

- [1] J.F. Ennever, A.F. McDonagh, W.T. Speck, Phototherapy for neonatal jaundice: optimal wavelengths of light, J. Pediatr. 103 (1983) 295–299.
- [2] P.S. Oh, K.S. Na, H. Hwang, H.S. Jeong, S. Lim, M.H. Sohn, H.J. Jeong, Effect of blue light emitting diodes on melanoma cells: involvement of apoptotic signaling, J. Photochem. Photobiol. B 142 (2015) 197–203.
- [3] P.S. Oh, H. Hwang, H.S. Jeong, J. Kwon, H.S. Kim, M. Kim, S. Lim, M.H. Sohn, H.J. Jeong, Blue light emitting diode induces apoptosis in lymphoid cells by stimulating autophagy, Int. J. Biochem. Cell. Biol. 70 (2016) 13–22.
- [4] I. Tanida, N. Minematsu-Ikeguchi, T. Ueno, E. Kominami, Lysosomal turnover, but not a cellular level, of endogenous LC3 is a marker for autophagy, Autophagy 1 (2005) 84–91.