

Lab resource: stem cell line

Episomal plasmid-based generation of an iPSC line from a 79-year-old individual carrying the *APOE4/4* genotype: i11001

Shadaan Zulficar^a, Barbara Fritz^b, Katja Nieweg^{a,*}^a Institute for Pharmacology und Clinical Pharmacy, Faculty of Pharmacy, Philipps University Marburg, Karl-von-Frisch-Str. 1, 35034 Marburg, Germany^b Centre for Human Genetics, Faculty of Medicine, Philipps University Marburg, Baldingerstraße, 35033 Marburg, Germany

ARTICLE INFO

Article history:

Received 14 September 2016

Accepted 22 September 2016

Available online 23 September 2016

ABSTRACT

In this study, lymphoblastoid cells derived from a 79-year old individual with a history of progressive presenile dementia, were used to generate iPSC cells, employing episomal plasmids expressing OCT4, SOX2, KLF4, LIN28, L-MYC and a p53 shRNA. The individual was homozygous for the *APOE4* allele. The resulting iPSC cells had a normal karyotype, retained the *APOE4/4* genotype, expressed pluripotency markers, were free of genomically integrated plasmids, and could be differentiated into cell type representatives from the three germ layers *in vitro*.

© 2016 The Authors. Published by Elsevier B.V. This is an open access article under the CC BY-NC-ND license (<http://creativecommons.org/licenses/by-nc-nd/4.0/>).

Resource table.

| | |
|---------------------------------|--|
| Name of Stem Cell line | i11001 |
| Institution | Philipps University Marburg |
| Person who created resource | Shadaan Zulficar |
| Contact person and email | Katja Nieweg, nieweg@staff.uni-marburg.de |
| Date archived/stock date | Nov 24, 2015 |
| Origin | EBV transformed human B lymphocytes (AG11001 from NIA Aging Cell Culture Repository) |
| Type of resource | Biological reagent: induced pluripotent stem cell, derived from human lymphoblastoid cells |
| Sub-type | Cell line |
| Key transcription factors | OCT4, SOX2, LIN28, L-MYC pCXLE-hOCT3/4-shp53-F (ID 27077) pCSLE-hSK (ID 27078) pCXLE-hUL (ID 27080) |
| Authentication | Identity and purity of cell line confirmed (as shown in Figure1) |
| Link to related literature | http://www.ncbi.nlm.nih.gov/pubmed/21556001 |
| Information in public databases | Not available |

Resource details

The generation of human iPSC line i11001 was carried out using E8 based media and feeder-free conditions. EBV-transformed B-lymphocyte cells from a 79-year old Alzheimer's Disease (AD) patient, generated by Coriell Institute For Medical Research, were employed. The patient had a history of progressive presenile dementia, with a parent and three siblings of the patient also being affected (iPSC line i10984).

The patient was homozygous for the *APOE4* allele. This was confirmed by a PCR-RFLP assay (Hixson and Vernier, 1990) carried out with the lymphoblastoid cell line (Fig. 1A).

The lymphoblastoid cells (LCs) were electroporated with OriP/EBNA-1 based episomal plasmids expressing reprogramming factors OCT4, SOX2, KLF4, L-MYC, LIN28 and a p53 shRNA (Okita et al., 2011). The *APOE4/4* genotype of the generated iPSC line was confirmed by sequencing (Fig. 1A). After several passages, a PCR for episomal vector backbone confirmed that exogenous episomal factors were not being expressed anymore (Fig. 1A).

i11001 iPSC colonies displayed a typical ES-like colony morphology (Fig. 1B). Immunocyto-chemical analysis revealed expression of transcription factors OCT4 and SOX2, and surface markers SSEA4 and TRA-1-81, characteristic of pluripotent stem cells (Fig. 1B). The iPSC cells exhibited a normal karyotype (46, XX) by G-Banding analysis (Fig. 1B). The cells also demonstrated differentiation capacity to three germ layers using an embryoid body-based directed differentiation protocol, followed by immunostaining with markers corresponding to the three germ layers i.e. ectoderm (PAX6), mesoderm (FOXC1) and endoderm (GATA4) (Fig. 1C).

Materials and methods

Cell culture

The LC line was cultured in RPMI-1640 Medium (Gibco™) supplemented with 15% FBS (MP Biomedicals), 1% Penicillin/Streptomycin (Invitrogen™) and 1% GlutaMAX (Invitrogen™) at 37 °C and 5% CO₂.

Derivation of the iPSC cell line and expansion

LCs were reprogrammed by electroporating 200,000 cells with OriP/EBNA-1-based episomal plasmids (the plasmids were a gift from Shinya

* Corresponding author.

E-mail address: nieweg@staff.uni-marburg.de (K. Nieweg).




Fig. 1. Generation and characterization of i11001 line. (A) Confirmation of APOE4/4 status of lymphoblastoid (source) line and i11001 iPS line by PCR-RFLP assay and Sanger sequencing, respectively. (Sequencing trace is reverse complement of sense strand). Product of OriP and EBNA PCR run with a positive control, shows episomal plasmids have been lost by the iPSCs. (B) i11001 cellular morphology resembles ES-like cell morphology. Immunofluorescence staining for pluripotency markers OCT4, SOX2, SSEA4, TR-A1-81 is positive. Karyotype is normal (46, XX). Scale bar: 100 μ m (C) EB formation was induced, followed by directed differentiation to, and immunofluorescence staining for markers of, the three germ layers: Ectoderm: PAX6, Mesoderm: FOXC1, Endoderm: GATA4. Scale bar: 100 μ m.

Table 1
Antibodies used for immunocytochemistry.

| | Antibody | Dilution | Company (Cat #) |
|-------------------------|-------------------------------------|----------|---|
| Pluripotency Markers | Rabbit anti-OCT4 | 1:300 | Cell Signalling Technology #9656 |
| | Rabbit anti-SOX2 | 1:300 | Cell Signalling Technology #9656 |
| | Mouse anti-SSEA4 | 1:300 | Cell Signalling Technology #9656 |
| | Mouse anti-TRA-1-81 | 1:300 | Cell Signalling Technology #9656 |
| Differentiation Markers | Rabbit anti-FOXC1 | 1:100 | Bethyl Laboratories Inc. #A303-520 A |
| | Rabbit anti-GATA4 | 1:200 | Bethyl Laboratories Inc. #A303-503 A |
| | Mouse anti-PAX6 | 1:100 | Developmental Studies Hybridoma Bank #AB 528427 |
| Secondary antibodies | CF488A Goat Anti-Mouse IgG (H + L) | 1:1000 | Biotium Inc. #20010 |
| | CF488A Goat Anti-Rabbit IgG (H + L) | 1:1000 | Biotium Inc. #20019 |

Yamanaka) encoding OCT3/4, SOX2, KLF4, L-MYC, LIN28, and shRNA of p53 (Okita et al., 2011) using a Neon transfection system (Invitrogen™). 1 μ g of each plasmid and 10 μ l tip of the Neon micropipettor kit were used for 2 pulses of 30 ms each at 1100 V. The cells were plated on Matrigel (Corning®) coated plates, in E8 medium without TGF- β , supplemented with 100 μ M sodium butyrate. After 15–17 days, colonies started to appear and when they reached a sufficient size, were manually picked and plated onto Matrigel (Corning®) coated plates, in complete E8 medium. Cells were passaged by exposing to DPBS (without Ca^{++} & Mg^{++}) for a few minutes, then scraping the cell aggregates with a scraper, followed by mechanical trituration and subculturing. Cells were frozen using CryoStor® (Sigma-Aldrich).

Polymerase chain reaction

PCR for OriP and EBNA-1 was carried out using primers as described in Cai et al. (2015). PCR reaction was performed using Kapa HiFi Polymerase (Kapa Biosystems) with cycling conditions: 94 °C for 30 s, 55 °C for 30 s, 72 °C for 1 min for 35 cycles. Episomal vector pCXLE-hUL was used as a positive control. The products were run on a 2% agarose gel.

Genotyping of APOE sequence

PCR of APOE exon 4 was carried out using the primers:

Forward: 5'CAGTCTCTCACACTCGTCTGGC3'; Reverse: 5'CTGCTCCTTCACCTCGTCCA3' PCR reaction was performed using Kapa HiFi Polymerase (Kapa Biosystems) with cycling conditions: 94 °C for 30 s, 61 °C for 30 s, 72 °C for 1 min for 30 cycles.

Sanger sequencing was then performed on the PCR product using the reverse primer by Eurofins Genomics.

APOE genotyping by RFLP assay

Genotyping for APOE status of the LCs was carried out by PCR amplification followed by *HhaI* digestion as described in Hixson and Vernier (1990). An APOE3/4 line DNA was run as a reference.

Karyotyping

Karyotyping of the iPSC cells was performed as described in Campos et al. (2009), followed by G-Banding. Metaphases were imaged using a Zeiss Axioskop, 60 \times , and analysed using Ikaros software (Metasystems). At least 30 metaphases were counted.

Embryoid body (EB) formation and directed differentiation

For EB formation, iPSC colonies were first dissociated to single cells using Accutase® (Sigma-Aldrich), counted and cultured, hanging as

single drops of 2000 cells/20 μ l on lids of petri dishes in E8 medium supplemented with 4 mg/ml polyvinyl alcohol.

For differentiation into germ layer precursors, small molecules were used as described in Frank et al. (2012). Briefly, for differentiation to neuroectoderm, 5 days of treatment with FTDA medium supplemented with 0.5 mM PD0325901, 15 mM SB431542 and 0.5 mM Dorsomorphin was applied. For mesodermal differentiation, treatment for 2 days with 10 ng/ml FGF2, 5 ng/ml Activin A, 10 ng/ml BMP4 (R&D), and 5 mM CHIR99021 was carried out. Endodermal differentiation was induced by treatment with 10 ng/ml FGF2, 100 ng/ml Activin A, 10 ng/ml BMP4, and 5 mM CHIR99021 for 1d, and then with 100 ng/ml Activin A for 3 days. Following this treatment, plated cells were immunostained.

Immuno cytochemistry

Expression of pluripotency and germ layer markers was verified using immunocytochemical staining on cells plated on coverslips. Cells were fixed in 4% PFA (15 min, room temperature (RT)), permeabilized (with the exception of surface markers) with 0.1 M Glycine and 0.25% triton in PBS (10 min, RT). This was followed by blocking in 30% normal goat serum (1 h, RT) and incubation with primary (overnight, 4 °C) and secondary (1 h, RT) antibodies (see Table 1) in 2% BSA/PBS. Nuclei were

stained using 1 μ g/ml DAPI (1 h, RT). Cells were imaged using a Leica DMI6000 B.

Acknowledgements

This work was supported by the Alzheimer Forschung Initiative e.V. (grant agreement 15032).

References

- Cai, J., Orlova, V.V., Cai, X., Eekhoff, E.M.W., Zhang, K., Pei, D., Pan, G., Mummery, C.L., Dijke, P., 2015. [Induced pluripotent stem cells to model human Fibrodysplasia ossificans progressiva](#). *Stem Cell Rep.* 5, 963–970 December 8.
- Campos, P.B., Sartore, R.C., Abdalla, S.N., Rehen, S.K., 2009. [Chromosomal spread preparation of human embryonic stem cells for karyotyping](#). *J. Vis. Exp.* 31, e1512.
- Frank, S., Zhang, M., Schoeler, H.R., Greber, B., 2012. [Small molecule-assisted, line-independent maintenance of human pluripotent stem cells in defined conditions](#). *PLoS One* 7 (7), e41958.
- Hixson, J.E., Vernier, D.T., 1990. [Restriction isotyping of human apolipoprotein E by gene amplification and cleavage with HhaI](#). *J. Lipid. Res.* 31, 545–548 Note on Methodology.
- Okita, K., Matsumura, Y., Sato, Y., Okada, A., Morizane, A., Okamoto, S., Hong, H., Nakagawa, M., Tanabe, K., Tezuka, K.-i., Shibata, T., Kunisada, T., Takahashi, M., Takahashi, J., Saji, H., Yamanaka, S., 2011. [A more efficient method to generate integration-free human iPS cells](#). *Nat. Methods* 8 (5), 409–412.