

RESEARCH NOTE

MYCOLOGY

EUCAST technical note on the EUCAST definitive document EDef 7.2: method for the determination of broth dilution minimum inhibitory concentrations of antifungal agents for yeasts EDef 7.2 (EUCAST-AFST)*

**Maiken C. Arendrup¹, Manuel Cuenca-Estrella²,
Cornelia Lass-Flörl³, William Hope⁴ and
the EUCAST-AFST***

1) Unit of Mycology, Department of Microbiological Surveillance and Research, Statens Serum Institute, Copenhagen, Denmark, 2) Mycology Reference Laboratory, National Center for Microbiology, Instituto de Salud Carlos III, Majadahonda, Spain, 3) Division of Hygiene and Microbiology, Innsbruck Medical University, Innsbruck, Austria and 4) The University of Manchester, Manchester Academic Health Science Centre, NIHR, Translational Research Facility in Respiratory Medicine, University Hospital of South Manchester NHS Foundation Trust, Manchester, UK

Abstract

The European Committee on Antimicrobial Susceptibility Testing-Subcommittee on Antifungal Susceptibility Testing (EUCAST-AFST) has revised the EDef 7.1 document on the method for the determination of broth dilution minimum inhibitory concentrations of antifungal agents for fermentative yeasts. Changes are: dimethylsulphoxide is now the recommended solvent for caspofungin, micafungin and fluconazole; the shelf-life of plates containing the echinocandins prepared from stock solutions in dimethylsulphoxide is extended to 6 months at -80°C ; testing of amphotericin and *Cryptococcus* has been incorporated; and minimum inhibitory concentration ranges for quality control strains and anidulafungin are included.

Keywords: antifungals, *Candida*, *Cryptococcus*, EUCAST, minimum inhibitory concentration, susceptibility testing, yeast

Original Submission: 21 February 2012; **Accepted:** 26 March 2012

Editor: E. Roilides

Article published online: 30 April 2012

Clin Microbiol Infect 2012; **18**: E246–E247

10.1111/j.1469-0691.2012.03880.x

Corresponding author: Maiken C. Arendrup, Unit of Mycology, Department of Microbiological Surveillance and Research, Statens Serum Institute, Artillerivej 5, DK-2300 Copenhagen, Denmark
E-mail: maiken@arendrup.dk

*EUCAST-AFST: MC Arendrup, Chairman (Denmark), WW Hope, Scientific secretary (UK), C Lass-Flörl, Steering Committee (Austria), M Cuenca-Estrella, Steering Committee (Spain), S Arikan (Turkey), F Barchiesi (Italy), J Bille (Switzerland), E Chrysanthou, (Sweden), A Groll (Germany), P Gaustad (Norway), H Järv (Estonia), N Klimko (Russia), O Lortholary (France), T Matos (Slovenia), C Moore (UK), A Schmalreck (Germany), A Velegraki (Greece), P Verweij (the Netherlands).

This Technical Note is based on the EUCAST definitive document EDef 7.2 (available at http://www.eucast.org/antifungal_susceptibility_testing/methods_of_antifungal_susceptibility_testing/susceptibility_testing_of_yeasts/).

The EUCAST Definitive Document EDef 7.1: method for the determination of broth dilution MICs of antifungal agents for fermentative yeasts was approved and published as a discussion document (7.1) in 2003 [1] and the definitive document was published in 2008 [2].

The European Committee on Antimicrobial Susceptibility Testing-Subcommittee on Antifungal Susceptibility Testing (EUCAST-AFST) has revised the document to reflect current knowledge regarding the testing of non-fermentative yeasts (most importantly *Cryptococcus*) and technical details regarding MIC determinations for echinocandins and amphotericin B. This Technical Note is based on the EUCAST definitive document EUCAST definitive document EDef 7.2: *Method for the determination of broth dilution minimum inhibitory concentrations of antifungal agents for yeasts* (available on the EUCAST website: <http://www.eucast.org>). In the full document the revised method is described in detail.

Main updates are as follows. (i) Dimethylsulphoxide (DMSO) is recommended as solvent for fluconazole and all three echinocandins and shelf-life is extended to 6 months. The motivation is that the commercially available fluconazole pure substance and three echinocandins are hydrophobic, that considerable variation in MICs has been demonstrated when water is used as solvent, particularly for caspofungin, and that susceptibility plates are stable at -80°C for up to 6 months when DMSO is used [3]. (ii) For the time being it is recommended that the RPMI 2% glucose is used also for amphotericin B. The non-synthetic broth Antibiotic Medium 3 (AM3), supplemented to a final concentration of 2% glucose, has been evaluated for detecting resistance to amphotericin B [4–7]. However, there is batch-to-batch variation in this medium, and also in the performance of the medium from different manufacturers. (iii) Susceptibility testing of *Cryptococcus* has been included. Main recommendations for the time being are that the EUCAST methodology is adopted and the plates are read when the optical density value is >0.2 . In cases with insufficient growth it is suggested that the test be repeated but that trays are incubated at 30°C [8].

(iv) The following quality control MIC ranges have been included for anidulafungin and *Candida krusei* ATCC 6258: MIC ≤ 0.06 mg/L and for *Candida parapsilosis* ATCC 22019: MIC range 0.25–1.0 mg/L.

Acknowledgements

William Hope is supported by a National Institute of Health (NIHR) Clinician Scientist Award.

Transparency Declaration

The authors do not have any potential conflicts of interests related particularly to this paper. Otherwise, MCA has received research grants and acted as speaker for Astellas, Gilead, MSD and Pfizer, and been a consultant for Gilead, MSD and Pcovery. MCE has received grant support from Astellas Pharma, bioMerieux, Gilead Sciences, Merck Sharp and Dohme, Pfizer, Schering Plough, Soria Melguizo SA, the European Union, the ALBAN program, the Spanish Agency for International Cooperation, the Spanish Ministry of Culture and Education, The Spanish Health Research Fund, The Instituto de Salud Carlos III, The Ramon Areces Foundation, The Mutua Madrileña Foundation. He has been an advisor/consultant to the Panamerican Health Organization, Gilead Sciences, Merck Sharp and Dohme, Pfizer and Schering Plough. He has been paid for talks on behalf of Gilead Sciences, Merck Sharp and Dohme, Pfizer and Schering Plough. CLF has research grants, acted as consultant and/or been on the speakers bureau, for Pfizer, Astellas, Gilead and Merck.

WWH has research grants, acted as consultant and/or been on the speakers bureau, for Pfizer, Astellas, Gilead, Merck, Vectura and F2G.

References

1. Subcommittee on Antifungal Susceptibility Testing (AFST) of the ESC-MID European Committee for Antimicrobial Susceptibility Testing (EUCAST), Rodriguez-Tudela JL, Barchiesi F et al. Method for the determination of minimum inhibitory concentration (MIC) by broth dilution of fermentative yeasts. *Clin Microbiol Infect* 2003; 9: i–viii.
2. Rodriguez-Tudela JL, Arendrup MC, Barchiesi F et al. EUCAST definitive document EDef 7.1: method for the determination of broth dilution MICs of antifungal agents for fermentative yeasts. *Clin Microbiol Infect* 2008; 14: 398–405.
3. Arendrup MC, Rodriguez-Tudela JL, Park S et al. Echinocandin susceptibility testing of *Candida* spp. using the EUCAST EDef 7.1 and CLSI M27-A3 standard procedures: analysis of the influence of Bovine Serum Albumin Supplementation, Storage Time and Drug Lots. *Antimicrob Agents Chemother* 2011; 55: 1580–1587.
4. Cuenca-Estrella M, Diaz-Guerra TM, Mellado E, Rodriguez-Tudela JL. Detection of resistance to amphotericin B in *Candida* isolates by using Iso-Sensitest broth. *Antimicrob Agents Chemother* 2001; 45: 2070–2074.
5. Lozano-Chiu M, Nelson PW, Lancaster M, Pfaller MA, Rex JH. Lot-to-lot variability of antibiotic medium 3 used for testing susceptibility of *Candida* isolates to amphotericin B. *J Clin Microbiol* 1997; 35: 270–272.
6. Rex JH, Cooper CR Jr, Merz WG, Galgiani JN, Anaissie EJ. Detection of amphotericin B-resistant *Candida* isolates in a broth-based system. *Antimicrob Agents Chemother* 1995; 39: 906–909.
7. Wanger A, Mills K, Nelson PW, Rex JH. Comparison of Etest and National Committee for Clinical Laboratory Standards broth macrodilution method for antifungal susceptibility testing: enhanced ability to detect amphotericin B-resistant *Candida* isolates. *Antimicrob Agents Chemother* 1995; 39: 2520–2522.
8. Zaragoza O, Mesa-Arango AC, Gomez-Lopez A, Bernal-Martinez L, Rodriguez-Tudela JL, Cuenca-Estrella M. Process analysis of variables for standardization of antifungal susceptibility testing of nonfermentative yeasts. *Antimicrob Agents Chemother* 2011; 55: 1563–1570.