



Fig.1 Comparison between Fluka MC simulation (red line), dose distribution evaluated with CERR using Double Gaussian model (green line) and dose predicted by standalone computational model (blue line) for a proton initial beam energy of 60MeV, initial σ of 0.4cm, at a normalized depth 0.9 (assuming Bragg peak position at 1)

The preliminary results show that the model is able to predict dose with better accuracy than DGA. So, this model could be a step forward in lateral dose prediction for current TPS and as a future application, a useful tool for online dose-evaluation and adaptive planning; next step will be the implementation of different media and geometries.

Keywords:

Proton-therapy, Treatment Planning System, Molière lateral dose distribution.

References:

- [1] Frühwirth R., Regler. M., Nuclear Instruments and Methods in Physics Research 369 (2000)
- [2] Parodi K., Mairani A., Sommerer F., J. Radiat. Res. J Radiat Res 2013 Jul;54 Suppl 1:i91-6, (2013)
- [3] Bellinzona V. E. et al., Physica Medica doi: 10.1016/j.ejmp.2015.05.004 (2015) [4] Soukup M. et al., Phys. Med. Biol. 50, 5089-5104 (2005)
- [5] Pedroni E. et al., Phys. Med and Biol. 50 541-61(2005) [6] Paganetti H., Phys Med Biol.; 57:R99-117 (2012)
- [7] Bellinzona E. V., Ciocca M., Embriaco A., Ferrari A., Fontana A., Mairani A., Parodi K., Rotondi, A., Sala P., Tessonier T., A model for the accurate computation of the lateral scattering of protons in water, (under review)
- [8] Bohlen T.T. et al., Nuclear Data Sheets 120, 211-214, (2014)
- [9] Ferrari A. et al., CERN-2005-10, INFN/TC05/11, SLAC-R-773, (2005)
- [10] Ulmer W., Rad. Phys. and Chem. 76 1089-107 (2007)
- [11] Deasy J.O. et al., Medical physics 30.5: 979-985(2003)
- [12] Schell S. and Wilkens J. J., Medical Physics, 37(10):5330-40, 2010.

19

A clinical protocol for Simultaneous Integrated Boost for proton treatment of Head and Neck carcinoma

M. F. Belosi, R. Malyapa (MD), A. Bolsi, A. J. Lomax, D. C. Weber (MD)¹

¹ Paul Scherrer Institut, Center for Proton Therapy (5232 Villigen)

Purpose: To define the clinical protocol for Head and Neck (H&N) patients treated with Proton Pencil Beam Scanning (PBS) at Paul Scherrer Institute (PSI, Villigen CH) with Simultaneous Integrated Boost (SIB). The focus was mainly to

define the optimization parameters and the normalization procedure to obtain a steep gradient between the elective and the boost volume, in order to reduce over-dosage in the fall off area of the boost.

Material and Methods: 7 patients, treated at PSI for different H&N carcinomas, were selected for this study. They were all originally treated with two sequential series at 2GyRBE; the 1st one up to 54GyRBE to the elective Planning Target Volume (PTV1), followed by a second series up to 72GyRBE for the boost target volume (PTV2). In this study, all have been re-planned using a SIB regimen of 1.8GyRBE (up to 54GyRBE) to PTV1 and 2.36GyRBE (up to 70.8GyRBE) to PTV2 for a total of 30 fractions. Dose constraints to the OARs were kept as for the sequential approach [1]. All plans were designed on the PSIPlan Treatment Planning System using Intensity Modulated Proton Therapy (IMPT) with non-coplanar fields. The prescription dose of 100% corresponded to the average dose to PTV1, while the dose to PTV2 was optimized prescribing a boosting factor based on the ratio of the two different prescribed dose levels (i.e. 54GyRBE and 70.8 GyRBE). Depending on the ratio of the PTV1 and PTV2 volumes the boosting factor was prescribed to PTV2 or to a smaller structure (GTV or PTV2-3mm). To guarantee dose homogeneity in the dose gradient region between PTV1 and PTV2, the optimized SIB plans were normalized such that 54GyRBE corresponded to the average dose to [PTV1-(Boost volume+3mm)]. This avoided over-dosages close to PTV2 being compensated by under-dosages at the PTV1 margin.

Results: The choice of the boosted volume, the boosting factor, the quality and conformity of the SIB plans were observed to be dependent on the volume ratio PTV2/PTV1. For the analyzed patients, this ranged from 28.8% to 61.7%. For differences in the PTV1 to PTV2 volume larger than 40%, either a PTV2-3mm or the GTV itself were selected for the boosting. Else no gradient could be observed. When compared to the nominal sequential plan, the SIB approach resulted in both a lower mean dose to the ring area (average value for SIB: 55.2 \pm 1.0 GyRBE; average value for sequential: 64.1 \pm 3.6 GyRBE), and 0.6 \pm 0.4% as average ring volume receiving the 95% of the dose prescribed to the boost (the same dose parameter was >10% for all the original plans). A similar sparing of OARs was obtained with both treatment schedules.

Conclusion: Planning H&N patients with SIB optimization resulted in dose distributions which guaranteed the PTV2 and PTV1 coverage and conformity whilst keeping dose to OARs within tolerance. Therefore this approach can be transferred to the clinical operation and has already been applied to a first patient.

Keywords: over-boosting, gradient, boosting factor

References:

- [1] Andrew Lauve, M. Morris, R. Schmidt-Ulrich et al. Simultaneous Integrated Boost intensity-modulated Radiotherapy for locally advanced Head-and-Neck squamous cell carcinomas: II-clinical results. Int J Rad Oncol Biol Phys 2004; 60 (2): 374-387.

20

Gamma Locator for Radionuclide Diagnostics Of Oncological Diseases

A. K. Berdnikova¹, V. N. Belyaev¹, A. I. Bolozdynya¹, V. A. Kantserov¹, V. V. Sosnovtsev¹, K. I. Zhukov²

¹ National Research Nuclear University MEPhI, Moscow, Russian Federation

² P.N. Lebedev Physical Institute of the Russian Academy of Sciences, Moscow, Russian Federation

Gamma locator is a handheld lightweight and compact gamma probe based on a scintillation crystal LaBr₃:Ce and silicone photomultiplier to be used for detection of gamma-radiation emitted by radionuclides such as ^{99m}Tc, ¹²⁵I, ¹¹¹In, ¹⁸F [1].

There are two main applications of gamma locator: intraoperative detection of sentinel lymph nodes and non-invasive scanning the surface of the body. In the first case, a radiotracer is injected into the patient preoperatively and

