

## Gene Section

### Review

# DIO2 (deiodinase, iodothyronine, type II)

Ana Luiza Maia, Simone Magagnin Wajner, Leonardo B Leiria

Thyroid Section, Endocrine Division, Hospital de Clinicas de Porto Alegre (HCPA), Universidade Federal do Rio Grande do Sul, Porto Alegre, RS, Brazil (ALM, SMW, LBL)

Published in Atlas Database: June 2010

Online updated version : <http://AtlasGeneticsOncology.org/Genes/DIO2ID44390ch14q31.html>

DOI: 10.4267/2042/44980

This work is licensed under a Creative Commons Attribution-Noncommercial-No Derivative Works 2.0 France Licence.  
© 2011 Atlas of Genetics and Cytogenetics in Oncology and Haematology

### Identity

**Other names:** 5DII; D2; SelY; TXDI2

**HGNC (Hugo):** DIO2

**Location:** 14q31.1

### DNA/RNA

#### Description

The Dio2 gene is composed of 3 exons comprising 14656 bp of the genomic DNA.

#### Transcription

The length of transcribed mRNA is about 6,8 kb and generates three variants of mRNA. Transcript variant 1 represents the longest transcript and encodes isoform a. Transcript variant 2 differs in the 5'UTR when compared to variant 1. Both variants 1 and 2 encode isoform a. Transcript variant 3 includes an alternate in-frame exon in the coding region, compared to variant 1. Variant 3 encodes isoform b, which is longer than isoform a.


#### Pseudogene

No pseudogene have been described.

### Protein

#### Description

The protein encoded by this gene belongs to the iodothyronine deiodinase family. This enzyme activates thyroid hormone by converting the prohormone thyroxine (T4) by outer ring deiodination to bioactive 3,3',5-triiodothyronine (T3). It is highly expressed in the thyroid, and may contribute significantly to the relative increase in thyroidal T3 production in patients with Graves' disease and thyroid adenomas. This protein contains selenocysteine (Sec) residues encoded by the UGA codon, which often signals the end of process of translation. The 3'UTR of Sec-containing genes have a common stem-loop structure, the sec insertion sequence (SECIS), which is necessary for the recognition of UGA as a Sec codon rather than a stop signal. Alternative splicing results in multiple transcript variants encoding different isoforms. Ubiquitination can also regulate proteins by transiently inactivating enzymatic function through conformational change in a dimeric enzyme, which can be reversed upon deubiquitination (post-translational).



**Organization of the Dio2 gene:** Yellow bars represent the coding region (exon) and red bars, the untranslated region.



**Schematic representation of D2 peptide structure (not on scale).** Isoform a (273 aa) and Isoform b (309 aa). In deep green transmembran domain (position 10-34). In yellow active site (position 133). In deep blue alternative sequence isoform b (position 74).

### Expression

Ohba et al. (2001) identified 2 alternatively spliced DIO2 transcripts that include intronic sequences between the 2 invariant DIO2 exons. These splice variants showed tissue-specific expression in brain, thyroid, liver, thymus, anterior pituitary gland and brown adipose tissue. In mesothelioma cell lysates, Curcio et al. (2001) determined that endogenous DIO2 gene had an apparent molecular mass of 31 kD. In normal tissues, D2 activity/mRNA ratio is variable, but the enzyme is expressed in rodents in the developing and adult testis, heart, muscle, thyroid, BAT, brain, pituitary, thymus, skin, spinal cord, placenta, liver and pancreas. In humans D2 is expressed in brain, BAT, heart, thyroid, muscle, placenta, skin and vascular smooth muscle cells.

### Localisation

Immuno location of the protein in cells showed D2 as an endoplasmic reticulum resident protein.

### Function

Type 2 deiodinase converts intracellular pro-hormone-3,3',5,5'-tetraiodothyronine (T4) into the active thyroid hormone 3,3',5-triiodothyronine (T3) thereby regulating intracellular levels of active T3 in target tissues.

#### Thermogenesis

The expression of D2 is increased in response to cold stimulation in brown adipocytes isolated from mice. Dio2 activation in the brown adipose tissue (BAT) of human newborns and rodents is known to play a role in adaptive energy expenditure during cold exposure.

#### Development

D2 activity is present in human placenta through all pregnancy, and is highly expressed during the first trimester. The level of activity is low in the non-pregnant uterus, but in pregnancy the level rises progressively to a maximum at gestation day 17 when it is increased threefold.

### Homology

Several homologues of Dio2 have been identified in Pan troglodytes and Macaca mulatta (100%). The chicken and mouse have similar domain structures with

human Dio2 (97%). Human Dio2 homology with D3 is expressed in *Sus scrofa*, *Equus caballus*, *Cricetus cricetus*, *Oryctolagus cuniculus*, *Pituophis deppei* (92% similarity) and limited domains with human D3 and D1.

## Mutations

### Note

No germinal or somatic mutations has been described. However, the polymorphism Thr92Ala in Dio2 gene is associated with increased risk of mental retardation, insulin resistance in type 2 diabetic patients, reduced glucose availability in obese women, symptomatic osteoarthritis, Graves' disease and arterial hypertension.

## Implicated in

### Various cancers

### Note

Although not completely understood, Dio2 gene expression and activity is altered in some tumors. It is under-expressed in papillary thyroid carcinomas (PTC). In follicular tumors, D2 activity is similar or elevated when compared to non tumoral tissues, and augmented in follicular adenomas. D2 is also highly expressed in medullary thyroid carcinoma. A higher expression of the Dio2 gene was also described in gliosarcoma, oligoastrocytoma, glioblastoma, oligodendroglioma and pituitary tumors. In contrast, meningioma does not express D2 activity. These differences might be related to the embryonic tumor origin. Mesothelioma expresses higher activity of D2, whereas osteosarcoma has diminished D2 activity.

### Insulin resistance

### Note

Dio2 polymorphism Thr92Ala interacts with a polymorphism in PPAR gamma 2 gene and is associated with insulin resistance in diabetic patients. This Dio2 polymorphism is associated with a ~20% lower rate of glucose disposal in obese women than in non-obese women. Although the association between those two genes occurs in patients with insulin resistance, these results are contradictory in non diabetic population.

## Hypothyroidism

### Note

Disruption in mouse Dio2 gene is associated with alterations in T4/T3 balance with elevated TSH levels, which demonstrates that the Dio2 gene is of critical importance in the feedback regulation of TSH secretion.

## Graves' disease

### Note

It is suggested that the Thr92Ala variant of the Dio2 gene is associated or might be in linkage disequilibrium with a functional DIO2 polymorphism which involves the development of Graves' disease in a Russian population.

## Mental retardation

### Note

A case control study in Chinese patients demonstrated that two allelic intronic SNPs (rs225010 (T/C) and rs225012 (A/G)) in the DIO2 gene could affect the amount of T3 available and in an iodine-deficient environment and partially determine on augmented risk of mental retardation. They found a positive association with mental retardation and the two intronic Dio2 polymorphisms but not with Dio2 Thr92Ala alone and concluded that the genetic variation in Dio2 determine the risk of development of mental retardation that could be due to alterations in the local amount of T3 available in the brain.

## Bone metabolism

### Note

Dio2 is expressed in human and mouse osteoblast cells. In patients with differentiated thyroid carcinoma, the Dio2 Thr92Ala polymorphism is associated with a decreased femoral neck bone mineral density and higher bone turnover independent of serum thyroid hormone levels.

## Cardiomyopathy and arterial hypertension

### Note

Dio2 gene expression is also markedly up-regulated in hearts of mice that develops hypothyroidism or eccentric hypertrophy after myocardial infarction. The Dio2 polymorphism Thr92Ala is also associated with increased risk for the development of hypertension.

## References

Campos-Barros A, Hoell T, Musa A, Sampaolo S, Stoltenburg G, Pinna G, Eravci M, Meinhold H, Baumgartner A. Phenolic and tyrosyl ring iodothyronine deiodination and thyroid hormone concentrations in the human central nervous system. *J Clin Endocrinol Metab.* 1996 Jun;81(6):2179-85

Croteau W, Davey JC, Galton VA, St Germain DL. Cloning of the mammalian type II iodothyronine deiodinase. A selenoprotein differentially expressed and regulated in human and rat brain and other tissues. *J Clin Invest.* 1996 Jul 15;98(2):405-17

Salvatore D, Tu H, Harney JW, Larsen PR. Type 2 iodothyronine deiodinase is highly expressed in human thyroid. *J Clin Invest.* 1996 Aug 15;98(4):962-8

Buettner C, Harney JW, Larsen PR. The 3'-untranslated region of human type 2 iodothyronine deiodinase mRNA contains a functional selenocysteine insertion sequence element. *J Biol Chem.* 1998 Dec 11;273(50):33374-8

Celi FS, Canettieri G, Yarnall DP, Burns DK, Andreoli M, Shuldiner AR, Centanni M. Genomic characterization of the coding region of the human type II 5'-deiodinase gene. *Mol Cell Endocrinol.* 1998 Jun 25;141(1-2):49-52

Araki O, Murakami M, Morimura T, Kamiya Y, Hosoi Y, Kato Y, Mori M. Assignment of type II iodothyronine deiodinase gene (DIO2) to human chromosome band 14q24.2-->q24.3 by in situ hybridization. *Cytogenet Cell Genet.* 1999;84(1-2):73-4

Bartha T, Kim SW, Salvatore D, Gereben B, Tu HM, Harney JW, Rudas P, Larsen PR. Characterization of the 5'-flanking and 5'-untranslated regions of the cyclic adenosine 3',5'-monophosphate-responsive human type 2 iodothyronine deiodinase gene. *Endocrinology.* 2000 Jan;141(1):229-37

Campos-Barros A, Amma LL, Faris JS, Shailam R, Kelley MW, Forrest D. Type 2 iodothyronine deiodinase expression in the cochlea before the onset of hearing. *Proc Natl Acad Sci U S A.* 2000 Feb 1;97(3):1287-92

Murakami M, Araki O, Morimura T, Hosoi Y, Mizuma H, Yamada M, Kurihara H, Ishiuchi S, Tamura M, Sasaki T, Mori M. Expression of type II iodothyronine deiodinase in brain tumors. *J Clin Endocrinol Metab.* 2000 Nov;85(11):4403-6

Curcio C, Baqui MM, Salvatore D, Rihn BH, Mohr S, Harney JW, Larsen PR, Bianco AC. The human type 2 iodothyronine deiodinase is a selenoprotein highly expressed in a mesothelioma cell line. *J Biol Chem.* 2001 Aug 10;276(32):30183-7

Mizuma H, Murakami M, Mori M. Thyroid hormone activation in human vascular smooth muscle cells: expression of type II iodothyronine deiodinase. *Circ Res.* 2001 Feb 16;88(3):313-8

Murakami M, Araki O, Hosoi Y, Kamiya Y, Morimura T, Ogiwara T, Mizuma H, Mori M. Expression and regulation of type II iodothyronine deiodinase in human thyroid gland. *Endocrinology.* 2001 Jul;142(7):2961-7

Murakami M, Kamiya Y, Morimura T, Araki O, Imamura M, Ogiwara T, Mizuma H, Mori M. Thyrotropin receptors in brown adipose tissue: thyrotropin stimulates type II iodothyronine deiodinase and uncoupling protein-1 in brown adipocytes. *Endocrinology.* 2001 Mar;142(3):1195-201

Ohba K, Yoshioka T, Muraki T. Identification of two novel splicing variants of human type II iodothyronine deiodinase mRNA. *Mol Cell Endocrinol.* 2001 Feb 14;172(1-2):169-75

Schneider MJ, Fiering SN, Pallud SE, Parlow AF, St Germain DL, Galton VA. Targeted disruption of the type 2 selenodeiodinase gene (DIO2) results in a phenotype of pituitary resistance to T4. *Mol Endocrinol.* 2001 Dec;15(12):2137-48

Zaninovich AA. [Thyroid hormones, obesity and brown adipose tissue thermogenesis]. *Medicina (B Aires).* 2001;61(5 Pt 1):597-602

Mentuccia D, Proietti-Pannunzi L, Tanner K, Bacci V, Pollin TI, Poehlman ET, Shuldiner AR, Celi FS. Association between a novel variant of the human type 2 deiodinase gene Thr92Ala and insulin resistance: evidence of interaction with the Trp64Arg variant of the beta-3-adrenergic receptor. *Diabetes.* 2002 Mar;51(3):880-3

Curcio-Morelli C, Zavacki AM, Christofollete M, Gereben B, de Freitas BC, Harney JW, Li Z, Wu G, Bianco AC.

- Deubiquitination of type 2 iodothyronine deiodinase by von Hippel-Lindau protein-interacting deubiquitinating enzymes regulates thyroid hormone activation. *J Clin Invest*. 2003 Jul;112(2):189-96
- Wagner MS, Morimoto R, Dora JM, Benneman A, Pavan R, Maia AL. Hypothyroidism induces type 2 iodothyronine deiodinase expression in mouse heart and testis. *J Mol Endocrinol*. 2003 Dec;31(3):541-50
- Chistiakov DA, Savost'yanov KV, Turakulov RI. Screening of SNPs at 18 positional candidate genes, located within the GD-1 locus on chromosome 14q23-q32, for susceptibility to Graves' disease: a TDT study. *Mol Genet Metab*. 2004 Nov;83(3):264-70
- Guo TW, Zhang FC, Yang MS, Gao XC, Bian L, Duan SW, Zheng ZJ, Gao JJ, Wang H, Li RL, Feng GY, St Clair D, He L. Positive association of the DIO2 (deiodinase type 2) gene with mental retardation in the iodine-deficient areas of China. *J Med Genet*. 2004 Aug;41(8):585-90
- Arnaldi LA, Borra RC, Maciel RM, Cerutti JM. Gene expression profiles reveal that DCN, DIO1, and DIO2 are underexpressed in benign and malignant thyroid tumors. *Thyroid*. 2005 Mar;15(3):210-21
- Bianco AC, Maia AL, da Silva WS, Christoffolete MA. Adaptive activation of thyroid hormone and energy expenditure. *Biosci Rep*. 2005 Jun-Aug;25(3-4):191-208
- Canani LH, Capp C, Dora JM, Meyer EL, Wagner MS, Harney JW, Larsen PR, Gross JL, Bianco AC, Maia AL. The type 2 deiodinase A/G (Thr92Ala) polymorphism is associated with decreased enzyme velocity and increased insulin resistance in patients with type 2 diabetes mellitus. *J Clin Endocrinol Metab*. 2005 Jun;90(6):3472-8
- Galton VA. The roles of the iodothyronine deiodinases in mammalian development. *Thyroid*. 2005 Aug;15(8):823-34
- Gouveia CH, Christoffolete MA, Zaitune CR, Dora JM, Harney JW, Maia AL, Bianco AC. Type 2 iodothyronine selenodeiodinase is expressed throughout the mouse skeleton and in the MC3T3-E1 mouse osteoblastic cell line during differentiation. *Endocrinology*. 2005 Jan;146(1):195-200
- Maia AL, Kim BW, Huang SA, Harney JW, Larsen PR. Type 2 iodothyronine deiodinase is the major source of plasma T3 in euthyroid humans. *J Clin Invest*. 2005 Sep;115(9):2524-33
- Morimura T, Tsunekawa K, Kasahara T, Seki K, Ogiwara T, Mori M, Murakami M. Expression of type 2 iodothyronine deiodinase in human osteoblast is stimulated by thyrotropin. *Endocrinology*. 2005 Apr;146(4):2077-84
- Canani LH, Leie MA, Machado WE, Capp C, Maia AL. Type 2 deiodinase Thr92Ala polymorphism is not associated with arterial hypertension in type 2 diabetes mellitus patients. *Hypertension*. 2007 Jun;49(6):e47; author reply e48
- Fiorito M, Torrente I, De Cosmo S, Guida V, Colosimo A, Prudente S, Flex E, Menghini R, Miccoli R, Penno G, Pellegrini F, Tassi V, Federici M, Trischitta V, Dallapiccola B. Interaction of DIO2 T92A and PPARgamma2 P12A polymorphisms in the modulation of metabolic syndrome. *Obesity (Silver Spring)*. 2007 Dec;15(12):2889-95
- Gumieniak O, Perlstein TS, Williams JS, Hopkins PN, Brown NJ, Raby BA, Williams GH. Ala92 type 2 deiodinase allele increases risk for the development of hypertension. *Hypertension*. 2007 Mar;49(3):461-6
- Maia AL, Dupuis J, Manning A, Liu C, Meigs JB, Cupples LA, Larsen PR, Fox CS. The type 2 deiodinase (DIO2) A/G polymorphism is not associated with glycemic traits: the Framingham Heart Study. *Thyroid*. 2007 Mar;17(3):199-202
- Sagar GD, Gereben B, Callebaut I, Mornon JP, Zeöld A, da Silva WS, Luongo C, Dentice M, Tente SM, Freitas BC, Harney JW, Zavacki AM, Bianco AC. Ubiquitination-induced conformational change within the deiodinase dimer is a switch regulating enzyme activity. *Mol Cell Biol*. 2007 Jul;27(13):4774-83
- Wajner SM, dos Santos Wagner M, Melo RC, Parreira GG, Chiarini-Garcia H, Bianco AC, Fekete C, Sanchez E, Lechan RM, Maia AL. Type 2 iodothyronine deiodinase is highly expressed in germ cells of adult rat testis. *J Endocrinol*. 2007 Jul;194(1):47-54
- Gereben B, Zavacki AM, Ribich S, Kim BW, Huang SA, Simonides WS, Zeöld A, Bianco AC. Cellular and molecular basis of deiodinase-regulated thyroid hormone signaling. *Endocr Rev*. 2008 Dec;29(7):898-938
- Maia AL, Hwang SJ, Levy D, Larson MG, Larsen PR, Fox CS. Lack of association between the type 2 deiodinase A/G polymorphism and hypertensive traits: the Framingham Heart Study. *Hypertension*. 2008 Apr;51(4):e22-3
- Meyer EL, Goemann IM, Dora JM, Wagner MS, Maia AL. Type 2 iodothyronine deiodinase is highly expressed in medullary thyroid carcinoma. *Mol Cell Endocrinol*. 2008 Jul 16;289(1-2):16-22
- Watanabe M, Yamamoto T, Mori C, Okada N, Yamazaki N, Kajimoto K, Kataoka M, Shinohara Y. Cold-induced changes in gene expression in brown adipose tissue: implications for the activation of thermogenesis. *Biol Pharm Bull*. 2008 May;31(5):775-84
- Williams AJ, Robson H, Kester MH, van Leeuwen JP, Shalet SM, Visser TJ, Williams GR. Iodothyronine deiodinase enzyme activities in bone. *Bone*. 2008 Jul;43(1):126-34
- Panicker V, Saravanan P, Vaidya B, Evans J, Hattersley AT, Frayling TM, Dayan CM. Common variation in the DIO2 gene predicts baseline psychological well-being and response to combination thyroxine plus triiodothyronine therapy in hypothyroid patients. *J Clin Endocrinol Metab*. 2009 May;94(5):1623-9
- Heemstra KA, Hofstijzer H, van der Deure WM, Peeters RP, Hamdy NA, Pereira A, Corssmit EP, Romijn JA, Visser TJ, Smit JW. The type 2 deiodinase Thr92Ala polymorphism is associated with increased bone turnover and decreased femoral neck bone mineral density. *J Bone Miner Res*. 2010 Jun;25(6):1385-91
- Wang YY, Morimoto S, Du CK, Lu QW, Zhan DY, Tsutsumi T, Ide T, Miwa Y, Takahashi-Yanaga F, Sasaguri T. Up-regulation of type 2 iodothyronine deiodinase in dilated cardiomyopathy. *Cardiovasc Res*. 2010 Sep 1;87(4):636-46

---

*This article should be referenced as such:*

Maia AL, Wajner SM, Leiria LB. DIO2 (deiodinase, iodothyronine, type II). *Atlas Genet Cytogenet Oncol Haematol*. 2011; 15(3):262-265.

---