

Engineering Conferences International ECI Digital Archives

Cell Culture Engineering XV

Proceedings

Spring 5-10-2016

Impact of harvest conditions on the glycosylation profile of a therapeutic antibody

Raghavan Venkat
MedImmune

Follow this and additional works at: http://dc.engconfintl.org/cellculture_xv



Part of the [Biomedical Engineering and Bioengineering Commons](#)

Recommended Citation

Raghavan Venkat, "Impact of harvest conditions on the glycosylation profile of a therapeutic antibody" in "Cell Culture Engineering XV", Robert Kiss, Genentech Sarah Harcum, Clemson University Jeff Chalmers, Ohio State University Eds, ECI Symposium Series, (2016). http://dc.engconfintl.org/cellculture_xv/92

This Abstract is brought to you for free and open access by the Proceedings at ECI Digital Archives. It has been accepted for inclusion in Cell Culture Engineering XV by an authorized administrator of ECI Digital Archives. For more information, please contact franco@bepress.com.

It is essential to maintain the glycosylation profile of a protein therapeutic as it can affect its efficacy (ADCC or CDC), pharmacokinetics/pharmacodynamics, safety, and/or immunogenicity. Pilot scale experiments for an antibody manufacturing process, which used low pH during harvest step to improve the robustness of a subsequent downstream filtration step, revealed the presence of certain peaks that were not present in the harvested material without such low pH treatment. Analytical characterization of these peaks identified them to be associated with M6, M7, M8, and M9 high mannose glycoforms. Further experiments, which looked into the individual effects of the feed pump and continuous centrifuge in the pilot scale harvest set up, concluded that the glycoforms were enriched in high mannose forms during the continuous centrifugation step. Intriguingly, the same was not observed when the experiments were done using a bench top centrifuge. This led to the hypothesis that the cells were exposed to varying shear rates across scales. The hypothesis was subsequently tested using the pilot scale centrifuge and a capillary shear device, in which the cells were exposed to different levels of shear stress. The data from the experiments confirmed that it was indeed the combination of high shear and low pH that resulted in the enrichment of high mannose glycoforms in the centrifuged material. Results from subsequent experiments also suggested that the enrichment potentially occurred due to the release of immature glycoforms associated with the endoplasmic reticulum (ER) and Golgi apparatus.