Frontiers in Molecular Neuroscience 2013 NAPR 2013

Calibration and functional analysis of three genetically encoded CI-/pH sensors

Mukhtarov M., Liguori L., Waseem T., Rocca F., Buldakova S., Arosio D., Bregestovski P. Kazan Federal University, 420008, Kremlevskaya 18, Kazan, Russia

Abstract

Monitoring of the intracellular concentrations of CI- and H+ requires sensitive probes that allow reliable quantitative measurements without perturbation of cell functioning. For these purposes the most promising are genetically encoded fluorescent biosensors, which have become powerful tools for non-invasive intracellular monitoring of ions, molecules and enzymatic activity. A ratiometric CFP/YFP-based construct with a relatively good sensitivity to CI- has been developed (Markova et al., 2008; Waseem et al., 2010). Recently, a combined Cl-/pH sensor (ClopHensor) opened the way for simultaneous ratiometric measurement of these two ions (Arosio et al., 2010). ClopHensor was obtained by fusion of a red-fluorescent protein (DsRedmonomer) to the E2GFP variant that contains a specific Cl--binding site. This construct possesses pKa = 6.8 for H+ and Kd in the 40-50 mM range for Cl- at physiological pH (\sim 7.3) As in the majority of cell types the intracellular CI- concentration ([CI-]i) is about 10 mM, the development of sensors with higher sensitivity is highly desirable. Here we report the intracellular calibration and functional characterization of ClopHensor and its two derivatives: the membrane targeting PalmPalm-ClopHensor and the H148G/V224L mutant with improved Claffinity, reduced pH dependence and pKa shifted to more alkaline values. For functional analysis, constructs were expressed in CHO cells and [Cl-]i was changed by using pipettes with different CI- concentrations during whole-cell recordings. Kd values for CI- measured at 33°C and pH ~ 7.3 were, respectively, 39 mM, 47 mM and 21 mM for ClopHensor, PalmPalm-ClopHensor and the H148G/V224L mutant. PalmPalm-ClopHensor resolved responses to activation of CI--selective glycine receptor channels better than did ClopHensor. Our observations indicate that these different ClopHensor constructs are promising tools for noninvasive measurement of [CI-]i in various living cells. © 2013 Mukhtarov, Liguori, Waseem, Rocca, Buldakova, Arosio and Bregestovski.

http://dx.doi.org/10.3389/fnmol.2013.00009

Keywords

Fluorescent biosensors, Intracellular chloride, Intracellular pH, Non-invasive monitoring, Optogenetics