## Solid-State NMR Study of Alkali Ion Dynamics in Layered Na<sub>2</sub>Ti<sub>3</sub>O<sub>7</sub> and Li<sub>2</sub>Ti<sub>3</sub>O<sub>7</sub>

Kai Volgmann<sup>1,4</sup>, Vanessa Werth<sup>1</sup>, Suliman Nakhal<sup>2</sup>, Martin Lerch<sup>2</sup>, Thomas Bredow<sup>3</sup>. Paul Heitians<sup>1,4</sup>

<sup>1</sup> Institut für Physikalische Chemie und Elektrochemie, Universität Hannover, Callinstr. 3-3a, 30167 Hannover, Germanu:<sup>2</sup> Institut für Chemie, Technische Universität Berlin, Straße des

17. Juni 135, 10623 Berlin, Germany; <sup>3</sup> Mulliken Center for Theoretical Chemistry, Universität

Bonn, Beringstr. 4-6, 53115 Bonn, Germany<sup>4</sup> Zentrum für Festköperchemie und Neue Materialien, Leibniz Universität Hannover, Callinstr. 3, 30167 Hannover, Germany

*E-mail:* volqmann@pci.uni-hannover.de

In  $Na_2Ti_3O_7$  having monoclinic layer structure the sodium can be replaced completely with lithium to achieve monoclinic  $Li_2Ti_3O_7$  [1, 2]. Both materials exhibit interesting properties such as zero-strain behavior upon intercalation and high cycling stability in case of potential battery application. However, there is still a lack of fundamental understanding of the ion diffusivity of both Na and Li in the corresponding host structure. Solid-state nuclear magnetic resonance (NMR) spectroscopy is used here to study the dynamics of both alkali ions in layered  $Na_2Ti_3O_7$  and  $Li_2Ti_3O_7$ . This includes activation energies for the ionic motion and jump rates on the microscopic scale from <sup>23</sup>Na and <sup>7</sup>Li NMR spin-lattice relaxation (SLR), <sup>7</sup>Li spin-alignment echo (SAE), and 2D <sup>6</sup>Li NMR exchange techniques. Moreover, the dimensionality of the ionic motion is investigated by frequency-dependent <sup>7</sup>Li NMR SLR. Structural details are studied using magic-angle spinning (MAS) <sup>23</sup>Na and <sup>7</sup>Li NMR spectroscopy. Results for the electric field gradient at the Na and Li site, respectively, are compared with those from theoretical calculations performed within this study. The dynamics are similar for both cations, and the frequency-dependence of the <sup>7</sup>Li NMR SLR rate indicates Li motion confined to two dimensions. Thus, these two materials may be regarded a model system for low-dimensional diffusion [2]. Another system studied by us also offering low-dimensional diffusion pathways both for Li and Na ions is mixed-alkali  $\text{Li}_x \text{Na}_{2-x} \text{Ti}_6 O_{13}$  [3].

- Synthesis, Structure, and Electrochemical Li-ion Intercalation Properties of Li<sub>2</sub>Ti<sub>3</sub>O<sub>7</sub> with [1]Na<sub>2</sub>Ti<sub>3</sub>O<sub>7</sub>-type Layered Structure, K. Chiba, N. Kijima, Y. Takahashi, Y. Idemoto, J. Akimoto, Solid State Ion. 178 (2008) 1725.
- Solid-State NMR Spectroscopy Study of Cation Dynamics in Layered Na<sub>2</sub>Ti<sub>3</sub>O<sub>7</sub> and Li<sub>2</sub>Ti<sub>3</sub>O<sub>7</sub>, [2]K. Volgmann, V. Werth, S. Nakhal, M. Lerch, T. Bredow, P. Heitjans, Z. Phys. Chem., in press, 2017 (DOI 10.1515/zpch-2016-0948).
- Solid-State NMR to Study Translational Li Ion Dynamics in Solids with Low-Dimensional Diffusion [3] Pathways, K. Volgmann, V. Epp, J. Langer, B. Stanje, J. Heine, S. Nakhal, M. Lerch, M. Wilkening, P. Heitjans, Z. Phys. Chem., submitted, 2017.

