Vol. 120, n. 1 (Supplement): 39, 2015

## Ultrastructural features of human sperm cells cryopreserved by different methods

Claudia Omes<sup>1</sup>, <u>Antonia Icaro Cornaglia<sup>2</sup></u>, Federica Riva<sup>2</sup>, Monica Savio<sup>3</sup>, Giuliano Mazzini<sup>4</sup>, Anna Lisa Marchetti<sup>5</sup>, Rossella Chiara Sanarica<sup>1</sup>, Carmine Tinelli<sup>6</sup>, Aurora Farina<sup>2</sup>, Marco Casasco<sup>2</sup>, Alberto Calligaro<sup>2</sup>, Rossella Elena Nappi<sup>7</sup>, Andrea Casasco<sup>2</sup>

<sup>1</sup>Research Center for Reproductive Medicine, Obstretrics and Gynecology Unit- Maternal Infant Dept., Fondazione IRCCS Policlinico San Matteo, Pavia, Italy - <sup>2</sup>Histology and Embryology Unit- Public Health, Experimental and Forensic Medicine Dept., University of Pavia, Pavia, Italy - <sup>3</sup>Immunology and General Pathology Unit- Molecular Medicine Dept., University of Pavia, Pavia, Italy - <sup>4</sup>IGM-CNR- Biology and Biotechnology Dept., University of Pavia, Pavia, Italy - <sup>4</sup>IGM-CNR- Biology and Biotechnology Dept., University of Pavia, Pavia, Italy - <sup>6</sup>Scientific Direction, Fondazione IRCCS Policlinico San Matteo, Pavia, Italy - <sup>6</sup>Scientific Direction, Clinical Epidemiology and Biometric Unit, Fondazione IRCCS Policlinico San Matteo, Pavia, Italy - <sup>7</sup>Clinical, Surgery, Diagnostic and Pediatric Sciences, University of Pavia, Pavia, Italy

Cryopreservation of human spermatozoa has been recognized as a key strategy for management of male fertility. Nevertheless, current protocols of sperm freezing are neither optimal nor standardized between different labs (1).

In this study we compare the ultrastructural features of human normospermic sperm samples (according to WHO parameters 2010) from 5 different freezing techniques in order to identify the best methods of cryopreservation. After informed consent, 21 normospermic patients (from the Medically Assisted Procreation PMA Center of the Fondazione IRCCS Policlinico San Matteo in Pavia) were recruited and both traditional and improved analysis of sperm quality were applied, in order to define critical steps of cryopreservation.

Cryopreservation of human spermatozoa has been related to decreased motility associated with impaired velocity and viability of sperm pre-freeze and post-thaw. For all applied methods there was a significant reduction of progressive and total motility (P) as a result of freezing. To investigate ultrastructural details, 5 additional cryopreserved samples by the best two freezing methods were analyzed with electron microscopy (TEM). Preliminary data showed the minimal differences between the protocols, with a large number of queues detached and large quantities of cytoplasmic debris after of the first protocol. Spermatozoa appear to be better preserved in the second analyzed method, despite both procedures induced deteriorations at ultrastructural level (2). Other non-routine analysis will be performed to determine whether the cooling time to +4°C may affect the procedure; Comet Assay (to assess the degree of sperm DNA fragmentation) (3) and flow cytometry (to study light scatters patterns and membrane integrity) (4) will be applied.

## References

- Omes C et al. (2013) Human spermatozoa cryopreservation: comparison of three different freezing protocols. Cryoletters, 34 (5), 535-543.
- [2] Barthelemy C. et al. (1990) Ultrastructural changes in membranes and acrosome of human sperm during cryopreservation. System Biology in Reproductive Medicine, 25, 29-40.
- [3] Simon L. et al. (2013) Sperm DNA Damage Measured by Comet Assay. Spermatogenesis 927, 137-146.
- [4] Ferrara F. et al. (1997) Automation of human sperm cell analysis by flow cytometry. Clinical Chemistry, 53(5), 801-807.

<sup>© 2015</sup> Firenze University Press

http://www.fupress.com/ijae