



POLITECNICO DI TORINO
Repository ISTITUZIONALE

Thermal diffusivity and conductivity in low-conducting materials: a new technique

Original

Thermal diffusivity and conductivity in low-conducting materials: a new technique / SPARAVIGNA A.; OMINI M; PASQUARELLI A; STRIGAZZI A. - In: INTERNATIONAL JOURNAL OF THERMOPHYSICS. - ISSN 0195-928X. - 13(1992), pp. 351-351.

Availability:

This version is available at: 11583/1405998 since:

Publisher:

Springer

Published

DOI:

Terms of use:

openAccess

This article is made available under terms and conditions as specified in the corresponding bibliographic description in the repository

Publisher copyright

(Article begins on next page)



Athens Authentication Point

Recognized as:

Sistema Bibliotecario del
Politecnico di Torino
(314-36-772)

1241 CILEA LN Italy
(704-88-862)

Open Consortium Italy
(589-89-716)

3960 CILEA-CIBER e-books
Italy (818-59-898)

Welcome!

To use the personalized
features of this site, please
log in or **register**.

If you have forgotten your
username or password, we
can help.

My SpringerLink

Marked Items

Alerts

Order History

Saved Items

All

Favorites



Content Types Subject Collections

English

Back to: All Search Results

Previous Result Next Result

Journal Article



Thermal diffusivity and conductivity in low-conducting materials: A new technique

Journal International Journal of Thermophysics
Publisher Springer Netherlands
ISSN 0195-928X (Print) 1572-9567 (Online)
Issue Volume 13, Number 2 / March, 1992
DOI 10.1007/BF00504442
Pages 351-360
Subject Collection Physics and Astronomy
SpringerLink Date Friday, September 30, 2005

Add to marked items

Add to shopping cart
Add to saved items
Request Permissions
Recommend this article



Disable Highlighting

Find

more options

**A. Sparavigna¹, M. Omini¹, A. Pasquarelli¹ and A. Strigazzi¹**

(1) Dipartimento di Fisica, Politecnico di Torino, CISM and INFN, Unità di Torino, Corso Duca degli Abruzzi 24, 10129 Torino, Italy

Received: 22 April 1991

Abstract A comparative method is presented, suitable to measure both thermal diffusivity and conductivity of low-conducting solids. The repeatability of the measurements of thermal conductivity is 3%, whereas for diffusivity is 6%. Data for some low-conducting materials are given, consistent with those reported in the literature.

Key words low-conducting materials - thermal conductivity - thermal diffusivity

References secured to subscribers.

 ...

- Within all content
 Within this journal
 Within this issue

Export this article

Export this article as RIS | Text

Text

PDF

The size of this document is 386 kilobytes. Although it may be a lengthier download, this is the most authoritative online format.

Open: Entire document

Ads by Google

Measurement services

Thermal Conductivity /
Diffusivity / CTE / CP
measurement services
www.tprl.com

Coating Services

Leading Supplier of
Plasma & HVOF Ceramic,
Metal and Carbide
Coatings
www.plasmathermalcoatings.com

Thermal Conductivity

of fluids. Precise
measurements according
to ASTM D 2717.
www.psl-systemtechnik.de

Selective gas sensors

For dissolved CO₂, N₂,
He, H₂ and D₂ low level
measurement
www.techtura.com

Conductivity Meter Search

Free Technical Search
Engine Search Thousands
of Catalogs Today
www.global-spec.com



[Frequently asked questions](#) | [General information on journals and books](#) | [Send us your feedback](#) | [Impressum](#)

© Springer. Part of Springer Science+Business Media

[Privacy](#), [Disclaimer](#), [Terms and Conditions](#), © [Copyright Information](#)

Remote Address: 130.192.10.82 • Server: mpweb20

HTTP User Agent: Mozilla/4.0 (compatible; MSIE 7.0; Windows NT 5.1; .NET CLR 2.0.50727; InfoPath.1)

