Kenyon College Digital Kenyon: Research, Scholarship, and Creative Exchange

Kenyon Summer Science Scholars Program

Summer Student Research Scholarship

Summer 2014

Kinetic Characterization of CHY4 and CHY5 in the Catabolism of Valine in Arabidopsis thaliana

Katherine C. Hoener

Follow this and additional works at: https://digital.kenyon.edu/summerscienceprogram Part of the <u>Chemistry Commons</u>


Recommended Citation

Hoener, Katherine C., "Kinetic Characterization of CHY4 and CHY5 in the Catabolism of Valine in Arabidopsis thaliana" (2014). *Kenyon Summer Science Scholars Program.* Paper 254. https://digital.kenyon.edu/summerscienceprogram/254

This Poster is brought to you for free and open access by the Summer Student Research Scholarship at Digital Kenyon: Research, Scholarship, and Creative Exchange. It has been accepted for inclusion in Kenyon Summer Science Scholars Program by an authorized administrator of Digital Kenyon: Research, Scholarship, and Creative Exchange. For more information, please contact noltj@kenyon.edu.

Kinetic Characterization of CHY4 and CHY5 in the Catabolism of Valine in Arabidopsis thaliana

Katherine C. Hoener '15 and Dr. Kerry Rouhier Kenyon College Department of Chemistry



Introduction


Results

CHY5

Valine is an essential branched chain amino acid that is used as an energy source during times of plant growth and extended periods of darkness when normal energy sources are not available.¹ The degradation of valine yields propionyl-CoA, which may be further metabolized through respiration. Unlike other organisms, in which valine degradation takes place exclusively in the mitochondria, plants possess enzymes involved in valine degradation in both the mitochondria and peroxisomes.² This study specifically investigates the enzymes that catalyze the conversion of β -hydroxyisobutyryl-CoA to β -hydroxyisobutyrate in an irreversible process. In the plant species A. thaliana, there are eight isoforms of the hydrolase enzymes that are believed to catalyze this reaction. Deemed CHY1-8, CHY1-3 have peroxisomal targeting sequences while CHY4-6 have mitochondrial leader sequences, and all are part of a larger enzyme family that encompass many functions.

To date, only one peroxisomal hydrolase, CHY1, has been classified as a HIBYL-CoA hydrolase.³ In order to confirm that valine degradation does in fact take place in the mitochondria, it is essential to classify each necessary enzyme with mitochondrial targeting sequences for the entire pathway. While previously none of the mitochondrial CHY enzymes have been classified, the present study investigates CHY4 and CHY5, both mitochondrial enzymes. Additionally, CHY4 was chosen for study as it is essential for plant germination. It is thought that in the absence of CHY4, a build up of the toxic intermediate methacryly-CoA proves fatal to the cell.

With previous work done by Scott Waters ('14) and Haley Dugan ('14) focused on expressing and purifying the CHY4 and CHY5 hydrolases, the present study seeks to classify CHY4 and CHY5 to determine if they are specific for this step in the valine degradation pathway, as HIBYL-CoA hydrolases, or if they also function in other pathways in the plant, such as the conversion of propionyl-CoA to acetyl-CoA.







Figure 2. Substrate specificity of CHY5 (1.57 ug) at 37°C in 0.2 M Tris-HCl pH 8, 1 mM DTNB buffer. CHY5 is most active with the HIBYL-CoA substrate.






Figure 3. CHY5 temperature assay in 0.2 M Tris-HCl pH 8, 1 mM DTNB buffer (1.57 ug CHY5, 0.14 mM S-HIBYL-CoA). CHY5 is most active at higher temperatures.




Figure 1. Standard curve of Absorbance vs. concentration of CoA ranging from 0.46-15 uM at 37°C in 0.2 M Tris-HCl pH 8, 1 mM DTNB (y=-11.9x+8.8).

Figure 4. Michaelis-Menten plot of rate vs. concentration of substrate ranging from 0.67-140 uM at 37°C in 0.2 M Tris-HCl pH 8, 1 mM DTNB buffer with 1.57 ug CHY5 (K_m : 9.5 +/- 0.8 uM, V_{max} : 3.60x10⁻² +/- 7x10⁻³ uM/s).

800





Figure 5. CHY5 pH assay in 0.2M Acetate-Mes-Tris-Glycine, 1 mM DTNB buffer at 37°C (1.57 ug CHY5, 0.14 mM S-HIBYL-CoA). CHY5 is most active at pH 8.





Conclusions

- CHY5 exhibited specificity for the HIBYL-CoA substrate and similar kinetic behavior to CHY1, strongly suggesting that it is in fact a HIBYL-CoA hydrolase that functions optimally at pH 8.
- CHY 4 displayed higher activity with isobutyryl-CoA, propionyl-CoA and acetyl-CoA, suggesting that this enzyme may have broader substrate specificity than other HIBYL-CoA hydrolases.
- Future work will re-investigate enzyme dependence on temperature with a modified assay procedure to overcome limitations of the current assay.

References and Acknowledgements

- 1. Harper, A., Miller, R., Block, K. "Branched Chain Amino Acid Metabolism." Annu Rev Nutr, 4 (1984): 409–454.

Figure 6. Substrate specificity of CHY4 (3.2 ug) at 37°C in 0.2 M Tris-HCl pH 8, 1 mM DTNB buffer. CHY4 is most active with isobutyryl-CoA.

Figure 7. Michaelis Michaelis-Menten plot of rate vs. concentration of substrate ranging from 13.4-696 uM at 37°C in 0.2 M Tris-HCl pH 8, 1 mM DTNB buffer with 1.57 ug CHY4 (K_m : 9.5 ± 0.8 uM, V_{max} : 3.21x10⁻² ± 6x10⁻⁴ uM/s).

2. Lucas, R., Filley, J., Erb, J., Graybill, E., Hawes, J. "Perixisomal Metabolism of Propionic Acid and Isobutyric Acid in Plants." J. Boil. Chem. 282.34 (2007): 24980-9. 3. Zolman BK, Monroe-Augustus M, Thompson B, Hawes JW, Krukenberg KA, Matsuda SPT, Bartel B. "CHY1, an Arabidopsis mutant with impaired b- oxidation, is defective in a peroxisomal b-hydroxyisobutyryl-CoA hydrolase." J Biol Chem (2001): 31037– 31046.

Many thanks to Dr. Rouhier for her advise and guidance throughout the summer. I would also like to thank the Kenyon College Summer Science program for their generous funding.