



Two novel synthetic peptoids exhibit rapid in vitro killing of methicillin-resistant *Staphylococcus pseudintermedius*

Damborg, Peter Panduro; Hansen, Paul Robert; Guardabassi, Luca

Publication date:
2011

Document version
Peer reviewed version

Citation for published version (APA):
Damborg, P. P., Hansen, P. R., & Guardabassi, L. (2011). *Two novel synthetic peptoids exhibit rapid in vitro killing of methicillin-resistant Staphylococcus pseudintermedius*. Poster session presented at ASM-ESCMID Conference on Methicillin-resistant Staphylococci in Animals, Washington DC, United States.

Two novel synthetic peptoids exhibit rapid *in vitro* killing of methicillin-resistant *Staphylococcus pseudintermedius*



Corresponding author: Peter Damborg
peda@life.ku.dk, +45 35332725

Peter Damborg¹, P. R. Hansen² & L. Guardabassi¹

Department of Veterinary Disease Biology (1) and Department of Basic Sciences and Environment (2), Faculty of Life Sciences, University of Copenhagen, 1870 Frederiksberg, Denmark.

Introduction

- Methicillin-resistant *S. pseudintermedius* (MRSP) are of increasing concern due to its recent spread amongst dog populations worldwide (Perreten et al., 2010).

- Limited treatment options are available for MRSP, thus alternative treatment is required to combat these multidrug-resistant microorganisms in the future.

- Peptoids are oligomers of N-substituted glycines and may mimic biologically active peptides by showing antimicrobial properties.

- Peptoids resemble peptides except that the side chains are appended to the nitrogen atom of the backbone rather than to the α -carbons. They are also characterised by being more stable *in vitro* towards heat, salt, pH fluctuations and organic solvents compared to peptides. *in vivo*, peptoids are stable to proteolysis whereas peptides are rapidly degraded.




Results

MIC determinations (Fig. 1)

- Low MIC's ranging from 1.56-6.25 μ M were observed for B1 and D2. MICs did not differ between MSSP and MRSP.
- MICs were normally distributed \rightarrow no obvious evidence of resistance in the 50 isolates tested.

Time kill kinetics (Fig 2a and 2b)

- Both B1 and D2 had a concentration-dependant antimicrobial effect on *S. pseudintermedius*.
- B1 acted more rapidly with complete killing at 4 * MIC in 30 min.
- D2 was slightly slower taking 2 hours to kill at 4 * MIC.



Objective

- To assess the *in vitro* efficacy of two newly developed peptoids against MRSP and methicillin-susceptible *S. pseudintermedius* isolates of canine origin

Materials and Methods

- Prior to the study, 10 newly developed peptoids were tested for their *in vitro* efficacy and hemolysis against various bacterial species. A lysine peptoid hybrid (B1) and a pure peptoid (D2) showing both a low minimum inhibitory concentration (MIC) and a low rate of hemolysis against *S. pseudintermedius* were selected for this study.

- B1 and D2 were tested for their *in vitro* activity against 50 clinical *S. pseudintermedius* isolates of canine origin using the broth microdilution method (CLSI, 2008). Ten of these isolates were MRSP and 40 isolates were MSSP, and all isolates had been collected at our diagnostic lab between 2007 and 2010.

- B1 and B2 were also tested in a time kill kinetic study against one representative *S. pseudintermedius* isolate. Survival of a standard inoculum (5×10^5) of this isolate was observed over time after inoculation without peptoid, at 1 * MIC and at 4 * MIC.

Conclusions

- Two novel peptoid compounds were shown to have a rapid concentration-dependant effect against *S. pseudintermedius* of canine origin, irrespective of antibiotic resistance phenotypes.
- The rapid killing resembles the pharmacodynamics of antiseptics but the mechanism of action is unknown.
- The next step will be to test the effect and toxicity of topical formulations of the two compounds *in vivo*, for example in a mice skin infection model.

References

- Clinical Laboratory Standards Institute (CLSI). 2008. Performance Standards for Antimicrobial Disk and Dilution Susceptibility Tests for Bacteria Isolated from Animals; Approved standard, third ed., M31-A3, vol. 28, no. 8. CLSI, Wayne, PA, USA.
- Perreten V, Kadlec K, Schwarz S, Grönlund Andersson U, Finn M, Greko C, Moodley A, Kania SA, Frank LA, Bemis DA, Franco A, Iurescia M, Battisti A, Duim B, Wagenaar JA, van Duijkeren E, Weese JS, Fitzgerald JR, Rossano A, Guardabassi L. 2010. Clonal spread of methicillin-resistant *Staphylococcus pseudintermedius* in Europe and North America: an international multicentre study. J. Antimicrob. Chemother. 65: 1145-54.

Acknowledgements

The study is funded by the Danish Centre for Antibiotic Research and Development (DanCARD). The centre is an interdisciplinary collaboration between 11 Danish research groups with the main objectives to develop new antibiotics and increase the efficacy of existing antibiotics (www.dancardproject.dk).