ORGANIC-INORGANIC HYBRID POLYMER COATINGS WITH CONTROLLED BIOFUNCTIONALITY

Svetlana A. Sukhishvili, Texas A&M University, Department of Materials Science & Engineering, USA svetlana@tamu.edu

Victoria Albright, Texas A&M University, Department of Materials Science & Engineering, USA Alexander Marin, Institute for Bioscience and Biotechnology Research, University of Maryland, USA Alexander K. Andrianov, Institute for Bioscience and Biotechnology Research, University of Maryland, USA

Key Words: polyphosphazenes; self-defensive coatings; layer-by-layer; antibacterial films

Fluorinated polyphosphazenes (FPs) offer important advantages as biocompatible coatings for coronary stents and other biomedical devices. Recently, a new class of FPs has been introduced, which integrates carboxylic^{1,2} or sulfonic acid³ and fluorinated moieties into a single macromolecular structure. Assemblies of such fluorinated polyelectrolytes with polyelectrolytes or charged small functional molecules can offer efficient modulation of




Figure 1 – Schematic representation of construction of LbL films of FPs and a polycation, as well as water contact angles of polyester surgical felt coated with 7layer FP-terminated FP/poly(ethylene imine) films, in which FPs were 86% fluorinated (left) or nonfluorinated (right).¹

hydrophobicity, improved biocompatibility, as well as biofunctionality, such as modulated drug release. Here, we have explored aqueous multilayer polyelectrolyte deposition as a convenient route to nanofabrication of layered coatings built from ionic FPs (iFPs) and polyelectrolytes^{1,2} or small molecule partners. The resulting laver-by-laver (LbL) assemblies displayed controlled film growth, modulated hydrophobicity, swelling, and protein adsorption characteristics. Hydrophobic interactions largely contributed to the formation of LbL films of iFPs with polycations, leading to linear growth and extremely low water uptake. As shown in neutron reflectometry (NR) studies, films of fluorinated polyphospazenes demonstrated superior layering and persistence of such layering in salt solution as compared to control

nonfluorinated polyphospha-zene/polycation films. Hydrophobicity-enhanced ionic pairing between iFP and linear polycations gave rise to large-amplitude oscillations in surface wettability as a function of capping layer. Importantly, hydrophobicity of iFP-capped LbL coatings could be further enhanced by using a highly porous polyester surgical felt rather than planar substrates for film deposition.²

Moreover, because of the unique combination of ionic and hydrophobic properties, iFPs enabled direct LbL assembly with cationic antibiotics at neutral pH – a feature not achievable with traditional synthetic or biological hydrophilic polymers. The amount of antibiotics included in the coatings could be precisely tuned by coating composition and thickness. Importantly, antibiotics could be retained within the coatings for at least 30 days, and released in response to a trigger that indicated the onset of bacterial colonization, such as locally induced acidification. In addition, the coatings have demonstrated extremely low hemolytic activity (<1%), and were nontoxic to fibroblasts. Taken together, the data suggest that iFPs are versatile building blocks for creating surface coatings with controllable interfacial adhesion, wettability, and programmable interactions with biological milieus, including bacteria-triggered release of antimicrobials for prevention of bacterial colonization of surfaces.

Acknowledgement: This work is supported by the National Science Foundation under Award DMR-1808483. References

1. V. Selin, V. Albright, J. F. Ankner, A. Marin, A. K. Andrianov, and S. A. Sukhishvili "Biocompatible Nanocoatings of Fluorinated Polyphosphazenes through Aqueous Assembly", ACS Appl. Mater. Interfaces 2018, 10, 9756–9764.

2. V. Albright, V. Selin, H. Hlushko, A. Palanisamy, A. Marin, A. K. Andrianov, and S.A. Sukhishvili "Fluorinated Polyphosphazene Coatings Using Aqueous Nano-Assembly of Polyphosphazene Polyelectrolytes" in Polyphosphazenes in Biomedicine, Engineering, and Pioneering Synthesis, A. Andrianov and H. Allcock, Eds., ACS Books, a Symposium Series #1298, 2018 Chapter 5, pp. 101-118. https://pubs.acs.org/doi/10.1021/bk-2018-1298.ch005

3. V. Albright, A. Marin, P. Kaner, S. Sukhishvili, and A. Andrianov "A New Family of Water-Soluble Sulfo-Fluoro Polyphosphazenes and their Assembly within Hemocompatible Nanocoatings", ACS Applied Bio Materials 2019, https://pubs.acs.org/doi/abs/10.1021/acsabm.9b00485