



Tansley insight

Optogenetics in plants

Author for correspondence
John M. Christie
Email: john.christie@glasgow.ac.uk

Matias D. Zurbriggen
Email: matias.zurbriggen@uni-duesseldorf.de

Received: 2 August 2020
Accepted: 18 September 2020

John M. Christie¹ and Matias D. Zurbriggen²

¹Institute of Molecular, Cell and Systems Biology, University of Glasgow, Glasgow, G12 8QQ, UK; ²Institute of Synthetic Biology and CEPLAS, University of Duesseldorf, Duesseldorf 40225, Germany

Contents

Summary	3108	V. Channelrhodopsin-based optogenetics in plants	3113
I. Introduction	3108	VI. Optical tuning to modulate photoreceptor sensitivity	3113
II. Plant photoreceptor contributions to the optogenetic toolkit	3109	VII. Conclusions and outlook	3113
III. Optogenetic control of gene expression in plant cells	3112	Acknowledgements	3113
IV. Optogenetic strategies to manipulate K ⁺ transport and plant physiology	3112	References	3114

New Phytologist (2021) **229**: 3108–3115
doi: 10.1111/nph.17008

Key words: light signalling, optogenetics, photoreceptors, plant optogenetics, synthetic biology.

Summary

The last two decades have witnessed the emergence of optogenetics; a field that has given researchers the ability to use light to control biological processes at high spatiotemporal and quantitative resolutions, in a reversible manner with minimal side-effects. Optogenetics has revolutionized the neurosciences, increased our understanding of cellular signalling and metabolic networks and resulted in variety of applications in biotechnology and biomedicine. However, implementing optogenetics in plants has been less straightforward, given their dependency on light for their life cycle. Here, we highlight some of the widely used technologies in microorganisms and animal systems derived from plant photoreceptor proteins and discuss strategies recently implemented to overcome the challenges for using optogenetics in plants.

I. Introduction

Optogenetics refers to the use of light-responsive proteins to control and dissect biological functions. Light-responsive proteins are abundant in nature and have been instrumental in generating an ever-growing suite of tools for optogenetic applications (Goglia & Toettcher, 2019). As these tools are genetically encoded, they can be introduced into any cell type of choice and used, together with optical methods, to artificially manipulate a specific signalling output. Light provides a powerful control tool as it can be delivered within defined spaces, and can elicit rapid, reversible and quantitative effects following photoreceptor activation with minimized toxicity.

The incentive for optogenetics came from the need to develop a technology that could dissect the complexity of the neural circuitry in animals (Deisseroth & Hegemann, 2017). Microbial opsins have contributed considerably in this regard, with these light-activated ion pumps and ion channels being used as first-generation optogenetic tools (Box 1) to introduce electrochemical changes across cellular membranes that can, in turn, excite or inhibit neural activity depending on their ion transport properties. Optogenetics now extends beyond the regulation of membrane potential changes. A second generation of optogenetic tools (Box 2) has been successful in placing a wide range of biological processes under optical control (Kolar & Weber, 2017; Krueger *et al.*, 2019). Many of these tools originate from plant photoreceptor proteins (Ziegler

Box 1 Optogenetics 1.0

Engineering of opsin-based light-gated ion channels into animal cells has revolutionized the neurosciences (Deisseroth & Hegemann, 2017). Optogenetics 1.0 began with the seminal work of Nagel, Bamberg, Hegemann, Boyden and Deisseroth, among others, in the early 2000s, which successfully used algal retinal-binding channel-rhodopsins (ChRs) to control membrane polarization and thus neuronal activity at high spatiotemporal resolution and reversibility. The following 15 yr witnessed an explosion of optogenetic tools and application, fuelled by a wide range of available ChRs and other opsins, including light-activated ion pumps (halorhodopsins, bacteriorhodopsins). The properties of -opsins have been engineered extensively, aided by structure-function analysis, to generate variants with faster or slower activation/deactivation kinetics, light adaptation/bistability, photosensitivity, colour tuning (wavelength dependency), ion selectivity (H^+ , K^+ , Na^+ , Ca^{2+} , Cl^-) and direction of ion transport (inwards vs outwards). Optogenetics 1.0 has led to major advances in our understanding of neuronal activity dynamics regulating cognition, perception and behaviour in health and disease that is likely to facilitate development of novel intervention therapies. Key advances in associated illumination technologies (optrodes) have also added to the various disciplines contributing to Optogenetics 1.0 (biophysics, biochemistry, neurobiology, biomedicine, material sciences and engineering). Engineering retinal cofactor biosynthesis into plant cells will enable broad application of opsin-based tools, opening up novel perspectives for plant developmental and physiological studies.

Box 2 Optogenetics 2.0

Optogenetics 2.0 deals with the engineering and implementation of photoreceptors for the control of cellular processes with light, going beyond ion transport control, and is being applied to a diversity of organisms ranging from bacteria and fungi to animals and plants (Kolar & Weber, 2017). For this purpose, light-sensitive proteins or domains are engineered to form synthetic photoswitches. They generally consist of two modules: a photosensory module coupled to an effector module of desired function. Key examples are given in the main text. To date, a multitude of photoreceptor sources from bacteria, fungi and plants have been incorporated into optogenetic devices. They vary in photobiology (wavelength dependency, kinetics of activation/deactivation, photosensitivity) and molecular mechanism of function. Their introduction into living cells has led to a myriad of applications, including optical control of gene expression, protein stability, subcellular localization of proteins and organelles, enzyme activity, kinase and receptor activity, and biohybrid materials. Multiplexing photoswitches permits us to achieve independent, orthogonal control over two or more processes simultaneously. Optogenetics 2.0 is still progressing at a fast pace, and current bottlenecks are being overcome. For example, if required, cells can be engineered to express the necessary enzymes required for chromophore production. Despite the advances made, a generalized application of advanced engineering strategies integrating mathematical models are required to optimize and customize photoswitches for their dynamic properties, functionalities and cell/organism of interest (Muller *et al.*, 2014a; Ruess *et al.*, 2015; Benzinger & Khamash, 2018; Golonka *et al.*, 2019). The web-based platform optobase (www.optobase.org) (Kolar *et al.*, 2018) provides a valuable resource in cataloguing the diversity of the optogenetic tools available, their different applications and organism of choice.

& Moglich, 2015; Glantz *et al.*, 2016). In this review, we highlight both the contributions made by plant photoreceptors to the optogenetic toolkit and how optogenetics is being used to study and manipulate cellular and physiological processes in plants.

II. Plant photoreceptor contributions to the optogenetic toolkit

Plants are a plentiful source of photoreceptors as light is pivotal for regulating their growth and development. PHYTOCHROMES (PHYA-E) absorb mainly in the red/far-red region (Legris *et al.*, 2019), whereas cryptochromes (CRY1, CRY2), PHOTOTROPINS (PHOT1, PHOT2) and members of the ZEITLUPE family (ZTL, FKF1, LKP2) respond to ultraviolet (UV)/blue wavelengths (Christie *et al.*, 2015). Plants also contain UV-RESISTANCE LOCUS 8 (UVR8) which specifically absorbs in the UV-B region (Christie *et al.*, 2012). Together, they absorb across much of the UV/visible spectrum (Fig. 1) and their photosensory properties exploited, singly or in combination, for a variety of optogenetic applications (Banerjee & Mitra, 2020).

Phytochromes bind a linear tetrapyrrole chromophore (phytochromobilin), synthesized from heme, that enables them to photoconvert between red and far-red absorbing forms known as Pr and Pfr, respectively (Rockwell & Lagarias, 2010). These properties make PHY particularly attractive for optogenetic applications as independent wavelengths can be used to toggle between active (Pfr) and inactive (Pr) states. Many photoreversible dimerization tools

have been created based on PHYs' light-regulated interaction with basic helix–loop–helix (bHLH) proteins known as phytochrome-interacting factors (PIFs) (Fig. 2a). For example, Shimizu-Sato *et al.* (2002) were the first to develop a red light-inducible gene expression system in yeast by fusing protein regions responsible for PHYB–PIF3 interactions to the DNA-binding domain and transactivation domain of GAL4. The utility of PHY-based tools in nonplant systems is limited by the need to supply the chromophore or to synthesize it *in vivo* through genetic manipulation. However, these have been engineered extensively for a variety of applications, including red/far-red control of gene expression, subcellular protein and organelle localization and signalling in mammalian cells (Levskaia *et al.*, 2009; Muller *et al.*, 2013a) and zebrafish (Beyer *et al.*, 2015; Buckley *et al.*, 2016) (Box 2).

Interaction modules derived from CRY2 and the transcription factor CRY2 INTERACTING bHLH 1 (CIB1) have also been used as optogenetic dimerizers (Fig. 2b) in mammalian cells to place protein translocation, transcription and Cre recombinase-mediated DNA recombination under blue light control (Kennedy *et al.*, 2010). CRYs bind flavin adenine nucleotide (FAD) as chromophore and dimerize following their photoactivation to initiate signalling (Wang *et al.*, 2016). The ability to also form homo-oligomers has been exploited to optically modulate protein function through clustering of cellular receptors (Bugaj *et al.*, 2015) or steric interference (Fig. 2c).




Fig. 1 Plant photoreceptors and their wavelength dependency.




Fig. 2 Optogenetic tool examples derived from plant photoreceptors that are widely used in microbial and mammalian cells. Their mode of action in controlling a protein of interest (POI) is illustrated. (a) Photoswitchable PHYB/PIF heterodimerization. (b) Photoreversible CRY2/CIB heterodimerization. (c) Photoreversible CRY2 dimerization/oligomerization. (d) Photoreversible LOV caging/uncaging. (e) Photoreversible UVR8/COP1 heterodimerization. (f) Photoreversible UVR8 monomerization/dimerization.

For instance, light-induced disruption of clathrin-mediated endocytosis could be achieved when the FAD-binding module of CRY2 was fused to the clathrin light chain (Taslami *et al.*, 2016). Conversely, CRY2 dimerization has been used to optically activate Raf1 kinase activity (Wend *et al.*, 2014). This process can be fine-tuned to increase or decrease Raf1 signalling by altering the clustering capacity of CRY2 by site-directed mutagenesis (Duan *et al.*, 2017).

Ultraviolet/blue light sensing by PHOTs and the ZTL family is mediated by a photosensory module known as the LOV domain which binds flavin mononucleotide (FMN) as chromophore (Christie *et al.*, 2015). LOV domains regulate protein activities by changing conformation upon irradiation. The photo-uncaging property of the LOV2-J α photosensory module from oat PHOT1 (Fig. 2d) has been used to develop a range of optogenetic tools (Losi *et al.*, 2018). Fusion of LOV2-J α to the small GTPase Rac1 represses its activity in darkness through steric interference (Wu *et al.*, 2009). Consequently, actin cytoskeletal dynamics can be controlled remotely upon Rac1 activation by blue light. Optogenetic regulation of Rac1 activity has also been achieved by controlling its plasma membrane recruitment using optogenetic dimerizers based on FKFI's light-regulated interaction with GIGANTEA (Yazawa *et al.*, 2009).

UVR8 is the most recently identified plant photoreceptor and is a homodimer that monomerizes upon UVB exposure to interact with its downstream signalling component CONSTITUTIVELY PHOTOMORPHOGENIC 1 (COP1) (Rizzini *et al.*, 2011). UVR8 requires no additional chromophore to sense UV-B and does this by a mechanism that involves tryptophan-mediated disruption of cross-dimer salt bridges (Christie *et al.*, 2012). A number of optogenetic applications have been created by hijacking UVR8-COP1 interactions (Fig. 2e), including UV-B-induced gene expression, which, when combined with red and blue light-regulated systems, can afford multichromatic control of transcription in mammalian cells (Muller *et al.*, 2013b). UVR8 has also been used to sequester protein targets at the endoplasmic reticulum and promote their trafficking through the secretory pathway to the plasma membrane following UV-B-dependent monomerization (Chen *et al.*, 2013) (Fig. 2f).

III. Optogenetic control of gene expression in plant cells

Plant optogenetics faces experimental constraints imposed by the fact that light is essential for shaping plant growth and development. This presents a formidable challenge as a plant's dependency on light would inadvertently activate optogenetic tools (Fig. 3a). That said, optogenetics has been applied successfully in plant cells and whole plants. The first demonstration employed a split transcription factor system to control gene expression in *Nicotiana benthamiana* and *Arabidopsis thaliana* protoplasts derived from PHYB-PIF interactions (Fig. 2a). Here, the N-terminal region of PHYB was fused to the VP16 transactivation domain from *Herpes simplex virus* in addition to a nuclear localization signal, whereas the N-terminus of PIF6 was fused to the E DNA-binding domain (Muller *et al.*, 2014b). Red light was used to induce the expression of a luciferase reporter, a microRNA or a phytohormone receptor in transformed protoplasts following reconstitution of the split


transcription factor complex and binding to a cognate minimal promoter sequence. This system exhibited a high dynamic range over dark-treated controls (*c.* 400-fold) that could be fine-tuned by varying the intensity of red light. Far-red light terminates PHYB-V16-E-PIF6 interactions, thereby switching the system to an off state. Moreover, supplemental far red repressed reporter activity in the presence of white light and offers a means to circumvent unwanted activation of the system under standard growth conditions (Muller *et al.*, 2014b; Ochoa-Fernandez *et al.*, 2016) with the caveat that endogenous PHY activity might also be affected.

A different engineering strategy was recently developed to help overcome this limitation. A system known as plant usable light-switch elements (PULSE) is insensitive to broad-spectrum white light but activated specifically by red light (Ochoa-Fernandez *et al.*, 2020). It combines two different photoswitches to control transcriptional initiation from a synthetic bipartite promoter: one engineered from the LOV-based transcription factor, EL222 (Nash *et al.*, 2011) (blue-off module), and the other based on the PHYB-PIF6 photoswitch (red-on module) discussed earlier. EL222 dimerizes and binds to its target DNA sequence in blue light (Nash *et al.*, 2011) and can inhibit transcriptional activation when fused to the EAR repressor domain designated SRDX (Ochoa-Fernandez *et al.*, 2020). Thus, PULSE enables effective optogenetic manipulation of plant gene expression under standard growth conditions as transcription is repressed in white light by the blue-off module and is activated only in monochromatic red light by the red-on module (Fig. 3b,c). PULSE has already proved successful in whole plants, as well as protoplasts, and has shown promising potential for a variety of applications, including light-dependent control of plant immunity.

Green/yellow light (500–600 nm) is an attractive wavelength for optogenetic control in plants as native photoreceptors absorb minimally within this region (Fig. 1). A green light-responsive system was recently engineered from the bacterial photoreceptor CarH to control gene expression. CarH uses a vitamin B12 derivative, adenosylcobalamin (AdoB12), as chromophore to regulate its function as transcription factor. CarH binds to the DNA operator CarO as a tetramer in darkness, leading to repression of the carotenogenic target genes, whereas green triggers a disassembly of the CarH tetramer to initiate gene expression. The utility of the CarH photoswitch has already proven useful for controlling reporter gene expression in *Arabidopsis* protoplasts (Chatelle *et al.*, 2018). However, a disadvantage of the CarH system for *in vivo* applications is that it relies on addition of the AdoB12 chromophore, as this is not produced in plants. It is also becoming more apparent that green wavelengths are effective in stimulating photoreceptor function in plants (Battle & Jones, 2020).

IV. Optogenetic strategies to manipulate K⁺ transport and plant physiology

Optogenetic tools are now being used to augment plant responses to light with the same wavelength dependency. Fusion of the LOV2-J α photoswitch to the miniature K⁺ channel Kcv from the *Chlorella* virus PBCV-1 resulted in a synthetic, blue light-induced K⁺ channel, designated BLINK1 (Cosentino *et al.*, 2015). BLINK1 is essentially closed in darkness and opens in response to blue light



and is capable of conducting both inward and outward currents depending on the K^+ concentration gradient across the membrane. This tool, along with its successor BLINK2 (Alberio *et al.*, 2018), has been used successfully to inhibit neural function in zebrafish, mice and rats (Cosentino *et al.*, 2015; Alberio *et al.*, 2018). Recent work by Papanatsiou *et al.* (2019) has shown that BLINK1 can improve stomatal movements in *Arabidopsis* when expressed in the guard cells surrounding the stomatal pore. PHOTs trigger blue light-induced stomatal opening by activating the plasma membrane H^+ -ATPase, the activity of which promotes K^+ uptake that is required to drive changes in guard cell volume and turgor. Guard cell-specific expression of BLINK1 works on top of this process by using the same wavelength dependency of PHOTs to enhance stomatal opening kinetics by adding another light-driven K^+ -uptake pathway. BLINK1 was also found to increase the rate of stomatal closure in darkness (Papanatsiou *et al.*, 2019), presumably by promoting K^+ efflux as a consequence of its slow recovery kinetics (Cosentino *et al.*, 2015). As a result, BLINK1-expressing plants produced more biomass, especially under fluctuating light conditions typical of outdoor growth, because their accelerated stomatal responses synchronize better with the changing light environment to enhance carbon assimilation while preserving water status (Papanatsiou *et al.*, 2019). Moreover, with this approach, unwanted stimulation of endogenous photoreceptor signalling becomes no longer an issue, as the same light quality is used to drive both stomatal opening and BLINK1 activation (Fig. 3d).

V. Channelrhodopsin-based optogenetics in plants

The algal retinal-binding protein, channelrhodopsin-2 (ChR2), is the most commonly used optogenetic tool in the neurosciences (Box 1). ChR2 allows the influx of Ca^{2+} and Na^+ when illuminated by blue light, resulting in membrane depolarization and neural activation. ChR2, like other algal and microbial opsins, is fully functional in mammalian cells where its all-trans-retinal chromophore is in plentiful supply. Supplementation with exogenous retinal has now been successful in fully reconstituting a functional derivative of ChR2 when expressed transiently in *N. benthamiana* or stably in transgenic *Arabidopsis*, leading to the first demonstration of ChR2-based optogenetics in plants. By studying ChR2-mediated membrane depolarization, Reyer *et al.* (2020) were able to show that the plasma membrane H^+ -ATPase is a major driver of membrane repolarization. This work gives a tantalizing glimpse of the potential for opsin-based tools (Box 1) in controlling plant ion transport. A requirement for cofactor addition, however, limits the implementation of these tools *in vivo*. Ongoing work towards engineering the biosynthesis of all-trans retinal in plants (G. Nagel, pers. comm.) will thereby facilitate the use of opsin-based tools for controlling ion fluxes, cellular metabolism and physiology by harnessing their full potential (Box 1).

VI. Optical tuning to modulate photoreceptor sensitivity

Optogenetic tools deactivate in darkness (from seconds to hours, depending on the photosensor) by a process called thermal

reversion. Slowing down thermal reversion (slow cycling) through protein engineering can increase photoproduct lifetime, thereby enhancing photosensitivity, whereas the opposite can be achieved when thermal reversion is accelerated (fast cycling). These principles can be used to fine-tune the photocycle kinetics and sensitivity of optogenetic tools (Ziegler & Moglich, 2015). This approach has also been applied to modify plant photoreceptor performance. For example, slow-cycle tuning of ZTL impacts the circadian period in *Arabidopsis* by altering the rate of clock-component degradation (Pudasaini *et al.*, 2017). A similar strategy has been used more recently to create slow-cycling variants of *Arabidopsis* PHOTs. Plants expressing these engineered photoreceptors exhibit more rapid and robust chloroplast movement responses and improved leaf positioning and expansion, leading to enhanced biomass accumulation under light-limiting conditions (Hart *et al.*, 2019). While this approach is technically not optogenetics *per se*, these findings demonstrate its feasibility to alter plant photoreceptor for benefits in growth. Changes in temperature can also impact the rate of thermal reversion and impact photoreceptor performance in plants (Jung *et al.*, 2016; Legris *et al.*, 2016; Burgie *et al.*, 2017; Fujii *et al.*, 2017).

VII. Conclusions and outlook


Plant optogenetics is still in its infancy and awaits extrapolation by incorporating additional approaches from the mammalian field, including optical control of kinase regulation, protein localization and protein stability (Box 2), and the use of other opsin-based tools, besides ChR2, to manipulate membrane potential changes and cellular signalling (Box 1). However, implementing such tools will be challenging as this may require additional protein engineering strategies to circumvent issues that could impede their functionality or provoke unwanted side-effects in plant cells. These include chromophore availability (e.g. a need for engineering cofactor production within plant cells), efficient expression and localization, as well as minimizing any interference between plant-derived tools and endogenous light signalling pathways. Given the high quantitative modulation, spatiotemporal resolution and the reversible control capabilities provided by optogenetic tools, we envisage that combinatorial approaches such as those used in PULSE will facilitate the targeted manipulation and study of biological processes in light-grown plants. Native and synthetic promoters (Ali & Kim, 2019) could also be used to drive tissue-specific and development regulation of optogenetic tools in plants, whereas multichromatic, independent and orthogonal control over the expression of several traits (multiplexed) will provide new ways to interrogate regulatory networks noninvasively, and with unmatched temporal and quantitative control. A generalized application of engineering strategies integrating mathematical models will also aid the optimization and customization of photoswitches for their dynamic properties, functionalities and cell/organism of choice.


Acknowledgements

We apologize to colleagues whose work could not be cited owing to space limitations and are grateful to Georg Nagel for

communicating unpublished data. This study was supported in part by Biotechnology and Biological Sciences Research Council (BBSRC, BB/M002128/1, BB/R001499/1 to JMC), the Gatsby Charitable Foundation (to JMC), the Deutsche Forschungsgemeinschaft (DFG) German Research Foundation (to MDZ) under Germany's Excellence Strategy (CEPLAS – EXC-2048/1 – Project ID 390686111), NEXT-Plant Plant GRK2466, Collaborative Research Center SFB1208 (Project-ID 267205415) and the European Commission – Research Executive Agency (H2020 Future and Emerging Technologies FET-Open Project ID 801041 CyGenTig to MDZ).

ORCID

John M. Christie  <https://orcid.org/0000-0002-9976-0055>

Matias D. Zurbriggen  <https://orcid.org/0000-0002-3523-2907>

References

- Alberio L, Locarno A, Saponaro A, Romano E, Bercier V, Albadri S, Simeoni F, Moleri S, Pelucchi S, Porro A *et al.* 2018. A light-gated potassium channel for sustained neuronal inhibition. *Nature Methods* 15: 969–976.
- Ali S, Kim WC. 2019. A fruitful decade using synthetic promoters in the improvement of transgenic plants. *Frontiers in Plant Science* 10: 1433.
- Banerjee S, Mitra D. 2020. Structural basis of design and engineering for advanced plant optogenetics. *Trends in Plant Science* 25: 35–65.
- Battle MW, Jones MA. 2020. Cryptochromes integrate green light signals into the circadian system. *Plant, Cell & Environment* 43: 16–27.
- Benzinger D, Khammash M. 2018. Pulsatile inputs achieve tunable attenuation of gene expression variability and graded multi-gene regulation. *Nature Communications* 9: 3521.
- Beyer HM, Juillot S, Herbst K, Samodelov SL, Muller K, Schamel WW, Romer W, Schafer E, Nagy F, Strahle U *et al.* 2015. Red light-regulated reversible nuclear localization of proteins in mammalian cells and zebrafish. *ACS Synthetic Biology* 4: 951–958.
- Buckley CE, Moore RE, Reade A, Goldberg AR, Weiner OD, Clarke JDW. 2016. Reversible optogenetic control of subcellular protein localization in a live vertebrate embryo. *Developmental Cell* 36: 117–126.
- Bugaj LJ, Spelke DP, Mesuda CK, Varedi M, Kane RS, Schaffer DV. 2015. Regulation of endogenous transmembrane receptors through optogenetic Cry2 clustering. *Nature Communications* 6: 6898.
- Burgie ES, Bussell AN, Lye SH, Wang T, Hu W, McLoughlin KE, Weber EL, Li H, Vierstra RD. 2017. Photosensing and thermosensing by phytochrome B require both proximal and distal allosteric features within the dimeric photoreceptor. *Scientific Reports* 7: 13648.
- Chatelle C, Ochoa-Fernandez R, Engesser R, Schneider N, Beyer HM, Jones AR, Timmer J, Zurbriggen MD, Weber W. 2018. A green-light-responsive system for the control of transgene expression in mammalian and plant cells. *ACS Synthetic Biology* 7: 1349–1358.
- Chen D, Gibson ES, Kennedy MJ. 2013. A light-triggered protein secretion system. *Journal of Cell Biology* 201: 631–640.
- Christie JM, Blackwood L, Petersen J, Sullivan S. 2015. Plant flavoprotein photoreceptors. *Plant and Cell Physiology* 56: 401–413.
- Christie JM, Gawthorne J, Young G, Fraser NJ, Roe AJ. 2012. LOV to BLUF: flavoprotein contributions to the optogenetic toolkit. *Molecular Plant* 5: 533–544.
- Cosentino C, Alberio L, Gazzarrini S, Aquila M, Romano E, Cerumenati S, Zuccolini P, Petersen J, Beltrame M, Van Etten JL *et al.* 2015. Optogenetics. Engineering of a light-gated potassium channel. *Science* 348: 707–710.
- Deisseroth K, Hegemann P. 2017. The form and function of channelrhodopsin. *Science* 357: ean5544.
- Duan L, Hope J, Ong Q, Lou HY, Kim N, McCarthy C, Acero V, Lin MZ, Cui B. 2017. Understanding CRY2 interactions for optical control of intracellular signaling. *Nature Communications* 8: 547.
- Fujii Y, Tanaka H, Konno N, Ogasawara Y, Hamashima N, Tamura S, Hasegawa S, Hayasaki Y, Okajima K, Kodama Y. 2017. Phototropin perceives temperature based on the lifetime of its photoactivated state. *Proceedings of the National Academy of Sciences, USA* 114: 9206–9211.
- Glantz ST, Carpenter EJ, Melkonian M, Gardner KH, Boyden ES, Wong GK, Chow BY. 2016. Functional and topological diversity of LOV domain photoreceptors. *Proceedings of the National Academy of Sciences, USA* 113: E1442–E1451.
- Goglia AG, Toettcher JE. 2019. A bright future: optogenetics to dissect the spatiotemporal control of cell behavior. *Current Opinion in Chemical Biology* 48: 106–113.
- Golonka D, Fischbach P, Jena SG, Kleeberg JRW, Essen LO, Toettcher JE, Zurbriggen MD, Moglich A. 2019. Deconstructing and repurposing the light-regulated interplay between Arabidopsis phytochromes and interacting factors. *Communication Biology* 2: 448.
- Hart JE, Sullivan S, Hermanowicz P, Petersen J, Diaz-Ramos LA, Hoey DJ, Labuz J, Christie JM. 2019. Engineering the phototropin photocycle improves photoreceptor performance and plant biomass production. *Proceedings of the National Academy of Sciences, USA* 116: 12550–12557.
- Jung JH, Domijan M, Klose C, Biswas S, Ezer D, Gao M, Khattak AK, Box MS, Charoensawan V, Cortijo S *et al.* 2016. Phytochromes function as thermosensors in Arabidopsis. *Science* 354: 886–889.
- Kennedy MJ, Hughes RM, Peteya LA, Schwartz JW, Ehlers MD, Tucker CL. 2010. Rapid blue-light-mediated induction of protein interactions in living cells. *Nature Methods* 7: 973–975.
- Kolar K, Knobloch C, Stork H, Znidaric M, Weber W. 2018. OptoBase: A web platform for molecular optogenetics. *ACS Synthetic Biology* 7: 1825–1828.
- Kolar K, Weber W. 2017. Synthetic biological approaches to optogenetically control cell signaling. *Current Opinion in Biotechnology* 47: 112–119.
- Krueger D, Izquierdo E, Viswanathan R, Hartmann J, Pallares Cartes C, De Renzi S. 2019. Principles and applications of optogenetics in developmental biology. *Development* 146: dev175067.
- Legris M, Ince YC, Fankhauser C. 2019. Molecular mechanisms underlying phytochrome-controlled morphogenesis in plants. *Nature Communications* 10: 5219.
- Legris M, Klose C, Burgie ES, Rojas CC, Neme M, Hiltbrunner A, Wigge PA, Schafer E, Vierstra RD, Casal JJ. 2016. Phytochrome B integrates light and temperature signals in Arabidopsis. *Science* 354: 897–900.
- Levkaya A, Weiner OD, Lim WA, Voigt CA. 2009. Spatiotemporal control of cell signalling using a light-switchable protein interaction. *Nature* 461: 997–1001.
- Losi A, Gardner KH, Moglich A. 2018. Blue-light receptors for optogenetics. *Chemical Reviews* 118: 10659–10709.
- Muller K, Engesser R, Metzger S, Schulz S, Kampf MM, Busacker M, Steinberg T, Tomakidi P, Ehrbar M, Nagy F *et al.* 2013a. A red/far-red light-responsive bistable toggle switch to control gene expression in mammalian cells. *Nucleic Acids Research* 41: e77.
- Muller K, Engesser R, Schulz S, Steinberg T, Tomakidi P, Weber CC, Ulm R, Timmer J, Zurbriggen MD, Weber W. 2013b. Multi-chromatic control of mammalian gene expression and signaling. *Nucleic Acids Research* 41: e124.
- Muller K, Engesser R, Timmer J, Zurbriggen MD, Weber W. 2014a. Orthogonal optogenetic triple-gene control in Mammalian cells. *ACS Synthetic Biology* 3: 796–801.
- Muller K, Siegel D, Rodriguez Jahnke F, Gerrer K, Wend S, Decker EL, Reski R, Weber W, Zurbriggen MD. 2014b. A red light-controlled synthetic gene expression switch for plant systems. *Molecular BioSystems* 10: 1679–1688.
- Nash AI, McNulty R, Shillito ME, Swartz TE, Bogomolni RA, Luecke H, Gardner KH. 2011. Structural basis of photosensitivity in a bacterial light-oxygen-voltage/helix-turn-helix (LOV-HTH) DNA-binding protein. *Proceedings of the National Academy of Sciences, USA* 108: 9449–9454.
- Ochoa-Fernandez R, Abel NB, Wieland FG, Schlegel J, Koch LA, Miller JB, Engesser R, Giuriani G, Brandl SM, Timmer J *et al.* 2020. Optogenetic control of gene expression in plants in the presence of ambient white light. *Nature Methods* 17: 717–725.
- Ochoa-Fernandez R, Samodelov SL, Brandl SM, Wehinger E, Muller K, Weber W, Zurbriggen MD. 2016. Optogenetics in plants: red/far-red light control of gene expression. *Methods in Molecular Biology* 1408: 125–139.

- Papanatsiou M, Petersen J, Henderson L, Wang Y, Christie JM, Blatt MR. 2019. Optogenetic manipulation of stomatal kinetics improves carbon assimilation, water use, and growth. *Science* 363: 1456–1459.
- Pudasaini A, Shim JS, Song YH, Shi H, Kiba T, Somers DE, Imaizumi T, Zoltowski BD. 2017. Kinetics of the LOV domain of ZEITLUPE determine its circadian function in Arabidopsis. *eLife* 6: e21646.
- Reyer A, Hassler M, Scherzer S, Huang S, Pedersen JT, Al-Rasheid KAS, Bamberg E, Palmgren M, Dreyer I, Nagel G *et al.* 2020. Channelrhodopsin-mediated optogenetics highlights a central role of depolarization-dependent plant proton pumps. *Proceedings of the National Academy of Sciences, USA* 117: 20920–20925.
- Rizzini L, Favory JJ, Cloix C, Faggionato D, O'Hara A, Kaiserli E, Baumeister R, Schafer E, Nagy F, Jenkins GI *et al.* 2011. Perception of UV-B by the Arabidopsis UVR8 protein. *Science* 332: 103–106.
- Rockwell NC, Lagarias JC. 2010. A brief history of phytochromes. *ChemPhysChem* 11: 1172–1180.
- Ruess J, Parise F, Miliadis-Argeitis A, Khammash M, Lygeros J. 2015. Iterative experiment design guides the characterization of a light-inducible gene expression circuit. *Proceedings of the National Academy of Sciences, USA* 112: 8148–8153.
- Shimizu-Sato S, Huq E, Tepperman JM, Quail PH. 2002. A light-switchable gene promoter system. *Nature Biotechnology* 20: 1041–1044.
- Taslami A, Zoltowski B, Miranda JG, Pathak GP, Hughes RM, Tucker CL. 2016. Optimized second-generation CRY2-CIB dimerizers and photoactivatable Cre recombinase. *Nature Chemical Biology* 12: 425–430.
- Wang Q, Zuo Z, Wang X, Gu L, Yoshizumi T, Yang Z, Yang L, Liu Q, Liu W, Han YJ *et al.* 2016. Photoactivation and inactivation of Arabidopsis cryptochrome 2. *Science* 354: 343–347.
- Wend S, Wagner HJ, Muller K, Zurbriggen MD, Weber W, Radziwill G. 2014. Optogenetic control of protein kinase activity in mammalian cells. *ACS Synthetic Biology* 3: 280–285.
- Wu YI, Frey D, Lungu OI, Jaehrig A, Schlichting I, Kuhlman B, Hahn KM. 2009. A genetically encoded photoactivatable Rac controls the motility of living cells. *Nature* 461: 104–108.
- Yazawa M, Sadaghiani AM, Hsueh B, Dolmetsch RE. 2009. Induction of protein-protein interactions in live cells using light. *Nature Biotechnology* 27: 941–945.
- Ziegler T, Moglich A. 2015. Photoreceptor engineering. *Frontiers in Molecular Biosciences* 2: 30.



About New Phytologist

- *New Phytologist* is an electronic (online-only) journal owned by the New Phytologist Foundation, a **not-for-profit organization** dedicated to the promotion of plant science, facilitating projects from symposia to free access for our Tansley reviews and Tansley insights.
- Regular papers, Letters, Viewpoints, Research reviews, Rapid reports and both Modelling/Theory and Methods papers are encouraged. We are committed to rapid processing, from online submission through to publication 'as ready' via *Early View* – our average time to decision is <26 days. There are **no page or colour charges** and a PDF version will be provided for each article.
- The journal is available online at Wiley Online Library. Visit **www.newphytologist.com** to search the articles and register for table of contents email alerts.
- If you have any questions, do get in touch with Central Office (np-centraloffice@lancaster.ac.uk) or, if it is more convenient, our USA Office (np-usaoffice@lancaster.ac.uk)
- For submission instructions, subscription and all the latest information visit **www.newphytologist.com**