University of Massachusetts Medical School eScholarship@UMMS

Open Access Articles

Open Access Publications by UMMS Authors

2016-05-24

Prevalent presence of periodic actin-spectrin-based membrane skeleton in a broad range of neuronal cell types and animal species

Jiang He Harvard University

Et al.

Let us know how access to this document benefits you.

Follow this and additional works at: https://escholarship.umassmed.edu/oapubs

Part of the Neuroscience and Neurobiology Commons

Repository Citation

He J, Zhou R, Wu Z, Carrasco MA, Kurshan PT, Farley JE, Simon DJ, Wang G, Han B, Hao J, Heller E, Freeman MR, Shen K, Maniatis T, Tessier-Lavigne M, Zhuang X. (2016). Prevalent presence of periodic actin-spectrin-based membrane skeleton in a broad range of neuronal cell types and animal species. Open Access Articles. https://doi.org/10.1073/pnas.1605707113. Retrieved from https://escholarship.umassmed.edu/oapubs/2832

This material is brought to you by eScholarship@UMMS. It has been accepted for inclusion in Open Access Articles by an authorized administrator of eScholarship@UMMS. For more information, please contact Lisa.Palmer@umassmed.edu.

Prevalent presence of periodic actin-spectrin-based membrane skeleton in a broad range of neuronal cell types and animal species

Jiang He^{a,b,c,1,2}, Ruobo Zhou^{a,b,c,2}, Zhuhao Wu^d, Monica A. Carrasco^{e,3}, Peri T. Kurshan^{f,g}, Jonathan E. Farley^{h,i}, David J. Simon^d, Guiping Wang^{a,b,c}, Boran Han^{a,b,c}, Junjie Hao^{a,b,c}, Evan Heller^{a,b,c}, Marc R. Freeman^{h,i}, Kang Shen^{f,g}, Tom Maniatis^e, Marc Tessier-Lavigne^d, and Xiaowei Zhuang^{a,b,c,4}

^aHoward Hughes Medical Institute, Harvard University, Cambridge, MA 02138; ^bDepartment of Chemistry and Chemical Biology, Harvard University, Cambridge, MA 02138; ^dLaboratory of Brain Development and Repair, The Rockefeller University, New York, NY 10065; ^eDepartment of Biochemistry and Molecular Biophysics, Columbia University Medical Center, New York, NY 10032; Howard Hughes Medical Institute, Stanford University, Stanford, CA 94305; ⁹Department of Biology, Stanford University, Stanford, CA 94305; ^hHoward Hughes Medical Institute, University of Massachusetts Medical School, Worcester, MA 01605; and ¹Department of Neurobiology, University of Massachusetts Medical School, Worcester, MA 01605

Contributed by Xiaowei Zhuang, April 8, 2016 (sent for review March 25, 2016; reviewed by Bo Huang and Feng Zhang)

Actin, spectrin, and associated molecules form a periodic, submembrane cytoskeleton in the axons of neurons. For a better understanding of this membrane-associated periodic skeleton (MPS), it is important to address how prevalent this structure is in different neuronal types, different subcellular compartments, and across different animal species. Here, we investigated the organization of spectrin in a variety of neuronal- and glial-cell types. We observed the presence of MPS in all of the tested neuronal types cultured from mouse central and peripheral nervous systems, including excitatory and inhibitory neurons from several brain regions, as well as sensory and motor neurons. Quantitative analyses show that MPS is preferentially formed in axons in all neuronal types tested here: Spectrin shows a long-range, periodic distribution throughout all axons but appears periodic only in a small fraction of dendrites, typically in the form of isolated patches in subregions of these dendrites. As in dendrites, we also observed patches of periodic spectrin structures in a small fraction of glial-cell processes in four types of glial cells cultured from rodent tissues. Interestingly, despite its strong presence in the axonal shaft, MPS is disrupted in most presynaptic boutons but is present in an appreciable fraction of dendritic spine necks, including some projecting from dendrites where such a periodic structure is not observed in the shaft. Finally, we found that spectrin is capable of adopting a similar periodic organization in neurons of a variety of animal species, including Caenorhabditis elegans, Drosophila, Gallus gallus, Mus musculus, and Homo sapiens.

actin | spectrin | cytoskeleton | neuron | STORM

A ctin is critically involved in the regulation of neuronal polar-ization, differentiation, and growth of neuronal processes, cargo trafficking, and plasticity of synapses (1-3). Spectrin is an actin-binding protein that is important for the development and stabilization of axons and maintenance of neuronal polarization (4-6). In *Caenorhabditis elegans*, spectrin is important for the stability and integrity of axons under mechanical stress (4, 6) and for mechanosensation (6), and spectrin depletion results in axon breakage during animal locomotion (4). In Drosophila, spectrin has been shown to be involved in axonal path finding (7) and stabilization of presynaptic terminals (8). In mice, spectrin null mutations are embryonically lethal, and neurons with spectrin knockdown display defects in axonal initial segment assembly (5, 9, 10).

Actin and spectrin form a 2D polygonal lattice structure underneath the membrane of erythrocytes (11). Recently, a novel form of actin-spectrin-based submembrane skeleton structure was discovered in neuronal axons (12) using superresolution STORM imaging (13, 14). This membrane-associated periodic skeleton (MPS) has been observed in both fixed and live cultured neurons (12, 15, 16) and in brain tissue sections (12). In this structure, short actin filaments are organized into repetitive, ring-like structures that wrap around the circumference of the axon with a periodicity of ~190 nm; adjacent actin rings are connected by spectrin tetramers, and actin short filaments in the rings are capped by adducin (12). This structure also appears to organize other associated molecules, such as ankyrin and sodium channels, into a periodic distribution in axons (12, 16). During neuronal development, MPS originates from the axonal region proximal to the soma and propagates to distal axonal terminals (16). At a relatively late stage during development, specific isoforms of ankyrin and spectrin molecules, ankyrin-G and BIV spectrin, are recruited to the axon initial segment (AIS) (17, 18), and these molecules are also assembled into the MPS structure, adopting a similar periodic distribution (16, 19). As in the AIS, this periodic structure is also present in the nodes of Ranvier (20). This periodic skeletal structure has been shown to preferentially form in axons compared with dendrites in primary neuronal cultures: actin and spectrin typically form a long-range, periodic lattice structure throughout the entire axonal shaft, except for the very distal region near the growth cone, in essentially

Significance

Actin, spectrin, and associated molecules form a submembrane periodic skeleton structure in neurons. In this study, we demonstrate that this membrane-associated periodic skeleton (MPS) is present in a broad range of neuronal cell types cultured from the central and peripheral nervous systems of rodents. The MPS structure is preferentially formed in axons compared with dendrites and is differentially regulated at the pre- and postsynaptic sites of neurons. Our data also suggest that MPS is conserved across a wide range of invertebrate and vertebrate animal species.

The authors declare no conflict of interest.

¹Present address: Institute of Medical Engineering and Sciences, Massachusetts Institute of Technology, Cambridge, MA 02139.

⁴To whom correspondence should be addressed. Email: zhuang@chemistry.harvard.edu.



Author contributions: J. He, R.Z., M.R.F., K.S., T.M., M.T.-L., and X.Z. designed research; J. He, R.Z., Z.W., M.A.C., P.T.K., J.E.F., D.J.S., G.W., J. Hao, and E.H. performed research; J. He, R.Z., G.W., and B.H. analyzed data; and J. He, R.Z., M.R.F., K.S., T.M., M.T.-L., and X.Z. wrote the paper

Reviewers: B.H., University of California, San Francisco; and F.Z., Massachusetts Institute of Technology.

Freely available online through the PNAS open access option.

²J. He and R.Z. contributed equally to this work.

³Present address: Department of Biomedical Science, Faculty of Health Sciences, University of Talca, Talca 3460000, Chile.

This article contains supporting information online at www.pnas.org/lookup/suppl/doi:10. 1073/pnas.1605707113/-/DCSupplemental

all observed axons. In contrast, such a periodic structure was observed in only a small fraction ($\sim 10-30\%$) of dendrites and typically appeared as short, isolated patches in portions of these dendrites (16, 20). The local concentration of spectrin is a key determinant for the preferential formation of MPS in axons: in wild-type neurons, β II spectrin is enriched in axons, and artificially increasing the concentration of β II spectrin through overexpression is sufficient to induce the formation of MPS in all dendrites (16). Ankyrin-B appears to be an important regulator of this structure: in ankyrin-B knockout mice, β II spectrin becomes evenly distributed between axons and dendrites, leading to the formation of the long-range MPS structure in all dendrites (16) without perturbing the MPS structure in axons (16, 21).

The ubiquitous expression of spectrin in the nervous systems of nearly all animal species (22) raises the questions of how widespread the MPS structure is in different nervous system cell types and distinct subcellular compartments and of how conserved this structure is across different animal species. A recent paper reports the presence of periodic actin structures in several nervous system cell types from rodents (23). Here we investigated these questions regarding the prevalence and conservation of the MPS structure by examining the distribution of spectrin in many different types of rodent neurons and glial cells and across a variety of organisms ranging from *C. elegans* to *Homo sapiens*. Furthermore, we examined the distribution of spectrin in presynaptic and postsynaptic compartments of axons and dendrites, respectively, to shed light on the relation between the MPS structure and synapses.


Results

The MPS Structure Is Present in Both Excitatory and Inhibitory Neurons Cultured from Mouse Cortex, Hippocampus, and Midbrain. We initially observed the membrane-associated, periodic actin–spectrin structure in axons of rodent hippocampal and cortical neurons (12). The neurons observed in this earlier study are most likely excitatory pyramidal neurons. To test whether our findings extend to other excitatory neuronal types, or inhibitory neurons in the mouse brain, we first cultured neurons from three different mouse brain regions—the cortex, hippocampus, and midbrain and distinguished excitatory and inhibitory neurons using immunofluorescence against vGlut1 and GAD2, respectively (Fig. S1). We then labeled β II spectrin in these neurons using immunofluorescence and imaged immunolabeled β II spectrin using 3D STORM imaging, a superresolution imaging method that uses stochastic switching and high-precision localization of single

Fig. 1. The MPS structure is present in both excitatory and inhibitory neurons from mouse cortex, midbrain, and hippocampus. Cultured neurons were immunostained with vGlut1 or GAD2 to label the excitatory and inhibitory neurons, respectively. (A) Reconstructed conventional image of a vGlut1+ cortical neuron shown together with 3D STORM images of βII spectrin in an axonal (A-1) and dendritic (A-2) region. The axonal (A-1) and dendritic (A-2) regions correspond to the regions indicated by arrows in A. Dotted lines in A-2 indicate patches of periodic pattern in the dendritic shaft. In 3D STORM images, localizations at different z values are depicted in different colors. (B, B-1, and B-2) Similar to A, A-1, and A-2 but for a GAD2+ cortical neuron. (C) Averaged autocorrelation functions calculated from multiple, randomly selected axonal (black) and dendritic (red) regions in neurons cultured from cortex (n = 18 for axons and n = 9 for dendrites of GAD2+ neurons; n = 9 for axons and n = 8 for dendrites of vGlut1+ neurons), midbrain (n = 8 for axons molecules to obtain subdiffraction-limit image resolution (13, 14). We distinguished axons from dendrites either by morphology, with the axon being a distinctly long process projecting from the soma, or by immunofluorescence, with dendrites being MAP2-positive. We observed highly periodic distributions of BII spectrin in the axons of nearly all (>90%) vGlut1+ excitatory neurons cultured from all three brain regions (Fig. 1A and Figs. S2A, S3A, and S44). Autocorrelation analysis indicated that the spacing was around 190 nm (Fig. 1C and Figs. S2A, S3A, and S4A), consistent with the value determined previously (12, 16). β II spectrin also adopted a periodic distribution in the axons of nearly all (90%) GAD2+ inhibitory neurons from mouse cortex, hippocampus, and midbrain (Fig. 1B and Figs. S2B, S3B, and S4B). We used autocorrelation analysis to systematically quantify the degree of periodicity for BII spectrin in axons of these neurons. The amplitude of the averaged autocorrelation function of many randomly selected segments of axons provided a measure of the degree of periodicity. As shown in Fig. 1C, both the spacing and the degree of periodicity in GAD2+ inhibitory neurons were quantitatively similar to those obtained for vGlut1+ excitatory neurons.

In contrast, β II spectrin exhibited a much less regular distribution in dendritic shafts from the same neurons: unlike in axons, where β II spectrin adopted a long-range, periodic pattern throughout nearly the entire axonal shaft, the β II spectrin distribution exhibited a periodic pattern only occasionally in a small fraction of dendritic shafts, often as short, isolated patches (Fig. 1*A*-2 and Figs. S2*A*-2 and S4 *A*-2 and *B*-2). This observation is similar to what we observed previously in mouse hippocampal neurons (16). The average amplitudes of the autocorrelation functions derived from many randomly selected segments of dendrites were much smaller than those observed for axons in both excitatory and inhibitory neurons from all three brain regions (Fig. 1*C*), quantitatively showing that the degree of periodicity in dendrites was much less than that in axons.

The MPS Structure Is Present in Multiple Specific Subtypes of Excitatory and Inhibitory Neurons in the Mouse Central Nervous System. The proteins vGlut1 and GAD2 are general markers for glutamatergic excitatory and GABAergic inhibitory neurons, respectively. Next, we selected some morphologically and functionally distinct subtypes of excitatory and inhibitory neurons from the mouse brain for further investigation. Specifically, for excitatory neurons, we imaged cerebellar granule cells, dopaminergic neurons from the midbrain, and pontine nuclei excitatory neurons. For inhibitory neurons, we imaged Golgi cells and Purkinje cells from cerebellum, as well as



and n = 7 for dendrites of GAD2+ neurons; n = 13 for axons and n = 11 for dendrites of vGlut1+ neurons), and hippocampus (n = 6 for axons and n = 8 for dendrites of GAD2+ neurons; n = 8 for axons and n = 9 for dendrites of vGlut1+ neurons), respectively. (Scale bar: 20 μ m for conventional images and 2 μ m for STORM images.)

Parvalbumin-positive cells from the cortex, midbrain, and hippocampus. Granule cells and Purkinje cells were chosen because of their distinct morphologies. Granule cells are among the smallest cells in the brain, and Purkinje cells are unique in their large size and elaborate dendritic branches (24). Dopaminergic neurons release the neurotransmitter dopamine and represent a type of excitatory neuron distinct from glutamatergic neurons (24). Parvalbumin cells belong to a specific subtype of GABAergic inhibitory neurons (24). In addition, we imaged olfactory neurons from the olfactory bulb.

Cultured dopaminergic neurons, olfactory neurons, Golgi cells, Purkinje cells, and Parvalbumin cells were identified using immunostaining against tyrosine hydroxylase (25), olfactory marker protein (26), metabotropic glutamate receptors (27), calbindin (28), and parvalbumin, respectively (Figs. S1 and S5), whereas granule cells cultured from cerebellum were distinguished by their lack of calbindin staining and relatively small soma size. We also used the general neuronal markers MAP2 and NeuN to exclude glial cells from our analysis. All neuronal subtypes were cultured for at least 10 d in vitro until a distinctly long axon had projected from the soma to ensure sufficient axon specification. The neurons were then fixed and immunostained with β II spectrin for STORM imaging. Axons were again distinguished from dendrites by either their morphology or lack of MAP2 staining. As shown in Fig. 2 A-I, all of these neuronal subtypes exhibited a periodic distribution of β II spectrin in their axons with a spacing of ~190 nm. Autocorrelation analyses showed that the degree of periodicity (i.e., autocorrelation amplitude) was similar among most of these excitatory and inhibitory neurons (Fig. 2 J and K). Two exceptions were Golgi and Purkinje cells, which showed less regularity (Fig. 2 H, I, and K). This is probably due to a relatively slow development process for these cultured cerebellar neurons as evidenced by less elaborate dendrites observed for these neurons in vitro (29, 30). In addition, all of these specific neuronal subtypes showed a substantially higher propensity for MPS formation in axons with only occasional presence of periodic patterns in dendrites; concordantly, quantitative autocorrelation analyses showed substantially lower average autocorrelation amplitudes in dendrites (Fig. 2 J and K). Together, our data

suggest that the MPS structure is prevalent in morphologically and functionally diverse excitatory and inhibitory neuronal types.

As many neurons are myelinated in the brain, we next examined whether myelination had any effect on the formation of MPS using a previously developed neuron–glia coculture system (31). Myelin sheaths were identified using an antibody against myelin basic protein (MBP). As shown by STORM images and autocorrelation analysis, β II spectrin was periodically distributed in both MBP-positive and MBP-negative axonal segments with a nearly identical degree of periodicity (Fig. S6), suggesting that myelination has little impact on the formation of the periodic membrane skeleton in vitro, consistent with the observation of the periodic structure in myelinated sciatic nerves (23).

The MPS Structure Is Present in Mouse Peripheral Sensory and Motor

Axons. The aforementioned results were all from neurons restricted to the mouse central nervous system. We next examined whether MPS was also present in peripheral neurons by studying two representative types: motor neurons and dorsal root ganglia (DRG) sensory neurons. Motor neurons were derived from mouse embryonic stem (mES) cells, and differentiated motor neurons were identified by Hb9 staining (32). BII spectrin was again detected by immunostaining and STORM imaging. The distribution of *βII* spectrin was highly periodic in the axons, but mostly nonperiodic in the dendrites of mES-derived motor neurons (Fig. 3 A and B). Different from motor neurons, DRG sensory neurons are pseudounipolar neurons that have axons but lack dendrites (24). We imaged immunostained BII spectrin in explants of DRG sensory neurons (Fig. 3C) and in dissociated DRG sensory neurons (Fig. 3D); the periodic distribution of β II spectrin was clearly evident in the axons of these neurons as well.

The Distribution of β II Spectrin in Presynaptic Boutons and Dendritic Spine Necks. We next extended our investigation to two neuronal compartments in axons and dendrites related to synapses: axonal boutons (presynaptic) and dendritic spines (postsynaptic). Spectrin is expressed in synapses and plays important roles in synapse stabilization and spine morphogenesis (8, 33, 34). To image β II




Fig. 2. MPS is present in multiple subtypes of inhibitory and excitatory neurons in the mouse central nervous system. (*A–I*) Representative STORM images of β II spectrin for a typical axonal region (*Left* panels in *A–I*) and averaged autocorrelation functions from multiple axons (*Right* panels *A–I*) for cultured cerebellar granule cells (*A*; *n* = 7 for averaged autocorrelation calculation), Pontine nuclei neurons (*B*; *n* = 7), dopaminergic neurons (*C*; *n* = 9), olfactory neurons (*D*; *n* = 16), Parvalbumin neurons from cortex (*E*; *n* = 9), Parvalbumin neurons from midbrain (*F*; *n* = 6), Parvalbumin neurons from hippocampus (*G*; *n* = 13), Golgi cells (*H*; *n* = 8), and Purkinje cells (*I*; *n* = 7) from cerebellum. (*J* and *K*) Averaged autocorrelation amplitudes of axons (blue) versus dendrites (red) for the neuronal subtypes we imaged in the mouse central nervous system. Error bars are SEM. (Scale bar: 2 µm.)

spectrin in presynaptic boutons, we transiently expressed a GFPtagged ßII spectrin (or GFP-tagged all spectrin) construct (16) in a sparse subpopulation of cultured hippocampal neurons, immununolabeled GFP-tagged spectrin with anti-GFP, and performed STORM imaging at day in vitro (DIV) 14. The presynaptic boutons were marked by immunofluorescence of Bassoon and imaged at the diffraction-limited resolution. Such transient expression was performed to allow for the spectrin labeling in only a sparse subset of neurons, thereby preventing the presynaptic spectrin signal from being obscured by postsynaptically distributed spectrin. Interestingly, despite the presence of MPS in nearly all axons, the periodic structure was disrupted at most (\sim 70%) of the presynaptic sites imaged (Fig. 4 A and B). Next, we imaged βII spectrin distribution in the dendritic spines of cultured hippocampal neurons at DIV 18, using immunofluorescence against endogenous BII spectrin, and observed the periodic structure of BII spectrin in a significant fraction (at least 25%) of the spine necks (Fig. 4C). Notably, the periodic pattern was observed even in some of the spine necks stemming from the shaft regions that exhibited an irregular β II spectrin distribution (Fig. 4*C*-2).

The Distribution of β II Spectrin in Glial Processes. We next examined whether the MPS structure also forms in glial-cell processes. Glial cells express β II spectrin, but at a lower level than neurons (35, 36). We cultured glial cells from mouse brain or rat sciatic nerve tissues, used specific immunomarkers to identify four different glial types, including astrocytes, microglia, NG2 glia, and Schwann cells, and imaged their endogenous β II spectrin distribution using immunolabeling and STORM. The β II spectrin distribution appeared to be largely irregular in most of the processes of all four glial types imaged, although patches of periodic β II spectrin patterns could be observed in a small fraction of these processes (Fig. 5A and Fig. S7), similar to the appearance of β II spectrin in dendrites. Autocorrelation analysis showed that the degree of periodicity of the β II spectrin distributions in the glial-




Fig. 3. MPS is present in mouse peripheral sensory and motor axons. (*A*) Representative STORM images of β II spectrin for a typical axonal region (*Left*) and averaged autocorrelation functions from multiple axons (*Right*, n = 12) for mES-derived motor neurons. (*B*) Same as in *A* but for dendrites in mES-derived motor neurons (n = 11). (*C* and *D*) Same as *A* but for axons in DRG neurons (*C*, explants, n = 13; *D*, dissociated DRG culture, n = 10). (Scale bar: 2 µm.)




Fig. 4. The distribution of BII spectrin in axonal boutons and dendritic spines. (A) Representative STORM images of *βII* spectrin in axons that contain presynaptic boutons marked by a presynaptic marker, Bassoon (magenta). Bll spectrin is labeled by expressing GFP-BII spectrin (or GFP-all spectrin) in a sparse subset of cultured neurons, and hence some bassoon-positive regions do not overlap with the GFP-positive axon. (Top and Middle) Examples of bassoonpositive presynaptic boutons (indicated by arrow) within which the periodic pattern of spectrin is disrupted. (Bottom) An example of bassoon-positive presynaptic bouton (indicated by arrow) within which the period pattern of spectrin is not disrupted. (B) Bar graph depicting relative fraction of boutons in which the periodic structure is disrupted or undisrupted. (C) The dendritic region of a cultured mouse hippocampal neuron immunostained with MAP2, a dendrite marker. (C-1 and C-2) STORM images of two selected regions showing periodic patterns of immunolabeled endogenous BII spectrin in spine necks (Left panels) with autocorrelation functions (Right panels) calculated from the boxed spine neck region (white boxes) and dendritic shaft (yellow boxes). (Scale bar: 10 µm for conventional images and 1 µm for STORM images.)

cell processes was similar among all four glial types and also similar to that observed in neuronal dendrites, but much smaller than that observed in neuronal axons (Fig. 5*B*).

The MPS Structure Was Observed in Neurons from Multiple Invertebrate and Vertebrate Animal Species. Spectrin tetramers are evolutionarily conserved with the lengths and major structural domains of the αII and BII spectrin homologs preserved across the metazoan kingdom (Fig. S8) (11, 22). Next, we examined whether MPS can also form in other animal species in addition to rodents. We first studied two additional vertebrate species, Gallus gallus (chicken), a nonmammalian vertebrate, and H. sapiens. We examined the endogenous BII spectrin distribution using immunofluorescence labeling and STORM imaging in cultured chicken neurons and human induced pluripotent stem (iPS) cell-derived motor neurons. BII spectrin exhibited a long-range, periodic pattern along the axonal shafts with a spacing of ~190 nm in both chicken neurons and human iPS-derived motor neurons (Fig. 6A and B). Next we extended our studies to two invertebrate systems: C. elegans and Drosophila. To this end, we created fluorescent fusion proteins of β spectrin for *C. elegans* (UNC-70) and Drosophila (d\beta-spectrin) and overexpressed these fusion constructs with neuronal-specific promotors in C. elegans and in glutamatergic neurons of Drosophila larva (SI Materials and Methods). We used structure illumination microscopy (SIM) (37) to image these neurons in vivo either in the whole animal (C. elegans) or in brain tissues (Drosophila), rather than in cultured neurons. In both C. elegans and Drosophila neurons, the expressed β-spectrin proteins adopted a periodic structure with a spacing of ~190–200 nm (Fig. 6 C and D), consistent with the conserved lengths




Fig. 5. Sparse presence of the periodic β II spectrin structure in glial cell processes. (*A*) An astrocyte stained for the astrocyte marker GFAP (green) and β II spectrin (red). (*A*-1 and *A*-2) STORM images of β II spectrin from two processes of that astrocyte (indicated by arrows) displaying relatively periodic (*A*-1) or irregular (*A*-2) spectrin distribution. (*Right* panels) Autocorrelation analysis of the boxed regions. (*B*) Averaged autocorrelation functions (blue) calculated from multiple randomly selected processes of astrocytes, microglia cells, NG2 glia, and Schwann cells (*n* = 10 for astrocytes, *n* = 9 for microglia, *n* = 8 for NG2 glia, *n* = 9 for Schwann cells). For comparison, the dotted gray and red curves show the averaged autocorrelation functions of β II spectrin distribution in axons and dendrites of cortical vGlut1+ neurons, respectively (reproduced from Fig. 1C). (Scale bar: 20 µm for conventional image and 1 µm for STORM images.)

and domain structures of spectrin tetramers across different species. These results indicate that the MPS structure can form not only in vertebrate neurons but also in invertebrate neurons. However, we note that, although overexpression of β II spectrin does not appear to perturb the formation of the MPS structure in rodent axons, such overexpression does promote MPS formation in dendrites of rodent neurons (16). Therefore, future experiments on imaging endogenous spectrin are needed to confirm whether the MPS structure is formed by endogenous proteins in *C. elegans* and *Drosophila*. Because spectrin is expressed in many neuronal and nonneuronal cells, the resulting high background makes such imaging experiments difficult using immunofluorescence. Using genome editing to introduce fluorescent tags to the endogenous spectrin in specific populations of neurons provides a promising solution.

Discussion

We have recently discovered that actin, spectrin, and associated molecules form a membrane-associated periodic skeletal structure in the axons of rodent neurons (12). Here, we show that this actin– spectrin-based MPS structure is prevalent in a wide spectrum of neuronal cell types cultured from the rodent central and peripheral nervous systems, including excitatory and inhibitory neurons in cortex, hippocampus, and midbrain, and many specific neuronal subtypes, such as granule cells, pontine nuclei neurons, dopaminergic neurons, olfactory neurons, Parvalbumin cells, Golgi cells, Purkinje cells, motor neurons, and sensory neurons. While we were finalizing this paper, an independent study recently reported the presence of periodic actin structure in several neuronal types in the rodent central and peripheral nervous systems (23), and the list of cell types reported in ref. 23 partially overlap with the neuronal types investigated here. Using specific neuronal markers to distinguish neuronal cell types, we further showed that, in all these neuronal types, the MPS structure preferentially forms in axons compared with dendrites. Interestingly, despite the prevalent presence of MPS in nearly all axons, we showed that this periodic structure is disrupted in most presynaptic boutons. This observation suggests the possibility that MPS may inhibit synapse formation, and thus needs to be disrupted for bouton formation, and that the small fraction of boutons that maintained the periodic structure might be relatively immature synapses. On the other hand, the periodic structure is formed in a significant fraction of dendritic spine necks, including ones that stemmed from the dendritic shafts lacking the MPS structure, suggesting a specific mechanism to promote the formation of this periodic structure in the spine neck. We also examined four types of rodent glial cells, including astrocytes, microglia, NG2 glia, and Schwann cells and observed isolated patches of periodic spectrin distributions in a small fraction of glial-cell processes, consistent with the observation of periodic actin and spectrin distributions in rodent oligodendrocytes (23). It is interesting to note that this periodic structure appeared to be absent in the majority of glial-cell processes in our images, and our quantitative analysis showed that the degree of periodicity for spectrin distributions in glial processes was similar to that observed in neuronal dendrites and much smaller than that observed in neuronal axons. Because the expression level of BII




Fig. 6. The observation of MPS in neurons across different animal species. (A) Representative STORM image of immunolabeled β II spectrin in the axon of a cultured chicken neuron (*Left*) and the autocorrelation function calculated from the boxed region in this image (*Right*). (Scale bar: 2 µm.) (*B*) Representative STORM image of immunolabeled β II spectrin in the axon of a human iPS-derived motor neuron (*Left*) and the autocorrelation function calculated from this image (*Right*). (Scale bar: 2 µm.) (*B*) representative SIORM image of β spectrin (UNC-70)-GFP in *C. elegans* neurons imaged directly in the animals (*Left*), and the corresponding autocorrelation functions of the boxed regions (*Right*). (Scale bar: 1 µm.) (*D*) Two representative SIM images of β spectrin-mMaple3 in *Drosophila* neurons imaged in *Drosophila* brain tissues (*Left*) and the corresponding autocorrelation functions of the boxed regions (*Right*). (Scale bar: 1 µm.)

spectrin in glial cells is similar to that in dendrites and lower than that in axons (35, 36), this observation is consistent with the notion that local concentrations of BII spectrin are a key determinant for MPS formation (16). Taken together, these results suggest that the MPS structure might have some tendency to self-assemble, but the axon is a specialized cellular compartment where the MPS structure forms with a much higher propensity. Whether this substantially higher propensity of MPS formation in neuronal axons observed in all examined cultured nervous system cell types can be extended to in vivo systems awaits future investigations. Finally, we found that spectrin is capable of adopting a periodic distribution in neurons from a wide variety of animal species across both invertebrates and vertebrates, ranging from C. elegans to H. sapiens, suggesting that the MPS structure is an evolutionarily conserved structure in neurons across different species. The high prevalence of this structure in the axons of diverse neuronal types and animal species suggests that the actinspectrin-based MPS is a key functional component of the axon. Indeed, such a periodic skeletal structure formed by actin rings connected by flexible spectrin tetramers provides a robust and flexible support for the axonal membrane, which could be particularly important for maintaining the integrity of axons under mechanical stress and can explain why spectrin depletion results in axon buckling and breakage during animal locomotion (4, 6). The ability of this periodic submembrane skeleton to organize membrane proteins, such as ion channels (12), into a periodic distribution might also impact the generation and/or propagation of action potentials as well as other signaling pathways at the axonal membrane.

- 1. Luo L (2002) Actin cytoskeleton regulation in neuronal morphogenesis and structural plasticity. Annu Rev Cell Dev Biol 18:601–635.
- 2. Cingolani LA, Goda Y (2008) Actin in action: The interplay between the actin cytoskeleton and synaptic efficacy. Nat Rev Neurosci 9(5):344–356.
- Stiess M, Bradke F (2011) Neuronal polarization: The cytoskeleton leads the way. Dev Neurobiol 71(6):430–444.
- Hammarlund M, Jorgensen EM, Bastiani MJ (2007) Axons break in animals lacking beta-spectrin. J Cell Biol 176(3):269–275.
- Galiano MR, et al. (2012) A distal axonal cytoskeleton forms an intra-axonal boundary that controls axon initial segment assembly. *Cell* 149(5):1125–1139.
- 6. Krieg M, Dunn AR, Goodman MB (2014) Mechanical control of the sense of touch by β -spectrin. Nat Cell Biol 16(3):224–233.
- Hülsmeier J, et al. (2007) Distinct functions of alpha-Spectrin and beta-Spectrin during axonal pathfinding. *Development* 134(4):713–722.
- Pielage J, Fetter RD, Davis GW (2005) Presynaptic spectrin is essential for synapse stabilization. Curr Biol 15(10):918–928.
- Tang Y, et al. (2003) Disruption of transforming growth factor-beta signaling in ELF beta-spectrin-deficient mice. Science 299(5606):574–577.
- Stankewich MC, et al. (2011) Cell organization, growth, and neural and cardiac development require αll-spectrin. J Cell Sci 124(Pt 23):3956–3966.
- Bennett V, Lorenzo DN (2013) Spectrin- and ankyrin-based membrane domains and the evolution of vertebrates. Curr Top Membr 72:1–37.
- 12. Xu K, Zhong G, Zhuang X (2013) Actin, spectrin, and associated proteins form a periodic cytoskeletal structure in axons. *Science* 339(6118):452–456.
- Rust MJ, Bates M, Zhuang X (2006) Sub-diffraction-limit imaging by stochastic optical reconstruction microscopy (STORM). Nat Methods 3(10):793–795.
- Huang B, Wang W, Bates M, Zhuang X (2008) Three-dimensional super-resolution imaging by stochastic optical reconstruction microscopy. Science 319(5864):810–813.
- Lukinavičius G, et al. (2014) Fluorogenic probes for live-cell imaging of the cytoskeleton. Nat Methods 11(7):731–733.
- Zhong G, et al. (2014) Developmental mechanism of the periodic membrane skeleton in axons. *eLife* 3:e04581.
- Jenkins SM, Bennett V (2001) Ankyrin-G coordinates assembly of the spectrin-based membrane skeleton, voltage-gated sodium channels, and L1 CAMs at Purkinje neuron initial segments. J Cell Biol 155(5):739–746.
- Yang Y, Ogawa Y, Hedstrom KL, Rasband MN (2007) betalV spectrin is recruited to axon initial segments and nodes of Ranvier by ankyrinG. J Cell Biol 176(4):509–519.
 Leterrier C, et al. (2015) Nanoscale architecture of the axon initial segment reveals an
- organized and robust scaffold. *Cell Reports* 13(12):2781–2793. 20. D'Este E, Kamin D, Göttfert F, El-Hady A, Hell SW (2015) STED nanoscopy reveals the
- ubiquity of subcortical cytoskeleton periodicity in living neurons. *Cell Reports* 10(8): 1246–1251. 21. Lorenzo DN, et al. (2014) A PIK3C3-ankyrin-B-dynactin pathway promotes axonal
- growth and multiographic transport. J Cell Biol 207(6):735–752.
- Bennett V, Baines AJ (2001) Spectrin and ankyrin-based pathways: Metazoan inventions for integrating cells into tissues. *Physiol Rev* 81(3):1353–1392.

Materials and Methods

All experimental procedures were performed in accordance with the *Guide* for the Care and Use of Laboratory Animals of the National Institutes of Health (38). The protocol was approved by the Institutional Animal Care and Use Committee of Harvard University.

Procedures for cell culture, antibodies, and immunofluorescence are described in *SI Materials and Methods*. STORM imaging was performed as previously described (12, 14). Briefly, labeled cells were imaged with continuous illumination of a 657-nm laser for imaging and 405 nm for activation of the photoswitchable dyes in an oblique-incidence geometry on an Olympus IX-71 inverted microscope. STORM images were generated from a sequence of about 20,000~40,000 image frames at a frame rate of ~60 Hz. SIM imaging (37) was performed on a Zeiss Elyra superresolution microscope at the Center for Biological Imaging at Harvard University. SIM imaging analyses were performed with the provided software for the Zeiss Elyra system. Autocorrelation analyses were performed as the difference between the first peak (at ~190 nm) and the average of the two first valleys of the autocorrelation curve. See *SI Materials and Methods* for more details.

ACKNOWLEDGMENTS. The preprint of this paper has been posted on bioRxiv (dx.doi.org/10.1101/045856) for discussion and we thank a wide community of colleagues for participating in the discussion. This work is supported in part by the Howard Hughes Medical Institute (X.Z., K.S., and M.R.F.); National Institutes of Health Grants NS089786 (to M.T.-L.), NS048392 (to K.S.), and NS053538 (to M.R.F.); and the ALS Association (T.M.). R.Z. is an HHMI Fellow of the Life Sciences Research Foundation. Z.W. is supported in part by a Bristol-Myers Squibb Postdoctoral Fellowship at the Rockefeller University. X.Z., K.S., and M.R.F. are Howard Hughes Medical Institute investigators.

- D'Este E, et al. (2016) Subcortical cytoskeleton periodicity throughout the nervous system. Sci Rep 6:22741.
- Kandel ER (2013) *Principles of Neural Science* (McGraw-Hill, New York), 5th Ed, 1709 pp.
 Daubner SC, Le T, Wang S (2011) Tyrosine hydroxylase and regulation of dopamine synthesis. *Arch Biochem Biophys* 508(1):1–12.
- Lee AC, He J, Ma M (2011) Olfactory marker protein is critical for functional maturation of olfactory sensory neurons and development of mother preference. *J Neurosci* 31(8):2974–2982.
- Neki A, et al. (1996) Metabotropic glutamate receptors mGluR2 and mGluR5 are expressed in two non-overlapping populations of Golgi cells in the rat cerebellum. *Neuroscience* 75(3):815–826.
- Whitney ER, Kemper TL, Rosene DL, Bauman ML, Blatt GJ (2008) Calbindin-D28k is a more reliable marker of human Purkinje cells than standard Nissl stains: A stereological experiment. J Neurosci Methods 168(1):42–47.
- Furuya S, Makino A, Hirabayashi Y (1998) An improved method for culturing cerebellar Purkinje cells with differentiated dendrites under a mixed monolayer setting. *Brain Res Brain Res Protoc* 3(2):192–198.
- Tabata T, et al. (2000) A reliable method for culture of dissociated mouse cerebellar cells enriched for Purkinje neurons. J Neurosci Methods 104(1):45–53.
- Gardner A, Jukkola P, Gu C (2012) Myelination of rodent hippocampal neurons in culture. Nat Protoc 7(10):1774–1782.
- Arber S, et al. (1999) Requirement for the homeobox gene Hb9 in the consolidation of motor neuron identity. *Neuron* 23(4):659–674.
- Pielage J, Fetter RD, Davis GW (2006) A postsynaptic spectrin scaffold defines active zone size, spacing, and efficacy at the Drosophila neuromuscular junction. J Cell Biol 175(3):491–503.
- 34. Gao Y, et al. (2011) β-III spectrin is critical for development of purkinje cell dendritic tree and spine morphogenesis. J Neurosci 31(46):16581–16590.
- Goodman SR, Lopresti LL, Riederer BM, Sikorski A, Zagon IS (1989) Brain spectrin(240/ 235A): A novel astrocyte specific spectrin isoform. *Brain Res Bull* 23(4-5):311–316.
- Susuki K, et al. (2011) Schwann cell spectrins modulate peripheral nerve myelination. Proc Natl Acad Sci USA 108(19):8009–8014.
- Gustafsson MG, et al. (2008) Three-dimensional resolution doubling in wide-field fluorescence microscopy by structured illumination. *Biophys J* 94(12):4957–4970.
- Committee for the Update of the Guide for the Care and Use of Laboratory Animals (2011) Guide for the Care and Use of Laboratory Animals (National Academies Press, Washington, DC), 8th Ed, p XXV.
- 39. He J, et al. (2013) Dual function of CD81 in influenza virus uncoating and budding. *PLoS Pathog* 9(10):e1003701.
- Jones SA, Shim SH, He J, Zhuang X (2011) Fast, three-dimensional super-resolution imaging of live cells. Nat Methods 8(6):499–508.
- Wang S, Moffitt JR, Dempsey GT, Xie XS, Zhuang X (2014) Characterization and development of photoactivatable fluorescent proteins for single-molecule-based superresolution imaging. *Proc Natl Acad Sci USA* 111(23):8452–8457.