

SUPPLEMENTARY INFORMATION

AFM cantilever with in-situ renewable mercury microelectrode

Peter Schön^{1*}, Joël Geerlings², Niels Tas², Edin Sarajlic^{2,3}

1) *Materials Science and Technology of Polymers, MESA+ Institute for Nanotechnology, University of Twente, Enschede, The Netherlands*

2) *Transducers Science and Technology, MESA+ Institute for Nanotechnology, University of Twente, Enschede, The Netherlands*

3) *SmartTip B.V., Enschede, The Netherlands*

*Corresponding author:

Peter M. Schön
Materials Science and Technology of Polymers
MESA+ Institute for Nanotechnology
University of Twente
Enschede,
The Netherlands

+31-53-4893170 (office)

+31 (0)53 489 3823 (fax)

E-mail: p.m.schon@utwente.nl

Experimental Details

Photographs of the experimental setup



Figure S1. Inverted optical microscope covered with plastic foil.



Figure S2. Experimental setup, on the right the syringe next to the pressure sensor and translation motor on the left.




Figure S3. Fountain pen probe connected to a polycarbonate block with fluidic and electric connections. The probe is placed above the ITO electrode with O-ring containing the electrolyte solution.



Figure S4. Syringe filled with mercury and positioning table to move the plunger of the syringe.