## TERMIS EU 2019

termis.

Tissue Engineering Therapies: From Concept to Clinical Translation & Commercialisation

27-31 May 2019 Rhodes, Greece

**Rodos Palace Hotel** 

Conference Chair: **Dimitrios I. Zeugolis, PhD** 

Conference Program Chair: Maria Chatzinikolaidou, PhD





Find us at:





🚺 @termis eu 2019

Organized by:











Organizing Secretariat:



T: +30 210 6833600 E: congress@convin.gr W: www.convin.gr Antibiotic loaded poly(\varepsilon-caprolactone) microspheres functionalized with poly(aspartic acid) as bone targeting delivery system to treat infection

S. Rotman<sup>1,2</sup>, M. Wang<sup>3</sup>, S. Yang<sup>3</sup>, D.W. Grijpma<sup>2</sup>, R.G. Richards<sup>1</sup>, F.T. Moriarty<sup>1</sup>, D. Eglin<sup>1</sup>, T.T. Tang<sup>3</sup>, O. Guillaume<sup>1</sup>

Presenting Author: Stijn Rotman, stijn.rotman@aofoundation.org

<sup>1</sup>AO Research Institute Davos, Davos, Switzerland <sup>2</sup>Department of Biomaterials Science and Technology, Faculty of Science and Technology, Technical Medical Centre, University of Twente, Enschede, The Netherlands <sup>3</sup>Shanghai Key Laboratory of Orthopedic Implants, Department of Orthopedic Surgery, Shanghai Ninth People's Hospital, Shanghai Jiao Tong University School of Medicine, Shanghai, China

**INTRODUCTION:** The recurrence rate of chronic osteomyelitis in adults is close to 30% [1]. Bacteria are known to migrate deeper into bone tissue through canaliculi and evade common systemic- and local-antibiotic therapies [2]. By fabricating antibiotic loaded poly( $\epsilon$ -caprolactone) (PCL) microspheres conjugated with the bone-chelator poly(aspartic acid) (PAA) we aim to prolong the microsphere residency near the site of infection, increasing bactericidal potential.

**METHODS:** Hydrophobic Gentamicin-dioctyl sulfosuccinate (Gen-AOT) loaded PCL microspheres were made by oil/water emulsion methods. In vitro antimicrobial properties were tested by zone of inhibition (ZOI) in a serial plate transfer test with S. aureus. In vivo antimicrobial efficacy of 1 mg of microspheres was tested in a femoral defect in rats (n = 5), infected with 2·10<sup>6</sup> colony forming units (CFU) of bioluminescent Xen-29 a week prior to treatment. In a 2<sup>nd</sup> study, the PCL microspheres underwent conjugation with PAA by carbodiimide chemistry. Interaction with bone mineral was assessed in the same model as above. IR780 iodide loaded PCL or PCL-PAA microspheres (1 mg) were injected in the bone defect and traced using an in vivo imaging system (IVIS Lumina III, Perkin Elmer).

**RESULTS & DISCUSSION:** ZOI of Gen-AOT loaded PCL was measurable for 5 days, while a ZOI for bactericidal collagen-sheets was visible for 3 days. The Gen-AOT loaded PCL microspheres caused an 81% reduction in CFU compared to untreated control. In vivo, a brighter signal was measured for PCL-PAA compared to PCL microspheres, validating the hypothesis that PAA-grafted PCL resides longer in bone as control PCL.

**CONCLUSIONS:** In the presented animal model, a monotreatment of 1 mg PCL microspheres caused an 81% CFU reduction. PCL-PAA microspheres enhance bone affinity by chelation of the PAA to bone mineral at the femoral defect. Further work is required to optimize the bone-targeted drug delivery system to bone.

**ACKNOWLEDGEMENTS:** Project funding was received from AO Trauma.

## **REFERENCES**

- [1] Tice AD et al. Am J Med. 2003; 114(9):723-728
- [2] Brady A et al. FEMS Immunol. Med. Microbiol. 2008; 52(1):13-22