



Validation of the ESC 0/1 algorithm in the EMS setting

K Wibring^{1,2}; J Herlitz³; M Lingman^{4,5} and A Bång¹

Institute of Health and Care Sciences, Sahlgrenska Academy, University of Gothenburg, Göteborg, Sweden.
Department of Prehospital Emergency Care, Region Halland, Sweden
The Prehospital Research Center Western Sweden, University of Borås, Borås, Sweden.
Halland Hospital, Region Halland, Sweden
Department of Molecular and Clinical Medicine/Cardiology, Institute of Medicine, Sahlgrenska Academy, University of Gothenburg, Göteborg, Sweden
Contact: kristoffer.wibring@regionhalland.se

Background

Chest pain is one of the most common reasons for contacting the emergency medical services (EMS). Troponins is a corner stone in ruling out high-risk conditions at the emergency department (ED) and the interest in using Troponins in the EMS setting is increasing. European Society of Cardiology (ESC) advocate the use of the 0/1 algorithm which allows rule-out after one Troponin test. However, no study has been validating the ESC 0/1

Methods

Observational cohort study including all patients with chest pain cared for by the EMS in the county of Halland, Sweden, during 2018. A blood sample was collected during the EMS mission and brought to the ED, where a high-sensitive Troponin T (hs-TnT) was analysed. The diagnosis at hospital discharge were retrieved from the hospital medical record.

algorithm in the EMS setting.

Aim


To examine the accuracy of the ESC 0/1 algorithm if being applied in the EMS setting using prehospital Troponin T.

Results

For AMI identification the sensitivity was 99.4 (CI 95% 96.8-100.0), specificity 18.0 (CI 95% 19.9-20.1), PPV 13.4 (CI 95% 11.6-15.4) and NPV 99.6 (CI 95% 97.7-). Giving a rule-out rate of 16%, with <0.1% patients with AMI wrongly ruled-out.

For high-risk conditions (AMI, pulmonary embolism, aortic dissection etc.) identification the sensitivity was 98.4 (Cl 95% 96.1-99.6), specificity 19.0 (Cl 95% 16.1-19.1), PPV 20.0 (Cl 95% 19.5-20.5) and NPV 98.3 (Cl 95% 95.7-99.4), with four false negatives. Giving a rule-out rate of 16%, with <0.2% patients with

In total 2 917 EMS missions were included, of which 1 501 had an EMS hs-TnT test result. The sensitivity, specificity, negative and positive predictive value (NPV, PPV) of the ESC 0/1 algorithms ability to identify acute myocardial infarction (AMI) and other high-risk diagnoses were calculated. Using a single hs-TnT test. Calculation were made twice. Once for patients with symptoms onset >3 hours ago and once regardless of time for symptom onset.



high-risk conditions wrongly ruled out.

Hs-cTnT 5th gen. (Roche)	5	14	4	52	10
--------------------------------	---	----	---	----	----

*symptom onset >3 hours ago

	Total number	True positive	False positive	True negative	False negative	Sens	Spec	PPV	NPV			
Debut >3 h och Tnt <5 ng/L, AMI	538	51	414	72	1	98	15	11	97			
Tnt <5 ng/L, AMI	1 501	169	1 0 9 2	239	1	99	18	13	100			
Debut >3 h och Tnt <5 ng/L, high-risk	538	76	389	70	3	96	15	16	96			
Tnt <5 ng/L, high-risk	1 501	252	1 009	236	4	98	19	20	98			
<u>ESC rule-in</u>												
	Total	True	False	True	False							
	number	positive	positive	negative	negative	Sens	Spec	PPV	NPV			
Tnt ≥52 ng/L, high-risk	1 501	89	78	1 167	167	35	94	53	87			
Tnt ≥52 ng/L, AMI	1 501	81	86	1 2 4 5	89	48	94	49	93			

Conclusions and clinical implications

By introducing new technology, enabling beside hs-Troponins analyses, and applying the ESC 0/1 algorithm in the EMS setting, high-risk conditions and especially AMI can be ruled out with high accuracy using only a single hs-Troponin test. Allowing referral of low-risk patients to self-care or alternative modes of transportation or destination. High-risk patients beneficial of increased triage level or transportation to hospitals providing PCI capabilities or cardiac care units could also be identified, allowing to by-pass the nearest hospital or the emergency department.

