



Tata, B. K., Harbulot, C., Csaba, Z., Peineau, S., Jacquier, S., & De Roux, N. (2017). Rabconnectin-3 is required for the morphological maturation of GnRH neurons and kisspeptin responsiveness. *Scientific Reports*, 7, [42463]. DOI: [10.1038/srep42463](https://doi.org/10.1038/srep42463)

Publisher's PDF, also known as Version of record

License (if available):

CC BY

Link to published version (if available):

[10.1038/srep42463](https://doi.org/10.1038/srep42463)

[Link to publication record in Explore Bristol Research](#)

PDF-document

This is the final published version of the article (version of record). It first appeared online via Scientific Reports at <https://doi.org/10.1038/srep42463>. Please refer to any applicable terms of use of the publisher.

University of Bristol - Explore Bristol Research

General rights

This document is made available in accordance with publisher policies. Please cite only the published version using the reference above. Full terms of use are available:
<http://www.bristol.ac.uk/pure/about/ebr-terms.html>

Rabconnectin-3 α is required for the morphological maturation of GnRH neurons and kisspeptin responsiveness

Brooke K. Tata¹, Carole Harbulot¹, Zsolt Csaba¹, Stéphane Peineau^{1,2}, Sandrine Jacquier¹, Nicolas de Roux^{1,3,*}.

Affiliations:

¹ Univ Paris Diderot, Sorbonne Paris Cité, U1141, Inserm, F- 75019, Paris, France.

² MRC Centre for Synaptic Plasticity; School of Physiology, Pharmacology, Neurosciences; University of Bristol; Bristol, UK

³ AP-HP, Laboratoire de Biochimie. Hôpital Robert Debré, Paris, F-75019, France.

Table S1. List of primary antibodies

Primary antibody	Host species	Source	Catalogue number	Concentration
Polyclonal anti- DMXL2	Rabbit	Sigma	HPA039375	1/500
Polyclonal anti-kisspeptin (mammalian)	Rabbit	Dr. Alain Caraty (INRA, Tours, France)	AC 564	1/10000
Polyclonal anti-GnRH	Guinea Pig	Dr. Greg Anderson (CNE, University of Otago)	GA04	1/15000
Polyclonal anti-cFos	Rabbit	Santa Cruz Biotechnology	SC-52	1/5000