



The motivation of our work is the spatial and temporal distribution analysis of methane around pockmarks and other CH₄ seeps.


Worldwide, the release of methane from sediments of lakes, coastal regions as well as ocean margins is observed. The gas release is often associated with specific features like pockmarks (morphological depressions at the seafloor), mud volcanoes, cold seeps and occurrence of gas hydrates. This gas plumes were observed by underwater camera systems and acoustic techniques.




Visual observation of the release of gas bubbles from the seafloor.



Acoustic blanking in surface sediments



Often the CH₄ concentrations around pockmarks are rather low.




Bottom Water Sampler Sauter et al. (2005)

The application of UWMS is a step towards a more detailed investigation of spatial and temporal variations of methane in aquatic systems.

Why under water mass spectrometry (UWMS)?

Compared to such semi-quantitative information, rather little is known about the concentration field of CH₄ as well as other gases around e.g., pockmarks. This is mainly due to the laborious sampling schemes and rather time consuming CH₄ analysis by gas chromatography.



Established method



Advantage of UWMS compared to established methods:

- NO sampling artifacts by e.g. de-pressurisation or sample warming due to in situ sampling
- Highest possible spatial and temporal resolution (up to 750 times higher than established techniques),
- Online and realtime measurements
- simultaneous measurements of major and trace gases

Underwater mass spectrometer:




Comparison UWMS versus gas chromatography

Mass spectrometer:


- Membrane to reduce water vapour
- Cryotrap to reduce water vapour
- Quadrupole based separation
- Varian turbo molecular with diaphragm backing pumps
- Embedded PC
- Online onboard visualization by ethernet connection

- + Detection limit of around 16nmol/l CH₄ by an implemented cryotrap
- + Results are comparable with established techniques like gas chromatography

Mode of operation



Ex situ



In situ in a frame including benthic chamber




In situ at sediment-water-transition-zone




Laboratory measurements


Outlook




Observation of a gas seep area (North Sea) in high resolution




Working area (Modified after Schroot et al. 2005)



Under water gas analyser, sampler and observing system



Transect of the observing system




Echosounding during transect




Video observation

Methane distribution visualized by ArcGIS



(Gentz et al. unpublished data)

Methane inventory calculated by Matlab



Results:

- A methane saturation of 23200 % was observed in 8 m water depth.
- The air sea exchange flux is calculated to $\sim 210 \pm 63 \mu\text{mol m}^{-2} \text{d}^{-1}$.
- Methane flux: 28.27 L min^{-1}
- Methane release: $35.3 \pm 17.65 \text{ t CH}_4 \text{ yr}^{-1}$

Entire interpolated inventory of methane (6.410.000 m³):

- $\sim 0.6 \text{ mol CH}_4$
- $\sim 1.000.000 \text{ m}^3$ (15.6 %) contain concentrations higher than 200 nmol L^{-1}
- 40 % of initial methane is dissolved above the pycnocline.

References:
Gentz T. & Schlüter M. (2012), Limnology and Oceanography: Methods 10 (2012): 317-328.
Mau, S. et al. (2015) Biogeosciences 12.18: 5261-5276.
Schlüter M. & Gentz T. (2008), American Society Mass Spectrometry 19: 1395-1402,
Short R.T. et al. (2006), Measurement Science and Technology 10: 1195-1201,

Gentz T. et al. (2014), Continental Shelf Research 72: 107-118.
Bell, R. J. et al. (2007), Environmental Science and Technology, 41: 8123-8128.
Sauter E.J., et al. (2005) Journal of Sea Research, Vol. 54: 204-210