Electronic Supplementary Information

Modulation of pK_a by cyclodextrins; subtle structural changes induce spectacularly different behaviors

Paweł M. Nowak^a, Michał Woźniakiewicz^a*, Mariusz P. Mitoraj^b, Magdalena Garnysz^a, Paweł Kościelniak^a

- a. Department of Analytical Chemistry, Faculty of Chemistry, Jagiellonian University, R. Ingardena 3, 30-060 Krakow, Poland.
- Department of Theoretical Chemistry, Faculty of Chemistry, Jagiellonian University, R. Ingardena 3, 30-060 Krakow, Poland

pH Phosphate buffer I	Buffer composition [mL]	
	H ₃ PO ₄ (100 mM)	NaH ₂ PO ₄ (100 mM)
2.50	18.52	4.80
3.50	1.83	4.98
Acetic buffer	CH ₃ COOH (500 mM)	CH ₃ COONa (500 mM)
4.50	14.16	10.00
5.50	1.42	10.00
Phosphate buffer II	NaH ₂ PO ₄ (100 mM)	Na ₂ HPO ₄ (100 mM)
6.00	3.59	4.70
7.00	1.01	13.29
8.00	0.12	16.25
Borate buffer	Na ₂ B ₄ O ₇ ·10H ₂ O (50 mM)	NaOH (1 M)
9.20	48.42	0.16

Table S-1. Composition and predicted pH of all BGEs (background electrolytes) prepared for experiments, calculated for 50 mL total volume and 100 mM ionic strength.

Above data were obtained from the PHoEBuS 1.3 software by Analis (Namur, Belgium).