

### King Saud University Journal of King Saud University – Science

www.ksu.edu.sa www.sciencedirect.com

# ORIGINAL ARTICLE EPIYA (or -like) motifs in mammalian proteins



Journal of King Saud University -Science

### Fatemeh Safari \*

Laboratory of Membrane and Cytoskeleton Dynamics, Institute of Molecular and Cellular Biosciences, The University of Tokyo, 1-1-1 Yayoi, Bunkyo-ku, Tokyo 113-0032, Japan

Received 21 July 2013; accepted 1 May 2014 Available online 14 May 2014

### **KEYWORDS**

EPIYA (or -like) motif; Mammalian proteins; SFKs; Bacterial effector proteins

Abstract Recently, two mammalian proteins with functional Glu-Pro-Ile-Tyr-Ala (EPIYA) or a similar sequence (EPIYA-like motif) have been identified. Pragmin contains a single EPIYA motif and p140Cap contains two EPIYA-like motifs (EPLYA/EGLYA). Searching the human proteome using the NCBI BLAST software revealed that EPIYA-like motifs such as EPLYA, ESIYE, EDLYA, ESIYA, EHIYD, EPIYD, ENIYE, EPVYA, EVVYA, and TPLYA (which are present in bacterial effectors and have roles in virulence) are present in mammalian proteins. Moreover, several of bacterial EPIYA-like motifs were duplicated in mammalian proteins. Since, only five mammalian proteins contain EPIYA motif, and there are no mammalian proteins with more than one copy of the EPIYA motif. We also showed that the most of EPIYA (or -like) motifs in mammalian proteins are predicted to be naturally ordered. Moreover, few mammalian proteins containing EPIYA (or -like) motif were tyrosine phosphorylated. On the other hand, most of the mammalian proteins contain "unfunctional EPIYA (or -like) motifs". Notably, regulation of Src family kinases' (SFKs) activity (including up and down regulation) originally is mediated by their own EPIYA like motifs (EPIYI and EPQYQ) and thereby, it seems that EPIYA (or -like) motifs are "potential regulators of SFKs". Therefore, it can explain why EPIYA (or -like) motifs are restricted sites with low copy number in most mammalian proteins whereas in bacterial effector proteins, they can be universal sites with repetitive copies. In this study, we identified new mammalian proteins which contain functional EPIYA (or -like) motifs and thereby, their possible roles in cell signal transduction pathways. In this respect, functional experiments will be necessary to confirm our findings.

© 2014 King Saud University. Production and hosting by Elsevier B.V. All rights reserved.

\* Current address: Laboratory of Molecular Medicine Cell Biology, Graduate School of Biological Sciences, Nara Institute of Science and Technology, 8916-5, Takayama, Ikoma, Nara 630-0192, Japan. Tel.: +81 743725470.

E-mail addresses: safarifatemeh@bs.naist.jp, safarif2003@gmail. com.

Peer review under responsibility of King Saud University.

ELSEVIER Production and hosting by Elsevier

#### 1. Introduction

The presence of the Glu-Pro-Ile-Tyr-Ala (EPIYA) motif or a sequence closely related to the EPIYA motif (EPIYA-like motif) in effector proteins that play roles in the virulence of pathogenic bacteria was analyzed (Table 1). Bacterial effector proteins enter mammalian cells through type III or type IV secretion systems (T3SS or T4SS), and these effector proteins are phosphorylated at tyrosine residues located at the EPIYA (or -like) motif by host kinases, which triggers an interaction

1018-3647 © 2014 King Saud University. Production and hosting by Elsevier B.V. All rights reserved. http://dx.doi.org/10.1016/j.jksus.2014.05.001

Pathogen	Effector protein	EPIYA (or -like) motifs	Tyrosine kinases involved	References
H. pylori	CagA	EPIYA	SFKs, Abl	Poppe et al. (2007), Tammer et al. (2007), Mueller et al. (2012)
A. phagocytophilum C. trachomatis H. ducreyi B. henselae EPEC; C. rodentium	AnkA Tarp LspA Bep Tir	ESIYE, EDLYA, ESIYA, EPIYA ENIYE EPIYG, EPVYA EPLYA, EVVYA, TPLYA, EPLYA VNPYA, EHIYD; EPIYD	SFKs, Abl SFKs, Abl, Syk SFKs SFKs (?) Fyn, Abl	Ijdo et al. (2007) Mehlitz et al. (2008) Deng et al. (2008) Schulein et al. (2005) Campellone and Leong (2005), Phillips et al. (2004), Swimm et al. (2004),

Table 1 List of bacterial effector proteins which can be tyrosine phosphorylated by host kinases at EPIYA (or -like) motif.

with host cell SH2 domain-containing proteins and manipulates the function of host cells for more effective infection and improved colonization (Selbach et al., 2009; Backert et al., 2010; Hayashi et al., 2013).

Among mammalian proteins, it was demonstrated that Pragmin EPIYA motif undergoes tyrosine phosphorylation at the EPIYA motif by Src family kinases (SFKs) or in response to EGF stimulation. Tyrosine phosphorylation at the EPIYA allows Pragmin to interact with the SH2 domain of Csk, and thereby sequesters Csk in the cytoplasm. Sequestration of Csk by Pragmin has a positive feedback on SFK activity (Safari et al., 2011). In another study, Repetto et al. (2013) showed that N-terminal of p140Cap (also known as SRC kinase signaling inhibitor 1) contains two EPIYA-like motifs (EPLYA and EGLYA) which can be tyrosine phosphorylated at EPLYA and EGLYA by c-Abl or in response to integrin-mediated adhesion and EGF stimulation. Upon tyrosine phosphorylation, p140Cap EPLYA and EGLYA sequences serve as binding sites for the SH2 domain of Csk (Repetto et al., 2013; Sharma et al., 2013). Moreover, C-terminal of p140Cap contains a proline rich sequence which interacts with Src SH3 domain and thereby, facilitates inhibition of SFK activity by Csk (Di Stefano et al., 2007). Taken together, it has recently been proposed that the mammalian EPIYA (or -like) motif might have been exploited by pathogenic bacteria (Safari et al., 2011).

In this study, we identified all of bacterial EPIYA-like motifs in mammalian proteins and by using PhosphoSite, we explored the functional EPIYA (or -like) motifs in mammalian proteins. Furthermore, by using RONN tool, we did structure disorder prediction for mammalian proteins containing EPI-YA (or -like) motifs. Our results provided new information about differences between mammalian and bacterial EPIYA (or -like) motifs. Furthermore, our findings may give insights into the role of EPIYA (or -like) motif in new mammalian proteins to regulate SFK activity as the target of EPIYA (or -like) motifs.



**Figure 1** (A) List of human EPIYA containing proteins. By using PhosphoSite, the functional EPIYA motifs containing proteins were found which are shown in red boxes. The logo plot was also generated with the WebLogo tool. (B) Structure disorder prediction of mammalian proteins which contain functional EPIYA motif calculated with the RONN tool. (See above-mentioned references for further information.)

A Protein name	EPIYA segments	Tyrosine phosphorylation	Interacting partner	References
p140Cap/SRC kinase signaling inhibitor 1 (ACC# Q9C0H9)	MVLVKGEGLYADPYGLL	+	-Y396: Csk	Repetto et al, 2013
Oxidative stress-induced growth inhibitor 1	YQSTRQEGLYAMGPLAG	_		1 ,
Short transient receptor potential channel 7,3,6	DPGIISEGLYAIAVVLS	_		
Zinc finger SWIM domain-containing protein 4	QQRALPEGLYAQDKVVR	_		
Zinc finger SWIM domain-containing protein 5	QQRLMPEGLYAQDKVCR	-		
Partitioning defective 3 homolog B (isoform b and c) (ACC# Q8TEW	(8) PERDHLEGLYAKVNKPY	+	-Y1000:Unknown	Gu et al, 2011
Citron Rho-interacting kinase isoform 1,2	VLVGTEEGLYALNVLKN	-		
Filamin-B isoform s 1,2, 3,4	<b>KFIPREEGLYAVDVTYD</b>	_		
Beta/gamma crystallin domain-containing protein 3	NHCVLEEGLYADLTSCG	_		
Striated muscle preferentially expressed protein kinase isoform	1 HVGSEDEGLYAVSAVNT	-		
WD repeat and FYVE domain-containing protein 3	AMASDVEGLYAAVKALV	-		
Malonyl-CoA-acyl carrier protein transacylase, (mitochondrial	<b>GAMEFAEGLYA</b> VKIRAE	-		
1.0     3.0       2.0     1.0       0.0     0.0		15 WebLogo 3.3		
R				
Protein name	EPIYA segments	Tyrosine phosphorylation	Interacting partner	References
p140Cap/SRC kinase signaling inhibitor 1 (ACC# Q9C0H9)	KIYRKEPLYAAFPGSH	+	-Y264: Csk	Repetto et al, 2013
Putative fidgetin-like protein 2	AGNLPEPLYAGNACGG	-		
Girdin isoform X1 (or Girdin iso 5)	AIDMPEPLYAQARNSR	-		
		0 3.3		
p140Cap/ SRC kinase signaling inhibitor 1	Partitionin	ng defective 3 homolog B (i	soforms b and c)	
$Frobality_{r}$	Probability of Disorder	Amm	MAA	
Residue position	0.0 50 100 150 200 250 30	0 350 400 450 500 550 600 650 700 75	30 800 850 900 950 1000 1050 110	0
		Residue positi	on	

Tyrosine residues in the EPIYA – like motifs

**Figure 2** (A and B) List of human EPIYA-like motifs (EGLYA and EPLYA) containing proteins. Functional EGLYA and EPLYA motifs in mammalian proteins are shown in red boxes. Alignment of sequence surrounding EGLYA and EPLYA motifs is displayed and structure disorder prediction of mammalian proteins which contain functional EPIYA-like motifs (EPLYA or EGLYA) is shown. (See above-mentioned references for further information.)

### 2. Material and methods

NCBI (National Center for Biotechnology Information, U.S. National Library of Medicine, www.ncbi.nlm.nih.gov) was used to obtain sequences of mammalian EPIYA (or -like) motif containing proteins.

PhosphoSitePlus was used for the detection of tyrosine phosphorylation at EPIYA (or -like) containing proteins (http://www.phosphosite.org) (Hornbeck et al., 2012).

RONN tool was used for prediction of structure disorder of EPIYA (or -like) motif in mammalian proteins, (www.strubi.ox.ac.uk/RONN) (Yang et al., 2009). WebLogo 3 was used to align and display protein sequences (http://code.google.com/p/weblogo/) (Crooks et al., 2004).

### 3. Results

3.1. Tyrosine phosphorylation and structure disorder prediction of EPIYA segments in mammalian proteins

It was assumed that flexible feature of bacterial EPIYA (or like) segments allows them to interact with huge number of SH2-domain containing proteins in a tyrosine-phosphorylation-dependent manner (Hayashi et al., 2012, 2013). Moreover,

EPIYA-lik motifs	e Protein name	Tyrosine phosphorylation	Interacting partner	References
	Retinal dehydrogenase 2 (ACC# 094788)	+	-Y333: Unknown	PhosphoSite: CST Curation Set: 5083, 2008
ESIYE	Carcinoembryonic antigen-related cell adhesion molecule 20	-		
	Trafficking protein particle complex subunit 8	-		
	Myotubularin-related protein 6 (ACC# Q9Y217)	+	-Y108: Unknown	Iliuk et al, 2010
	Myotubularin-related protein 8	-		
EDLYA	Early endosome antigen 1	-		
	Pre-mRNA-processing-splicing factor 8	-		
	L Piezo-type mechanosensitive ion channel component 2	-		
ESIYA	Solute carrier family 25 member 38	-		
	Protein unc-119 homolog A	-		
EHIYD	Occludin isoforms (ACC# Q16625)	+	-Y287: Unknown	Rikova et al, 2007
Linib	Clathrin interactor 1	-		
	[	-		
EPIYD	Nuclear receptor ROR-beta	-		
	LCUB and sushi domain-containing protein 3	-		
	Torsin-1A precursor	-		
	Zinc finger protein 568	-		
	Zinc finger protein 674	-		
ENIYE	Zinc finger protein 699	-		
	Sodium/hydrogen exchanger 9 (ACC# Q8IVB4)	+	-Y622: Unknown	St-Germain et al, 2009
	MICAL-like protein 2 (ACC# Q8IY33)	+	-Y59: Unknown	Rikova et al, 2007
	Putative ATP-dependent RNA helicase DHX30	-		
	LTestis- and ovary-specific PAZ domain-containing protein 1	-		
	Dalmdalnhin (regulation)			Cure et al. 2008
EPVVA	Tringertite motif containing protein (TDIM161)	+	-Y262: Unknown	Guo et al, 2008
	The GTP as activating protain 27 isoforms (ACCU occurrence)			D.1 1 2007
	Kilo OTPase-activating protein 27 Isoformis (ACC#Q62UM4)	+	-Y228: Unknown	Rikova et al, 2007
	r			
EVVYA ·	l Probable ubiquitin carboxyl-terminal hydrolase FAF-X	-		
	[39S ribosomal protein L3, mitochondrial	-		
TDIVA	Zinc finger MIZ domain-containing protein 2	-		
ITLIA .	Bile salt export pump	-		
	Xin actin-binding repeat-containing protein 2	-		
	Zinc finger protein 469	-		

**Table 2** List of Human EPIYA-like motif containing proteins (EPIYA-like motifs which are present in bacterial effectors and have roles in virulence). Mammalian proteins with functional EPIYA-like motifs are shown in red boxes.

it was shown that five mammalian proteins contain EPIYA motif (Safari et al., 2011). Therefore, the question was whether EPIYA segments in mammalian proteins show order feature which is different from flexible feature of bacterial EPIYA segment and thereby, by using Phosphosite, we first explored tyrosine phosphorylation of EPIYA motif mammalian containing proteins. Our results showed that two out five mammalian proteins containing EPIYA motif were tyrosine phosphorylated (Fig. 1A). We then used RONN tool to predict disorder structure of EPIYA segments of mammalian proteins which were tyrosine phosphorylated. Our results indicated that Pragmin EPIYA segment revealed flexible feature which is similar to that of bacterial EPIYA segment (Fig. 1B). This result was not unexpected because of similar sequences distal between Pragmin EPIYA motif (EPIYAE-SAKR) and CagA-B EPIYA motif (EPIYAQVAKK).

# 3.2. Identification of functional EPIYA-like motifs in mammalian proteins

p140Cap contains two EPIYA-like motifs (EPLYA and EGLYA) (Repetto et al., 2013). This result indicates that EPI-YA-like motif in mammalian proteins has a function similar to that of EPIYA motif. Moreover, it raises a possibility that the EPIYA-like motifs (which are present in bacterial effectors and have roles in virulence) are present in mammalian proteins. To examine this possibility, first, we searched EPIYA-like motifs (EPLYA or EGLYA) in the human proteome by using NCBI BLAST software and our results showed that 11 mammalian proteins contain single EGLYA sequences and two mammalian proteins contain single EPLYA sequences (Fig 2A and B). However, only two out 14 mammalian proteins containing EGLYA or EPLYA were tyrosine phosphorylated. Disorder structure prediction for EPIYA-like segments of two mammalian proteins indicated that p140Cap has rigid feature and the other protein (Partitioning defective 3 homolog B) showed flexible feature. Interestingly, sequences distal Partitioning defective 3 homolog B EGLYA (EGLYAKVNKP) showed homology with CagA-A EPIYA motif (EPIYAKVNKK) (Fig. 2A). In this respect, it was shown that EPIYA like motifs such as ESIYE, EDLYA, ESIYA, EPIYG, VNPYA, EHIYD, EPIYD, ENIYE, EPVYA, EVVYA and TPLYA, which have roles in virulence, can be found in bacterial effector proteins (Selbach et al., 2009; Backert et al., 2010; Hayashi et al., 2013). We then investigated whether mammalian proteins contain such EPIYA-like motifs. Our results revealed that except EPIYG and VNPYA, all bacterial EPIYA-like motifs are present in mammalian proteins (Table 2), and bacterial EPIYA-like motif (ENIYE) is more common in mammalian proteins (Table 3). Among mammalian proteins containing

**Table 3** Frequencies of the 12 types of bacterial EPIYA-like motifs in mammalian proteins. <sup>\*</sup>EGLYA sequence is present in p140Cap (i.e. EGLYA sequence is not present in bacterial effector proteins which were discussed in this study).

EPIYA-like motif	Frequencies
EPLYA	3
ESIYE	3
EDLYA	5
ESIYA	1
EPIYG	0
VNPYA	0
EHIYD	3
EPIYD	2
ENIYE	8
EPVYA	3
EVVYA	1
TPLYA	5
*EGLYA	11
Total	45

bacterial EPIYA (or -like) motifs, only seven proteins were tyrosine phosphorylated. Disorder structure prediction of seven mammalian proteins containing bacterial EPIYA-like motifs showed that four proteins are naturally disordered (Fig. 3).

### 3.3. Identification of mammalian proteins with duplication of bacterial EPIYA like motifs

It was assumed that the presence of numerous EPIYA (or like) copies in bacterial effector proteins with disordered features allowed them to interact with several SH2-domaincontaining proteins (Hayashi et al., 2012, 2013). The presence of two EPIYA like motifs (EPLYA and EGLYA) in p140Cap, raises a possibility that the mammalian proteins may contain duplication of EPIYA like motifs and we wished to know the structure disorder prediction of them. So, we explored human proteome to find all proteins with duplication of bacterial EPIYA-like motifs and we found that several bacterial EPIYA-like motifs are duplicated in mammalian proteins (eight proteins) and only two of them were tyrosine phosphorylated (Fig. 4A). However, none of them are predicted to be disordered in EPIYA-like segments (Fig. 4B).



Figure 3 Structure disorder prediction of mammalian proteins which contain functional EPIYA-like motifs.

A Protein name	EPIYA segments	Tyrosine phosphorylation	Interacting partner	References
p140Cap/ SRC kinase signaling inhibitor 1 (ACC# Q9C0H9)	EPLYA / EGLYA	+/+	-Y264/ -Y396: Csk	Repetto et al, 2013
Myotubularin-related protein 8	EDLYA / EDNYA	-/-		
Piezo-type mechanosensitive ion channel component	EDLYA / EDIYA	-/-		
CUB and sushi domain-containing protein 3	EPIYD/ EPNYD	-/-		
Bile salt export pump	TPLYA / TPSYA	-/-		
Xin actin-binding repeat-containing protein 2	TPLYA / QPLYA	-/-		
Striated muscle preferentially expressed protein kinase isoforms	EGLYA / EGEYA	-/-		
WD repeat and FYVE domain-containing protein 3 isoform (ACC# Q8IZQ1)	FGLYA/ EGLYA	+/-	Y1327/-: Unknown	Huttlin et al, 2010



Figure 4 (A) List of human proteins which contain double EPIYA-like motif. (B) Structure disorder prediction of mammalian proteins with double EPIYA like motifs. (See above-mentioned references for further information.)

<b>A</b> Protein name	<b>EPIYA</b> segments	Tyrosine phosphorylation	Interacting partner	References
<b>c-Src</b> (ACC# P12931)	<b>AVSSEEPIYIVTEYM</b>	+	-Y338: Shp1	Barker et al, 1995; Somani et al, 1997
Fyn (ACC# P06241)	<b>AVVSEEPIYIVTEYM</b>	+	-Y339: Unknown	Moritz et al, 2010
c-Yes (ACC# P07947)	<b>AVVSEEPIYIVTEFM</b>	+	-Y345: Unknown	Oppermann et al, 2009
Lyn (ACC# P07948)	<b>VVTREEPIYI</b> ITEYM	+	-Y316: Unknown	Wolf-Yadlin et al, 2007
Blk (ACC# P51451)	<b>AVVTKEPIYIVTEYM</b>	+	-Y309: Unknown	PhosphoSite: CST Curation Set: 2062, 2006
Hck (ACC# P08631)	<b>AVVTKEPIYIITEFM</b>	+	-Y330: Unknown	Oppermann et al, 2009
Lck (ACC# P06239)	AVVTQ <b>EPIYI</b> ITEYM	-	-Y313	_
c- Fgr (ACC# P09769)	AVVSEEPIYIVTEFM	-	-Y331	-
4.0				



5

10 WebLogo 3.3

B Protein name	EPIYA segments	Tyrosine phosphorylation	Corresponding tyrosine	References
c-Src (ACC# P12931)	<b>YFTSTEPOYOPGENL</b>	+	-Y530	Bjelfman et al, 1990
Fyn (ACC# P06241)	YFTATEPOYOPGENL	+	-Y531	Peters et al, 1990
c-Yes (ACC# P07947)	YFTATEPOYOPGENL	+	-Y537	Park et al, 1993
Lyn (ACC# P07948)	FYTATEGOYOOOP	+	-Y508	Hirao et al, 1997
Blk (ACC# P51451)	FYTATERQYELQP	+	-Y501	PhosphoSite: CST Curation Set: 1789, 2006
Hck (ACC# P08631)	FYTATESQYQQP	+	-Y522	Ziegler et al, 1989
Lck (ACC# P06239)	FFTATEGQYQPQP	+	-Y505	Marth et al, 1988
c- Fgr (ACC# P09769)	YFTSAEPQYQPGDQT	+	-Y523	PhosphoSite: CST Curation Set: 492, 2004
4.0	FFTAT <mark>ERQYO</mark> BS			

Figure 5 (A and B) The presence of two EPIYA-like motifs (EPIYI and EPQYQ) in SFK members. (See above-mentioned references for further information.)

# 3.4. The presence of functional EPIYA-like motifs in SFK members and regulation of SFK activity via their own EPIYA-like motifs

In the current study, we focused on the functional EPIYA (or like) motifs (that have virulence roles in bacterial effector proteins) in mammalian proteins. However, it was shown that p140Cap contains EGLYA sequence that is not present in bacterial effector proteins. Therefore, it suggests that the non-bacterial EPIYA-like motifs are also present in mammalian proteins. In this respect, we found that members of Src family kinases contain duplication of functional EPIYA-like motifs. It was shown that N-terminal SH2 domain of SHP-1 interacts with phosphorylated tyrosine residue (or-Y338) which is located at EPIYA-like motif," EPIYI", of c-Src and thereby, up regulates Src activity in platelets and lymphocytes (Smart et al., 1981; Somani et al., 1997). Moreover, Y338 is conserved in all SFK members, although, the function of this tyrosine residue in the other SFK members is still unknown (Fig. 5A). Of note, the activity of SFKs was down regulated by Csk via EPIYA-like motif "EPQYQ" (Okada et al., 1988a,b; Okada et al., 1991). On the other hand, Csk is the responsible kinase which phosphorylates tyrosine residue (or-Y530 in c-Src) that is located at EPIYA-like motif "EPQYQ" (Fig. 5B). These results suggest that EPIYA (or -like) motifs can originally play regulatory roles for SFK members.

#### 4. Conclusion

The sequence EPxYAxV (where x is any amino acid) is significantly underrepresented in mammalian proteomes (Selbach et al., 2009). However, the number of mammalian proteins containing EPIYA-like motifs is not as low as it was proposed before (Xu et al., 2010; this paper). This difference may be due to explanation of "EPIYA-like motif". In the current study, we identified 50 mammalian proteins which contain bacterial EPIYA (or -like) motifs. Our results showed that most of mammalian proteins containing EPIYA (or -like) motifs were not tyrosine phosphorylated at EPIYA (or -like) motifs and thereby, they are "unfunctional EPIYA (or -like) motifs". Since, bacterial effector proteins containing same EPIYA (or -like) motifs were tyrosine phosphorylated. In this regard, we found that out of 11 functional EPIYA (or -like) motifs in mammalian proteins, six proteins were predicted to be naturally unfolded. Moreover, out of 39 proteins containing unfunctional EPIYA (or -like) motif, 32 proteins were predicted to be ordered at EPIYA (or -like) segments. Therefore, in mammalian proteins containing EPIYA (or -like) motif, tyrosine phosphorylation at unfolded EPIYA (or -like) segments is more frequent than ordered EPIYA (or -like) segments. Although, bacterial effector proteins were continuously tyrosine phosphorylated by host tyrosine kinases which seems to be due to unfolded nature of bacterial EPIYA (or -like) motifs. Our results also raise a possibility that EPIYA (or like) motifs can be potential regulators of SFK members in the cell machinery. On the other hand, SFK activity can be directly regulated via their own EPIYA-like motifs (EPIYI and EPQYQ) or via indirect feedback by functional EPIYA (or -like) motifs containing proteins (including mammalian proteins and pathogenic effector proteins). During long-term co-evolution with hosts, pathogenic bacteria enable to employ

host key regulators such as EPIYA (or -like) to better colonization. These key regulators were expanded via recombination (Furuta et al., 2011; Majazki et al., 2013). However, in the mammalian hosts, such key regulators are mostly unfunctional with low copy number. Furthermore, it seems that intramolecular interactions in mammalian proteins containing EPIYA (or -like) allowed the EPIYA (or -like) segments to be as unaccessed sites. Notably, SFKs have pleiotropic functions in the regulation of fundamental cellular processes. SFKs are highly conserved across the animal kingdom, and regulation of SFKs is a prerequisite for evolution of multicellular animals (Segawa et al., 2006; Okada, 2012). Our findings showed that via tyrosine phosphorylation at EPIYA (or -like) motifs, mammalian proteins can interact with Csk (Safari et al., 2011; Repetto et al., 2013) and SHP-1 (Somani et al., 1997). In this study, several mammalian proteins with functional EPIYA (or -like) motifs were identified for which interacting partner(s) are still unknown. So, it would be interesting to gain a better understanding of SFK activity via EPIYA (or -like) motifs. Moreover, elucidation of the structure of mammalian proteins containing EPIYA (or -like) motif proteins will provide more information on the specificity of mammalian EPIYA (or -like) segments. It is known that SFKs are activated in some human cancers, including colon and breast cancer (Summy and Gallick, 2003). In this respect, the role of Pragmin EPIYA motif in activation of SFKs should be explored. Notably, SFKs are candidates for the treatment of human cancers (Kopetz et al., 2007; Hollande et al., 2010). Therefore, the specificity of EPIYA (or -like) segments to regulation of SFK activity should be further explored.

### Acknowledgment

We would like to thank the Reviewers for their valuable comments and suggestions.

#### References

- Backert, S., Tegtmeyer, N., Selbach, M., 2010. The versatility of helicobacter pylori CagA effector protein functions: the master key hypothesis. Helicobacter 15, 163–176.
- Barker, S.C., Kassel, D.B., Weigl, D., Huang, X., Luther, M.A., Knight, W.B., 1995. Characterization of pp60c-src tyrosine kinase activities using a continuous assay: autoactivation of the enzyme is an intermolecular autophosphorylation process. Biochemistry 34, 14843–14851.
- Bjelfman, C., Meyerson, G., Cartwright, C.A., Mellström, K., Hammerling, U., Påhlman, S., 1990. Early activation of endogenous pp60src kinase activity during neuronal differentiation of cultured human neuroblastoma cells. Mol. Cell. Biol. 10, 361–370.
- Campellone, K.G., Leong, J.M., 2005. Nck-independent actin assembly is mediated by two phosphorylated tyrosines within enteropathogenic *Escherichia coli* Tir. Mol. Microbiol. 56, 416–432.
- Crooks, G.E., Hon, G., Chandonia, J.M., Brenner, S.E., 2004. WebLogo: a sequence logo generator. Genome Res 14, 1188–1190.
- Deng, W., Vallance, B.A., Li, Y., Puente, J.L., Finlay, B.B., 2003. Citrobacter rodentium translocated intimin receptor (Tir) is an essential virulence factor needed for actin condensation, intestinal colonization and colonic hyperplasia in mice. Mol. Microbiol. 48, 95–115.
- Deng, K., Mock, J.R., Greenberg, S., van Oers, N.S., Hansen, E.J., 2008. *Haemophilus ducreyi* LspA proteins are tyrosine phosphorylated by macrophage-encoded protein tyrosine kinases. Infect Immun. 76, 4692–4702.

- Di Stefano, P., Damiano, L., Cabodi, S., Aramu, S., Tordella, L., Praduroux, A., Piva, R., Cavallo, F., Forni, G., Silengo, L., Tarone, G., Turco, E., Defilippi, P., 2007. P140Cap protein suppresses tumour cell properties, regulating Csk and Src kinase activity. EMBO J. 26, 2843–2855.
- Furuta, Y., Yahara, K., Hatakeyama, M., Kobayashi, K.S., 2011. Evolution of cagA oncogene of *Helicobacter pylori* through recombination. PLoS One 6, e23499.
- Gu, T.L., Deng, X., Huang, F., Tucker, M., Crosby, K., Rimkunas, V., Wang, Y., Deng, G., Zhu, L., Tan, Z., Hu, Y., Wu, C., Nardone, J., MacNeill, J., Ren, J., Reeves, C., Innocenti, G., Norris, B., Yuan, J., Yu, J., Haack, H., Shen, B., Peng, C., Li, H., Zhou, X., Liu, X., Rush, J., Comb, M.J., 2011. Survey of tyrosine kinase signaling reveals ROS kinase fusions in human cholangiocarcinoma. PLoS One 6, e15640.
- Guo, A., Villén, J., Kornhauser, J., Lee, K.A., Stokes, M.P., Rikova, K., Possemato, A., Nardone, J., Innocenti, G., Wetzel, R., Wang, Y., MacNeill, J., Mitchell, J., Gygi, S.P., Rush, J., Polakiewicz, R.D., Comb, M.J., 2008. Signaling networks assembled by oncogenic EGFR and c-Met. Proc. Natl. Acad. Sci. U.S.A. 105, 692– 697.
- Hayashi, T., Senda, M., Morohashi, H., Higashi, H., Horio, M., Kashiba, Y., Nagase, L., Sasaya, D., Shimizu, T., Venugopalan, N., Kumeta, H., Noda, N., Inagaki, F., Senda, T., Hatakeyama, M., 2012. Tertiary structure-function analysis reveals the pathogenic signaling potentiation mechanism of Helicobacter pylori oncogenic effector CagA. Cell Host Microbe 12, 20–33.
- Hayashi, T., Morohashi, H., Hatakeyama, M., 2013. Bacterial EPIYA effectors – where do they come from? What are they? Where are they? Where are they going? Cell. Microbiol. 15, 377–385.
- Hirao, A., Hamaguchi, I., Suda, T., Yamaguchi, N., 1997. Translocation of the Csk homologous kinase (Chk/Hyl) controls activity of CD36-anchored Lyn tyrosine kinase in thrombin-stimulated platelets. EMBO J. 16, 2342–2351.
- Hollande, F., Pannequin, J., Joubert, D., 2010. The long road to colorectal cancer therapy: searching for the right signals. Drug Resist. Updat. 13, 44–56.
- Hornbeck, P.V., Kornhauser, J.M., Tkachev, S., Zhang, B., Skrzypek, E., Murray, B., Latham, V., Sullivan, M., 2012. PhosphoSitePlus: a comprehensive resource for investigating the structure and function of experimentally determined post-translational modifications in man and mouse. Nucleic Acids Res. 40, D261–D270.
- Huttlin, E.L., Jedrychowski, M.P., Elias, J.E., Goswami, T., Rad, R., Beausoleil, S.A., Villén, J., Haas, W., Sowa, M.E., Gygi, S.P., 2010. A tissue-specific atlas of mouse protein phosphorylation and expression. Cell 143, 1174–1189.
- Ijdo, J.W., Carlson, A.C., Kennedy, E.L., 2007. Anaplasma phagocytophilum AnkA is tyrosine-phosphorylated at EPIYA motifs and recruits SHP-1 during early infection. Cell. Microbiol. 9, 1284– 1296.
- Iliuk, A.B., Martin, V.A., Alicie, B.M., Geahlen, R.L., Tao, W.A., 2010. In-depth analyses of kinase-dependent tyrosine phosphoproteomes based on metal ion-functionalized soluble nanopolymers. Mol. Cell. Proteomics 9, 2162–2172.
- Kopetz, S., Shah, A.N., Gallick, G.E., 2007. Src continues aging: current and future clinical directions. Clin. Cancer Res. 13, 7232– 7236.
- Luo, W., Slebos, R.J., Hill, S., Li, M., Brábek, J., Amanchy, R., Chaerkady, R., Pandey, A., Ham, A.J., Hanks, S.K., 2008. Global impact of oncogenic Src on a phosphotyrosine proteome. J. Proteome Res. 7, 3447–3460.
- Majazki, J., Wüppenhorst, N., Hartelt, K., Birtles, R., von Loewenich, F.D., 2013. Anaplasma phagocytophilum strains from voles and shrews exhibit specific ankA gene sequences. BMC Vet. Res. 9, 235. http://dx.doi.org/10.1186/1746-6148-9-235.
- Marth, J.D., Cooper, J.A., King, C.S., Ziegler, S.F., Tinker, D.A., Overell, R.W., Krebs, E.G., Perlmutter, R.M., 1988. Neoplastic

transformation induced by an activated lymphocyte-specific protein tyrosine kinase (pp56lck). Mol. Cell. Biol. 8, 540–550.

- Mehlitz, A., Bänhart, S., Hess, S., Selbach, M., Meyer, T.F., 2008. Complex kinase requirements for Chlamydia trachomatis Tarp phosphorylation. FEMS Microbiol. Lett. 289, 233–240.
- Moritz, A., Li, Y., Guo, A., Villén, J., Wang, Y., MacNeill, J., Kornhauser, J., Sprott, K., Zhou, J., Possemato, A., Ren, J.M., Hornbeck, P., Cantley, L.C., Gygi, S.P., Rush, J., Comb, M.J., 2010. Akt-RSK-S6 kinase signaling networks activated by oncogenic receptor tyrosine kinases. Sci. Signal. 3, ra64. http://dx.doi.org/ 10.1126/scisignal.2000998.
- Mueller, D., Tegtmeyer, N., Brandt, S., Yamaoka, Y., De Poire, E., Sgouras, D., Wessler, S., Torres, J., Smolka, A., Backert, S., 2012. C-Src and c-Abl kinases control hierarchic phosphorylation and function of the CagA effector protein in Western and East Asian *Helicobacter pylori* strains. Clin. Invest. 122, 1553–1566.
- Okada, M., Nakagawa, H., 1988a. Protein tyrosine kinase in rat brain: neonatal rat brain expresses two types of pp60c-src and a novel protein tyrosine kinase. J. Biochem. 104, 297–305.
- Okada, M., Nakagawa, H., 1988b. Identification of a novel protein tyrosine kinase that phosphorylates pp60c-src and regulates its activity in neonatal rat brain. Biochem. Biophys. Res. Commun. 154, 796–802.
- Okada, M., Nada, S., Yamanashi, Y., Yamamoto, T., Nakagawa, H., 1991. CSK: a protein-tyrosine kinase involved in regulation of src family kinases. J. Biol. Chem. 266, 24249–24252.
- Okada, M., 2012. Regulation of the Src family kinases by Csk. Int. J. Biol. Sci. 8, 1385–1397.
- Oppermann, F.S., Gnad, F., Olsen, J.V., Hornberger, R., Greff, Z., Kéri, G., Mann, M., Daub, H., 2009. Large-scale proteomics analysis of the human kinome. Mol. Cell. Proteomics 8, 1751–1764.
- Park, J., Meisler, A.I., Cartwright, C.A., 1993. C-Yes tyrosine kinase activity in human colon carcinoma. Oncogene 8, 2627–2635.
- Peters, D.J., McGrew, B.R., Perron, D.C., Liptak, L.M., Laudano, A.P., 1990. In vivo phosphorylation and membrane association of the fyn proto-oncogene product in IM-9 human lymphoblasts. Oncogene 5, 1313–1319.
- Phillips, N., Hayward, R.D., Koronakis, V., 2004. Phosphorylation of the enteropathogenic *E. coli* receptor by the Src-family kinase c-Fyn triggers actin pedestal formation. Nat. Cell Biol. 6, 618–625.
- Poppe, M., Feller, S.M., Römer, G., Wessler, S., 2007. Phosphorylation of *Helicobacter pylori* CagA by c-Abl leads to cell motility. Oncogene 26, 3462–3472.
- Repetto, D., Aramu, S., Boeri Erba, E., Sharma, N., Grasso, S., Russo, I., N., Jensen, O., Cabodi, S., Turco, E., Di Stefano, P., Defilippi, P., 2013. Mapping of p140Cap phosphorylation sites: the EPLYA and EGLYA motifs have a key role in tyrosine phosphorylation and Csk binding, and are substrates of the Abl kinase. PLoS One 8, 1–12.
- Rikova, K., Guo, A., Zeng, Q., Possemato, A., Yu, J., Haack, H., Nardone, J., Lee, K., Reeves, C., Li, Y., Hu, Y., Tan, Z., Stokes, M., Sullivan, L., Mitchell, J., Wetzel, R., Macneill, J., Ren, J.M., Yuan, J., Bakalarski, C.E., Villen, J., Kornhauser, J.M., Smith, B., Li, D., Zhou, X., Gygi, S.P., Gu, T.L., Polakiewicz, R.D., Rush, J., Comb, M.J., 2007. Global survey of phosphotyrosine signaling identifies oncogenic kinases in lung cancer. Cell 131, 1190–1203.
- Safari, F., Murata-Kamiya, N., Saito, Y., Hatakeyamaa, M., 2011. Mammalian Pragmin regulates Src family kinases via the Glu-Pro-Ile-Tyr-Ala (EPIYA) motif that is exploited by bacterial effectors. Proc. Natl. Acad. Sci. U.S.A. 108, 14938–14943.
- Schulein, R., Guye, P., Rhomberg, T.A., Schmid, M.C., Schröder, G., Vergunst, A.C., Carena, I., Dehio, C., 2005. A bipartite signal mediates the transfer of type IV secretion substrates of Bartonella henselae into human cells. Proc. Natl. Acad. Sci. U.S.A. 102, 856– 861.
- Segawa, Y., Suga, H., Iwabe, N., Oneyama, C., Akagi, T., Miyata, T., Okada, M., 2006. Functional development of Src tyrosine kinases

during evolution from a unicellular ancestor to multicellular animals. Proc. Natl. Acad. Sci. U.S.A. 103, 12021–12026.

- Selbach, M., Ernst Paul, F., Brandt, S., Guye, P., Daumke, O., Backert, S., Dehio, C., Mann, M., 2009. Host cell interactome of tyrosine-phosphorylated bacterial proteins. Cell Host Microbes 5, 397–403.
- Sharma, N., Repetto, D., Aramu, S., Grasso, S., Russo, I., Fiorentino, A., Mello-Grand, M., Cabodi, S., Singh, V., Chiorino, G., Turco, E., Di Stefano, P., Defilippi, P., 2013. Identification of two regions in the p140Cap adaptor protein that retain the ability to suppress tumor cell properties. Am. J. Cancer Res. 3, 290–301.
- Smart, J.E., Oppermann, H., Czernilofsky, A.P., Purchio, A.F., Erikson, R.L., Bishop, J.M., 1981. Characterization of sites for tyrosine phosphorylation in the transforming protein of Rous sarcoma virus (pp6Ovsrc) and its normal cellular homologue (pp6oc rc). Proc. Natl. Acad. Sci. U.S.A. 78, 6013–6017.
- Somani, A.K., Bignon, J.S., Mills, G.B., Siminovitch, K.A., Branch, D.R., 1997. Src kinase activity is regulated by the SHP-1 proteintyrosine phosphatase. J. Biol. Chem. 272, 21113–21119.
- St-Germain, J.R., Taylor, P., Tong, J., Jin, L.L., Nikolic, A., Stewart, I.I., Ewing, R.M., Dharsee, M., Li, Z., Trudel, S., Moran, M.F., 2009. Multiple myeloma phosphotyrosine proteomic profile associated with FGFR3 expression, ligand activation, and drug inhibition. Proc. Natl. Acad. Sci. U.S.A. 106, 20127–20132.

- Summy, J.M., Gallick, G.E., 2003. Src family kinases in tumor progression and metastasis. Cancer Metastasis Rev. 22, 337–358.
- Swimm, A., Bommarius, B., Li, Y., Cheng, D., Reeves, P., Sherman, M., Veach, D., Bornmann, W., Kalman, D., 2004. Enteropathogenic *E. coli* use redundant tyrosine kinases to form actin pedestals. Mol. Biol. Cell. 15, 3520–3529.
- Tammer, I., Brandt, S., Hartig, R., König, W., Backert, S., 2007. Activation of Abl by *Helicobacter pylori*: a novel kinase for CagA and crucial mediator of host cell scattering. Gastroenterology 132, 1309–1319.
- Wolf-Yadlin, A., Hautaniemi, S., Lauffenburger, D.A., White, F.M., 2007. Multiple reaction monitoring for robust quantitative proteomic analysis of cellular signaling networks. Proc. Natl. Acad. Sci. U.S.A. 104, 5860–5865.
- Xu, S., zhang, C., Gao, J., Xu, D., 2010. Effector prediction in hostpathogen interaction based on a Markov model of a ubiquitous EPIYA motif. BMC Genomics 11 (Suppl. 3), S1.
- Yang, Z.R., Thomson, R., McNeil, P., Esnouf, R.M., 2009. RONN: the bio-basis function neural network technique applied to the detection of natively disordered regions in proteins. Bioinformatics 21, 3369–3376.
- Ziegler, S.F., Levin, S.D., Perlmutter, R.M., 1989. Transformation of NIH 3T3 fibroblasts by an activated form of p59hck. Mol. Cell Biol. 9, 2724–2727.